



In your pursuit to advance science, every experiment matters. There's no time to start again, no need to wonder if the products you chose will set you back or propel you forward. With our comprehensive portfolio of thermal cyclers, PCR plastics, reagents, and service plans, you know exactly what you're getting—quality made to support your every win, big and small, at every step.



Contents

Which instrument his your needs?	4
ProFlex PCR System	6
VeritiPro Thermal Cycler	8
SimpliAmp Thermal Cycler	10
MiniAmp Thermal Cycler	12
Automated Thermal Cycler	14
Thermal Cycler Fleet Control Software	16
Thermal cycler services and support	16
MicroAmp PCR plastics	18
SuperScript IV Reverse Transcriptase	22
PCR enzymes	24
PCR primers and dNTPs	26
Nucleic acid separation and analysis	27
Specialty enzymes	31
Molecular biology education	31

Which instrument fits your needs?

Providing precise, consistent results, Applied Biosystems[™] thermal cyclers are designed for every challenge, application, and budget. Be confident in your results with instruments designed for reliability, accuracy, and intuitive programming. Reduce time spent on PCR optimizations with Applied Biosystems[™] VeriFlex[™] Block temperature control technology. Manage your instruments remotely with secure access.





	Ultimate flexibility and throughput	Ultimate performance
	ProFlex PCR System	VeritiPro Thermal Cycler
Max. sample throughput	480,000 reactions	384 reactions
Max. block ramp rate	6.0°C/sec	6.0°C/sec
Block formats (temperature optimization)	 3 x 32-well, 0.2 mL (2-zone VeriFlex Block) 96-well, 0.2 mL (6-zone VeriFlex Block) 2 x 96-well, 0.2 mL 2 x flat block 2 x 384-well, 0.02 mL 	 96-well, 0.2 mL (6-zone VeriFlex Block) 384-well, 0.02 mL

For Research Use Only. Not for use in diagnostic procedures.













Elegantly simple and precise	Routine PCR	Designed for easy robotic integration
SimpliAmp Thermal Cycler	MiniAmp Thermal Cycler	Automated Thermal Cycler
96 reactions	96 reactions	384 reactions
4.0°C/sec	3.0°C/sec	3.5°C/sec
96-well, 0.2 mL (3-zone VeriFlex Block)	• 96-well, 0.2 mL	 96-well, 0.2 mL, compatible with full- or semi-skirted plates 384-well, 0.02 mL

Looking for a thermal cycler for *in vitro* diagnostic use? Learn more about the Applied Biosystems™ Veriti™ Dx Thermal Cycler at **thermofisher.com/veritidx**



Don't forget reagents—choose from industry-leading enzymes such as Invitrogen™ SuperScript™ IV reverse transcriptases and Invitrogen™ Platinum™ SuperFi™ II reagents (see pages 22–25).



Interested in private-label thermal cyclers or PCR plastics? To find out more, go to thermofisher.com/oem

Ultimate flexibility and throughput

ProFlex PCR System

The Applied Biosystems[™] ProFlex[™] PCR System combines flexible configuration and control features to fit how you work today and tomorrow with the reliability you've come to expect from Applied Biosystems[™] products. Interchangeable block formats allow you to maximize your throughput or run independent experiments concurrently.

The ProFlex PCR System is cloud-enabled, giving you the freedom to design and securely upload your methods, monitor runs, and check instrument availability from any mobile device or desktop computer with the Thermo Fisher[™] Connect Platform.

- Multi-user accessible—run three experiments at once
- Flexible block configuration—accepts five different block formats for optimization and throughput
- Cloud-enabled—conveniently access your instrument anytime, anywhere with the Connect Platform
- Fleet control-compatible-manage multiple instruments, users, and methods with Applied Biosystems™ Thermal Cycler Fleet Control Software



Five interchangeable block options

The ProFlex PCR System has five different blocks that can be changed with the flip of a switch, including a 3 x 32-well block. This allows up to three experiments to be run simultaneously and completely independently of each other.

Dual 96-well and dual 384-well blocks are available for high-throughput needs. A dual flat block is also available to support Applied Biosystems[™] OpenArray[™] plate technology for genotyping analysis on the Applied Biosystems™ QuantStudio™ 12K Flex Real-Time PCR System as well as our sealed-chip technology on the QuantStudio™ 3D Digital PCR System.



3 x 32-well



96-well



Dual 96-well



Dual 384-well



Dual flat

Specifications

Block format	3 x 32-well, 0.2 mL, 2-zone VeriFlex Block independent control	96-well, 0.2 mL, 6-zone VeriFlex Block	2 x 96-well, 0.2 mL	2 x flat block for chips and arrays	2 x 384-well, 0.02 mL
Features	Run three experiments at once or at different times	Perform complete optimization work with full 96-well VeriFlex Block	High throughput in 96-well format	High throughput capability: 8 x 3,072 OpenArray Plate* or 24 x 20K Chip**	High throughput in 384-well format
Max. block ramp rate	6.0°0	C/sec	3.0°C/sec	1.6°C/sec	3.0°C/sec
Max. sample ramp rate	4.4°(C/sec	1.6°C/sec	NA	1.6°C/sec
Temperature accuracy			±0.25°C (35–99.9°C)		
Temperature range			0-100.0°C		
Temperature uniformity		<0.5	5°C (30 sec after reaching 95°C)		
Dimensions (H x W x D)		27.2 x	27.2 x 33.0 x 56.5 cm (10.6 x 13 x 22 in.)		
Weight	18.8 kg	g (41 lb)		20.4 kg (45 lb)	
PCR volume range	10-8	10-80 μL		33 nL	5-20 µL
Instrument memory			USB, onboard		
Display interface			8.4-inch color TFT LCD		
Power		100-	-240 V, 50-60 Hz, max. 95	50 VA	
VeriFlex Block	2 temperature zones per block (5°C zone-to-zone)	6 temperature zones, 25°C range (5°C zone-to-zone)	NA		
Data connectivity [†]		Clou	id or mobile via Ethernet or Wi-Fi		

^{*} OpenArray Plate is compatible with the QuantStudio 12K Flex Real-Time PCR System.

Ordering information

Product	Complete system Cat. No.	Block only Cat. No.	Instrument + 5-year* Rapid Exchange warranty Cat. No.
ProFlex 96-Well PCR System	4484075	4483637	A27934
ProFlex 3 x 32-Well PCR System	4484073	4483638	A28986
ProFlex Dual 96-Well PCR System	4484076	4484071	A27937
ProFlex Dual Flat PCR System	4484078	4484074	A27931
ProFlex Dual 384-Well PCR System	4484077	4484072	A30229

Recommended plastics					
3 x 32-well block	Cat. No.	96-well block	Cat. No.	384-well block	Cat. No.
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32811	MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354	MicroAmp Optical 384-Well Reaction Plate	4343370
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589	MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32811	MicroAmp Optical 384-Well Reaction Plate with Barcode	4309849
MicroAmp Reaction Tube with Cap, 0.2 mL	N8010540	MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589	MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode	4483285
MicroAmp 32-Well Clear Adhesive Film	A32812	MicroAmp Clear Adhesive Film	4306311	MicroAmp Clear Adhesive Film	4306311
		MicroAmp 32-Well Clear Adhesive Film	A32812		

^{*} The system comes with a 2-year standard warranty; 3 years of an extended warranty may be added for a total of 5 years of service.



Did you know?

The ProFlex PCR System, Applied Biosystems™ VeritiPro™ Thermal Cycler, and the Applied Biosystems™ SimpliAmp™ Thermal Cycler feature VeriFlex temperature control technology, which enables more precise and efficient PCR optimization.



Find out more at thermofisher.com/veriflextechnology

^{**} The Applied Biosystems™ QuantStudio™ 3D Digital PCR 20K Chip is compatible with the QuantStudio 3D Digital PCR System.

[†] The Instrument Connect app, available at Apple™ and Google™ app stores, can be used to monitor your instrument.

The Connect Platform can be used to create and share protocols, schedule an instrument, and start and monitor runs remotely.

Ultimate performance

VeritiPro Thermal Cycler

The VeritiPro Thermal Cycler delivers proven reliability with advanced temperature control technology and connectivity. Take advantage of next-level PCR optimization with the precision offered by VeriFlex temperature control technology. Connect to the cloud-enabled VeritiPro Thermal Cycler remotely. Conveniently design and securely upload your methods, monitor runs, and check instrument availability from any mobile device or desktop computer with the Connect Platform.

- Innovative design—ramp rate of 6.0°C/sec, quiet fan, and ergonomic soft-close lid
- VeriFlex Block—precise control over six independent temperature zones for PCR optimization
- Fleet control-compatible—manage multiple instruments, users, and methods with Thermal Cycler Fleet Control Software






Precise control with VeriFlex temperature control technology

VeriFlex Blocks are constructed of segmented metal blocks with separate heating and cooling elements below each, enabling:

- More precise control over PCR optimization—each block can be set with up to six specific temperatures
- Precise incubation—use the six temperature zones of the VeriFlex Block to do enzyme studies, restriction digests, or any other process that requires precise temperature control







Tired of water baths?

Incubate samples at up to six different temperatures simultaneously for enzyme activation studies, restriction digests, or sequencing library preps with the VeritiPro Thermal Cycler.

Specifications

Block format	96-well, 0.2 mL, 6-zone VeriFlex Block	384-well, 0.02 mL, isothermal		
Max. block ramp rate	6.0°C/sec	5.0°C/sec		
Max. sample ramp rate	4.4°C/sec	3.5°C/sec		
Temperature accuracy	±0.25°C (35–99.9°C)		
Temperature range	0–10	0.0°C		
Temperature uniformity	<0.5°C (30 sec at	ter reaching 95°C)		
Dimensions (H x W x D)	21.7 x 24.5 x 46.5 cm (8.5 x 9.6 x 18.3 in.)			
Weight	12.0 kg (26.5 lb)			
PCR volume range	10–100 μL	5–20 μL		
Instrument memory	USB port and 16 GB onboard memory; onboard capacity >1,000 protocols			
Display interface	8-inch color TFT LCD			
Power	100-240 V, 50-60 Hz, max. 700 W			
VeriFlex Block range	30°C range across block, 6 temperature zones (up to 10°C zone-to-zone)	NA		
Data connectivity*	Cloud or mobile via Ethernet or Wi-Fi			

^{*} The Instrument Connect app, available at Apple and Google app stores, can be used to monitor your instrument.

Ordering information

Product	Cat. No.	Instrument + 5-year* Rapid Exchange warranty Cat. No.
VeritiPro Thermal Cycler, 384 well	A48140	A50326
VeritiPro Thermal Cycler, 96 well	A48141	A48765
Recommended plastics		
384-well block		Cat. No.
MicroAmp Optical 384-Well Reaction Plate		4343370
MicroAmp Optical 384-Well Reaction Plate with Barcode		4309849
MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode		4483285
MicroAmp Clear Adhesive Film		4306311
96-well block		Cat. No.
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode		4483354
MicroAmp Optical 96-Well Reaction Plate		N8010560
MicroAmp Clear Adhesive Film	4306311	
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate		A32811
MicroAmp 32-Well Clear Adhesive Film		A32812
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589	

^{*} The system comes with a 2-year standard warranty; 3 years of an extended warranty may be added for a total of 5 years of service.

The Connect Platform can be used to create and share protocols, schedule an instrument, and start and monitor runs remotely.

Elegantly simple and precise

SimpliAmp Thermal Cycler

The SimpliAmp Thermal Cycler is an easy-to-use, compact, and accurate thermal cycler designed to fit every lab's essential PCR workflow. Features like a responsive color touchscreen and VeriFlex temperature control technology enable simple, accurate optimization. Plus, the SimpliAmp Thermal Cycler is cloud-enabled, giving you the freedom to design and securely upload your methods, monitor runs, and check instrument availability from any mobile device or desktop computer with the Connect Platform.

- Intuitive interface—large, easy-to-use color touchscreen for easy programming and quick status checks
- VeriFlex Block—three independent temperature zones for PCR optimization
- Cloud-enabled—conveniently access your instrument anytime, anywhere with the Connect Platform
- Compact design—helps save bench space
- Fleet control-compatible—manage multiple instruments, users, and methods with Thermal Cycler Fleet Control Software





The status dial displays the current block temperature and elapsed run time.



Simple, intuitive user interface

The SimpliAmp Thermal Cycler has an 8-inch color touchscreen, making navigation of the intuitive menu options fast and efficient.



Specifications

Block format	96-well, 0.2 mL, 3-zone VeriFlex Block		
	Enabled to run Fast chemistry	Auto restart (after power outages)	
Features	Controllable ramp rate	Edit program during experiment	
	Program overwrite protection	One-touch incubation	
Max. block ramp rate*	4°C	C/sec	
Max. sample ramp rate*	3°0	C/sec	
Temperature accuracy	±0.25°C (35-99.9°C)		
Temperature range	0–100.0°C		
Temperature uniformity	<0.5°C (30 sec after reaching 95°C)		
Temperature calibration	Calibrated to standards traceable to the National Institute of Standards and Technology (NIST)		
Dimensions (H x W x D)	21.0 x 24.0 x 46.0 cm (8.3 x 9.5 x 18.1 in.)		
Weight	8.3 kg (18.3 lb)		
PCR volume range	10–100 μL		
nstrument memory	2,000 MB onboard memory (capacity for >1,000 protocols); USB port for additional external storage		
Display interface	8-inch color TFT LCD		
Power	100-240 V, 50-60 Hz, max. 600 W		
VeriFlex Block	3 temperature zones, 20°C range (10°C zone-to-zone)		
Data connectivity**	Cloud or mobile via Ethernet or Wi-Fi		

^{*} At a reaction volume of 1 µL.

The Connect Platform can be used to create and share protocols, schedule an instrument, and start and monitor runs remotely.

Ordering information

Product	Cat. No.	Instrument + 5-year* warranty with Rapid Exchange plan Cat. No.
SimpliAmp Thermal Cycler	A24811	A27603
High-Power USB Wi-Fi Module	A26774	NA
Recommended plastics		
96-well block	Cat. No.	
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354	
MicroAmp Optical 96-Well Reaction Plate	N8010560	
MicroAmp Clear Adhesive Film	4306311	
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32811	
MicroAmp 32-Well Clear Adhesive Film	A32812	
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589	



Did you know?

ProFlex, VeritiPro, SimpliAmp, and Applied Biosystems[™] MiniAmp[™] Thermal Cyclers have thermal simulation modes that make the transition from other thermal cyclers simple, accurate, and efficient. A library of modes that mimic the ramp rates of other instruments is available on each instrument.

^{**} The Instrument Connect app, available at Apple and Google app stores, can be used to monitor your instrument.

Routine PCR

MiniAmp Thermal Cycler

The MiniAmp Thermal Cycler delivers the reliability you've come to expect from Applied Biosystems™ technology for routine PCR at every lab bench.

- Compact design—fits in everyone's lab space at just 7.5 in. (19 cm) wide
- Cloud-enabled—conveniently access your instrument anytime, anywhere with the Connect Platform
- Fleet control-compatible—manage multiple instruments, users, and methods with Thermal Cycler Fleet Control Software

Secure remote access with cloud-enabled instruments

The ProFlex, VeritiPro, SimpliAmp, and MiniAmp Thermal Cyclers are cloud-enabled instruments, allowing you secure, private access with an account on the Connect Platform. Anywhere, anytime with any mobile device or desktop computer you can:

- Design and share protocols
- · Schedule an instrument
- Start or stop a run
- Check run status







Did you know?

Most Applied Biosystems thermal cyclers come with a two-year standard warranty and a starter kit that includes Applied Biosystems™ PCR plate and tube samples and all the tools you need. Packages that include the instrument, an extended warranty, and consumables are available. Visit the instrument webpage to view available packages.

Specifications

	MiniAmp Thermal Cycler	MiniAmp Plus Thermal Cycler	
Block format	96-well, 0.2 mL isothermal block	96-well, 0.2 mL, 3-zone VeriFlex Block	
Max. block ramp rate*	3.0°C/sec	3.5°C/sec	
Max. sample ramp rate*	2.2°C/sec	2.7°C/sec	
Temperature accuracy	±0.25°C (35–99.9°C)	
Temperature range	0-10	00.0°C	
Temperature uniformity	<0.5°C (30 sec after reaching 95°C)		
Temperature calibration	Calibrated to standards traceable to the NIST		
Dimensions (H x W x D)	20 x 19 x 39 cm (7.9 x 7.5 x 15.4 in.)		
Weight	5.9 kg (13.0 lb)		
PCR volume range	10–1	00 μL	
Instrument memory	2,000 MB onboard memory (capacity for >1,000 protocols); USB port for additional external storage		
Display interface	5-inch color TFT LCD		
Power	100-240 V, 50-60 Hz, max. 500 W		
VeriFlex Block	NA	3 temperature zones, 20°C range (10°C zone-to-zone)	
Data connectivity**	Cloud or mobile via Ethernet or Wi-Fi		

 $^{^{\}star}$ At a reaction volume of 1 $\mu L.$

Ordering information

		Instrument + 5-year warranty with
Product	Cat. No.	Rapid Exchange plan Cat. No.
MiniAmp Thermal Cycler	A37834	A38081
MiniAmp Plus Thermal Cycler	A37835	A38077
High-Power USB Wi-Fi Module	A26774	NA
Recommended plastics		
96-well block	Cat. No.	
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354	
MicroAmp Optical 96-Well Reaction Plate	N8010560	
MicroAmp Clear Adhesive Film	4306311	
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32811	
MicroAmp 32-Well Clear Adhesive Film	A32812	
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589	



 $^{^{\}star\star} \text{ The Instrument Connect app, available at Apple and Google app stores, can be used to monitor your instrument.}$

The Connect Platform can be used to create and share protocols, schedule an instrument, and start and monitor runs remotely.

Designed for easy robotic integration

Automated Thermal Cycler

The Applied Biosystems™ Automated Thermal Cycler offers the flexibility, reliability, and performance needed in a complete PCR automation system. The small, easy-to-integrate format of the Automated Thermal Cycler (ATC) enables hands-free PCR results.

- Flexible modular design and small footprint—helps save space on deck
- Automated lid—easy, hands-free operation with a liquid handler or plate stacker
- Free desktop software—for PCR optimization prior to robotic integration
- Plug-and-play drivers and Standardization in Lab Automation (SiLA) compatibility—for easy integration on your liquid handler of choice
- Fleet control-compatible—manage multiple instruments, users, and methods with Thermal Cycler Fleet Control Software





Designed for any stage of your workflow automation journey



Stand-alone

Optimize assays before robotic integration with our direct software.



Plug-and-play drivers

Ask about available drivers for leading robotic platforms.



SiLA-compatible

Maximum robotic platform integration flexibility with SiLA rapid integration coding.



Find out more at thermofisher.com/atc

Specifications

Block format	96-well, 0.2 mL*	384-well, 0.02 mL				
PCR volume range	10–100 μL for full-skirted plates; 20–100 μL for semi-skirted plates	5–20 μL for full-skirted plates				
	Predrilled mounting and/or alignment points at each co	orner of the chassis				
Hardware integration features	3-side and top-plate access					
	Available in 3-connector configurations					
	Application programming interfaces (APIs) available for	integration with robotics systems				
Software features	SiLA Rapid Integration software-standardized program	nming access**				
	Free software available for instrument demonstration and stand-alone operation					
Block module dimensions (H x W x D)	13.3 x 17.9 x 31.7 cm (5.2 x 7.0 x 12.5 in.)					
Control module dimensions (H x W x D)	7.0 x 25.7 x 33.1 cm	(2.8 x 10.1 x 13.0 in.)				
Temperature accuracy	±0.25°C (3	5.0-99.9°C)				
Max. block ramp rate [†]	3.5°C/sec	2.8°C/sec				
Max. sample ramp rate [†]	1.8°C/sec	1.6°C/sec				
Temperature range	4–105°C (no condensation risk	with subambient temperatures)				
Temperature uniformity	≤0.50°C (20 sec a	Ifter reaching 95°C)				
Temperature calibration	Calibrated to standard	s traceable to the NIST				
Service options	2-year standard warranty includes Rapid Exchange ser	vice plan				
Cervice options	NIST-traceable temperature probe equipment available					
Weight	9.4 kg (20.7 lb) total (block module 6.0 k	kg (13.2 lb), control module 3.4 kg (7.5 lb))				
Power	100–240 V, 50–6	60 Hz, max. 600 W				
Flexible ramp rates	Program your own ramp rates, or use preprogrammed simulation modes					
Data connectivity	L	AN				

^{*} Compatible with full- or semi-skirted plates. ATC semi-skirted adaptor required for use with semi-skirted 96-well plates.

ATC ships with the full-skirted adaptor installed, which is required for use with full-skirted 96-well plates.

Ordering information

Product	96-well Cat. No.	384-well Cat. No.	3-year extended warranty* with Rapid Exchange plan Cat. No.
Automated Thermal Cycler System, laptop, 1 m cable	A31486	A33977	
Automated Thermal Cycler System, laptop, 3 m cable	A31487	A33978	
Automated Thermal Cycler System, laptop, 10 cm cable	A31488	A33979	ZGEXSCATC3Y/ZGEXSCATC3843Y
Automated Thermal Cycler System, 1 m cable	A31489	A33980	ZGEXSCATO31/ZGEXSCATO36431
Automated Thermal Cycler System, 3 m cable	A31490	A33981	
Automated Thermal Cycler System, 10 cm cable	A31491	A33982	
Automated Thermal Cycler Semi-Skirted Adaptor	A33044	NA	NA
Automated Thermal Cycler Full-Skirted Adaptor	A33045	NA	NA

^{*} The Automated Thermal Cycler comes standard with a 27-month warranty.

Recommended plastics							
96-well block	Cat. No.	384-well block	Cat. No.				
MicroAmp EnduraPlate Optical 96-Well Full-Skirted Plate with Barcode, Clear*	A31728	MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode	4483273				
MicroAmp EnduraPlate Optical 96-Well Multicolor Reaction Plate with Barcode, semi-skirted**	4483356	MicroAmp EnduraPlate Optical 384-Well Multicolor Reaction Plate with Barcode	4483317				
MicroAmp Clear Adhesive Film	4306311	MicroAmp Clear Adhesive Film	4306311				

^{*} Available in multiple colors.

^{**} Requires semi-skirted plate adaptor (Cat. No. A33044); included with instrument purchase.



The Applied Biosystems[™] MicroAmp[™] EnduraPlate[™] plate (Cat. No. A31728) was specifically designed for the Automated Thermal Cycler and has three-way barcoding. Learn about the amplification uniformity benefits at thermofisher.com/automatedpcrplate

^{**} sila-standard.org

 $[\]dagger$ At a reaction volume of 1 $\mu\text{L}.$

Manage multiple users, methods, and instruments with a single intuitive interface

Thermal Cycler Fleet Control Software


The Thermal Cycler Fleet Control Software is as intuitive and easy to use as our thermal cyclers are. With this tool, you'll be able to view and control all Applied Biosystems™ thermal cyclers in your facility.

Secure—comes preinstalled on a server that sits securely behind your firewall

Flexible—create custom permissions and rules for users, methods, and instruments

Powerful—quickly view your thermal cycler inventory and availability right from the dashboard

Find out more at thermofisher.com/fleetcontrol



Convenient—control your thermal cyclers, reports, methods, and firmware upgrades from your office, another lab computer, or any other networked device

Science can't wait on instrument downtime

Thermal cycler services and support

You can't afford downtime in your lab. Built on more than 35 years of service expertise, our superior service solutions for Applied Biosystems™ instruments and applications help keep your lab up and running. More than 1,400 trained professionals make up the industry's largest network, ready to assist you when you need it.

Rapid Exchange

When fast turnaround time is critical and asset tag consistency is not required, choose our Rapid Exchange plan. If our Remote Service Center engineers can't repair your instrument through telephone support, we will ship you a factory-certified refurbished replacement instrument within one business day—and it's yours to keep. AB Repair Center (ABRC) Support Plus

1	•	Receive a replacement unit (your existing service coverage now applies to this instrument)
2		Package the broken instrument in the same box
_		. 45.449
3		Ship it back to us free of charge

Off-site service plans at a glance

	Rapid Exchange	AB Repair Center (ABRC) Support Plus
Plan designed for	Labs that operate under critical timelines	Regulated labs that must maintain asset tag consistency to comply with regulated protocols
Repair time	Replacement shipped in 1 business day*	Loaner instrument shipped in 1 business day**
Off-site repair service including shipping, parts, and labor	✓	✓
Telephone and email access to our Remote Service Center	✓	✓

^{*} Subject to regional availability.

^{**} Your instrument repaired and returned in ~3 weeks. Service plan without loaner instrument also available.

Our AB Repair Center Support Plus plan is the best choice for customers in regulated environments that must maintain asset tag consistency to comply with regulated protocols and need their instrument returned after repair. With this mail-in option, we provide a loaner instrument for you to use while factory-trained engineers maintain and repair your instrument at our Remote Service Center. All instruments undergo certification performance tests, temperature verification, instrument and optical calibration, ground continuity, and firmware upgrades, when applicable.

Just mail us your instrument and we'll repair and ship it back to you within three weeks. All parts, labor, and shipping charges are included. Planned maintenance, temperature verification, and loaner instrument service* can be added to any off-site service plan.

Check your instrument's coverage status at thermofisher.com/manage

Digital service innovations

Our instrument service plans include digital innovations that help keep your instruments and your lab running smoothly. With pioneering on-demand tools and capabilities such as remote support using mobile augmented-reality (AR) technology, instrument-driven support, and on-demand instrument training, we're constantly looking ahead so your lab never falls behind.

Explore our digital service innovations at thermofisher.com/innovations

Temperature verification services

Multiplex dynamic temperature verification (MDTV) services are designed to deliver fast, accurate temperature readings and documentation to assess the accuracy and uniformity of your PCR systems. MDTV is available for Applied Biosystems thermal cyclers.**

Have confidence in the accuracy of your results, learn more about the MDTV service at **thermofisher.com/mdtv**

Education services

It can be difficult to prepare yourself for what's next while you're focused on the work you have now. Our interactive professional training courses make it easier. We offer a combination of virtual and in-person classroom instruction, and hands-on learning in your lab to match your schedule, budget, and learning preferences. Whichever course style you choose, you'll learn from one of our 300 highly skilled application scientists who are available to lead sessions online, at your location, or at one of our 12 training centers located worldwide.

Explore courses at thermofisher.com/educationservices

Technical support


If you have questions about product selection or use, assay or experimental design, data analysis, or troubleshooting, contact our team of technical support scientists or access our online product and application support tools.

How to reach us

To find your local support or technical support team, go to

thermofisher.com/contactus

For product FAQs, protocols, training courses, and webinars, go to thermofisher.com/technicalresources



An example arrangement of five Applied Biosystems instruments receiving the MDTV service simultaneously.



Find out more at thermofisher.com/thermalcyclerservice

 $^{^{\}star}$ Loaner instruments are not available in all regions. Contact your services and support representative for availability.

^{**} MDTV service not currently offered on the 60-well Applied Biosystems™ GeneAmp™ PCR System 9700.

High-performance PCR plastics for optimal PCR results

MicroAmp PCR plastics

Applied Biosystems™ PCR plastics have been designed and validated to work with our thermal cyclers for more than 25 years. That's why they are engineer approved to enable optimal PCR performance.

Applied Biosystems[™] MicroAmp[™] PCR plastics are:


- Validated on Applied Biosystems thermal cyclers for optimal fit and performance
- Designed for optimal heat transfer with thin-walled polypropylene wells
- Designed to reduce cross-contamination with raised well rims for effective sealing

Applied Biosystems[™] MicroAmp[™] TriFlex 3 x 32-Well PCR Reaction Plate and adhesive film

- One plate that can be used as separated, individual 32-well plates or as a full 96-well plate
- Easy-to-tear dual side tabs
- · Alphanumeric labeling on each plate segment
- Precut Applied Biosystems[™] MicroAmp[™] 32-Well Clear Adhesive Film



Unique, high-performance features of Applied Biosystems™ MicroAmp™ EnduraPlate™ plastic consumables





Options for every format and throughput need

Choose from tubes, tube strips, plates, adhesive film, and accessories for any throughput need. MicroAmp EnduraPlate plastics offer a solution for experiments that require special handling, such as automated or high-throughput workflows, and an even greater degree of durability for use with multi-instrument experiments.



MicroAmp 8-tube strip with attached domed caps; optical flat caps are also available for use with real-time PCR (qPCR).



The Applied Biosystems™ MicroAmp™ 8-Tube Strip with attached optical or domed caps offers a combination of features to help prevent cross-contamination, pipetting errors, and sample identification errors in your PCR and real-time PCR applications.

- Attached caps that open and close independently of each other
- Etched A-H letter labeling for individual tubes and caps
- Dual side tabs for strip labeling
- Graduated 20 µL measuring markers on every tube
- Available exclusively for Applied Biosystems thermal cyclers to enable optimal PCR results



Find out more at **thermofisher.com/pcrplastics**



Did you know?

Proper plate sealing helps reduce evaporation and well-to-well contamination.

- 1. Remove the backing of the Applied Biosystems™ adhesive film.
- 2. Align the adhesive film so as to cover all wells while placing on the plate.
- 3. Rub the flat edge of the applicator along the long edge (length) of the plate, then along the short edge (width). Finally, rub the applicator between all the wells and around the outside edges of the plate using small back-and-forth motions to form a complete seal.



Which PCR plastic fits your needs?

Find the PCR plastic format with the throughput and features for your application

Use for:	Small-scale experiments with a few samples	Routine experiments	Automation	Laboratory use
	Single tubes, strips, caps, adhesive film, and accessories	MicroAmp optical microplates	MicroAmp EnduraPlate optical microplates	MicroAmp EnduraPlate optical microplates GPLE*
Formats	Single tubes	• 32-well	• 96-well	• 96-well
	Single tubes with caps	48-well Fast	96-well Fast	96-well Fast
	8-strip tubes with caps	• 96-well	• 384-well	• 384-well
	• 12-strip caps	96-well Fast		
	Adhesive film/plate seal	• 384-well		
DNA-, RNase-, PCR inhibitor-free	Yes	Yes	Yes	Yes
ANSI/SBS standard dimension color	Clear, or mixed packs containing red, orange, blue, and green	Clear	Single-color packs (red, blue, green, yellow, or clear) and 5-plate sampler (1 of each color)	Clear
Instrument compatibility	Use our plastics selection tool	Use our plastics selection tool	Use our plastics selection tool	Use our plastics selection tool
Barcode available	No	Yes (1 or 2 sides)	Yes (3 sides)	Yes (3 sides)
Multiple applications	No	No	Yes	Yes
Optical compatibility	Yes (applicable for optical version)	Yes	Yes	Yes
Use	Research use only	Research use only	Research use only	For laboratory use**

^{*}General purpose laboratory equipment (GPLE); for laboratory use.




Did you know?

Need high-quality PCR plastics for non–Applied Biosystems instruments?

Visit <u>thermofisher.com/thermoscientificplastics</u> for a wide range of Thermo Scientific™ PCR plastics.

Custom and OEM plastics for PCR and qPCR are available.

Learn more at thermofisher.com/oemplastics



^{**} Lot-based contamination test with Certificate of Analysis.

Quickly find the plastics and accessories you need for your instrument

		3 x				96-well Fast	384-we	II		Genetic
Applied Biosystems MicroAmp Product	Cat. No.	32-well ProFlex	ProFlex, SimpliAmp, VeritiPro, Veriti, MiniAmp Plus, MiniAmp	2720	9700	Veriti	ProFlex, VeritiPro	0026	310	3130, 3130xl, 3500, 3500xL, 3730, 3730xl
96-well, 0.2 mL reaction plates	10010500 1010010									
Optical 96-Well plate	N8010560, 4316813		•	•	•					•
Optical 96-Well plate with barcode	4306737, 4326659		•	•	•					•
Optical 96-Well plate with barcode and optical caps	403012		•	•	•					
Optical 96-Well plate with barcode and optical adhesive films	4314320		•	•	•					
EnduraPlate Optical 96-Well Clear plate	A36924		•	•	•					•
EnduraPlate Optical 96-Well Clear plate with barcode*	4483354, 4483352		•	•	•					•
TriFlex 3 x 32-Well PCR plate	A32810, A32811	•**	•	•	•					
96-well, 0.1 mL reaction plates										
Fast Optical 96-Well plate	4346907								•	•
Fast Optical 96-Well plate with barcode	4346906, 4366932					•			•	•
EnduraPlate Optical 96-Well Fast Clear plate	A36930					•			•	•
EnduraPlate Optical 96-Well Fast Clear plate with barcode*	4483485, 4483494					•			•	•
384-well reaction plates										
Optical 384-Well plate	4343370				П		•	•		•
Optical 384-Well plate with barcode	4309849, 4326270, 4343814						•	•		•
EnduraPlate Optical 384-Well Clear plate	A36931						•	•		•
EnduraPlate Optical 384-Well Clear plate with barcode*	4483285, 4483273						•	•		•
Strip tubes and caps										
Fast 8-Tube Strip, 0.1 mL	4358293					•				
Optical 8-Tube Strip with Attached Optical Caps, 0.2 mL	A30588	•	•	•	•					
8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589	•	•	•	•					
8-Tube Strip, 0.2 mL*	N8010580	•	•	•	•				•	
Optical 8-Tube Strip, 0.2 mL	4316567	•	•	•	•					
8-Cap Strip*, clear	N8010535, N8011535	•	•	•	•	•				
Optical 8-Cap Strip	4323032	•	•	•	•	•				
12-Cap Strip*	N8010534, N8011534		•	•	•	•				
Single tubes										
Fast Reaction Tube with Cap, 0.1 mL	4358297, 4358293 N8010540,					•				
Reaction Tube with Cap, 0.2 mL*	N8010612, N8011540	•	•	•	•					
Reaction Tube without Cap, 0.2 mL*	N8010533, N8011533	•	•	•	•					
Optical Tube without Cap, 0.2 mL	N8010933	•	•	•	•					
Seals and covers	1000011									
Clear Adhesive Film	4306311		•		•	•	•	•		
Optical Adhesive Film	4360954, 4311971		•		•	•	•	•		
96-Well Full Plate Cover	N8010550			•	•				_	
32-Well Clear Adhesive Film	A32812	•**	•	•	•					
Accessories	404000									
Splash-Free 96-Well Base	4312063		•	•	•	•				
96-Well Support Base	4379590		•	•	•	•				•
96-Well Base * Multiple colors are available.	N8010531		•	•	•					

^{*} Multiple colors are available.

Note: Experiments using one or two 8-tube strips with attached caps require blank tube strips to balance lid pressure on the block or the use of the Applied Biosystems™ MicroAmp™ 96-Well Tray/Retainer Set (Cat. No. 4381850)—bottom part of tray only. For use with 96-well block of ProFlex, SimpliAmp, and VeritiPro thermal cyclers.

^{**} Do not use Applied Biosystems™ MicroAmp™ 3 x 32-Well Retainer (Cat. No. 4481669).

Superior cDNA synthesis for any application

SuperScript IV Reverse Transcriptase

With more than 50,000 citations, reviews, and publications, Invitrogen™ SuperScript™ reverse transcriptases are among the most trusted and widely used products for cDNA synthesis. Invitrogen™ SuperScript™ IV Reverse Transcriptase is the latest enzyme in the portfolio, engineered to deliver superior cDNA synthesis performance with even the most challenging RNA samples.

- Super-efficient—up to 100x higher cDNA yields than with other reverse transcriptase enzymes
- Super-sensitive—C_t values reduced by as many as eight cycles for RT-qPCR
- Super-robust—high thermostability and processivity for superior cDNA synthesis, even from degraded or inhibitor-containing RNA samples
- Super-fast—10 min cDNA synthesis protocol



Select the right format for your application

We offer a comprehensive portfolio of enzymes and kits within the SuperScript IV family to suit your research needs. Start with the selection guide below to find the best format for common cDNA synthesis applications.

Would you like to have	the ability to optimize reaction components and conditions?	a complete kit with all cDNA synthesis reaction components?	ultimate convenience and minimal pipetting steps for RT-PCR?	the ability to skip the RNA extraction process and synthesize cDNA directly from mammalian cell lysate?
	•			
Product format	Stand-alone enzyme	First-strand cDNA synthesis kit	One-Step RT-PCR kit	Direct RT kit
Recommended product	SuperScript IV Reverse Transcriptase	SuperScript IV First-Strand Synthesis System		SuperScript IV CellsDirect cDNA Synthesis Kit
Applications	RT-PCR, RT-qPCR, sequence detection, gene expression analysis, transcript variant detection, cloning, cDNA library construction, RACE, RNA-Seq	RT-PCR, RT-qPCR, sequence detection, gene expression analysis, transcript variant detection, cloning, cDNA library construction, RACE, RNA-Seq	RT-PCR, sequence detection, cloning, genotyping, high-throughput analysis	RT-PCR, RT-qPCR
Input total RNA	1 pg-5 μg	1 pg-5 μg	0.01 pg–1 μg	1–10,000 cells
Optimal reaction temperature	50-55°C	50-55°C	50-55°C	50-55°C
Reaction time	10 min	10 min	10 min	10 min
cDNA synthesis with challenging or degraded RNA	Yes	Yes	Yes	Yes

Order Invitrogen™ Oligo(dT) or random primers for cDNA synthesis at thermofisher.com/rtprimers



Find out more at thermofisher.com/ssiv

Complete kit with flexible priming options

SuperScript IV First-Strand Synthesis System

The Invitrogen™ SuperScript™ IV First-Strand Synthesis System is optimized for synthesis of first-strand cDNA from purified poly(A)⁺ or total RNA. The kit contains all components needed for reverse transcription, plus an additional control gene and primers, providing the flexibility to customize the reaction setup to fit the needs of your application.

Find out more at thermofisher.com/ssiv-firststrand



Conducting a single-cell low-RNA input study?

Invitrogen[™] SuperScript[™] IV Single Cell/Low Input cDNA PreAmp Kit combines Invitrogen[™] SuperScript[™] IV and Invitrogen[™] Platinum[™] SuperFi[™] enzymes along with template switching technology, which allows generation of high-quality full-length cDNA with high yields, sensitivity, and accuracy for RNA sequencing or qPCR. Find out more at thermofisher.com/ssivpreamp

Enabling faster, more efficient RT-PCR

SuperScript IV One-Step RT-PCR System

Even with challenging RNA samples, you can get more efficient results faster and easier than with any other RT-PCR reagent. The Invitrogen™ SuperScript™ IV One-Step RT-PCR System combines high-processivity SuperScript IV Reverse Transcriptase and high-fidelity Invitrogen™ Platinum™ SuperFi™ DNA Polymerase to provide superior one-step RT-PCR performance.

Find out more at

thermofisher.com/ssiv-onestep

- Two-phase hot-start activation mechanism—for high specificity, improved yields, and easy room temperature setup
- Superior sensitivity and speed—down to 0.01 pg of RNA, target length up to 13.8 kb, and the fastest one-step RT-PCR protocol
- Reliable target detection—even in RNA samples with suboptimal purity
- Fast and easy gDNA removal—for superior accuracy and confidence in your results

Doing RT-qPCR?

Invitrogen™ SuperScript™ IV VILO™ Master Mix offers exceptional cDNA synthesis for RT-qPCR applications, while maintaining superior linearity across the broadest range of input RNA. Learn more when you flip the brochure to the real-time PCR side, or visit thermofisher.com/4vilo

Ordering information

Product	Quantity	Cat. No.
	2,000 units	18090010
SuperScript IV Reverse Transcriptase	10,000 units	18090050
aperScript IV Reverse Transcriptase aperScript IV First-Strand Synthesis System aperScript IV First-Strand Synthesis System with ezDNase Enzyme aperScript IV One-Step RT-PCR System aperScript IV One-Step RT-PCR System with ezDNase Enzyme aperScript IV One-Step RT-PCR System with ezDNase Enzyme aperScript IV VILO Master Mix aperScript IV CellsDirect cDNA Synthesis Kit	4 x 10,000 units	18090200
Cupar Carint IV First Strand Synthasia System	50 reactions	18091050
uperScript IV First-Strand Synthesis System uperScript IV First-Strand Synthesis System with ezDNase Enzyme	200 reactions	18091200
Cupar Carint IV First Strand Synthasia System with az DNasa Enzyma	50 reactions	18091150
Superscript in First-Strand Synthesis System with ezdinase Enzyme	200 reactions	18091300
Cupar Carint IV One Stee DT DCD System	25 reactions	12594025
Superscript in One-Step ni-Fon System	100 reactions	12594100
SuperSeriet IV One Step DT DCD System with a DNage Enzyme	50 reactions	12595025
Superscript in One-Step ni-Pon System with ezonase Enzyme	100 reactions	12595100
SuperScript IV VILO Master Mix	50 reactions	11756050
SuperScript IV CellsDirect cDNA Synthesis Kit	50 reactions	11750150
SuperScript IV Single Cell/Low Input cDNA PreAmp Kit	48 reactions	11752048

DNA polymerases for PCR

A comprehensive portfolio of PCR enzymes and master mixes is available with the high performance and consistency you need. Start with the selection guide below to find the enzyme best suited for your PCR applications.





PCR type	High-fidelity PCR	Hot-start PCR	Hot-start PCR	Direct PCR
Do you need	accurate DNA sequences?	cleaner bands or to detect low-abundance targets?	a chemical hot start?	direct amplification without genomic DNA purification?
	Platinum SuperFi II	Platinum II <i>Tag</i> Hot-Start	AmpliTag Gold 360	Platinum Direct PCR
Recommended DNA polymerase	DNA Polymerase	DNA Polymerase	DNA Polymerase	Universal Master Mix
Applications				
Cloning and subcloning	•			
Site-directed mutagenesis	•			
GC-rich templates	•	•	•	•
Sanger sequencing templates	•	•	•	•
High-throughput PCR	•	•		•
Long-range PCR (20-40 kb)	•			
Genotyping	•	•	•	•
Amplification of samples with suboptimal purity	•	•		•
Colony PCR	•	•	•	•
Multiplex PCR	•	•	•	•
Fast PCR	•	•		•
Technical specifications				
Fidelity compared to <i>Taq</i> DNA polymerase	>300x	1x	1x	1x
Target length	Up to 40 kb*	Up to 5 kb	Up to 5 kb	Up to 8 kb
Hot-start modification	Antibody-mediated	Antibody-mediated	Chemical modification	Antibody-mediated
Speed	15-30 sec/kb	15 sec/kb	1 min/kb	20 sec/kb
Universal primer annealing	Yes	Yes	NA	Yes
Inhibitor tolerance	Yes	Yes	NA	Yes
Blunt or 3'-A end	Blunt	3´-A	3´-A	3′-A
Compatibility with TaqMan probes	NA	Yes	Yes	NA
Formats				
Master mix	Colorless/green**	Colorless/green**	Colorless	Green**
Stand-alone enzyme	Colorless	Colorless/green [†]	Colorless	NA

^{*} Depends on complexity of DNA templates; may require additional optimization of reaction conditions and primer design.



Did you know?

Assembled PCR reactions that include Invitrogen™ Platinum™ SuperFi™ II or Invitrogen™ Platinum™ II *Taq* Hot-Start DNA Polymerase are stable for 24 hours at room temperature, enabling high-throughput applications with automation or long waits between runs.

^{**} Direct gel loading with green buffer options.

[†] Green buffer available as separate item for use with stand-alone enzyme.

Ultimate accuracy and robustness

Platinum SuperFi II DNA Polymerase

Platinum SuperFi II DNA Polymerase is a hot-start, engineered, proofreading DNA polymerase for PCR applications requiring the highest sequence accuracy. Its fidelity is >300x that of *Taq* DNA polymerase and its buffer is specially formulated for a universal annealing temperature of 60°C.

- Exceptional fidelity—>300x more accurate than Taq DNA polymerase
- **Simplified workflow**—no need for primer melting temperature (T_m) calculation; enables cocycling of multiple PCR targets
- Increased PCR success—robust amplification of GC-rich targets, DNA of suboptimal purity, and long sequences
- Invitrogen™ Platinum™ hot-start technology—enables superior specificity, sensitivity, and yields; allows for room temperature reaction setup and benchtop stability

Find out more at thermofisher.com/platinumsuperfi



Enhanced PCR specificity and yield with fast cycling

Platinum II Taq Hot-Start DNA Polymerase

Platinum II *Taq* Hot-Start DNA Polymerase is an enzyme engineered for simplicity, specificity, and speed. Its universal primer annealing feature simplifies optimization and allows cocycling of different PCR assays.

- Universal primer annealing at 60°C—reduces tedious optimization steps and enables cocycling of different PCR assays
- Engineered Taq DNA polymerase—allows 4x faster synthesis and successful amplification even in the presence of inhibitors
- Platinum hot-start technology—enables superior specificity, sensitivity, and yields; allows for room temperature reaction setup and benchtop stability

Find out more at thermofisher.com/platinumiitaq

Ordering information

Product	Quantity*	Cat. No.
Platinum SuperFi II DNA Polymerase	100 reactions	12361010
Platinum SuperFi II PCR Master Mix	100 reactions	12368010
Platinum SuperFi II Green PCR Master Mix	100 reactions	12369010
Platinum II Taq Hot-Start DNA Polymerase	100 reactions	14966001
Platinum II Hot-Start PCR Master Mix (2X)	50 reactions	14000012
Platinum II Hot-Start Green PCR Master Mix (2X)	50 reactions	14001012

^{*} Additional product sizes available.

Everything you need for reliable PCR

PCR primers and dNTPs

Primers

We offer high-quality Invitrogen™ custom DNA oligos for successful PCR. Our oligos synthesis services feature over 25 years of customer service and popular primer design tools, as well as an informative, easy-to-use ordering portal to meet your needs. Invitrogen custom DNA oligos offer:

- Purity that you require: desalting, cartridge, HPLC, or PAGE
- Scales that suit your research needs:
 25 nmol, 50 nmol, 200 nmol, 1 µmol, and
 10 umol
- Product format you prefer: tube or plate
- Affordable price
- · Enhanced bulk ordering

Find out more at thermofisher.com/primers



- Gain confidence—Primer3-based design algorithm
- Speed up—design primers for PCR or capillary electrophoresis (CE), up to 50 genes at the same time
- Work smarter-recognizes .txt and .fasta file types
- Order with ease—seamlessly integrates with the Invitrogen™ ordering portal
- Store your data—ability to save your projects

Try OligoPerfect Primer Designer at

thermofisher.com/oligoperfect-designer

Complex and large-scale oligos

Easily customize your oligos based on the unique scale, modification, and purification requirements for your research.

- Milligram to >15-gram delivered yields
- Unmodified and modified DNA or RNA up to 100 bases
- ISO 9001– and ISO 13485–controlled process
- Custom formulations and packaging solutions
- OEM collaboration and commercial supply support available
- Dedicated technical project team

Our experienced technical team can help you with your scale-up from step one to launch.

Learn more at thermofisher.com/largescaleoligos

dNTPs

Our dNTPs have been extensively tested and verified for use in a wide variety of molecular biology applications, including highly sensitive techniques such as RT-qPCR and next-generation sequencing.

Learn more at thermofisher.com/dntp





Did you know?

We offer additional manufacturing services, including OEM scale-up and support?

Leverage our expertise and infrastructure. We can manufacture products with your labels and packaging.

Explore OEM and commercialization options for oligo synthesis at thermofisher.com/oem.

Simplify nucleic acid electrophoresis

Nucleic acid separation and analysis

E-Gel precast agarose gels

Using precast agarose gels can simplify the nucleic acid electrophoresis workflow. Invitrogen™ E-Gel™ precast gels are ready-to-use, prestained agarose gels packaged in a disposable cassette. There are no gels to pour, buffers to make, staining or destaining steps to perform, or gel boxes to assemble. Just load your samples and start the run.

- Ready to use—precast with agarose and DNA stain packaged inside a disposable cassette
- Time-saving-ideal for customers looking to streamline workflows
- Versatile—designed for rapid analysis of PCR products, plasmid preparations, and restriction digests to improve genotyping workflows



Sample analysis in three simple steps—load, run, and analyze.



Approximate time to complete an E-Gel electrophoresis system workflow and a traditional DNA electrophoresis workflow

Ordering information

Product	Quantity	Cat. No.
E-Gel Agarose Gels with SYBR Safe DNA Gel Stain, 1%	10 gels	A42100
E-Gel Agarose Gels with SYBR Safe DNA Gel Stain, 2%	10 gels	A42135
E-Gel EX Agarose Gels, 2%	10 gels	G401002
E-Gel 48 Agarose Gels with SYBR Safe DNA Gel Stain, 2%	8 gels	G820802
E-Gel 96 Agarose Gels with SYBR Safe DNA Gel Stain, 2%	8 gels	G720802
E-Gel CloneWell II Agarose Gels with SYBR Safe DNA Gel Stain, 0.8%	10 gels	G661818



Find out more at thermofisher.com/egel



Did you know?

E-Gel precast gels are available in a variety of formats for routine and high-throughput applications, with different stains and agarose percentages (0.8%, 1.2%, 2%, and 4%). To find the right gel for your needs, see the selection guide at thermofisher.com/egelselection

Electrophoresis elevated to empower

E-Gel Power Snap Plus Electrophoresis System

The Invitrogen™ E-Gel™ Power Snap Plus Electrophoresis System combines the convenience of quick real-time nucleic acid analysis with high-quality visualization in both high- and low-throughput Invitrogen™ E-Gel™ products. It offers:

- Fast and safe workflow—uses Invitrogen™ E-Gel™ technology, Invitrogen™ SYBR™ Safe DNA Gel Stain, and a blue-light transilluminator
- Convenience—gel-running and documentation on the same small device
- Easy analysis—available on both the device and on the Thermo Fisher™
 Connect Platform with Invitrogen™ iBright™ image analysis software



Specifications

Invitrogen™ E-Gel™ Power Snap Plus Electrophoresis Device		
Dimensions (L x W x H)	29.1 x 22.0 x 10.1 cm	
Weight	2.15 kg	
Touchscreen LCD display	3.5 in. TFT module with capacitive touch	
LED light	Blue LED center wavelength: 465 nm; full width, half max: 20 nm	
Viewing surface dimensions	15.4 x 11.2 cm	
Amber filter dimensions	16.0 x 13.0 cm	
LED life span	50,000 hours	
LED specification	Array of 39 high-power LEDs with emission wavelength at 465 ± 10 nm	


Invitrogen™ E-Gel™ Power Snap Plus Camera		
Dimensions (L x W x H) 23.3 x 21.0 x 22.8 cm		
Weight	2 kg	
Touchscreen LCD display	8 in. TFT module with capacitive touch	
Camera type	Color, complementary metal-oxide semiconductor (CMOS)	
Gel-image resolution	4,208 x 3,120 pixels (13 MP), 8 bits	
Dynamic range	69 dB	
Image output	TIF, JPG, PNG	
Lens f-number	2.8	
Internal memory	64 GB SD card	

Ordering information

Product	Cat. No.*	Cat. No.**
E-Gel Power Snap Plus Electrophoresis Device	G9101	G9110
E-Gel Power Snap Plus Camera	G9201	G9200
E-Gel Power Snap Plus Electrophoresis System	G9301	G9311
E-Gel Power Snap Plus Electrophoresis System with E-Gel 48 1% agarose gels	G9341	G9331
E-Gel Power Snap Plus Electrophoresis System with E-Gel 48 2% agarose gels	G9342	G9332
E-Gel Power Snap Plus Electrophoresis System with E-Gel 96 1% agarose gels	G9391	G9381
E-Gel Power Snap Plus Electrophoresis System with E-Gel 96 2% agarose gels	G9392	G9382

^{*} North America/Europe/Middle East/Africa

^{**} Asia Pacific/Japan/Latin America/Greater China



E-Gel precast agarose gels

E-Gel precast agarose gels in 48- and 96-well with the option of 11- and 12-well gel analysis can be used with the Invitrogen E-Gel Power Snap Plus Electrophoresis System, and are offered with easy-to-use SYBR Safe stain and at percentages suitable for your lab's needs.

Find out more at thermofisher.com/powersnapplus

Fully integrated routine electrophoresis

E-Gel Power Snap Electrophoresis System

The Invitrogen™ E-Gel™ Power Snap Electrophoresis System combines rapid, real-time, low-throughput nucleic acid analysis with high-resolution image capture for ultimate convenience. The integrated design helps reduce workflow time and accelerate discovery.

Features of the E-Gel Power Snap Electrophoresis System include:

- Faster analysis—go from sample loading to image capture in as few as 15 minutes
- **Simple operation**—intuitive user interface with large touchscreen and integrated operating system
- Safer workflow—minimize handling of hazardous chemicals when used with E-Gel precast gel cassettes



Specifications

Invitrogen™ E-Gel™ Power Snap Electrophoresis Device		
Dimensions (L x W x H)	24 x 13 x 7 cm	
Weight	1 kg	
Touchscreen LCD display	7.7 x 4.4 cm	
LED light Blue LED; center wavelength: 468 full width, half max.: 20 nm		
Viewing surface dimensions	9 x 11 cm	
Amber filter dimensions	8.6 x 10.5 cm	
LED life span	50,000 hours	
LED specification	Array of 12 high-power LEDs with emission wavelength at 465 ± 10 nm	

Invitrogen™ E-Gel™ Power Snap Camera		
Dimensions (L x W x H)	26 x 13 x 15 cm	
Weight	1 kg	
Touchscreen LCD display	11.5 x 8.6 cm	
Camera type	Color, CMOS	
Gel-image resolution	1,600 x 1,944 pixels (3 MP), 8 bits	
Dynamic range	68 dB	
Image output	TIF (grayscale) and JPG (color) file formats	
Lens f-number	2.8	
Internal memory	32 GB	

Ordering information

Product	Cat. No.
E-Gel Power Snap Electrophoresis Device	G8100
E-Gel Power Snap Electrophoresis Camera	G8200
E-Gel Power Snap Electrophoresis System	G8300
E-Gel Power Snap Electrophoresis System Starter Kit, E-Gel EX gels (1%)	G8341ST
E-Gel Power Snap Electrophoresis System Starter Kit, E-Gel EX gels (2%)	G8342ST
E-Gel Power Snap Electrophoresis Device Starter Kit, E-Gel SizeSelect II gels (2%)	G8162ST



To learn more, go to thermofisher.com/powersnap



Did you know?

Samples can be run up to twice as fast with E-Gel precast gels compared to conventional handcast gels. E-Gel precast gels are available in a variety of formats: from low to high throughput, with SYBR Safe stain, and with agarose percentages suitable for either general-purpose or high-resolution separations.

Pour your own agarose gels

Stains and reagents for electrophoresis

Pour-your-own electrophoresis reagents

Choosing the right tools for nucleic acid electrophoresis can significantly improve and accelerate results, enabling you to address downstream applications sooner.

DNA stains

Detection of nucleic acid samples in gels can be improved using fluorescent dyes that are safer and/or more sensitive than ethidium bromide. The Invitrogen™ SYBR™ Safe, SYBR™ Green I, and SYBR™ Gold stains are less toxic and/or more sensitive with lower background fluorescence compared to conventional ethidium bromide stain.

Find out more at thermofisher.com/stains

UltraPure reagents for electrophoresis

Invitrogen™ UltraPure™ reagents are specifically formulated to meet your nucleic acid analysis and purification needs. UltraPure agarose and reagents are made from highly purified biochemicals for maximum reliability and superior performance.

Find out more at thermofisher.com/ultrapure







Ordering information

Product	Quantity	Cat. No.
SYBR Safe DNA Gel Stain	400 μL	S33102
SYBR Gold Nucleic Acid Gel Stain	500 μL	S11494
UltraPure Agarose	100 g	16500100
Tracklt 1 Kb Plus DNA Ladder	100 applications	10488085
UltraPure TAE Buffer, 10X	4 L	15558026
UltraPure DNase/RNase-Free Distilled Water	10 x 500 mL	10977023
UltraPure Ethidium Bromide, 10 mg/mL	10 mL	15585011



Specialty enzymes for molecular diagnostics development and commercial supply

Whether you need flexibility in assay design or custom functional testing, we can help accelerate your nucleic acid-based assay development with our innovative DNA-free and lyo-ready enzymes.

DNA-free PCR enzymes*

We've pioneered the use of single-use systems (SUS) to manufacture DNA-free PCR enzymes. The advantages include:

- Manufactured in a completely closed system using dedicated or single-use equipment
- Verified free of contaminating DNA from host, human operator, and environment
- Produced in an ISO 13485–certified facility for high standards of quality
- * Invitrogen™ Platinum™ *Taq* DNA Polymerase, DNA-free, is available as a catalog product.

Learn more at thermofisher.com/dna-free

Lyo-ready DNA polymerases and reverse transcriptases

From Platinum *Taq* DNA polymerases to SuperScript reverse transcriptases, we offer the largest selection of lyophilization-compatible enzymes, providing:

- The same functional enzyme performance as with the conventional format
- Tailor-made solutions for your specific applications, including custom formulation
- Higher confidence in results with low residual DNA contamination (human and bacterial)

Learn more at thermofisher.com/lyo-ready

Molecular biology education

Have you ever had questions about PCR but didn't want to ask? Find answers in our online education hub for molecular biology that features free technical content, fun videos, educational webinars, and application notes. Learn about PCR, reverse transcription, and more—at your own pace and on your own time—to amplify your research to the next level.

Find out more at thermofisher.com/molbioschool





Educational webinars



<u>Technical</u> <u>videos</u>



Application notes



Mobile app: PCR Quest lab game


Test your PCR knowledge with our lab game, PCR Quest, where you travel from lab to lab crushing the world's toughest diseases. Download the game at thermofisher.com/pcrquest





Find out more at thermofisher.com/pcrworkflow

For Research Use Only. Not for use in diagnostic procedures. © 2019–2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a trademark of Roche Molecular Systems, Inc., used under permission and license. AmpliTaq Gold is a trademark of Roche Molecular Systems, Inc., used under permission and license. Apple is a trademark of Apple Inc. Google is a trademark of Google Inc. COL35052 1222



Real-time PCR solutions

More efficiency to take your research further

Discover more efficient real-time PCR (qPCR)

The family of most-cited Applied Biosystems[™] qPCR solutions delivers the caliber of results you need, faster and easier than you could have imagined. Free yourself from extra steps, preventable rework, and other inefficiencies so you can focus on making meaningful discoveries.

Experience the difference the efficient workflow of Applied Biosystems qPCR solutions can make to your research success:

- Connected instruments incorporate smart technologies to enable hands-free operation and remote support, while a suite of cloud-based analysis software provides exceptional flexibility.
- Predesigned Applied Biosystems[™] TaqMan[™] Assays provide everything you need, including Applied Biosystems[™] TaqMan[™] probe and PCR primer sets, formulated to work right out of the box. No additional design, optimization, or post-run analysis steps are required.
- Complement your assay with an Applied Biosystems[™] TaqMan[™]
 master mix and experience the gold standard in qPCR
 performance. We also offer Applied Biosystems[™] SYBR[™] dye
 chemistry with advanced features, like dye tracking, that help
 prevent pipetting errors.

More choices. More flexibility. More efficiency. Experience more from your qPCR.

Contents

Which qPCR instrument fits your needs? Real-time PCR applications Real-time PCR software	4 6 7
QuantStudio 3 and 5 Real-Time PCR Systems QuantStudio 6 Pro and 7 Pro Real-Time PCR Systems QuantStudio 12K Flex Real-Time PCR System Orbitor RS2 Microplate Mover	8 10 12 14
TaqMan Array Plates with RFID	15
MicroAmp qPCR plastics	16
Which qPCR plastic fits your needs?	18
TaqMan and SYBR Green chemistries for qPCR	20
TaqMan Assays to accelerate your research	21
TaqMan and SYBR Green master mixes	23
SuperScript IV VILO Master Mix	25
SuperScript IV CellsDirect cDNA Synthesis Kit	26
nstrument Connect app	26
Digital PCR (dPCR) QuantStudio Absolute Q Digital PCR System nstrument services and support	27 28 30

Which qPCR instrument fits your needs?

QuantStudio real-time PCR systems

	Ultimate simplicity	Total control
	QuantStudio 3 system	QuantStudio 5 system
		The second of th
Colors	4 colors	5 or 6 colors (21 filter combinations)
Available formats*	96-well (0.1 mL block) 96-well (0.2 mL block)	96-well (0.1 mL block) 96-well (0.2 mL block) 384-well
Dimensions (H x W x D)	40 x 27 x 50 cm	40 x 27 x 50 cm
Block change	Fixed	Fixed
VeriFlex temperature control	Yes, 3 zones	Yes, 6 zones (96-well blocks only)
Throughput	Medium	Medium
21 CFR Part 11-enabled	Security	Security, auditing, e-signature package
Touchscreen	Yes, interactive	Yes, interactive
Cloud-enabled	Yes	Yes
Automation compatible	NA	NA
Key applications	 Gene expression miRNA profiling SNP genotyping Copy number variation Protein thermal shift High-resolution melt Pathogen detection 	 Gene expression miRNA profiling SNP genotyping Copy number variation Protein thermal shift High-resolution melt Pathogen detection
* Some blocks may not be available		I

 $^{^{\}ast}$ Some blocks may not be available.

Some instruments are also available in a diagnostic format. Learn more at $\underline{\text{thermofisher.com/qsdx}}$

For information on other Applied Biosystems[™] PCR instruments, please visit **thermofisher.com/qpcrinstruments**

Smart and connected	Smarter productivity	Maximum productivity
QuantStudio 6 Pro system*	QuantStudio 7 Pro system*	QuantStudio 12K Flex system
5 colors	6 colors (21 filter combinations)	6 colors (21 filter combinations)
96-well (0.1 mL block) 96-well	96-well (0.1 mL block) 96-well	96-well (0.1 mL block) 96-well
(0.2 mL block)	(0.2 mL block)	(0.2 mL block)
384-well	384-well TaqMan Array Card (384-well microfluidic card)	384-well TaqMan Array Card (384-well microfluidic card) OpenArray plates (3,072 through-holes)
54.7 x 33.8 x 52.5 cm	54.7 x 33.8 x 52.5 cm	73.8 x 50.1 x 66.1 cm
Interchangeable; no tools required	Interchangeable; no tools required	Interchangeable; no tools required
Yes, 3 zones (96-well blocks only)	Yes, 6 zones (96-well blocks only)	NA
High	High	Very high
Upgradeable	Security, auditing, e-signature package	Optional security, auditing, e-signature package
Yes, interactive	Yes, interactive	Yes
Yes	Yes	No
Upgradeable	Yes	Yes
Gene expression	Gene expression	Gene expression
 miRNA profiling 	miRNA profiling	miRNA profiling
 SNP genotyping 	SNP genotyping	SNP genotyping
 Copy number variation 	Copy number variation	Copy number variation
 Protein thermal shift 	Protein thermal shift	Protein thermal shift
High-resolution melt	High-resolution melt	High-resolution melt
Pathogen detection	Pathogen detection	Pathogen detection
 Pharmacogenomics 	Pharmacogenomics	Pharmacogenomics
		Growing menu of qualified solutions



Real-time PCR applications

Real-time PCR is used for sensitive, specific detection and quantification of nucleic acid targets. We have developed powerful assay design algorithms, optimized master mixes, intuitive data analysis software, and flexible instrumentation to help harness the power of qPCR across a rich and diverse set of applications. Explore our robust solutions for your qPCR-based research.

Infectious disease research

See our growing catalog of sensitive, specific qPCR probe and primer sets to study microbe biology and pathogenicity, host response, and disease pathology.

Food pathogen detection

Detect multiple bacteria in the same run, including Salmonella, Campylobacter, E. coli O157:H7, and Listeria monocytogenes.

Waterborne pathogen detection

Designed to detect and monitor waterborne pathogens in recreational and drinking water supplies.

Pharmaceutical analytics

Designed to detect mycoplasmas, viruses, and residual host cell contamination for pharmaceutical, cosmetics, and personal care product manufacturing.

Qualified solutions

A growing menu of new and valuable content for you to use on your high-throughput qPCR systems, including pharmacogenomics, vaginal microbiota, and *CFTR* mutation analysis.

Stem cell research

Solutions for analyzing stem cells, determining stemness, and studying gene regulation and translation in stem cells.

Pharmacogenomics research

Predesigned TaqMan Assays for more than 175 ADME and CYP targets, including >95% of ADME core markers and a warfarin metabolism panel.

Oncology and genetic disease research

Enabling robust, reliable detection and quantitation of markers for cancer and genetic diseases.

Plant sciences and agricultural biotechnology

Instruments, reagents, and kits designed for plant researchers that enable remarkable agricultural discoveries—from improved crops that feed more people to sustainable biofuels.

Other key applications include gene expression, genotyping, and sequencing.

Real-time PCR software

A suite of analysis modules for the QuantStudio instrument family is available on the Thermo Fisher™ Connect Platform, enabling users to access and analyze their data anytime and anywhere. Just create and sign in to your account to use the Connect Platform for secure, cloud-based data storage, scientific analysis apps, and peer collaboration tools. Use the asset management apps and web tools to schedule time on your lab's instruments via your mobile device.

thermofisher.com/connect

Applied Biosystems qPCR analysis modules are a set of modules that provide an online toolkit for the analysis of qPCR data. The software takes advantage of the Connect Platform to provide highly versatile analysis tools that are flexible, fast, and easy to use and facilitate functional understanding of qPCR and related data. These modules allow users to combine over 100 genotyping, expression, or qPCR experiments into a project and analyze the data within minutes. Featured modules include:

- Design and Analysis
- High-Resolution Melt (HRM) Analysis
- Presence/Absence Analysis
- Standard Curve
- Genotyping
- hPSC Scorecard Analysis
- Relative Quantification

Design and Analysis

The Design and Analysis app offers the ability to create, edit, and analyze qPCR instrument files.

HRM

High-Resolution Melt (HRM) Analysis

This HRM Analysis app is designed for post-PCR analysis to identify variation in nucleic acid sequences. The method is based on detecting small differences in PCR melting (dissociation) curves. It is enabled by high-brightness, dsDNA-binding dyes used in conjunction with qPCR instrumentation that has precise temperature ramp control, advanced data capture capabilities, and access to software designed specifically for HRM analysis.

PA Presence/Absence Analysis

The Presence/Absence Analysis app analyzes and interprets qPCR data or post-read data to determine if a specific target sequence is present in a sample or not. It provides an easy-to-view presence/absence result in a plate grid

Standard Curve

view.

The Standard Curve app offers reliable quantification of unknown quantities of genes and enables importing of standard curves from other experiments, providing analysis flexibility.

Genotyping

The Genotyping app includes improved visuals and integrated traces of allelic discrimination plots to allow thorough QC of SNP assays to accurately reflect the true signals versus background noise.

hPSC Scorecard Analysis The hPSC Scorecard Analysis

The hPSC Scorecard Analysis app scores gene expression profiles generated using the Applied Biosystems™ TaqMan™ hPSC Scorecard™ Panel compared to a reference set of well-characterized pluripotent stem cell (PSC) lines. It provides raw data as well as box plots, heat maps, and correlation plots for export.

Relative Quantification

The Relative Quantification app allows fast and powerful gene expression analysis with enhanced visual capabilities for relative quantification, including integrated correlation and volcano and cluster analysis with the ability to drill down to amplification plots.



Intuitive and easy to use for all levels of experience

QuantStudio 3 and 5 Real-Time PCR Systems

The Applied Biosystems™ QuantStudio™ 3 and QuantStudio™ 5 Real-Time PCR Systems provide intuitive touchscreen usability, allowing you to stay connected to your data easily. They're designed for both new and experienced users who need simple and affordable instruments without compromising performance or quality.

Get a premium instrument at an affordable price

Access, analyze, and share data anytime and anywhere with the Connect Platform—cloud-enabled services that support remote monitoring of your runs in real time, analyze sophisticated datasets in minutes, securely store data, and share results online with colleagues across institutions and around the world

Obtain results you can trust—detect differences in target quantity as small as 1.5-fold in singleplex reactions, and obtain 10 logarithmic units of linear dynamic range

Multiplex with ease—up to six excitation and six emission filters offer 21 different color combinations, allowing a broad range of detection chemistries and maximum multiplexing

Helps save valuable time—QuantStudio 3 and 5 systems offer three or six independent temperature zones for flexibility to run multiple experiments simultaneously; fast thermal cycling enables results in less than 30 minutes

Get up and running quickly—instruments are factory-calibrated for accuracy, quick installation, and immediate use; preoptimized protocol templates help minimize training for new users, and the included SmartStart Orientation provides basic qPCR training and setup for both the Connect Platform and the Instrument Management tool

Maximize benchtop space—the compact instrument can be configured stand-alone or with a computer









Specifications

	QuantStudio 3 system	QuantStudio 5 system
Sample capacity (wells)	96	96 or 384
Reaction volume	0.1 mL block: 10-30 μL 0.2 mL block: 10-100 μL	96-well, 0.1 mL block: 10–30 μL 96-well, 0.2 mL block: 10–100 μL 384-well block: 5–20 μL
Footprint (H x W x D)	40 x 27	x 50 cm
Excitation source	Bright w	hite LED
Optical detection	4 coupled filters	96-well: 6 decoupled filters 384-well: 5 coupled filters
Excitation/detection range	450–600 nm/500–640 nm	96-well: 450–680 nm/500–730 nm 384-well: 450–650 nm/500–700 nm
Multiplexing	Up to 4 targets	96-well: up to 6 targets 384-well: up to 5 targets
2D barcode reading	Opti	ional
Heating/cooling method	Pel	tier
Temperature zone function	3 VeriFlex zones	96-well: 6 VeriFlex zones 384-well: NA
Max. block ramp rate	0.2 mL block: 6.5°C/sec 0.1 mL block: 9.0°C/sec	0.2 mL block: 6.5°C/sec 0.1 mL block: 9.0°C/sec 384-well block: 6.0°C/sec
Average sample ramp rate	3.66°	C/sec
Temperature uniformity	0.4	1°C
Temperature accuracy	0.2	5°C
Run time	<30 min runs	96-well block: <30 min runs 384-well block: <35 min runs
Dye compatibility (name)	FAM/SYBR Green, VIC/JOE/HEX/TET, ABY/NED/TAMRA/Cy3, JUN, ROX/Texas Red	FAM/SYBR Green, VIC/JOE/HEX/TET, ABY/NED/TAMRA/Cy3, JUN, ROX/Texas Red, Mustang Purple, Cy5/LIZ, Cy5.5
Chemistry capabilities	Fast/st	andard
Features to assist with 21 CFR Part 11 compliance	No	Yes, with no additional fees
Detection sensitivity	1 c	ору
Sensitivity	Detect differences as small as 1.5-fold in	n target quantities in singleplex reactions

Ordering information

Product		Instrument + 1-year extended warranty with AB Assurance Cat. No.*
QuantStudio 3 Real-Time PCR System (96-well, 0.1 mL block)**	A28136	A33777
QuantStudio 3 Real-Time PCR System (96-well, 0.2 mL block)**	A28137	A33779
QuantStudio 5 Real-Time PCR System (96-well, 0.1 mL block)**	A28138	A33619
QuantStudio 5 Real-Time PCR System (96-well, 0.2 mL block)**	A28139	A33624
QuantStudio 5 Real-Time PCR System (384-well block)**	A28140	A33628

^{*} Includes SmartStart Orientation.

Recommended plastics

96-well block	Cat. No.	96-well block Fast	Cat. No.	384-well block	Cat. No.
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354	MicroAmp EnduraPlate Optical 96-Well Fast Clear Reaction Plate with Barcode	4483485	MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode	4483285
MicroAmp Optical 96-Well Reaction Plate with Barcode	4306737	MicroAmp Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL	4346906	MicroAmp Optical 384-Well Reaction Plate with Barcode	4309849
MicroAmp Optical 8-Tube Strip with Attached Optical Caps, 0.2 mL	A30588	MicroAmp Fast 8-Tube Strip	4358293	MicroAmp Optical 384-Well Reaction Plate	4343370
MicroAmp Optical Adhesive Film	4360954	MicroAmp Optical 8-Cap Strip	4323032	MicroAmp Optical Adhesive Film	4360954
		MicroAmp Optical Adhesive Film	4360954		

 $^{^{\}star\star}\, \text{Does not include computer. Additional Cat. Nos. are available that include laptop or desktop computer.}$

Smarter productivity and improved workflow

QuantStudio 6 Pro and 7 Pro Real-Time PCR Systems

The Applied Biosystems™ QuantStudio™ 6 Pro and 7 Pro Real-Time PCR Systems are designed to save space on the lab bench, increase format flexibility, and support quick and simple block changes. These systems also deliver a smart workflow experience with advanced features that are designed for efficiency and productivity.

Personalized

- · Automatically load your settings and plate setup
- · Simply log in with facial authentication, no passwords required
- SmartStart Orientation, which includes one day of on-site training that covers basic instrument operation and maintenance, and a choice of hands-on application training (including reagents) on either absolute quantification, gene expression, or genotyping

Efficient

- Minimize hands-on time with voice commands facilitating hands-free operations
- Eliminate manual steps to obtain plate layout, protocol, and assay information on the instrument when using Applied Biosystems™ TaqMan™ Array Plates with RFID (see <u>page 15</u> for more information)
- Maximize uptime with the push-button Smart Help feature to access technical support scientists for faster, more efficient troubleshooting and resolution
- Reduce downtime using the Smart Remote Support collaboration tool with real-time video and audio, and desktop support to guide you through repairs remotely

Simple

- Streamlined workflow directly from touchscreen
- Simple, tool-free block changes



* Automation is available for the QuantStudio 7 Pro system and is available for the QuantStudio 6 Pro system after upgrade.

Productive

- Access data anytime and anywhere with cloud-enabled services from the Connect Platform
- Built-in help videos for commonly used activities
- Increased ergonomics due to larger touchscreen, improved screen angle, and optional Bluetooth[™] keyboard

The QuantStudio 7 Pro system also features automation compatibility and a security, auditing, and e-signature (SAE) package.

Specifications

	QuantStudio 6 Pro system	QuantStudio 7 Pro system				
Sample capacity (wells)	96, 384	96, 384, TaqMan Array Card				
Reaction volume	0.1 mL block: 10–30 μL 0.2 mL block: 10–100 μL 384-well block: 5–20 μL	0.1 mL block: 10–30 μL 0.2 mL block: 10–100 μL 384-well block: 5–20 μL TaqMan Array Card: ~1 μL				
Footprint (H x W x D)	54.7 x 33.8	3 x 52.5 cm				
Excitation source	Bright w	hite LED				
Filter or color combinations	5	21				
Multiplexing	5 targets	6 targets				
Excitation/emission range	450–680 nm	/500–730 nm				
2D barcode reading	Via USB connection					
Heating/cooling method	Peltier					
Temperature zone function	3 independent temperature zones	6 independent temperature zones				
Maximum ramp rate	6.5°C	6.5°C/sec				
Average sample ramp rate	3.66°	C/sec				
Temperature uniformity	0.4	4°C				
Temperature range	4–99	9.9°C				
Temperature accuracy		5°C				
Run time	Less than 30 min	Less than 30 min				
Dye compatibility (name)	FAM/SYBR Green, VIC/JOE/TET/HEX, TAMRA/NED/ABY, ROX/	JUN/Texas Red, LIZ/Cy5/Mustang Purple, CY5.5/Alexa Fluor dye				
Features to assist with 21 CFR Part 11 compliance	Upgradeable	Yes				
Detection sensitivity		ору				
Sensitivity	Detect differences as small as 1.5-fold in	n target quantities in singleplex reactions				

Ordering information

Product	Cat. No.	Instrument + extended warranty package Cat. No.*
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.1 mL block	A43160	A47200
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.1 mL block, laptop	A43167	A47202
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.1 mL block, desktop	A43181	A47204
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.2 mL block	A43159	A44288
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.2 mL block, laptop	A43166	A44290
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.2 mL block, desktop	A43180	A44292
QuantStudio 6 Pro Real-Time PCR System, 384-well	A43161	A45582
QuantStudio 6 Pro Real-Time PCR System, 384-well, laptop	A43168	A45584
QuantStudio 6 Pro Real-Time PCR System, 384-well, desktop	A43182	A45586
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.1 mL block	A43163	A47201
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.1 mL block, laptop	A43170	A47203
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.1 mL block, desktop	A43184	A47205
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.2 mL block	A43162	A44289
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.2 mL block, laptop	A43169	A44291
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.2 mL block, desktop	A43183	A44293
QuantStudio 7 Pro Real-Time PCR System, 384-well	A43164	A4558383
QuantStudio 7 Pro Real-Time PCR System, 384-well, laptop	A43171	A455855585
QuantStudio 7 Pro Real-Time PCR System, 384-well, desktop	A43185	A4558745587
QuantStudio 7 Pro Real-Time PCR System, TaqMan Array Card block	A43165	A47651
QuantStudio 7 Pro Real-Time PCR System, TaqMan Array Card block, laptop	A43172	A47653
QuantStudio 7 Pro Real-Time PCR System, TaqMan Array Card block, desktop	A43186	A47655

^{*} The extended warranty package includes the instrument, SmartStart Orientation, and a 1-year AB Assurance Plan with 1 planned maintenance (PM) visit. Packages are not available in all countries. Contact your sales representative for more information.

Recommended consumables can be found on <u>page 15</u> (RFID-tagged plates), <u>page 16</u> (plastics), <u>page 21</u> (assays), <u>page 23</u> (master mixes), and <u>page 25</u> (reverse transcriptases).

Maximum productivity with minimum effort

QuantStudio 12K Flex Real-Time PCR System

This one instrument enables multiple users to conduct a wide range of experiments, from low- to high-throughput sample processing and for virtually any PCR application, such as:

- Drug discovery
- Pharmacogenomics research
- MicroRNA profiling
- · Agriculture molecular testing
- CFTR mutation analysis
- · Vaginal microbiota research

Miniaturization at a lower cost

Applied Biosystems™ OpenArray™ technology is a broadly applicable nanoliter fluidics platform for low-volume solution-phase reactions, and enables lower reagent and assay costs, and rapid parallel processing.

Superior throughput

The Applied Biosystems™ QuantStudio™ 12K Flex Real-Time PCR System can simultaneously run up to four 3,072-reaction Applied Biosystems™ OpenArray™ plates in about 4 hours. Up to 110,000 data points can be produced in an 8-hour day with the addition of a single Applied Biosystems™ ProFlex™ PCR System.

Outstanding flexibility

Easily switch between five available thermal cycling blocks: OpenArray plate, TaqMan Array Card, 384-well, and standard or Fast 96-well blocks. Load the thermal cycling block in less than 1 minute—no tools required.



Increased data integrity and quality control—with results you can trust


The integrated sample tracking and SAE module assist you in supporting 21 CFR Part 11 compliance.

Application spotlight

Pharmacogenomics research

Pharmacogenomics is the study of drug efficacy based on a subject's unique genomic composition. The QuantStudio 12K Flex system and OpenArray technology provide a simple, cost-effective, and fast workflow for the analysis of mutations and copy number variants associated with drug metabolism enzyme (DME) genes.

OpenArray workflow



Find out more at thermofisher.com/quantstudio12k

As low as 2-fold change for singleplex reaction

Specifications

	QuantStudio 12K Flex system				
Block change design	96-well (10-100	uL reactions)			
3 3	Fast 96-well (15-3	1 /			
	384-well (5–20	μL reactions)			
	TaqMan Array Card (~1 µL reactions)				
	OpenArray plate (33 nL reactions)				
Block change design	Block change from front in less than 1	min; no tools or service call required			
Excitation source	Enhanced OptiFlex s	system, white LED			
Instrument control	Instrument touchscreen, networked	computer, or attached computer			
Detection channels	Decoupled: 6 emission, 6 excitation (9)	96-well/Fast, 384-well, TaqMan blocks)			
Detection channels	Coupled: 4 emission, 4 ex	citation (OpenArray block)			
21 CFR Part 11 enablement	Optional software module				
Dimensions (H x W x D)	73.8 x 50.1 x 66.1 cm				
Weight	70 kg				
Remote monitoring	Available to monitor up to 15 netwo	orked instruments simultaneously			
	96-well, 96-well Fast, 384-well,				
	TaqMan Array Card blocks	OpenArray block			
Detection channels	Decoupled: 6 emission, 6 excitation	Coupled: 4 emission, 4 excitation			
Well-to-well variability	±0.25°C	±0.75°C			
Max. block ramp rate	3.0°C/sec (384-well)	3.0°C/sec			
Run time	30 min expected (Fast 96-well block) 2 hr (gene expression)				
	35 min (384-well block, using Fast master mix)	4 hr (genotyping)			
Demonstrated sensitivity	То 1 сору	To 1 copy			
Dynamic range	To 9 logarithmic units	To 7 logarithmic units			

As low as 1.5-fold change for singleplex reaction

Ordering information

Resolution

Product	Cat. No.	Instrument + 1-year extended warranty with AB Complete Cat. No.*
QuantStudio 12K Flex Real-Time PCR System, OpenArray block with AccuFill System, desktop configuration	4471090	4480621
QuantStudio 12K Flex Real-Time PCR System with OpenArray block without AccuFill System, desktop configuration	4472380	NA
QuantStudio 12K Flex Real-Time PCR System with TaqMan Array Card block, desktop configuration	4471089	4480622
QuantStudio 12K Flex Real-Time PCR System 384-well block, desktop configuration	4471134	4480623
QuantStudio 12K Flex Real-Time PCR System 96-well Fast block, desktop configuration	4471088	4480625
QuantStudio 12K Flex Real-Time PCR System 96-Well block, desktop configuration	4471087	4480631
96-Well Block Upgrade Kit	4453543	NA
96-Well Fast Block Upgrade Kit	4453544	NA
384-Well Block Upgrade Kit	4453545	NA
TaqMan Array Card Block Upgrade Kit	4453546	NA
OpenArray Block with AccuFill System	4471067	NA
QuantStudio 12K Flex AccuFill Upgrade Kit (For existing AccuFill system users)	4471022	NA
<u> </u>	·	·

^{*} Includes SmartStart Orientation.

Recommended plastics

96-well block	Cat. No.	96-well block Fast	Cat. No.	384-well block	Cat. No.
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354	MicroAmp EnduraPlate Optical 96-Well Fast Clear Reaction Plate with Barcode	4483485	MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode	4483285
MicroAmp Optical 96-Well Reaction Plate with Barcode	4306737	MicroAmp Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL	4346906	MicroAmp Optical 384-Well Reaction Plate with Barcode	4309849
MicroAmp Optical 8-Tube Strip with Attached Optical Caps, 0.2 mL	A30588	MicroAmp Fast 8-Tube Strip	4358293	MicroAmp Optical 384-Well Reaction Plate	4343370
MicroAmp Optical Adhesive Film	4360954	MicroAmp Optical 8-Cap Strip	4323032	MicroAmp Optical Adhesive Film	4360954
		MicroAmp Optical Adhesive Film	4360954		

Multiple colors are available for most Cat. Nos.

The perfect solution to automate a variety of workflows

Orbitor RS2 Microplate Mover

The easy-to-use, flexible, and reliable Thermo Scientific™ Orbitor™ RS2 Microplate Mover is a collaborative benchtop mover that provides superior reliability and improved process efficiency. Its innovative bidirectional telescoping arm, coupled with its expansive 360° workspace, provides exceptional reach and precision, making the Orbitor RS2 Microplate Mover your trusted laboratory productivity partner.

Designed with you in mind

- Integrated barcode reader enables sample tracking, barcode transmission, and inventory management
- Plate detection in the gripper helps to eliminate labware handling errors and reduces the risk of lost samples
- Variable force settings make it compatible with many plate types for a variety of applications



Specifications

<u>'</u>					
General specifications Axis/joint specifications					
Power requirements	110-230 V, automatic switching	Axis Range Velocity			Acceleration
Operating environment	4–40°C, noncondensing	Z (height)	575 mm	750 mm/sec	2,250 mm/sec ²
Safety compliances	CE, CSA	Radius (reach)	±406 mm	200 mm/sec	3,600 mm/sec ²
Weight	25 kg	Theta (rotation)	Infinite	225 deg/sec	675 deg/sec ²
Motion control	Closed-loop servo, brushless AC motor	Gripper stroke	51 mm		
Control	Directly from PC, requires Windows [™] 7 or 10 operating system				
Storage options					
Random access hotel	8 plates per hotel, maximum capacity of 24 plates, best for 384-well plates				
Sequential access stack	40 plates per stack, maximum capacity of 80 plates, best for 384	4-well and 96-well (0).2 mL and 0	.1 mL) plates	

Ordering information

Product	Cat. No.
Orbitor RS2 Microplate Mover Bundle	A43780
Orbitor RS2 Microplate Storage, Sequential Access Stacks (3-pack)	A43781
Orbitor RS2 Microplate Storage, Random Access Hotels (3-pack)	A43782
Dual QuantStudio System-Orbitor RS2 Microplate Mover Integration Kit	A43779

One simple step gets you to your data faster

TaqMan Array Plates with RFID


TaqMan Array Plates with RFID and QuantStudio 6 Pro and 7 Pro systems create an improved user experience with smart features that eliminate the manual steps currently necessary to set up the plate file for your qPCR instrument run. Just insert the RFID plate, and the QuantStudio 6 Pro or 7 Pro system captures the plate setup information, including:

- Plate layout—assay IDs in each well
- Plate type
- · Expiration date
- Cat. No.


- Lot. No.
- Reaction volume per well
- Passive reference dye
- Thermal protocol



Conventional workflow



Smart qPCR workflow



Ordering information

Description	No. of assays	No. of samples	Cat. No. (96-well 0.1 mL)	Cat. No. (96-well 0.2 mL)
TaqMan Array Plates with RFID	NA	NA	A43	3823
TaqMan Array, 96-well Plate; format 8	8	12	4413263	4413266
TaqMan Array, 96-well Plate; format 16/16 Plus	16	6	4413261/4413262	4413264/4413265
TaqMan Array, 96-well Plate; format 32/32 Plus	32	3	4413259/4413260	4391528/4391529
TaqMan Array, 96-well Plate; format 48/48 Plus	48	2	4413257/4413258	4391526/4391527
TaqMan Array, 96-well Plate; format 96/96 Plus	96	1	4413255/4413256	4391524/4391525

High-performance qPCR plastics for optimal qPCR results

MicroAmp qPCR plastics


Applied Biosystems™ PCR plastics have been designed and validated to work with our thermal cyclers for more than 25 years. That's why they are engineer approved to enable optimal PCR performance.

Applied Biosystems[™] MicroAmp[™] qPCR plastics are:

- Validated on Applied Biosystems thermal cyclers for optimal fit and performance
- Designed to perform on all Applied Biosystems qPCR instruments
- Designed for optimal heat transfer with thin-walled polypropylene wells
- Designed to reduce cross-contamination with raised well rims for effective sealing



Unique, high-performance features of Applied Biosystems™ MicroAmp™ EnduraPlate™ plastic consumables





Options for every format and throughput need

Choose from tubes, tube strips, plates, adhesive film, and accessories for any throughput need. MicroAmp EnduraPlate plastics offer a solution for experiments that require special handling, such as automated or high-throughput workflows, and an even greater degree of durability for use with multi-instrument experiments.





The Applied Biosystems™ MicroAmp™ 8-Tube Strip with attached optical caps* is optimally designed for precise real-time PCR with lid and tube labeling, dual end tabs, and 20 µL graduation marks on each tube to prevent pipetting errors. The 8-tube strips fit in all 0.2 mL Applied Biosystems real-time PCR instruments.

 * MicroAmp 8-Tube Strip with attached domed caps is also available for PCR.



Find out more at **thermofisher.com/findplastics**



Did you know?

Proper plate sealing helps reduce evaporation and well-to-well contamination.

- 1. Remove the backing of the Applied Biosystems™ adhesive film.
- 2. Align the adhesive film so as to cover all wells while placing on the plate.
- 3. Rub the flat edge of the applicator along the long edge (length) of the plate, then along the short edge (width). Finally, rub the applicator between all the wells and around the outside edges of the plate using small back-and-forth motions to form a complete seal.



Which qPCR plastic fits your needs?

Find the plastic format with the throughput and features for your application

	Small-scale experiments with a few samples	Routine experiments	Automation	Laboratory use
Use for:	Single tubes, strips, caps, adhesive film, and accessories	MicroAmp optical microplates	MicroAmp EnduraPlate optical microplates	MicroAmp EnduraPlate optical microplates GPLE*
Formats	Single tubes	48-well Fast	• 96-well	• 96-well
	Single tubes with caps	• 96-well	96-well Fast	96-well Fast
	8-strip tubes with caps	96-well Fast	• 384-well	• 384-well
		• 384-well		
DNA-, RNase-, PCR inhibitor-free	Yes	Yes	Yes	Yes
ANSI/SBS standard dimension color	Clear	Clear	Single-color packs (red, blue, green, yellow, or clear) and 5-plate sampler (1 of each color)	Clear
Instrument compatibility	Use our plastics selection tool	Use our plastics selection tool	Use our plastics selection tool	Use our plastics selection tool
Barcode available	No	Yes (1 or 2 sides)	Yes (3 sides)	Yes (3 sides)
Multiple applications	No	No	Yes	Yes
Optical compatibility	Yes (applicable for optical version)	Yes	Yes	Yes
Use	Research use only	Research use only	Research use only	For laboratory use**

^{*} General purpose laboratory equipment (GPLE); for laboratory use.




Did you know?

Need high-quality PCR plastics for non–Applied Biosystems instruments?

Visit <u>thermofisher.com/thermoscientificplastics</u> for a wide range of Thermo Scientific™ PCR plastics.

Custom and OEM plastics for PCR and qPCR are available. Learn more at thermofisher.com/oemplastics



^{**} Lot-based contamination test with Certificate of Analysis.

Quickly find the plastics and accessories you need for your instrument

		40 !!				20	_		204 !!
		48-well			96-well 0.2 mL	96	-wel	I 0.1 mL	384-well
Applied Biosystems MicroAmp Product	Cat. No.	StepOne	7000	7300, 7500	QuantStudio 3, 5, 6 Flex, 7 Flex, 12K Flex, 6 Pro, 7 Pro; ViiA 7; 7900HT	StepOnePlus	7500	QuantStudio 3, 5, 6 Flex, 7 Flex, 12K Flex, 6 Pro, 7 Pro; ViiA 7; 7900HT	QuantStudio 5, 6 Flex, 7 Flex, 12K Flex, 6 Pro, 7 Pro; ViiA 7; 7900HT
96-well, 0.2 mL reaction plates									
Optical 96-Well plate	N8010560, 4316813		•	•	•				
Optical 96-Well plate with barcode	4306737, 4326659		•	•	•				
Optical 96-Well plate with barcode and optical caps	403012		•	•	•				
Optical 96-Well plate with barcode and optical adhesive films	4314320		•	•	•				
EnduraPlate Optical 96-Well Clear plate	A36924			•**	•				
EnduraPlate Optical 96-Well Clear plate with barcode*	4483354, 4483352			•**	•				
96-well, 0.1 mL reaction plates									
Fast Optical 96-Well plate	4346907					•	•	•	
Fast Optical 96-Well plate with Barcode	4346906, 4366932					•	•	•	
EnduraPlate Optical 96-Well Fast Clear plate	A36930					•	•	•	
EnduraPlate Optical 96-Well Fast Clear plate with barcode*	4483485, 4483494					•	•	•	
384-well reaction plates									
Optical 384-Well plate	4343370								•
Optical 384-Well plate with barcode	4309849, 4326270, 4343814								•
EnduraPlate Optical 384-Well Clear plate	A36931								•
EnduraPlate Optical 384-Well Clear plate with barcode*	4483285, 4483273								•
48-well reaction plates									
Fast Optical 48-Well Plate	4375816	•							
Strip tubes and caps									
Fast 8-Tube Strip, 0.1 mL	4358293	•				•	•	•	
Optical 8-Tube Strip with Attached Optical Caps, 0.2 mL	A30588		•	•	•				
Optical 8-Tube Strip, 0.2 mL	4316567		•	•	•				
Optical 8-Cap Strip	4323032	•	•	•	•	•	•	•	
Single tubes and caps									
Fast Reaction Tube with Cap, 0.1 mL	4358297	•				•		•	
Optical Tube without Cap, 0.2 mL	N8010933		•	•					
Seals and covers						1			
Optical Adhesive Film	4360954, 4311971		•	•	•	•	•	•	•
48-Well Optical Adhesive Film	4375323	•							
Reaction trays	1								
96-Well Tray/Retainer Set	403081		•						
Fast 48-Well Tray	4375282	•							
96-Well Tray for VeriFlex Block	4379983					•			
Accessories	1010000								
Splash-Free 96-Well Base	4312063		•	•	•	•	•	•	
96-Well Support Base	4379590		•	•	•	•	•	•	

^{*} Multiple colors are available.

Note: Experiments using one or two 8-tube strips with attached caps require blank tube strips to balance lid pressure on the block or the use of the Applied Biosystems™ MicroAmp™ 96-Well Tray/Retainer Set (Cat. No. 4381850)—bottom part of tray only. For use with 96-well block of Applied Biosystems™ 7000, 7300, 7500, and ViiA 7 systems, and QuantStudio 3/5/6/7/12K instruments.

^{**} Requires use of proper adapter, Cat. No. A24820.

TaqMan and SYBR Green chemistries for real-time PCR

We offer two types of chemistries to detect PCR products using real-time PCR instruments:

- Applied Biosystems[™] TaqMan[™] chemistry (also known as fluorogenic 5' nuclease chemistry)
- Applied Biosystems[™] SYBR[™] Green I dye chemistry

With more than 21 million predesigned assays, including 2.8 million gene expression assays for more than 32 species, 5 assay formats, and >200,000 publications, TaqMan Assays represent the most trusted and comprehensive collection of qPCR assays available.

Find master mixes, reagents, and kits to power your applications.

	TaqMan chemistry—based detection	SYBR Green dye-based detection
Chemistry overview	Uses a fluorogenic probe to enable detection of a specific PCR product as it accumulates during PCR cycles	Uses SYBR Green I, or similar dye that binds to double-stranded DNA to detect the PCR product as it accumulates during PCR

	TaqMan Assays and reagents	SYBR Green reagents
Specificity	High	Low
Sensitivity—low copy number	High	Variable*
Reproducibility	High	Variable*
Multiplexing	Yes	No
Predesigned assays	Yes	No
Custom assays	Yes	No
User design and optimization	No	Yes
Cost	High	Low*
Gene expression quantitation	High	Low
DNA quantitation	Yes	Yes (pathogen detection)
ChIP	Yes	Yes
SNP genotyping	Yes	No
MicroRNA	Yes	No
Copy number	Yes	No
Somatic mutation detection	Yes	No
Pathway analysis	Yes	No

^{*} Depends on template quality, and primer design and optimization.

Chemistry	Application	Starting material	Recommended master mix
	Gene expression, DNA quantification	cDNA, gDNA	TaqMan Fast Advanced Master Mix
5' nuclease (TaqMan Assays)	Gene expression, RNA virus quantification	RNA	TaqMan Fast Virus 1-Step Master Mix
	Genotyping	gDNA	TaqPath ProAmp Master Mix
SYBR Green dye	Gene expression	cDNA	PowerTrack SYBR Green Master Mix



TaqMan Assays to accelerate your research

Comprehensive, high-quality solutions for genetic analysis

TaqMan Assays are the industry-leading choice for 5' nuclease qPCR assays. They are cited in more publications than any other qPCR assay product and are considered the gold standard for quantitative genomic analysis. Backed by a performance guarantee,* TaqMan Assays are consistently chosen as a proven solution to reliably provide fast and accurate results.

- Specificity—advanced primer/probe sequence selection criteria plus minor groove binder (MGB) probe enhancement deliver the specificity and reproducibility you need for confidence that your results are generated from amplification of the intended target and not from nonspecific dye binding or amplification of closely related genes or pseudogenes
- Sensitivity—the nonfluorescent quencher (NFQ) on Applied Biosystems™ TaqMan™ probes minimizes background, and intelligent PCR primer and probe design maximizes amplification efficiency; get better sensitivity and accuracy to reliably detect targets present at 10 or fewer copies
- Reproducibility—accurately reproduce results from well to well, day to day, and lab to lab, even across manufacturing lots
- **Proven technology**—referenced by more than 200,000 publications to date

Guaranteed to perform for all your research needs*



TaqMan Assays cover a broad range of research for a wide variety of applications:

wide variety of applications:				
Expression	Genetic variation			
Gene expression	Single-nucleotide polymorphism			
 MicroRNA (miRNA) 	(SNP) genotyping			
expression	Drug metabolism enzyme genotyping			
 Long noncoding RNA 	Copy number variation (CNV)			
(IncRNA) expression	Rare somatic mutation detection			
 Fusion transcript detection 				
Protein expression				

^{*} Terms and conditions apply. See full details of the quarantee at thermofisher.com/tagmanguarantee

Applied Biosystems[™] TaqMan[™] Gene Expression Assays

- Designed to detect virtually any gene product, with more than 2.8 million predesigned assays
- Best-coverage assays available to detect the highest number of transcript variants possible
- Available for 32 species and some microbial pathogens
- Multiple format options—single tubes, 96-well and 384-well plates, TaqMan Array Cards, and OpenArray plates
- Can't find a predesigned assay to meet your needs?
 Use our Custom Assay Design Tool to design an assay based on your specific target sequence

Applied Biosystems[™] TaqMan[™] Advanced miRNA Assays

- Detect targets with as few as 60 copies of input microRNA (miRNA) in the cDNA synthesis reaction
- One universal reverse transcription (RT) step for all miRNA assays
- cDNA can be archived for future miRNA studies
- Identify only mature miRNA and distinguish related highly homologous miRNAs with gold-standard TaqMan probe specificity
- Detect and quantify mature miRNA from as little as 1 pg of total RNA or 2 µL of purified plasma or serum

Applied Biosystems[™] TaqMan[™] Noncoding RNA Assays

- For reliable detection and quantitation of noncoding transcripts longer than 200 nucleotides
- Designed to detect only noncoding transcript targets

Applied Biosystems[™] TaqMan[™] Fusion Assays

- Quantify the expression level of gene fusions using qPCR
- Orthogonal validation method for confirming next-generation sequencing (NGS) results

Applied Biosystems[™] TagMan[™] SNP Genotyping Assays

- Comprise the world's largest predesigned collection, with more than 17 million assays
- Rigorous design pipeline with >90 parameters for optimal primer-probe combinations
- Functional quality testing on at least 20 gDNA samples for each assay
- Flexible formats to accommodate any number of targets and samples

Applied Biosystems[™] TaqMan[™] Drug Metabolism Genotyping Assays

- Detect polymorphisms in high-value drug metabolism enzyme (DME) pharmacogenetics markers
- Cover specific SNP alleles, multinucleotide polymorphisms (MNPs), and insertions and deletions (indels)
- Targets derived from public databases, consortia, and published articles

Applied Biosystems[™] TaqMan[™] Copy Number Assays

- Evaluate copy number of genomic DNA targets
- Easy to interpret—Applied Biosystems[™] CopyCaller[™]
 Software provides the calculated copy number and predicted copy number, along with confidence value and z-score quality metrics
- Fast and simple—setup to primary analysis in 3–4 hours

Applied Biosystems[™] TaqMan[™] Mutation Detection Assays

- Powered by competitive allele-specific Applied Biosystems[™]
 TaqMan[™] PCR (castPCR[™]) technology
- Detect and measure somatic mutations in genes associated with cancer research
- Detect rare amounts of mutated DNA in a sample that contains large amounts of normal, wild-type DNA
- Compatible with different sample types, such as cell lines; formalin-fixed, paraffin-embedded (FFPE) tissue; and fresh-frozen tissue samples

Custom TagMan Probes

Design your own TaqMan primers and probes by choosing from dual-labeled Applied Biosystems™ TaqMan™ MGB, QSY™, or TAMRA™ custom probes. TaqMan MGB probes include an MGB moiety at the 3′ end that increases the melting temperature (T_m) of the probe and stabilizes probe—target hybrids. TaqMan MGB probes can be significantly shorter than traditional probes, providing better sequence discrimination and flexibility to accommodate more targets. TaqMan QSY probes are available for multiplexing three or more targets and are a seamless replacement for BHQ™ probes.

Flexible formats

A variety of formats to meet your system and throughput needs

TaqMan Assay formats

Configurations to fit your research goals

Are you analyzing hundreds (or thousands) of samples, and expression from a handful of genes? Or does your research involve a few samples that need to be analyzed for a long list of mRNA targets? No matter your experimental design, there's a TaqMan Assay format* that fits your research needs.



Single tubes

- Low entry price
- Flexible for small projects
- Run on any qPCR instrument



96- or 384-well plates

- Optimal for medium to large projects
- Balance flexibility with streamlined reaction setup
- Run on 96- or 384-well qPCR instrument



384-well TaqMan Array Cards

- Low cost per reaction
- Optimal for medium to large projects
- Run on QuantStudio 7 Pro, 7 Flex, and 12K Flex systems**

^{*} Also available on OpenArray plates for the QuantStudio 12K Flex system.

^{**} Also available on Applied Biosystems™ ViiA™ 7 and 7900HT systems.

TaqMan and SYBR Green master mixes

TaqMan and Applied Biosystems™ SYBR™ Green master mixes have been designed to give you the reliable results you need in run times as short as 40 minutes. Master mixes are optimized for specific applications and contain everything needed for reliable qPCR—buffer, dNTPs, passive reference dye, thermostable hot-start DNA polymerase, and other components. Simply add your sample and assay, and start your reaction.

Applied Biosystems[™] TaqMan[™] Fast Advanced Master Mix

Accurate, dependable gene expression quantification in less time

- Wide linear dynamic range even on fast cycling instruments (<40 min run time)
- 72-hour pre-PCR benchtop stability enables flexibility in automated workflows
- Compatible with single and multiplex assays

Applied Biosystems[™] TaqMan[™] Fast Virus 1-Step Master Mix

Sensitive 1-step quantification of RNA viruses and transcripts, even with challenging samples

- Designed for reliable, high-sensitivity 1-step RT-qPCR to enhance virus and transcript quantification
- Formulated to handle common RT-PCR inhibitors found in blood, stool, and other difficult samples
- Single-tube 4X master mix format enables use of more sample for enhanced sensitivity
- Works in singleplex and multiplex, and with exogenous or endogenous internal controls



Applied Biosystems[™] TaqPath[™] ProAmp[™] and ProAmp[™] Multiplex Master Mixes*

Accurate SNP and CNV genotyping calls in both purified and crude samples

- Reproducible genotyping results even in the presence of inhibitors
- Two formulations that enable single- to four-target detection per reaction
- Compatible with buccal swab and blood crude lysates prepared with the Applied Biosystems[™] DNA Extract All Reagents Kit
- 72-hour pre-PCR benchtop stability enables flexibility in automated workflows

Applied Biosystems[™] TaqPath[™] BactoPure[™] Microbial Detection Master Mix*

Less background DNA for lower-level microbial detection

- Low-level detection of bacterial, fungal, mammalian, and viral DNA, as well as of antibiotic resistance markers
- Tolerance of PCR inhibitors typically found in both purified samples and crude lysates in biopharmaceutical, molecular diagnostic, and research applications
- Wide dynamic range** (up to 8 orders of magnitude) enabling accurate detection from both low- and high-concentration samples

Explore Applied Biosystems™ QualTrak™ qPCR solutions for biopharma discovery at thermofisher.com/qpcr/biopharma

^{**} Dynamic range is a property of both the assay and template concentration in the sample, as well as the formulation of the master mix; thus, individual results may vary.



Learn more about master mixes and how to use our selection tool at **thermofisher.com/qpcrmastermixes**

^{*} TaqPath ProAmp and TaqPath BactoPure master mixes are General Purpose Reagents for Laboratory Use.

Applied Biosystems[™] PowerTrack[™] SYBR[™] Green Master Mix

SYBR Green dye-based gene expression quantification with unrivaled specificity

- Helps reduce pipetting errors with a two-color tracking dye system
- Formulated for superior specificity and tight reproducibility of C_t values over a broad dynamic range
- Broad primer T_m and primer concentration compatibility allows flexibility in qPCR reaction setup with minimal optimization
- Compatible with Invitrogen™ SuperScript™ IV VILO™ Master Mix for fast, reproducible reverse transcription results

Applied Biosystems[™] PowerUp[™] SYBR[™] Green Master Mix

Exceptional performance for even the most challenging qPCR applications

- High specificity with proprietary *Taq* DNA polymerase with dual-lock, hot-start mechanism
- Room-temperature stability of preassembled qPCR reactions for up to 72 hours
- Compatibility with standard or fast cycling for results in less than 50 minutes

Ordering information (more sizes available)

Product	Quantity	Cat. No.
	1 x 1 mL	4444556
	1 x 5 mL	4444557
TaqMan Fast Advanced Master Mix	1 x 50 mL	4444558
	1 x 1 mL (no UNG)	A44359
	1 x 5 mL (no UNG)	A44360
	1 x 1 mL	4444432
TaqMan Fast Virus 1-Step Master Mix	5 x 1 mL	4444434
	1 x 10 mL	4444436
	1 x 1 mL	A30865
TaqPath ProAmp Master Mix	1 x 10 mL	A30866
	1 x 50 mL	A30867
	1 x 1 mL	A30868
TaqPath ProAmp Multiplex Master Mix	1 x 10 mL	A30869
	1 x 50 mL	A30870
	1 x 1 mL	A52699
TaqPath BactoPure Microbial Detection Master Mix	1 x 5 mL	A52700
	1 x 50 mL	A52702
	1 mL	A46012
PowerTrack SYBR Green Master Mix	5 mL	A46109
	50 mL	A46113
	1 mL	A25741
PowerUp SYBR Green Master Mix	5 mL	A25742
	50 mL	A25743

Superior cDNA synthesis performance in RT-qPCR applications

SuperScript IV VILO Master Mix

SuperScript IV VILO Master Mix is a first-strand cDNA synthesis reaction mix for two-step RT-qPCR. The master mix format elevates the trusted VILO technology (variable input, linear output) to the next level by combining further optimized buffer conditions with the highly processive and thermostable Invitrogen™ SuperScript™ IV Reverse Transcriptase. The master mix offers exceptional performance features while maintaining superior linearity across the broadest range of input RNA.

- Super-efficient—C_t values earlier by an average of 2 cycles compared to other reverse transcription reagents, in a 10 min reaction
- Super-strong—reliable results even with degraded or inhibitor-containing RNA samples
- Super-reliable—improved RT-qPCR data reproducibility due to single-tube master mix format
- Super-safe—integrated, easy, and RNA-friendly genomic DNA removal

Find out more at thermofisher.com/4vilo



Doing other cDNA synthesis applications?

The SuperScript IV Reverse Transcriptase portfolio of products is engineered to offer superior cDNA synthesis performance with even the most challenging RNA samples. The portfolio includes a stand-alone enzyme, first-strand cDNA synthesis kit, and one-step RT-PCR kit.

thermofisher.com/ssiv

Ordering information

Product	Quantity	Cat. No.
SuperScript IV VILO Master Mix	50 reactions	11756050
Superscript iv vico iviaster iviix	500 reactions	11756500
CuparCarint IV//II O Master Mix with a DNIggs Faryma	50 reactions	11766050
SuperScript IV VILO Master Mix with ezDNase Enzyme	500 reactions	11766500



Did you know?

SuperScript IV VILO Master Mix is available in a format with the novel dsDNA-specific Invitrogen™ ezDNase™ Enzyme, which offers convenient and fast genomic DNA removal from RNA samples to help ensure high confidence in RT-qPCR results.

Why isolate RNA when you don't have to?

SuperScript IV CellsDirect cDNA Synthesis Kit

The Invitrogen™ SuperScript™ IV CellsDirect™ cDNA Synthesis Kit is designed to synthesize first-strand cDNA directly from mammalian cell lysates without first isolating the RNA. With lysis and reverse transcription (RT) performed in the same tube, the resulting first-strand cDNA is ready to use in many downstream applications such as PCR and qPCR.



Key benefits of this easy-to-use format

- Time-saving—save up to 2 hours in overall workflow time
- Superior sensitivity—easy detection of low-abundance targets
- Direct convenience—no need to isolate RNA
- qPCR compatibility—superior performance with PowerTrack SYBR Green Master Mix

Find out more at thermofisher.com/cellsdirect

Ordering information

Product	Cat. No.
SuperScript IV CellsDirect cDNA Synthesis Kit	
50 reactions	11750150
500 reactions	11750350
SuperScript IV CellsDirect Lysis Reagents	
500 reactions	11750550

Connect with your instrument and achieve lab-life balance

Instrument Connect app

Thermo Fisher Scientific has pioneered a way to connect you to your instruments, giving you real-time updates on your run and access to data as they come up. We connect you to your cloud-enabled instruments and benchtop devices anytime and anywhere using the Instrument Connect mobile app.

The Instrument Connect remote monitoring app allows

you to stay connected to any of our cloud-enabled instruments, including the QuantStudio 3, 5, 6 Pro, and 7 Pro qPCR instruments, as well as endpoint PCR devices including Applied Biosystems™ ProFlex™, SimpliAmp™, and MiniAmp™ thermal cyclers.

With the Instrument Connect app, you can:

- Check the availability of your cloud- and network-connected device
- Monitor run progress
- View amplification plots in real time (available for QuantStudio 3, 5, 6 Pro, and 7 Pro qPCR instruments)
- View plots and filter by sample or target in real time
- Schedule an instrument

It's easy to get started.

Just download the
Instrument Connect
app from the Apple™
App Store or Google
Play™ Store and log in
to your account on the
Connect Platform. You
can view your connected
instruments, monitor ren

can view your connected instruments, monitor remaining time in your run, and view your amplification plots in real time.

Learn more about the Instrument Connect app at thermofisher.com/connect



Discover powerfully simple digital PCR


Digital PCR (dPCR) is quickly becoming a reliable, complementary technology to qPCR for precise nucleic acid quantification because it can deliver a high level of precision, sensitivity, and robust performance. dPCR overcomes common limitations of qPCR, such as the need for routine standard curves, low precision when measuring rare targets, and lack of sensitivity in high-background conditions.

Experience the difference that the efficient workflow of Applied Biosystems dPCR solutions can make to your research success:

- Single instrument, single plate, simple workflow—eliminate cumbersome reaction preparation steps from your dPCR workflow. The Applied Biosystems™ QuantStudio™ Absolute Q™ Digital PCR System requires only one hands-on step that takes less than five minutes to complete with minimal technical skill.
- Fast time to answer—avoid wasting precious time generating microchambers, moving plates, and reading samples. In 90 minutes, reagent digitization, thermal cycling, and data collection are integrated into a single system with no manual transfer steps required.
- Industry-leading choice—the Applied Biosystems brand is relied on for genetic analysis, offering robust assays powered by our unique design algorithms and years of experience. Our easy-to-use predesigned assay technology makes dPCR simple with Applied Biosystems™ Absolute Q™ digital PCR assays.

dPCR that works how you need it, when you need it

QuantStudio Absolute Q Digital PCR System



The higher precision, sensitivity, and absolute nature of dPCR is ideal for research in:

- Rare-target detection, such as somatic mutation detection in oncology research
- Precise quantification of viral targets
- · Pathogen detection and load determination
- Generation of references and standards
- Copy number variation

Simple workflow

The simplified dPCR workflow requires one instrument, one consumable, and one hands-on step that takes less than 5 minutes to complete with minimal technical skill.

Fast turnaround time

In 90 minutes, you can generate highly precise dPCR data in a fraction of the time required for typical dPCR workflows.

Industry-leading consistency

MAP technology enables exceptional consistency in total microchambers analyzed per reaction.

Maximum reagent efficiency

More than 95% of the input sample is analyzed per reaction compared to the 25–60% range from other dPCR platforms.

Plate flexibility

Each 16-sample dPCR plate may be run one column at a time to preserve unused wells for future runs, minimizing waste.


Absolute quantification

The >20,000 microchambers per sample generate data expressed as copies/µL, enabling high precision and sensitivity without the need for a standard curve.

Application spotlight

Rare cancer mutation quantification in liquid biopsy

Liquid biopsies use cell-free DNA (cfDNA) that is derived from both normal and cancerous cells that have undergone apoptosis or necrosis and have released their DNA contents into their environment. The small fraction of cfDNA molecules originating from the tumor are called circulating tumor DNA (ctDNA), and these show cancer mutations unique to the tumor. dPCR is an established method for the fast, reliable, and accurate quantification of these cancer mutations that occur in low abundance. dPCR offers outstanding specificity and precision, making it a cost-effective and rapid research method for routine monitoring of ctDNA from liquid biopsy samples.



Find out more at thermofisher.com/absoluteq

Specifications

QuantStudio Absolute Q Digital PCR System	
Time to results	<2 hr
Colors	4 (endpoint detection)
Illumination	Rax, blue, phosphor green high-power LED
Detection channels	FAM, VIC, and HEX
Reader size (H x W x D)	54.0 x 60.0 x 62.0 cm
Weight	60 kg
QuantStudio Absolute Q Digital PCR MAP	
Samples per plate	16
Targets per plate	4
Microchambers per sample	20,480
Loading volume	9 μL
Sample waste/dead volume	<5%
Performance	
Dynamic range	5 logarithmic units
Precision at 95% confidence interval	±10%
Compatible chemistries	Absolute Q Digital PCR Assays and TaqMan Assays

Ordering information

Product	Cat. No.
QuantStudio Absolute Q Digital PCR System, desktop, extended warranty, SmartStart Orientation	A53267
QuantStudio Absolute Q Digital PCR System, desktop	A52864
QuantStudio Absolute Q MAP16 Plate Kit and Master Mix	A53301
QuantStudio Absolute Q MAP16 Plate Kit	A52865
Absolute Q DNA Digital PCR Master Mix (5X)	A52490
QuantStudio Absolute Q Isolation Buffer	A52730

Keep working wonders. We've got your back.

Count on superior services and support to keep your real-time and digital PCR instruments up and running



Service plans help protect your investment

All new instruments are protected by a standard factory warranty. This warranty covers all costs for travel, labor, and parts for repairs. Extended coverage service plans are also available at the time of instrument purchase. Services at additional charge include:

- Service plans with guaranteed response times*
- Priority technical support—priority phone access to our trained and certified technical and instrument support specialists
- 24-hour response time—expedited access to remote service engineers within 24 hours

Check your instrument's coverage status at

thermofisher.com/manage

Digital service innovations

Our instrument service plans include digital innovations that help keep your instruments and your lab running smoothly. With pioneering on-demand tools and capabilities such as remote support using mobile augmented-reality (AR) technology, instrument-driven support, and on-demand instrument training, we're constantly looking ahead so your lab doesn't fall behind.

Explore our digital service innovations

at thermofisher.com/innovations

Instrument hardware qualification (IQ/OQ/PQ)

Our manufacturer-trained and certified field service engineers will conduct and document comprehensive tests, including software and hardware compatibility matrixing, component verification, and site requirements, to verify performance and provide reliable, audit-style documentation to meet your regulatory requirements. Installation qualification (IQ), operational qualification (OQ), and performance qualification (PQ) or instrument performance verification (IPV) services are recommended at installation, after moving the instrument, after software or hardware upgrades, and after planned maintenance or critical repairs.

Contact an instrument qualifications specialist at thermofisher.com/igogpg

Analytical validation consulting services

Why handle your analytical validation (AV) in-house, when using our AV consulting services can speed up the launch process by up to 75% and cut costs by up to half? We'll consult with you to develop and optimize your validation workflow, while providing data analysis support and template documentation to fully maximize your instrument and reagent investment.

Get a complete list of assays and panels we support, or contact an AV consulting services specialist at thermofisher.com/av

^{*} Guaranteed response times are dependent on type of service plan and geographic region.

Service plans at a glance

		Repair Center service plan*		
	AB Platinum	AB Assurance**	AB Maintenance Plus	AB Repair Center Support Plus
Response time	Next business day [†]	2 business days [†]	3 business days‡	
Planned maintenance	•	•	•	
Access to technical support (Mon-Fri, standard business hours)	24/7/365	•	•	•
Parts, labor, and travel	•	•	10% discount	•
Digital remote support	•	•	•	•
Qualification service	•	Available as add-on	Available as add-on	
Field application scientist (FAS) consultation	•	Available as add-on	Available as add-on	
Loaner instrument issued during repair (Repair Center plan only)				Available as add-on

^{*} Repair Center service plan available for QuantStudio 3 and 5 qPCR systems only.

Education services

It can be difficult to prepare yourself for what's next while you're focused on the work you have now. Our professional, interactive training courses make it easier.

We offer a combination of virtual and in-person classroom instruction, and hands-on learning in your lab to match your schedule, budget, and learning preferences. Whichever course style you choose, you'll learn from one of our 300 highly skilled application scientists who are available to lead sessions online, at your location, or at one of our 12 training centers located worldwide.

Explore courses at <u>thermofisher.com/educationservices</u>

Technical support

If you have questions about product selection or use, assay or experimental design, data analysis, or troubleshooting, contact our team of technical support scientists or access our online product and application support tools.

How to reach us

To find your local support or technical support team, go to thermofisher.com/contactus


For product FAQs, protocols, training courses, and webinars, go to thermofisher.com/technicalresources



^{**} AB Assurance service plan offered for QuantStudio Absolute Q Digital PCR System.

[†] Availability limited in some geographic areas.

[‡] After receipt of purchase order.





Find out more at **thermofisher.com/qpcr**

For Research Use Only. Not for use in diagnostic procedures. © 2019–2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a trademark of Roche Molecular Systems, Inc., used under permission and license. AmpliTaq Gold is a trademark of Roche Molecular Systems, Inc., used under permission and license. Cy is a registered trademark of GE Healthcare. Apple is a trademark of Apple Inc. Google Play is a trademark of Google Inc. Bluetooth is a trademark of Bluetooth SIG, Inc. Windows is a trademark of Microsoft Corp. BHQ is a trademark of Biosearch Technologies Inc. COL35052 1222