

# Introduction to FT-IR Sample Handling

**Thermo**  
SCIENTIFIC

# Introduction


*Why is it important to know about different methods of sample handling?*

Certain techniques of sample handling are more effective than others for specific sample types. In order to obtain the best quality spectrum from your sample, it is important to know which handling technique works best for your sample type. Acquiring the best spectral data possible will give you more confidence in your results.

## Sampling Techniques

*Transmission*

### How does it work?



Conceptual diagram of the beam path through a transmission sample

The transmission technique does not require a separate accessory. The user simply places a sample directly into the infrared (IR) beam. As the IR beam passes through the sample, the transmitted energy is measured and a spectrum is generated. However, the analyst must often prepare the sample into a pellet, mull, film, etc. before the transmission measurement can be made. This requires expertise and can be time consuming.

### What types of samples can you analyze?

Excellent quality spectra can be obtained for many types of samples using transmission. The transmission technique can be used alone or in conjunction with accessories such as microscopes and liquid or gas cells to analyze:

- Organic powders in pellet or mull form
- Thermoplastic powders
- Soluble polymers
- Thin polymer films
- Regular-shaped polymers (with preparation)
- Irregular-shaped polymers (with preparation)
- Dark polymer films (not carbon-filled)
- Liquids (free-flowing or viscous)
- Gases (high concentrations to trace amounts)


### What are the advantages of transmission?

- **Economical** – cells and mounts are generally inexpensive
- **Well established** – most traditional form of sample measurement
- **Excellent spectral information** – ideal for qualitative measurements
- **Great for quantitative work** – many standard operating procedures are based on transmission

## Attenuated Total Reflection (ATR)

### How does it work?

An ATR accessory operates by measuring the changes that occur in an internally reflected IR beam when the beam comes into contact with a sample. An IR beam is directed onto an optically dense crystal with a high refractive index at a certain angle. This internal reflectance creates an evanescent wave that extends beyond the surface of the crystal into the sample held in contact with the crystal. In regions of the IR spectrum where the sample absorbs energy, the evanescent wave will be attenuated. The attenuated beam returns to the crystal, then exits the opposite end of the crystal and is directed to the detector in the IR spectrometer. The detector records the attenuated IR beam as an interferogram signal, which can then be used to generate an IR spectrum.



### What types of samples can you analyze?

ATR is ideal for strongly absorbing or thick samples which often produce intense peaks when measured by transmission. ATR works well for these samples because the intensity of the evanescent waves decays exponentially with distance from the surface of the ATR crystal, making the technique generally insensitive to sample thickness.

Other solids that are a good fit for ATR include homogeneous solid samples, the surface layer of a multi-layered solid or the coating on a solid. Even irregular-shaped, hard solids can be analyzed using a hard ATR crystal material such as diamond. Ideal solids include:

- Laminates
- Paints
- Plastics
- Rubbers
- Coatings
- Natural powders
- Solids that can be ground into powder

In addition, ATR is often the preferred method for liquid analysis because it simply requires a drop of liquid to be placed on the crystal. ATR can be used to analyze:

- Free-flowing aqueous solutions
- Coatings
- Viscous liquids
- Biological materials


### What are the advantages of ATR?

- **Minimal sample preparation** – place the sample on the crystal and collect data
- **Fast and easy cleanup** – simply remove the sample and clean the surface of the crystal
- **Analysis of samples in their natural states** – no need to heat, press into pellets, or grind in order to collect spectra
- **Excellent for thick or strongly absorbing samples** – ideal for difficult samples like black rubber

## Diffuse Reflectance (DRIFTS)

### How does it work?

When an IR beam is focused onto a fine particulate material, the incident beam can interact with the particle in one of several ways. First, radiation can be reflected off the top surface of the particle without penetrating the particle. Second, the light can undergo multiple reflections off particle surfaces without penetrating into the particle. True diffuse reflectance results from the penetration of the incident radiation into one or more sample particles and subsequent scatter from the sample matrix.



The IR beam interacting with a sample in a diffuse reflectance experiment

A DRIFTS accessory operates by directing the IR energy into a sample cup filled with a mixture of the sample and an IR transparent matrix (such as KBr). The IR radiation interacts with the particles and then reflects off their surfaces, causing the light to diffuse, or scatter, as it moves throughout the sample. The output mirror then directs this scattered energy to the detector in the spectrometer. The detector records the altered IR beam as an interferogram signal, which can then be used to generate a spectrum. Typically, a background is collected with the DRIFTS accessory in place and the cup filled with just the IR matrix. Excellent quantitative and qualitative data can be collected with proper sample preparation. However, transmission and ATR techniques are preferable to diffuse reflectance for quantitative data due to pathlength.

### What types of samples can you analyze?

DRIFTS is commonly used for the analysis of both organic and inorganic samples that can be ground into a fine powder (less than 10 microns) and mixed in a powder matrix such as potassium bromide (KBr). Typical sample types include:

- Soft powders and powder mixtures
- Hard polymers
- Rigid polymers

The DRIFTS technique can also be used with silicon carbide paper for the analysis of large intractable surfaces. Silicon carbide paper can be used to rub off a small amount of a variety of samples for analysis. This technique is a viable alternative to traditional sampling techniques for:

- Paint and varnish surfaces
- Tablets
- Rigid polymers

### What are the advantages of diffuse reflectance?

- **Little to no sample preparation** – just place in the sample cup
- **Fast and easy cleanup** – dump the cup and blow or rinse clean
- **No need for pressed KBr pellets or messy mulls** – samples can be run neat or diluted with KBr powder

# True Specular Reflectance/Reflection-Absorption

## How does it work?

True specular reflectance is a surface measurement technique that works on the principle of reflective efficiencies. This principle states that every sample has a refractive index that varies with the frequency of light to which it is exposed. Instead of examining the energy that passes through the sample, true specular reflectance measures the energy that is reflected off the surface of a sample, or its refractive index. By examining the frequency bands in which the rate of change in the refractive index is high, users can make assumptions regarding the absorbency of the sample. The true specular reflectance technique provides excellent qualitative data.




Diagram of the interaction of the beam using true specular reflectance




Diagram of the interaction of the beam using reflection-absorption

Reflection-absorption works on the same principle, but due to sample properties, some of the energy passes through the surface layer, is absorbed into the bulk of the sample, and then reflects off a substrate below the surface layer. A combination of true specular reflectance and reflection-absorption can occur when criteria for both techniques are met. If a qualitative comparison to transmission spectra is desired, users can apply the Kramers-Kronig correction to the data to remove the effects of dispersion.

## What types of samples can you analyze?

Specular reflectance is commonly used for the analysis of both organic and inorganic samples having large, flat, reflective surfaces. Reflection-absorption can occur when one of the above criteria is compromised and the sample has a reflective substrate present just below the surface. This type of analysis is commonly used for:

- Metallic surfaces
- Silicon wafers
- Thin films on reflective substrates
- Laminated materials on metals

## What are the advantages of specular reflectance?

- **Sensitivity to monolayer samples** – can detect Angstrom thick coatings on metal substrates
- **Nondestructive analysis** – no contact or sample damage during analysis
- **Wide range of accessories available** – can utilize main spectrometer and microscope accessories depending on the size of the sample and the thickness of surface layer

## RATINGS

■ Excellent
 ■ Good
 ■ Adequate

## KEY

\* capable of being ground

† incapable of being ground

Transmission  
 Transmission – IR Microscope  
 Diffuse Reflectance  
 Diffuse Reflectance – Si Carb Sampler  
 Diffuse Reflectance – IR Microscope  
 ATR  
 ATR – IR microscope  
 Specular Reflectance  
 Specular Reflectance – IR Microscope

	Transmission	Transmission – IR Microscope	Diffuse Reflectance	Diffuse Reflectance – Si Carb Sampler	Diffuse Reflectance – IR Microscope	ATR	ATR – IR microscope	Specular Reflectance	Specular Reflectance – IR Microscope
<b>Solids</b>	Thermoplastic Polymers (can be melted)	Excellent	Excellent		Adequate		Good	Good	
	Thermoplastic Polymers (can't be melted)				Adequate		Good	Good	
	Soluble Polymers	Excellent					Good		
	Thin Polymer Films	Excellent	Excellent				Good	Good	
	Thick Polymer Films						Excellent		Adequate
	Flat, Smooth Polymers *	Excellent	Excellent				Good	Good	Adequate
	Flat, Smooth Polymers †						Good	Good	Adequate
	Irregularly Shaped Polymers *	Good	Excellent				Excellent	Excellent	
	Irregularly Shaped Polymers †						Good	Good	
	Thin, Dark Polymer Films	Good	Excellent				Excellent	Excellent	
	Thick, Dark Polymer Films		Excellent				Excellent	Excellent	
	Layered Polymer Films		Excellent				Adequate	Adequate	
	Thin Polymer Film on Reflective Substrates				Adequate		Good		Excellent
	Thick Polymer Film on Reflective Substrates				Adequate		Good		Excellent
	Organic Powders	Excellent	Excellent	Good		Good	Good		
	Adhesives	Adequate	Good				Excellent	Excellent	
	Rubber		Good		Adequate		Excellent	Excellent	
	Thin Fibers	Adequate	Excellent				Good	Excellent	
	Thick Fibers	Adequate	Excellent				Good	Excellent	
	Surface Analysis	Adequate		Good			Excellent	Excellent	Good
<b>Liquids</b>	Free-Flowing Aqueous Solutions	Excellent					Excellent		
	Other Free-Flowing Liquids	Excellent					Excellent		
	Viscous Liquids	Excellent					Excellent		
<b>Gas</b>	Gases (ppb to 100% concentration)	Excellent							

# Index of Sample Types

**Powders** – organic and inorganic solids that can be ground into a powder (2–5 micron particle size); Examples: chemicals, pharmaceuticals, crystalline materials, pigments, fibers, polymers and powders

**Thermoplastic Polymers** – polymers that can be pressed into free-standing thin films

**Soluble Polymers** – polymers that can be dissolved in a solvent or cast as a thin film

**Thin Polymer Films** – free-standing polymer films that are not thermoplastic or soluble and are less than 50 microns thick

**Thick Polymer Films** – free-standing polymer films that are not thermoplastic or soluble and are more than 50 microns thick

**Regularly Shaped Polymers** – polymers, films, and plaques that are hard or soft with a smooth surface, capable of being ground, not thermoplastic or soluble and regularly shaped

**Regularly Shaped Polymers** – polymers, films, and plaques that are hard or soft with a smooth surface, incapable of being ground, not thermoplastic or soluble and regularly shaped

**Irregularly Shaped Polymers** – polymers that are hard or soft with a rough or uneven surface, capable of being ground, not thermoplastic, or soluble and irregularly shaped; Examples: formed polymers, polymer beads and pellets

**Irregularly Shaped Polymers** – polymers that are hard or soft with a rough or uneven surface, incapable of being ground, not thermoplastic or soluble and irregularly shaped

**Thin, Dark Polymers** – carbon-filled polymers high in inorganic content that are not thermoplastic or soluble and less than 10 microns thick, such as carbon black

**Thick, Dark Polymers** – carbon-filled polymers high in inorganic content that are not thermoplastic or soluble and more than 10 microns thick

**Layered Polymer Films** – polymers that contain two or more layers or thin or thick films; Examples: layered paints and packaging materials

**Thin Polymer Film on Reflective Substrate** – polymer film on any kind of surface that reflects IR energy (usually metal) that is less than 15 microns thick; Examples: lubricants on hard disk media and layers on silicon wafers

**Thick Polymer Film on Reflective Substrate** – polymer film on any kind of surface that reflects IR energy (usually metal) that is more than 15 microns thick; Examples: coatings on containers (such as soda cans)

**Adhesives** – solid adhesives like tapes and solid glues

**Rubbers** – irregular-shaped rubber items that are not thermoplastic or soluble; Examples: o-rings, gaskets, and fittings

**Thin Fibers** – thin and bundled fibers

**Thick Fibers** – thick and bundled fibers

**Surface Analysis** – for qualitative analysis of the outermost layer of any solid or film

**Free-Flowing Aqueous Solutions** – liquids that contain any amount of water; Examples: inks, dyes, solvents, and paints

**Other Free-Flowing Liquids** – liquids that do not contain water

**Viscous liquids** – thick liquids, pastes, and emulsions; Examples: polyols, greases, and heavy oils

**Gases (ppb to 100% concentration)** – any sample that is a gas at room temperature or slightly above room temperature

## Smart Accessories

### *Design elements to consider when choosing an accessory*

We appreciate the need for efficiency and reproducibility in today's lab. We conducted extensive research regarding the features and benefits our users desire in their sampling accessories. This exhaustive process led to the development of unique Thermo Scientific™ Smart Accessories™. These accessories offer the following value-added features:

- **Permanently aligned optics** – needs no optical adjustment, so results are reproducible and quantifiable
- **Rugged design** – protects optics from daily use and dust
- **Automatic and fast purge** – achieves purge up to three times faster than a standard accessory
- **Accessory recognition** – automatically identifies the accessory as soon as it is snapped in place and records its serial number in a non-editable history file
- **Experiment setup** – automatically sets up your experiment parameters so you can start sampling immediately
- **Accessory performance checks** – tests the accessory to ensure that it is performing optimally
- **Spectral quality checks** – examines data as it is collected and rates the quality of the spectra you have collected, offering suggested improvements when necessary
- **Multi-media tutorials and on-line help** – answers any questions you may have while conducting your analysis

We offer a comprehensive line of both Smart Accessories and standard sampling accessories to meet the needs of your laboratory. In addition, we offer a selection of microscope objectives to facilitate in the analysis of small samples.



<b>Africa</b>	+27 11 822 4120
<b>Australia</b>	+61 3 9757 4300
<b>Austria</b>	+43 1 333 50 34 0
<b>Belgium</b>	+32 53 73 42 41
<b>Canada</b>	+1 800 530 8447
<b>China</b>	+86 10 8419 3588
<b>Denmark</b>	+45 70 23 62 60
<b>Europe-Other</b>	+43 1 333 50 34 0
<b>Finland/Norway/Sweden</b>	+46 8 556 468 00
<b>France</b>	+33 1 60 92 48 00
<b>Germany</b>	+49 6103 408 1014
<b>India</b>	+91 22 6742 9434
<b>Italy</b>	+39 02 950 591
<b>Japan</b>	+81 45 453 9100
<b>Latin America</b>	+1 561 688 8700
<b>Middle East</b>	+43 1 333 50 34 0
<b>Netherlands</b>	+31 76 579 55 55
<b>New Zealand</b>	+64 9 980 6700
<b>Russia/CIS</b>	+43 1 333 50 34 0
<b>Spain</b>	+34 914 845 965
<b>Switzerland</b>	+41 61 716 77 00
<b>UK</b>	+44 1442 233555
<b>USA</b>	+1 800 532 4752

[www.thermoscientific.com](http://www.thermoscientific.com)

©2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

BR50557\_E 06/13M

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific