

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Chromatography for Foods and Beverages Adulteration and Authentication Applications Notebook

Analytical Methods for Effective and Efficient Confirmation

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Adulteration and Authentication

Introduction

Guaranteeing the authenticity of a product is integral to quality procedures, essential in the compliance of food labelling regulations and important to consumers when paying for premium products.

Adulteration refers to the addition of substances to foods or substitution of components with cheaper alternatives, usually for economic gain (also called food fraud). Adulterants can lead to serious health issues. Food authenticity and adulteration testing encompasses multiple approaches and techniques that are constantly evolving to meet emerging challenges.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Analytical Technologies



High-Performance Liquid Chromatography

Thermo Scientific™ Vanquish™ UHPLC System and Thermo Scientific™ Dionex™ UltiMate™ 3000 UHPLC+ systems offer excellent chromatographic performance, operational simplicity and unrivaled flexibility. Choose from a wide range of standard and unique specialty detectors to extend your laboratory's analytical capabilities.




Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

The next generation in UHPLC innovations

The Vanquish system takes high-end UHPLC to a new level, offering more resolution while meeting the throughput demands of modern laboratories.

The system delivers better separations, more results and easier interaction, simultaneously, without compromise.



The Vanquish UHPLC System

Delivering the new standard in UHPLC

- More powerful separations with 1500 bar of pump pressure at flow rates up to 5 mL/min
- Industry-leading flow and gradient precision
- Excellent injections up to 100 μ L in 0.01 μ L increments
- Automated workflows with barcode reading for simplified setup and tracking
- Maximum sample capacity with up to 23 well plates, or 8832 samples
- More confident separations with a wide temperature range of 5 $^{\circ}$ C to 120 $^{\circ}$ C for two thermostating modes and active column pre-heating for improved precision
- UV detection with linear response up to 3000 mAu and noise levels as low as 3 μ Au
- Thermo Scientific™ LightPipe™ technology assures lowest peak dispersion with UV detection
- Available Vanquish Charged Aerosol detector for quantification of non-chromophoric compounds



Vanquish Diode Array Detector with LightPipe technology

UHPLC Portfolio

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)


[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)





UltiMate 3000 UHPLC+ Systems

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Best-in-class HPLC systems for all your chromatography needs

UltiMate 3000 UHPLC+ Systems provide excellent chromatographic performance while maintaining easy, reliable operation. The basic and standard analytical systems offer ultra HPLC (UHPLC) compatibility across all modules, ensuring maximum performance for all users and all laboratories.

Covering flow rates from 20 nL/min to 10 mL/min with an industry-leading range of pumping, sampling, and detection modules, UltiMate 3000 UHPLC+ Systems provide solutions from nano to semipreparative, from conventional LC to UHPLC.

Superior chromatographic performance

- UHPLC design philosophy throughout nano, standard analytical, and rapid separation liquid chromatography (RSLC)
- 620 bar (9,000 psi) and 100 Hz data rate set a new benchmark for basic and standard analytical systems
- RSLC systems go up to 1000 bar and data rates up to 200 Hz
- ×2 Dual System for increased productivity solutions in routine analysis
- Fully UHPLC compatible advanced chromatographic techniques
- Thermo Scientific™ Dionex™ Viper™ and nanoViper™ fingertight fittings—the first truly universal, fingertight fitting system even at UHPLC pressures

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



UltiMate 3000 UHPLC+ Systems

We are uniquely focused on making UHPLC technology available to all users, all laboratories, and for all analytes.



Rapid Separation LC Systems

The extended flowpressure footprint of the RSLC system provides the performance for ultrafast high-resolution and conventional LC applications.



RSLCnano Systems

The Rapid Separation nano LC System (RSLCnano) provides the power for high resolution and fast chromatography in nano, capillary, and micro LC.



Standard LC Systems

Choose from a wide variety of standard LC systems for demanding LC applications at nano, capillary, micro, analytical, and semipreparative flow rates.



Basic LC Systems

UltiMate 3000 Basic LC Systems are UHPLC compatible and provide reliable, high performance solutions to fit your bench space and your budget.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice


Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Charged Aerosol Detection

Charged Aerosol Detection provides near universal detection independent of chemical structure for non- or semi-volatile analytes with HPLC and UHPLC. Thermo Scientific™ Dionex™ Corona™ Veo™ and Vanquish Charged Aerosol detectors are ideally suited as a primary detector for any laboratory, while providing complementary data to UV or MS methods. No other LC detector available today can match the performance of a Corona Veo detector.

- High sensitivity – single-digit nanogram on column
- Consistent response – independent of chemical structure
- Wide dynamic range – to four orders of magnitude or greater
- Simple to use – easy to integrate with any HPLC/UHPLC system

Charged aerosol detectors give the simplicity, reproducibility and performance required for a full range of applications from basic research to manufacturing QC/QA. With charged aerosol detection you get predictable responses to measure analytes in direct proportion to their relative amounts for quantitation without actual standards.

This detector offers the flexibility to use reversed-phase gradients, as well as normal phase and HILIC modes of separation on any LC system. And, in many cases eliminates the need for derivatization or sample pre-treatment to provide real dilute-and-shoot simplicity.

Advanced Detection Capabilities



Corona Veo Charged Aerosol Detector



Vanquish system with Charged Aerosol Detector

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice


Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Advanced Detection Capabilities

CoulArray Multi-electrode Array Detector

The Thermo Scientific™ Dionex™ CoulArray™ Multi-electrode Array detector is the only practical multi-channel electrochemical detection system that allows you to measure multiple analytes simultaneously, including those that are chromatographically unresolved. The CoulArray detector delivers the widest dynamic range of any available electrochemical detector with unmatched selectivity for detection of trace components in complex matrixes, even when used with aggressive gradients.

- Measures analytes from femtomole to micromole levels
- Greatly simplify sample preparation and eliminate interferences
- Simultaneously analyze multiple analytes in very complex samples
- Easily produce qualitative information for compound identification

Multiple system configurations offer 4, 8, 12, or 16 channels that can be upgraded anytime. The unique data acquisition and processing software uses automatic signal ranging and a unique patented baseline correction algorithms to provide identification and quantitation of single or multiple analytes and powerful 3D data for quick sample fingerprint confirmation with integration to pattern recognition platforms.

With the power of coulometric array technology, the CoulArray detector can give you the qualitative data of a optical PDA with 1,000 fold greater sensitivity to profile the characteristic qualities of products, determine integrity, identify adulteration and even evaluate competitors' products.



CoulArray Multi-electrode Array Detector

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)


[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Advanced Detection Capabilities

RefractoMax 521 Refractive Index Detector


The Thermo Scientific RefractoMax 521 Refractive Index Detector from ERC Inc. This detector, in combination with the UltiMate 3000 system, is the right choice for the isocratic analysis of sugars, polymers, and fatty acids. It features fast baseline stabilization and excellent reproducibility, combined with high sensitivity. The RefractoMax 521 is fully controlled by Thermo Scientific™ Dionex™ Chromeleon™ Chromatography Data System Software (CDS), and can also operate in stand-alone mode.

- The detector is highly sensitive and applicable universally. It provides very stable baselines with a drift of 0.2 μ RIU/h and a noise specification of 2.5 nRIU or less
- The optical bench, thermostatically regulated from 30 °C to 55 °C, and the superior signal-to-noise ratio ensure highly precise measurement results

- The extended flow rate range from 1 mL/min up to 10 mL/min and the operating range of 1.00 to 1.75 RIU enable the use of this detector for a wide range of applications
- Applications include the analysis of all compounds with low UV-Vis activity, such as alcohols, mono- and polysaccharides, esters, fatty acids, or polymers
- An Auto Set-up function automates purging, equilibration, autozero, and the control baseline stability and noise
- Operation with Chromeleon CDS makes the detector easy to use and ensures maximum productivity in instrument control, data processing, and reporting of results



RefractoMax 521 Refractive Index Detector



Advanced Detection Capabilities

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

UltiMate 3000 Diode Array and Multiple-Wavelength Detectors

The Thermo Scientific Dionex UltiMate DAD 3000 detector is a high-resolution, 1024-element diode array detector (DAD) available in Rapid Separation (200 Hz) and Standard (100 Hz) versions. It operates with Chromeleon CDS software to provide a variety of spectra views, including 3-D plotting and automated chromatogram handling. The high resolution and low-noise performance of the DAD-3000 family makes it ideal for the most sensitive and accurate library searches and peak purity analyses.

The detector is also available as a multiple wavelength detector (MWD) in Standard (100 Hz) and Rapid Separation (200 Hz) versions.

- Data collection at up to 200 Hz using a maximum of eight single-wavelength data channels and one 3-D field (3-D only with DAD-3000 (RS)) for best support of ultrafast separations
- Standard versions operate at up to 100 Hz data collection rate for optimum support of 62 MPa (9000 psi) UltiMate 3000 Standard systems
- Accurate compound confirmation with a 1024-element, high resolution photodiode array
- Flexibility in both UV and Vis applications with 190–800 nm wavelength range
- Low-noise over the full spectral range using deuterium and tungsten lamps
- Fast and accurate wavelength verification using a built-in holmium oxide filter

- The detector can be upgraded with the UltiMate PCM 3000 for accurate monitoring pH gradients
- Excellent reliability and reproducibility with low baseline drift (typically < 500 μ AU/h)
- Simplified routine maintenance with front access to pre-aligned cells and lamps
- ID chips on flow cells and lamps for identification and life-span monitoring
- Chromeleon CDS software for full control and flexible data handling
- Front-panel display for easy monitoring of detector status to maximize uptime
- Flow cells for semi-micro, semi-analytical, analytical, and semi-preparative applications
- Flow cells available in stainless steel and biocompatible versions



UltiMate 3000 DAD-3000 Diode Array Detector

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Advanced Detection Capabilities

UltiMate 3000 Electrochemical Detector

Electrochemical detection delivers high sensitivity for neurotransmitter analysis, simplicity and robustness for pharmaceutical or clinical diagnostics, and the selectivity for the characterization of complex samples such as natural products, biological tissues and fluids. For today's researcher, there is a continuing need for detecting vanishingly small quantities of analyte and often in complex samples. Because electrochemical detection measures only compounds that can undergo oxidation or reduction it is both highly sensitive and very selective.


The Thermo Scientific Dionex UltiMate 3000 Electrochemical Detector, designed by the pioneers of coulometric electrochemical detection, delivers state-of-the-art sensor technologies complete with an entire range of high performance and ultra-high performance LC systems optimized for electrochemical detection. The UltiMate 3000 ECD-3000RS takes electrochemical detection to the next level with UHPLC compatibility, total system integration, and selection of detection mode, all with unprecedented operational simplicity.

Features include:

- Detection Modes – choose from DC and PAD for optimum analyte response
- Choice of sensors – both coulometric and amperometric sensors to meet the demands of any application
- UHPLC compatibility – ultralow peak dispersion and high data acquisition rates for conventional or fast, high resolution chromatography
- Modularity – easily expandable to multiple independent sensors for unrivaled flexibility
- Autoranging – simultaneously measure both low and high levels of analytes without losing data
- SmartChip™ technology – easy operation with automatic sensor recognition, event logging and electrode protection



UltiMate 3000 Electrochemical Detector



Advanced Detection Capabilities

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

UltiMate 3000 Fluorescence Detector

The Thermo Scientific Dionex UltiMate 3000 FLD-3000 is a high-sensitivity fluorescence detector series for UltiMate 3000 HPLC systems. It is available in Rapid Separation (RS) and Standard (SD) versions. The optics of the FLD-3000 series provide maximum stray-light suppression for best detection sensitivity. Operated with the Chromeleon CDS software, the detector provides automated qualification, various tools for method development, and instrument wellness monitoring for ease of use, maximum uptime, and the highest degree of regulatory compliance.

- Data collection at up to 200 Hz for optimal support of even the fastest UHPLC separations (FLD-3400RS)
- Standard detectors operate at up to 100 Hz data rate for optimum support of 62 MPa (9,000 psi) UltiMate 3000 standard systems
- Lowest limits of detection with a Raman signal-to-noise ratio (S/N): > 550 ASTM (> 2100 using dark signal as noise reference)

- Unsurpassed reproducibility with active flow cell temperature control for stable fluorophore activity independent of changes in ambient temperature
- Long-life xenon flash lamp for highest sensitivity and long-term operation without the need for frequent lamp changing
- Optional second photomultiplier (PMT) for unique Dual-PMT operation, offering an extended wavelength range up to 900 nm without sacrificing sensitivity in the standard wavelength range
- Two-dimensional (2D) or three dimensional (3D) excitation, emission, or synchro scans to provide the highest degree of flexibility for method development or routine sample characterization
- Innovative Variable Emission Filter for real-time compound-related sensitivity optimization (FLD-3400RS only)
- Large front-panel display for easy monitoring of the detector status
- Two flow-cell sizes for easy optimization to application requirements: the 8 μ L flow cell is ideal for trace analysis, and the 2 μ L flow cell offers best peak resolution with narrow-bore HPLC and UHPLC columns



Ultimate 3000 Fluorescence Detector

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



UltiMate 3000 Variable Wavelength Detectors

The Thermo Scientific Dionex UltiMate 3000 VWD-3000 is a variable wavelength detector (VWD) series for industry leading UV-Vis detection. The forward optics design and wide range of available flow cells ensure optimal performance over a flow rate range of five orders of magnitude. Automated qualification, performance optimization, and instrument wellness monitoring deliver maximum uptime, simplify work-flow, and give you full confidence in your analytical results. The detector is available in a standard 100 Hz (VWD-3100) and a 200 Hz Rapid Separation version (VWD-3400RS) for the most challenging UHPLC applications.

High-Performance UV-Vis Detection

- The VWD-3400RS variant provides data collection rates of up to 200 Hz for optimal support of today's and tomorrow's UHPLC separations
- The VWD-3100 standard detector operates at up to 100 Hz data rate for optimum support of 62 MPa (9000 psi) UltiMate 3000 Standard systems
- Superior detection of trace analytes with low noise ($< -2.0 \mu\text{AU}$) and drift ($< 100 \mu\text{AU/h}$)
- The detector's large linearity range of up to 2.5 AU is ideal for applications with widely varying analyte concentrations
- Up to four absorption channels (VWD-3400RS) and spectral scans support effective method development
- Active temperature control of optics and electronics for data acquisition independent of ambient conditions

Advanced Detection Capabilities

- Front panel access for quick and easy lamps and flow cells changes
- Automated qualification monitoring for full regulatory compliance
- Large front panel display for monitoring the detector status even from a distance
- Maximize uptime using predictive performance-based on monitoring the life cycle of detector lamps
- The detector can be upgraded with the Thermo Scientific Dionex pH/Conductivity Monitor (PCM-3000) for accurate and precise pH- and conductivity monitoring
- Unique 45 nL ultra-low dispersion UV monitor for dispersion-free UV detection in LC/MS



UltiMate 3000 VWD-3400 Variable Wavelength Detector.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Analytical Technologies



Ion Chromatography

Thermo Scientific Dionex IC systems have led the analytical instrument industry for over 30 years with solutions that represent state-of-the-art technological advancements and patented technologies.

IC and RFIC Systems

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)


[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Innovative Ion Chromatography Solutions

Our High-Pressure™ Ion Chromatography (HPIC™) systems include the Thermo Scientific Dionex ICS-5000+ HPIC system, which is optimized for flexibility, modularity, and ease-of-use, combining the highest chromatographic resolution with convenience. In addition, the Thermo Scientific Dionex ICS-4000 Capillary HPIC system is the world's first commercially available dedicated capillary high-pressure Reagent-Free™ (RFIC™) IC system. The Dionex ICS-4000 system is always ready for the next analysis, delivering high-pressure IC on demand.

Reagent-Free IC systems eliminate daily tasks of eluent and regenerant preparation in turn saving time, preventing errors, and increasing convenience. RFIC-EG systems use electrolytic technologies to generate eluent on demand from deionized water, and to suppress the eluent back to

pure water to deliver unmatched sensitivity. RFIC-ER systems are designed to use carbonate, carbonate/ bicarbonate, or MSA eluents for isocratic separations.

At the heart of our ion chromatography portfolio is a unique set of column chemistries that provide high selectivities and efficiencies with excellent peak shape and resolution. Thermo Scientific™ Dionex™ IonPac™ chromatography columns address a variety of chromatographic separation modes including ion exchange, ion exclusion, reversed-phase ion pairing, and ion suppression. Our column chemistries are designed to solve specific applications, and we offer a variety of selectivities and capacities for simple and complex samples. Additionally, our Dionex IonPac column line is available in standard bore, microbore and capillary formats for the ultimate application flexibility.



Thermo Scientific Dionex IC instrument family

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements


Tryptophan Supplements

Olive Oil

References




Analytical Technologies



Mass Spectrometry

We provide advanced integrated IC/MS and LC/MS solutions with superior ease-of-use and modest price and space requirements. UltiMate 3000 System Wellness technology and automatic MS calibration allow continuous operation with minimal maintenance. The Dionex ion chromatography family automatically removes mobile phase ions for effort-free transition to MS detection.



Mass Spectrometry Instruments

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Single-Point Control and Automation

We provide advanced integrated IC/MS and LC/MS solutions with superior ease-of-use and modest price and space requirements. UltiMate 3000 System Wellness technology and automatic MS calibration allow continuous operation with minimal maintenance. The Dionex ion chromatography family automatically remove mobile phase ions for effort-free transition to MS detection.

- Thermo Scientific™ MSQ Plus™ mass spectrometer, the smallest and most sensitive single quadrupole on the market for LC and IC
- Self-cleaning ion source for low maintenance operation

- Chromeleon CDS software for single-point method setup, instrument control, and data management compatible with existing IC and LC methods
- The complete system includes the MSQ Plus mass spectrometer, PC data system, electrospray ionization (ESI) and atmospheric pressure chemical ionization (APCI) probe inlets, and vacuum system

Now, you no longer need two software packages to operate your LC/MS system. Chromeleon CDS software provides single-software method setup and instrument control; powerful UV, conductivity, and MS data analysis; and fully integrated reporting.



MSQ Plus Mass Spectrometer

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Analytical Technologies



Chromatography Data Systems

Tackle chromatography management challenges with the world's most complete chromatography software. Whether your needs are simple or complex or your scope is a single instrument, a global enterprise, or anything in between – the combination of Chromeleon CDS' scalable architecture and unparalleled ease-of use, makes your job easy and enjoyable with one Chromatography Data System for the entire lab.




Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

The Fastest Way from Samples to Results

The 7.2 release of Chromeleon Chromatography Data System software is the first CDS that combines separation (GC/IC/LC) and Mass Spectrometry (MS) in an enterprise (client/server) environment. By extending Chromeleon 7.2 CDS beyond chromatography into MS, lab technicians can now streamline their chromatography and MS quantitation workflows with a single software package. MS support in Chromeleon 7.2 CDS is focused on routine and quantitative workflows, which provides access to rich quantitative data processing and automation capabilities — ultimately boosting your overall lab productivity and increasing the quality of your analytical results.



Chromeleon CDS Software

- Enjoy a modern, intuitive user interface designed around the principle of operational simplicity
- Streamline laboratory processes and eliminate errors with eWorkflows™, which enable anyone to perform a complete analysis perfectly with just a few clicks
- Access your instruments, data, and eWorkflows instantly in the Chromeleon Console
- Locate and collate results quickly and easily using powerful built-in database query features
- Interpret multiple chromatograms at a glance using MiniPlots
- Find everything you need to view, analyze, and report data in the Chromatography Studio
- Accelerate analyses and learn more from your data through dynamic, interactive displays
- Deliver customized reports using the built-in Excel® compatible spreadsheet

Excel is a registered trademark of Microsoft Corporation.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Analytical Technologies



Process Analytical Systems

Thermo Scientific Dionex process analytical systems provide timely results by moving chromatography-based measurements on-line.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Process Analytical Systems and Software

Improved Process Monitoring with On-line Chromatography IC and LC Systems

Information from the Thermo Scientific Dionex Integral process analyzer can help reduce process variability, improve efficiency, and reduce downtime. These systems provide comprehensive, precise, accurate information faster than is possible with laboratory-based results. From the lab to the factory floor, your plant's performance will benefit from the information provided by on-line LC.

- Characterize your samples completely with multicomponent analysis
- Reduce sample collection time and resources with automated multipoint sampling
- Improve your process control with more timely results
- See more analytes with unique detection capabilities
- The Thermo Scientific Integral Migration Path approach lets you choose the systems that best meets your needs



Integral process analyzer

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Analytical Technologies



Automated Sample Preparation

Solvent extractions that normally require labor-intensive steps are automated or performed in minutes, with reduced solvent consumption and reduced sample handling using the Thermo Scientific™ Dionex™ ASE™ Accelerated Solvent Extractor system or Thermo Scientific™ Dionex™ AutoTrace™ 280 Solid-Phase Extraction instrument.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Accelerated Solvent Extractor System

Complete Extractions in Less Time Using Less Solvent

Thermo Scientific Dionex ASE systems extract solid and semisolid samples using common solvents at elevated temperature and pressure. The Dionex ASE 150 and 350 systems feature pH-hardened pathways with Dionium™ components to support extraction of acidic or alkaline matrices, and combine pretreatment, solvent extraction, and cleanup into one step. Dionium is zirconium that has undergone a proprietary

hardening process that makes it inert to chemical attack by acids and bases at elevated temperatures.

Dionex ASE systems are dramatically faster than Soxhlet, sonication, and other extraction methods, and require significantly less solvent and labor. Accelerated solvent extraction methods are accepted and established in the environmental, pharmaceutical, foods, polymers and consumer product industries. Accelerated solvent extraction methods are accepted and used by government agencies worldwide.



Dionex ASE 150/350 and Dionex AutoTrace 280 SPE instruments

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Adulteration and Authentication

Part 1 Targeted Methods

The measurement of identified analytes as indicators of possible adulteration and to verify product authenticity.



Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Targeted Approaches

Adulteration and Authentication



Coffee Adulteration

The less costly Robusta varietal is commonly used in decaffeinated coffee and as an adulterant in the preferred and more costly Arabica varietal.



Targeted Approaches

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Coffee Adulteration

Organic acid profiles of coffee are used to determine Robusta content in Arabica coffee. Additionally, organic acids are of interest because they are primarily responsible for the coffee's acidity. The chromatograms shown here compare the differences in decaffeinated and caffeinated

brewed coffee samples with the characteristic higher quinate (quinic acid) concentration present in the caffeinated coffee. Lower quinate and lower malate concentrations than normal would be indicative of adulteration by a Robusta varietal coffee.

Column: Thermo Scientific™ Dionex™ IonSwift™ MAX-100G, MAX-100, 0.25 × 250 mm
 Flow: 12 µL/min
 Column Temp.: 30 °C
 Injection Volume: 0.4 µL
 Eluent Source: Dionex EGC-KOH capillary cartridge
 Gradient: 0.1 mM KOH (-10–4 min), 0.1–2 mM (4–6 min), 2–15 mM (6–12 min), 15–35 mM (12–16 min), 65 mM (17–30 min)
 Detection: Suppressed conductivity, Dionex ACES 300 suppressor, recycle mode
 Sample Prep.: 1:50 dilution




Figure 7-1. Determination of anions and organic acids in a brewed decaffeinated coffee sample more typical of the lower cost Robusta coffee.

Column: Dionex IonSwift MAX-100G, MAX-100, 0.25 × 250 mm
 Flow: 12 µL/min
 Column Temp.: 30 °C
 Injecton Volume: 0.4 µL
 Eluent Source: Dionex EGC-KOH capillary cartridge
 Gradient: 0.1 mM KOH (-10–4 min), 0.1–2 mM (4–6 min), 2–15 mM (6–12 min), 15–35 mM (12–16 min), 65 mM (17–30 min)
 Detection: Suppressed conductivity, Dionex ACES 300 suppressor, recycle mode
 Sample Prep.: 1:50 dilution




Figure 7-2. Determination of anions and organic acids in a brewed caffeinated coffee sample typical of the more costly Arabica coffee.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Targeted Approaches

Adulteration and Authentication



Controlled Substances

The identification of controlled substances in food samples is challenging due to the wide range of drugs available and the complex nature and diversity of the sample being analyzed.



Controlled Substances: Cannabinoids

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Marijuana is the most common illegal drug in the United States, and each year U.S. law enforcement agencies seize more than two million pounds of marijuana in various forms. Seized evidence submitted to forensic laboratories is screened for marijuana by microscopic inspection and simple chemical tests such as the Duquenois-Levine test. Presumptive positive results are confirmed by using gas chromatography-mass spectrometry (GC/MS) to positively identify cannabinoids including Δ^9 -tetrahydrocannabinol (THC, the main psychoactive component), cannabinol (the main degradation product of THC) and cannabidiol.

This traditional approach works fairly well for leaf marijuana, hashish, hash oil, and residue collected from smoking paraphernalia.

An alternative method to positively identify marijuana cannabinoids in complex food matrices is to use UHPLC/MS. UHPLC/MS offers a threefold benefit compared to GC/MS; simpler sample preparation, no derivatization, and less instrument clean up time. Application Note 433 demonstrates how a working forensic laboratory uses UHPLC/MS to analyze baked goods for three cannabinoids of forensic importance.

Column:	hermo Scientific™ Hypersil GOLD™ PFP (perfluorinated phenyl) 1.9 μ m, 100 \times 2.1 mm			
Flow:	1 mL/min			
Column Temp:	45 °C			
Injection Volume:	2 μ L partial loop injection, 25 μ L loop size			
Mobile Phase:	A: Water with 0.06 % acetic acid B: Acetonitrile (MeCN) with 0.06% acetic acid C: Methanol with 0.06% acetic acid			
Gradient:	T (min)	A%	B%	C%
	0.00	95.0	0.0	5.0
	1.00	60.0	32.5	7.5
	2.00	50.0	40.0	10.0
	5.00	45.0	45.0	10.0
	6.00	25.0	60.0	15.0
	6.50	5.0	0.0	95.0
	7.50	5.0	0.0	95.0
	7.51	95.0	0.0	5.0
	8.00	95.0	0.0	5.0
Syringe Speed:	8 μ L/sec			
Flush Speed:	1 100 μ L/sec			
Flush Volume:	400 μ L			
Wash Volume:	100 μ L			
Flush/Wash Source:	Bottle with methanol			




Figure 7-3. Extracted ion chromatograms (m/z 310.5–311.5 and 314.5–315.5) of cannabinoid standards: (A) extracts from brownie, (B) extracts from cookie, and a standard (C) by UHPLC/MS.

Controlled Substances: Cannabinoids



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

MS analysis was carried out on a MSQ Plus single quadrupole LC/MS detector with Thermo Scientific™ Xcalibur™ 2.05 software.




Figure 7-4. MS spectra of cannabinoid standards: cannabidiol (A), THC (B), cannabinol (C), eluted at 4.1 min, 5.1 min, and 5.4 min respectively, and extracts from brownie (D) and cookie (E), eluted at 5.1 min.



Targeted Approaches

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Controlled Substances: LSD

Lysergic acid diethylamide (LSD) is a controlled substance in forensic chemistry that is notorious for being difficult to identify. Its myriad evidentiary forms include paper tabs, eye drops, sugar cubes and small sugary candies such as sweet tarts, valentine hearts, or mints. Because it is such a potent hallucinogen, typical street doses require only 40 to 120 µg of LSD. The small personal-use amounts seized by state and local law enforcement often lack sufficient drug to allow both forensic analysis by traditional means and archiving of some of the evidence for follow-up testing. Most forensic laboratories confirm the presence of LSD by using gas chromatography with mass spectrometry (GC/MS).

An alternative method to positively identify LSD in complex food matrices is to use ultra high performance liquid chromatography with mass spectrometric detection (UHPLC/MS). UHPLC/MS offers a threefold benefit compared to GC/MS; simpler sample preparation, no derivatization, and less time wasted baking out or cleaning the instrument. Application Note 432 demonstrates how a working forensic laboratory uses UHPLC/MS to analyze sugar candies for LSD.





Targeted Approaches

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Column: Hypersil GOLD PFP (perfluorinated phenyl)
 1.9 μm , 100 x 2.1 mm
 Flow: 1 mL/min
 Column Temp.: 45 $^{\circ}\text{C}$
 Injection Volume: 2 μL partial loop injection, 25 μL loop size
 Mobile Phase: A: Water with 0.06 % acetic acid
 B: Acetonitrile (MeCN) with 0.06% acetic acid
 C: Methanol with 0.06% acetic acid

Syringe Speed: 8 $\mu\text{L}/\text{sec}$
 Flush Speed: 100 $\mu\text{L}/\text{sec}$
 Flush Volume: 400 μL
 Wash Volume: 100 μL
 Flush/Wash Source: Bottle with methanol

Gradient:

T (min)	A%	B%	C%
0.00	95.0	0.0	5.0
1.00	95.0	0.0	5.0
1.50	90.0	5.0	5.0
2.70	70.0	10.0	20.0
3.00	5.0	15.0	80.0
7.00	5.0	0.0	95.0
7.10	95.0	0.0	5.0
8.00	95.0	0.0	5.0




Figure 7-5. Extracted ion chromatograms at $m/z = 324 \pm 0.5$ amu obtained by UHPLC/MS on a Hypersil GOLD PFP column: (A) LSD standard; (B) methanol extract of LSD-doped candy heart; (C) methanol extract of LSD-doped sugar cube.

Controlled Substances: LSD

Mass Spectrometer Conditions

Ionization: Electrospray (ESI)
 Polarity: Positive
 Probe Temperature: 500 $^{\circ}\text{C}$
 Cone Voltage: 90 V
 Scan Mode: Full scan with mass range of $m/z = 125-425$ amu
 ESI Voltage: 3.5 kV
 Scan Time: 0.2 s




Figure 7-6. MS spectra of LSD obtained by UHPLC/MS on a Hypersil GOLD PFP column: (A) LSD standard; (B) methanol extract of LSD-doped candy heart; (C) methanol extract of LSD-doped sugar cube.



Controlled Substances: Psychotropics

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

The identification of controlled substances in food matrices is a challenge for forensic laboratories. Classical techniques, such as color tests and thin-layer chromatography, do not provide molecular structural information and cannot be used as a principal means of identification. Infrared spectroscopy (IR) and GC/MS can identify controlled substances, but suffer from several shortcomings. To ameliorate these issues, extensive wet chemistry preparation methods have been developed; unfortunately, they are time-consuming and often require greater amounts of the controlled substances than are present in the evidence.

Given the increasing emphasis on instrumental methods of analysis and the limitations of the traditional instruments in the field of forensics, another solution is necessary. LC/MS holds several advantages over the traditional methods of analysis. For example, psilocybin does not decompose at the lower temperatures used in HPLC. The low concentrations of psilocybin and psilocin are not an issue due to increased sensitivity of this technique. Many matrix components that interfere in GC/MS methods do not interfere in HPLC methods because of greater differences in analyte solubility as compared to analyte volatility.

Application Note 428 describes a method that can be used to separate and confirm psilocybin and psilocin in mushroom and chocolate samples using UHPLC/MS.




Figure 7-7. Methanol extract of chocolates containing psilocybin and psilocin.

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Targeted Approaches

Adulteration and Authentication



Fruit Juice Adulteration

Fruit juice adulteration presents an economic and regulatory problem. The United States orange juice industry estimates that orange juice sales gross more than one billion dollars annually. The most common forms of adulteration include simple dilution and blending of inexpensive and synthetically produced juices into the more expensive ones.

[Learn more about US Fruit Juice Adulteration Regulations](#)

[Learn more about the History of Fruit and Vegetable Juice US Regulations](#)



Targeted Approaches

Fruit Juice Adulteration: Organic Acids

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

The concentration of organic acids can be used to assess whether an expensive juice has been illegally adulterated with a cheaper juice. Because organic acid profiles are distinct to each type of fruit juice, evidence of tampering can be evaluated by comparing the known juice fingerprint to that of the suspected adulterated juice.

As shown in these examples, the major organic acids in grape juice are malate, tartrate, and citrate, whereas in orange juice the primary analyte is citric acid. Both lactic and acetic acids are typically low in these juice samples – elevated levels of lactic and acetic acids may be caused by microbiological spoilage, so it is important to monitor the concentrations of these organic acids as a measure of product quality.

Malate is the major constituent of apple juice, whereas in cranberry juice cocktail, there are high concentrations of succinate and quinate, which provide cranberry's tart taste.

The presence of galacturonate in cranberry juice cocktail and grape juice can be attributed to the degradation of pectins in the skins of fruit. Freshly squeezed juices generally show lower levels of galacturonate.

Peaks:

1. Quinate	55 ppm	9. Malate	116.0
2. Fluoride	1.3	10. Malonate/Tartrate	190.0
3. Lactate	46.2	11. Maleate	4.0
4. Galacturonate	60	12. Sulfate	22.1
5. Chloride	1.2	13. Oxalate	19.4
6. Nitrate	0.6	14. Phosphate	27.0
7. Glutarate	0.7	15. Citrate	80.0
8. Succinate	1.2	16. Isocitrate	1.8




Figure 7-8. Determination of anions and organic acids in grape juice.





Targeted Approaches

Fruit Juice Adulteration: Organic Acids

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peaks:			
1. Quinate	4.0 ppm	9. Succinate/Malate	161.0
2. Lactate	9.0	10. Malonate	3.0
3. Acetate	0.5	11. Sulfate	4.4
4. Propionate	1.9	12. Oxalate	13.8
5. Formate	3.0	13. Phosphate	23.0
6. Chloride	2.0	14. Citrate	400
7. Nitrate	0.1	15. Isocitrate	4.0
8. Glutarate	13.0		




Figure 7-9. Determination of anions and organic acids in orange juice.

Peaks:			
1. Quinate	34.5 mg/L	10. Malate	250.1
2. Lactate	14.8	11. Malonate/Tartrate	0.8
3. Glycolate	0.6	12. Maleate	1.6
4. Formate	0.5	13. Sulfate	1.7
5. Pyruvate	1.2	14. Oxalate	0.5
6. Galacturonate	1.5	15. Phosphate	12.6
7. Chloride	0.4	16. Citrate	6.6
8. Nitrate	0.8	17. Isocitrate	0.4
9. Succinate	9.0	18. <i>Cis</i> -aconitate	0.4




Figure 7-10. Determination of anions and organic acids in apple juice.



Fruit Juice Adulteration: Organic Acids

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peaks:

1. Quinate	210	mg/L	11. Succinate	257
2. Fluoride	<0.1		12. Unknown	-
3. Lactate/Acetate	10		13. Sulfate	10.3
4. Glycolate	2.6		14. Oxalate	14.8
5. Formate	3.7		15. Phosphate	1.8
6. Pyruvate	2.1		16. Unknown	-
7. Unknown	-		17. Citrate	163
8. Galacturonate	16.9		18. Isocitrate	1.0
9. Chloride	2.3		19. <i>Trans</i> -aconitate	2.7
10. Nitrate	<0.1		20. Unknown	-




Figure 7-11. Determination of anions and organic acids in cranberry juice cocktail.



Did You Know?

There are three types of adulteration:

- Replacement (complete or partial replacement of a food ingredient with a less expensive substitute). Examples include addition of melamine to milk, and water/citric acid to lemon juice. This form of adulteration also includes false declaration of geographical, species, botanical or varietal origin.
- Addition (addition of small amounts of non-authentic substance to mask inferior quality ingredient). E.g., addition of Sudan red dyes to enhance the color of poor quality paprika.
- Removal (Removal or intentional omission of an authentic and valuable constituent). An example is filtering of poor quality honey to remove pollen or other residue in order to make it difficult to trace and identify the source of the honey.



Targeted Approaches

Fruit Juice Adulteration: Oligosaccharides

Oligosaccharides are routinely determined in food and beverage products for a variety of purposes, including quality control, verifying food-labeling claims, establishing product authenticity, and monitoring fermentation processes. Oligosaccharide profiles obtained using HPAE-PAD with the Thermo Scientific™ Dionex™ CarboPac™ PA100 column can be used to establish the “fingerprint” of food samples. Suspect samples can be analyzed and compared to the known profiles. This technique is useful in detecting adulteration and in quality control. For example, oligosaccharide profiles are used to detect adulteration of natural fruit juices, establish the geographic origin of molasses, and analyze for polysaccharides in hydrolyzed glucose syrup.



Column: Dionex CarboPac PA-100
 Flow: 1 mL/min
 Injection Volume: 25 μ L
 Eluent: Sodium hydroxide/Sodium acetate gradient
 Detection: Pulsed amperometry, PAD, gold electrode




Figure 7-12. Oligosaccharide profiles of pure and adulterated orange juice. Oligosaccharide composition profiles can be used to detect the adulteration of natural fruit juices by inexpensive sweeteners such as beet sugar.

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Targeted Approaches

Fruit Juice Adulteration: Profiling

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

One adulterant currently in use is partially inverted sucrose, wherein about one-half of the sucrose has been hydrolyzed to glucose and fructose. This ratio of approximately 1:2 (glucose: fructose: sucrose) closely matches the ratio found in orange juice. When cane sugar is the source of inverted sucrose, Stable Isotope Ratio Analysis (SIRA) can be used to identify adulterated juices because the ratio of ^{13}C to ^{12}C is different for sugars in orange juice and cane sugar. Beets, on the other hand, produce sugar via a metabolic pathway different from cane and similar to many fruits,

so that the ratio of ^{13}C to ^{12}C is about the same for sugars in orange juice and beet sugar. This fact renders SIRA inadequate for detecting adulteration by beet sugar.

Recently, investigators using high performance anion exchange chromatography with pulsed amperometric detection (HPAE-PAD) have discovered several components in beet medium invert sugar (BMIS) that are not present in orange juice – such as late eluting components as well as raffinose as markers of orange juice adulteration.




Figure 7-13. Chromatograms showing (A) unadulterated orange juice and (B) orange juice adulterated with medium invert sugar. Note the late-eluting fingerprint between 50 and 60 minutes.



Targeted Approaches

Fruit Juice Adulteration: Profiling

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)




Figure 7-14. Chromatograms showing (A) unadulterated orange juice (note the lack of any peaks eluting at 20 minutes) and (B) orange juice adulterated to 12% with medium invert beet sugar. Adulteration can be detected by the presence of raffinose at 20 minutes.



Targeted Approaches

Fruit Juice Adulteration: Profiling

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Pomegranate is extensively cultivated worldwide and has become a high-value crop for juice production. The retail market now contains numerous pomegranate-related products such as juices, smoothies, flavored waters, and sports and energy drinks. Due to the high demand for pomegranates outstripping the supply, adulteration of pomegranate juice (PJ) has become widespread. Manufacturers have attempted to extend the limited supply of PJ by blending with filler ingredients such as cane sugar, corn syrup sweeteners, and lower-quality juices containing sorbitol, malic acid, and sucrose (e.g., grape, apple, and blackberry). To establish an authentication criterion, an International Multidimensional Authenticity Specifications algorithm was developed based on the analysis of commercial juice samples from 23 manufacturers in the United States, Iran, Turkey, Azerbaijan, Syria, India, and China. There is universal agreement that the anthocyanin profile in PJ consists of a constant group of six anthocyanins, regardless of the origin. However, the anthocyanin concentrations can vary depending on the geographic source of the PJ. The anthocyanin profile is one of several chemical analyses that are required to determine the authenticity of PJ. Additional chemical profiling methods include measuring other polyphenols (i.e., ellagitannins), monosaccharides (e.g., fructose and glucose), organic acids, amino acids, and potassium in PJ samples. Determinations of monosaccharides, organic acids, and punicalagins in fruit juices have been previously described in AN 82, 143, and CAN 106, respectively




Figure 7-15. Typical profile of anthocyanins in pomegranate wildberry juice.

Download Application Note 82: Analysis of Fruit Juice Adulterated with Medium Invert Sugar from Beets

Download Application Note 143: Determination of Organic Acids in Fruit Juices

Download Customer Application Note 106: Determination of the Punicalagins Found in Pomegranate by High Performance Liquid Chromatography



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Grape juice contains four of the six anthocyanins present in PJ but at much lower concentrations. Simulated adulterated PJ was prepared by combining PJ and grape juice, then diluting 1:5 in mobile phase A prior to analysis. The adulterated juice shows all of the signature anthocyanins and several other late-eluting peaks not characteristic of PJ. The anthocyanin content of the adulterated juice is also lower than that of PJ, as expected.



Fruit Juice Adulteration: Profiling


Column: Acclaim RSLC 120, C18, 2.2 μ m, Analytical (2.1 \times 150 mm)
 Flow: 0.475 mL/min
 Temperature: 30 $^{\circ}$ C
 Injection Volume: 0.5 μ L

Eluent: A: 9% CH₃CN, 10% formic acid
 B: 36% CH₃CN, 10% formic acid
 Gradient: 0.0–0.9 min, 100% A
 0.9–8.00 min to 28.5% B
 Hold for 2 min at 28.5% B

Detection: Absorbance, vis at 540 nm
 Sample: 5A: Adulterated PJ
 5B: Grape Juice

Peaks: 1.Dp3,5 33.50 μ g/mL*
 2. Cy3,5 61.00
 3. Dp3 26.00
 4. Pg3,5 5.15
 5. Cy3 73.00
 6. Pg3 2.10

*Calculated Concentrations



Peaks: 1.Dp3,5 —
 2. Cy3,5 4.95 μ g/mL*
 3. Dp3 33.95
 4. Pg3,5 7.25
 5. Cy3 110.70

*Calculated Concentrations




Figure 7-16. Separation of anthocyanins in simulated adulterated pomegranate juice (5A) overlaid with a separation of anthocyanins in grape juice (5B).



Targeted Approaches

Fruit Juice Adulteration: Profiling

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Anthocyanins are water-soluble plant pigments widely present in fruits, vegetables, and flowers. Anthocyanins have gained considerable interest in the scientific community and consumer market due to their anti-inflammatory action and strong antioxidant and radiation-protection properties. Bilberries are known to have a high anthocyanin content and are, therefore, one of the most expensive botanical ingredients in the health food industry. The high price of the extract makes bilberries more susceptible to adulteration.

Analytical characterization and quantification methods for anthocyanins specific to bilberries are therefore required. This study demonstrates a rapid separation liquid chromatography-mass spectrometry (RSLC-MS) method using the MSQ Plus mass spectrometer for the determination of anthocyanins in bilberry extract and pomegranate juice based on previously developed LC separations.



Chromatographic Conditions

System: UltiMate 3000 RSLC System
 Column: Acclaim RSLC 120 C18, 2.1 × 100 mm, 2.2 μm
 Flow: 0.5 mL/min
 Temperature: 40 °C
 Injection Volume: 5 μL
 Mobile Phase: A) CH₃CN, B) DI water, C) 20% formic acid
 Gradient: C held constant at 10%, A from 0% to 8% from 11 to 42 min, held for 13 min, return to 0% in 5 min
 Sample: NIST Bilberry Extract Standard Reference Material

Mass Spectrometric Conditions

System: MSQ Plus single quadrupole mass spectrometer
 Ionization interface: ESI
 Probe Temperature: 500 °C
 Needle Voltage: 2000 V
 Nebulizer Gas: Nitrogen at 80 psi
 Detection Mode: SIM
 Refer to chromatogram for SIM Acquisitions




Figure 7-17. Determination of 15 anthocyanins in bilberry extract by LC-MS.



Targeted Approaches

Fruit Juice Adulteration: Profiling

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

The primary mission of the Office of Dietary Supplements (ODS) at the National Institutes of Health (NIH) is to promote the quality, safety, and efficacy of dietary supplements. To accomplish this mission, authentic reference materials that closely match the matrix components of the dietary supplements are needed.

Standard Reference Materials (SRMs) for vaccinium (e.g., cranberries, blueberries, and bilberries) are being developed at the National Institute of Standards and Technologies (NIST) in collaboration with the NIH-ODS to evaluate these types of dietary supplements. Several SRMs with certified values for organic acids are currently available from NIST to aid dietary supplement and juice manufacturers in their analytical method development and QA/QC operations. Further work to certify anthocyanins and anthocyanidins is under way.

The popularity of cranberries is primarily associated with their reported ability to combat urinary tract infection and their high organic acid, anthocyanin, and polyphenol content. The primary health benefits of cranberries may be derived from their antioxidant properties, and their anthocyanin profiles are important in determining good quality and authenticity of the products. Cranberry extracts are prone to adulteration with lower value products in order to offer consumers competitive prices, therefore making it important to ensure the quality, safety and efficacy of these products. Organic acids present in cranberries impart flavor, and their specific ratios can be used to detect adulteration.

Organic acids in cranberry extract were measured using suppressed conductivity detection. Ratios of quinate, malate, and citrate can be used to identify potential adulteration.

Anthocyanin content were measured by 14 participating laboratories and results were submitted to the NIST. An overall consensus of the average reported concentration was observed, indicating validity of the data and the method.

Column:	Dionex IonPac AG11-HC, AS11-HC, 2 mm						
Flow:	0.38 mL/min						
Temperature:	30 °C						
Injection Volume:	5 µL						
Eluent:	1 mM KOH from -8 to 8 min, 1—30 mM from 8 to 20 min, 30—60 mM from 20 to 30 min, 60 mM from 30—45 min						
Eluent Source:	Dionex EGC II KOH with Dionex CR-ATC						
Detection:	Dionex ASRS 300 suppressor. 2mm, recycle mode, 57 mA						
Sample Prep:	0.1 g cranberry extract/160 mL DI water (centrifugation at 25 °C for 15 min for each 40 mL aliquot)						
Peaks:	<table> <tr> <td>1. Quinate</td> <td>11.7 mg/g</td> </tr> <tr> <td>2. Malate</td> <td>5.5</td> </tr> <tr> <td>3. Citrate</td> <td>21.7</td> </tr> </table>	1. Quinate	11.7 mg/g	2. Malate	5.5	3. Citrate	21.7
1. Quinate	11.7 mg/g						
2. Malate	5.5						
3. Citrate	21.7						




Figure 7-18 Separation of organic acids in cranberry powder.



Targeted Approaches

Fruit Juice Adulteration: Profiling

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Column: Thermo Scientific™ Accucore™ C18, 2.6 μm, analytical, 2.1 × 150 mm
 Flow: 0.65 mL/min
 Temperature: 42 °C
 Injection Volume: 2.0 μL
 Eluent: A: 10% Formic acid
 B: 10% Formic acid, 22.5% methanol, 22.5% MeCN
 Gradient: 10% B, 0–15 min
 Detection: Absorbance, Vis 520 nm
 Sample: 25 mg/mL cranberry extract
 Peaks: *Conc. in mg/g

1. Dp3Gal	0.021
2. Dp3Glu	0.037
3. Dp3Ara	-
4. Cy3Gal	0.090
5. Cy3Glu	0.041
6. Cy3Ara	0.157
7. Peo3Gal	0.176
8. Peo3Glu	0.032
9. Peo3Ara	0.170
10. Mal3Gal	0.088

*Calculated Concentration




Figure 7-19. Separation of anthocyanins in cranberry powder.

Table 7-1. Comparison of the experimentally determined anthocyanin values to the average values determined by the NIST collaborative study.

Analyte	Experimental Values (mg/g) n = 3	Average Values Reported by Collaborative Study (mg/g)
Dp3Gal	0.021 ± 0.01	0.019 ± 0.03
Dp3Glu	0.097 ± 0.05	0.057 ± 0.06
Cy3Gal	0.090 ± 0.04	0.130 ± 0.04
Cy3Glu	0.041 ± 0.07	0.059 ± 0.08
Cy3Ara	0.157 ± 0.03	0.200 ± 0.06
Peo3Gal	0.176 ± 0.02	0.047 ± 0.05
Peo3Ara	0.170 ± 0.02	0.150 ± 0.05
Mal3Gal	0.088 ± 0.05	0.060 ± 0.09



Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Targeted Approaches

Adulteration and Authentication



Melamine Adulteration

In 2008 investigations of death and health problems of babies in China revealed that some baby foods (milk powder) were contaminated by melamine. Some manufacturers illegally used melamine as an adulterant to increase the apparent protein content.



Targeted Approaches

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Standard methods enacted by the Chinese government for determining melamine in raw milk and dairy products included HPLC-UV, LC-MS, and GC-MS methods. However, the high cost of operation and maintenance of GC/LC-MS systems as well as the labor intensive derivatization that GC-MS requires limits their use in the milk product factories. The HPLC-UV method therefore is presently the popular choice for most factories.



Melamine in Milk

HPLC-UV Detection

Guard Column: Acclaim 120 C18, 5 μ m, 4.3 \times 10 mm
 Anal. Column: Acclaim 120 C18, 5 μ m, 4.6 \times 250 mm
 Flow: 1.0 mL/min
 Column Temp.: 40 $^{\circ}$ C
 Injection Volume: 20 μ L
 Eluents: Mix of 10 mM citric acid and
 10 mM sodium 1-octane sulfonate (pH 3)
 CH₃CN (92 : 8, v/v)
 Detection: UV at 240 nm

Traces: A. 0.2 μ g/mL
 B. 0.5 μ g/mL
 C. 2.0 μ g/mL
 D. 20 μ g/mL
 E. 25 μ g/mL
 F. 50 μ g/mL
 G. 100 μ g/mL
 Peak: 1. Melamine





Figure 7-20. Analysis of standards.



Targeted Approaches

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Melamine in Milk

HPLC-UV Detection

Guard Column: Acclaim 120 C18, 5 μ m, 4.3 \times 10 mm
 Anal. Column: Acclaim 120 C18, 5 μ m, 4.6 \times 250 mm
 Flow: 1.0 mL/min
 Column Temp.: 40 $^{\circ}$ C
 Injection Volume: 20 μ L
 Eluents: Mix of 10 mM citric acid and
 10 mM sodium 1-octane sulfonate (pH 3)
 CH_3CN (92 : 8, v/v)
 Detection: UV at 240 nm
 Peak: 1. Melamine




Figure 7-21. Analysis of milk powder samples.

HPLC-Charged Aerosol Detection

HPLC-Charged Aerosol Detection Parameters

Column: Amino, 4.6 \times 250, 5 μ m
 Flow: 1.0 mL/min
 Column Temp.: 30 $^{\circ}$ C
 Injection Volume: 25 μ L
 Mobile Phase: 80/10/10 MeCN/MeOH/ H_2O

Sample Preparation

A sample of 0.5 mL of milk was placed in a glass tube and 2.0 mL of 80% acetonitrile added to help precipitate any proteins. The tube was mixed well and centrifuge. The supernatant was then poured onto previously condition SPE columns as described below.

Solid Phase Extraction Protocol

1. Condition SPE* Column: 2.0 mL MeOH, 2.0 mL 0.1% TFA,
2. 0.5 mL milk + 2.0 mL 80% ACN, mix, centrifuge
3. Pour onto column
4. Wash: 2.0 mL 0.1% TFA,
5. 2.0 mL MeOH, vacuum dry for 30 seconds
6. Elute: 1.0 mL 5% NH_3 in MeOH, Dry for 30 seconds
7. Repeat Step 6
8. Evaporate at 40 $^{\circ}$ C
9. Reconstitute with 0.50 mL 50% CAN
10. Analyze




Figure 7-22. Analysis of milk spiked with melamine.



Targeted Approaches

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

A few years ago contaminated food led to numerous pet animal deaths and health problems. The toxins were found to be melamine and cyanuric acid, which when present together form an insoluble crystal, causing kidney failure. Since contaminated wheat gluten, rice protein concentrate, and corn gluten used in animal feed can be also used in human food (for example, bread, pasta, and baby food), it is crucial to monitor the presence of melamine and cyanuric acid in raw materials and animal tissue.




Figure 7-23. Structures of melamine and cyanuric acid and full scan MS spectra.

Melamine and Cyanuric Acid

Column:	Acclaim Mixed-Mode WAX-1, 5 μ m, 2.1 \times 150 mm	Mass Spectrometric Conditions
Flow:	0.25 mL/min	Ionization Interface: Electrospray Ionization (ESI)
Temperature:	20 $^{\circ}$ C	Detection Mode: Selected Ion Monitoring (SIM)
Injection Volume:	5 μ L	Probe Temp.: 500 C
Eluent:	90% Acetonitrile 10% Ammonium acetate buffer, 20 mM, pH 4	Needle Voltage: 2000 V
Detection:	MS	Cone Voltage: 50 V
		Dwell Time: 0.5 s for each SIM




Figure 7-24. SIM chromatogram of an FDA control sample (pork) spiked with melamine, cyanuric acid, and internal standards.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Targeted Approaches

Adulteration and Authentication



Curry powder

Mixtures of Complex Materials

Understanding the dynamics of a mixture of complex materials poses difficult analytical problems. The most common analytical approach for solving this issue is to assay for a single, hopefully unique, analyte for each of the constituents in the mixture. However, such marker peaks are not always seen by the analytical approach being used, or if found are not always known.



Targeted Approaches

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Mixtures of Complex Materials

Deconvolution of a Spice Mix

Gradient HPLC with coulometric electrochemical array detection is showing promise for distinguishing components contained in complex mixtures. In Application Brief 168 this approach was used to measure the contribution of cumin and coriander to a curry powder.




Figure 7-25. Curry powder mixture constituents (A) cumin and (B) coriander.




Figure 7-26. Curry powder blend.




Figure 7-27. Contribution of single components to a complex mixture based on marker peaks. Peak 2: abundance of cumin in mix. Peak 3: abundance of coriander in the mix. Only a single channel is presented for clarity.

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)


[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Adulteration and Authentication

Part 2 Non-targeted Methods

The measurement of patterns of both known and unknown analytes in a sample to authenticate product and to highlight possible adulteration. This is also referred to as metabolite profiling or metabolomic approaches.




Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Metabolomics

Adulteration and Authentication



Screening Herbs for Authenticity and Blending with Adulterants

Herb adulteration can include both inadvertent misidentification of the herb and deliberate blending of the herb with cheaper or more available substitutes. Although the adulteration of herbs is a frequent phenomenon, there are few simple methods available for the screening of large numbers of commercial batches of product.

Screening Herbs: Chemical Fingerprinting

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



The challenge to herb screening arises from the complexity and variability of genuine material combined with unrelenting conduct of adulteration. Herb and fruit variety, growing region, season, ripeness, and processing methods all contribute to the variability of the authentic product, making unambiguous characterization difficult.

Currently, one of the most reliable and applicable authentication methods is based on analytical chemical fingerprinting (non-targeted metabolomic) techniques. Gradient HPLC with coulometric electrochemical array detection is particularly suitable for generating information rich metabolite fingerprints of endogenous electroactive metabolites termed the “redoxome”.

The analytes that form the redoxome are compounds that influence color, flavor, nutritional, value, stability, and aroma. Such fingerprints can

be interrogated using pattern recognition and unsupervised statistical programs such as principal component analysis (PCA) to evaluate the authenticity (or possible adulteration) of an unknown sample by comparing its chromatogram with a compiled population of authenticated reference samples in the database. Poster Note 70534 illustrates the use of this approach to examine deliberate adulteration of herbs.

Authentic herbs and deliberate blends were extracted and analyzed using a gradient HPLC with coulometric electrochemical array detection method. The resulting chromatograms, containing patterns of both known and unknown analytes, were then used to form a database. An herb extract sample can then be analyzed and its metabolite fingerprint checked against the database to confirm authenticity or identify adulteration.



Oregano



Marjoram



Thyme

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Screening Herbs: Reference Samples

Column: Acclaim 120, C18, 3 × 150 mm, 3 μm
 Flow: 0.65 mL/min
 Injection Volume: 10 or 20 μL
 Mobile Phase: A: 20 mM Monobasic sodium phosphate, 3% acetonitrile, 0.2% tetrahydrofuran, pH 3.35
 B: 20 mM Monobasic sodium phosphate, 50% acetonitrile, 10% Tetrahydrofuran, pH 3.45
 C: 90% Methanol
 Gradient: 0–2 min: 2% B/3% C.; 30 min: 97% B/3% C.; 45 min: 97% B/3% C
 Curve 7
 EC Parameters: 16-channel array from 0 to +900 mV in +60 mV increments
 Data Station: Chromeleon CDS 6.8 SR9 and CoulArray software 3.1; Pirouette® software V4.5




Figure 7-28. EC array chromatograms of A) dried pure oregano, B) dried pure marjoram and C) dried pure thyme herbs (low sensitivity presented for clarity).

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



For easy visualization, metabolite patterns for each sample can be plotted in 3D space using principal component analysis. The more similar samples analytical fingerprints are, the closer they reside in the 3D plot.

In this way the location of an unknown sample can readily indicate authenticity (fingerprint similarity), level of adulteration, or whether a sample is totally different from samples in the database. As discussed in Application Note 1064 this approach can also be used to authenticate or assess adulteration of many different samples such as botanicals, wine, beer, tea, coffee, honey, essential oils, olive oil, and juices.



Screening Herbs: Principal Component Analysis




Figure 7-29. Principal component analysis of oregano-marjoram and oregano-thyme blends showing how the position in 3D space shifts according to the parent herb and the amount of blending taking place (shown as % in the figure).

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Metabolomics

Adulteration and Authentication



Screening Beer for Origin and Style

Beer is the most widely consumed alcoholic beverage in the world and the third most popular drink after water and tea. It is typically brewed from four basic ingredients: water, a starch source (e.g., malted barley), brewer's yeast, and a flavoring agent such as hops. Many varieties of beer result from differences in these ingredients, the additives used and the brewing process followed.

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



As discussed in Application Note 1065, gradient HPLC with coulometric electrochemical array detection could readily differentiate different beer styles, country of origin, and whether samples were normal or light.



Screening Beer: Principal Component Analysis




Figure 7-30. Various beer samples differentiated using principal component analysis.

Did You Know?

- Beer is the second most popular beverage in the world, coming in behind tea.
- After consuming a good quantity of a potent brew they called aul, or ale, the Vikings would head fearlessly into battle, often without armor or even shirts. In fact, “berserk” means “bare shirt” in Norse, and eventually took on the meaning of their wild battles.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Metabolomics

Adulteration and Authentication



Screening for Fraudulent Wine

Wine fraud is a major problem dating back to Roman times. Fraud includes adulteration and packaging cheap, poor quality wine with labels of more expensive, better products.

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Wine: Adulteration and Mislabeling

Wine adulteration is usually done the addition of cheaper products such as water, juices, or cheaper wines, but can also include the addition of harmful chemicals and sweeteners (e.g., ethylene glycol and lead acetate) to compensate color or flavor.

Conformation of fraud is not easy, and many analytical approaches have to be used. As discussed in Application Note 1064, one approach that can be used to study differences between grape varieties, blending, vintage, mislabeling, and effects of growing region is Gradient HPLC with coulometric electrochemical array detection. This approach can target potential marker peaks and can differentiate samples based on their metabolic fingerprints.

Table 7-2. Targeted analysis of some of the more abundant analytes measured in different wine samples. Wine 1: Cabernet Sauvignon, Argentina; Wine 2: Cabernet Sauvignon, South Africa; Wine 3: Cabernet Sauvignon, US; Wine 4: Cabernet Sauvignon, Chile; Wine 5: Hearty Burgundy, US.

Compound	Wine #1 Cabernet Sauvignon, Argentina (mg/L)	Wine #2 Cabernet Sauvignon, So. Africa (mg/L)	Wine #3 Cabernet Sauvignon, U.S. (mg/L)	Wine #4 Cabernet Sauvignon, Chile (mg/L)	Wine #5 Hearty Burgundy, U.S. (mg/L)
Apigenin	16	17.5	9.5	13	41
Caffeic Acid	8	13	5	17	3
Catechin Hydrate	37	26	26.5	24	22
Ellagic Acid Dihydrate	52	133	84	94	100
Epicatechin	19	15	16.5	11	4
Ferulic Acid	1	1	2	3	2
Gallic Acid	57	33.5	37	35	29.5
Isorhamnetin	6	5.5	2.5	6.5	2
Kaempferol	0.5	0.5	0.5	1	1
Myricetin	11	11	5	8	1.5
p-Coumaric Acid	8.5	16	2.5	14.5	3.5
Quercetin Dihydrate	13.5	15.5	3	14	4
cis-Resveratrol	1	1.5	0.5	2	0.5
trans-Resveratrol	2.5	2	1	2.5	1.5
Sinapic Acid	2	2	2	2	2
Syringic Acid	19	9.5	9	12	7
Vanillic Acid	6.5	4.5	2.5	8	4




Figure 7-31. Initial study showing principal component analysis of wines and the ability to distinguish between wine types.

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Metabolomics

Adulteration and Authentication



Fruit Juice Adulteration

Fruit juice adulteration presents an economic and regulatory problem. The most common forms of adulteration include simple dilution and blending of inexpensive and synthetically produced juices into the more expensive ones. Adulteration can also include dilution followed by addition of peel and/or pulp wash.

[Learn more about US Fruit Juice Adulteration Regulations](#)

[Learn more about the History of Fruit and Vegetable Juice US Regulations](#)

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Orange Juice Adulteration

The United States orange juice industry estimates that orange juice sales gross more than one billion dollars annually. While it is possible to identify orange juice adulteration through the measurement of some marker peak, or changes in levels of targeted analytes (see carbohydrate and organic ions earlier in the chapter), sometimes this is not possible.

As discussed in Application Note 1064 one approach that is showing promise to distinguish between orange juice varieties, juice dilution, the inclusion of peel or pulp wash or the effect of growing region is gradient HPLC with coulometric electrochemical array detection. A database of sample analyte fingerprints can be used to assess whether adulteration of pure orange juice has taken place.




Figure 7-32. Measurement of orange juice adulteration by blending with other juices (A) or by addition of orange peel or pulp wash (B) using the Spectro-Electro Array and PCA. GF – grapefruit; OJ – orange juice; OJ10%GF – orange juice blended with 10% grapefruit juice. POOL – equal blend of several orange juice samples.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Metabolomics

Adulteration and Authentication



Various types of tea leaves

Tea Adulteration

Tea is another product that can readily be adulterated, including the addition of iron or leather flakes, kaolin, and turmeric. Spent tea leaves can be darkened with coal tar and added to normal tea leaves to bulk up the amount for sale. Non-tea leaves may be used. This is especially true for green tea where paddy husks and “lie tea” can be added, as the latter closely resembles tea leaves.



Metabolomics

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

As discussed in Application Note 1064, one approach that is showing promise for authenticating tea, distinguishing between tea varieties, and with the potential to identify adulteration, is gradient HPLC with coulometric electrochemical array detection. A database of sample analyte fingerprints can be used to assess whether adulteration has taken place.

Table 7-3. Targeted analysis of some of the more abundant analytes measured in different types of tea.

Compound	Green Tea (mg/g)	Black Tea (mg/g)	White Tea (mg/g)
Catechin Hydrate	3.7	3.0	8.1
Epicatechin	50.8	9.3	39.8
Epicatechin Gallate	65.3	40.6	95.9
Epigallocatechin	49.2	2.5	32.3
Epigallocatechin Gallate	180	31.3	211
Gallocatechin	18.8	3.2	22.0
Gallocatechin Gallate	5.9	7.0	3.0

Did You Know?

All true tea (white, green, black, and oolong) comes from the same plant, the *Camellia sinensis*. It takes a minimum of three to five years for a tea bush to grow and be ready for harvesting.

Tea Adulteration




Figure 7-33. Initial study showing principal component analysis of tea and the ability to distinguish between tea types.



Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



Metabolomics

Adulteration and Authentication



Adulteration of Botanical Supplements

It is estimated that there are more than 50,000 dietary supplements currently available to the consumer. More than half of the U.S. adult population (53%–55%) consume dietary supplements. Although the most common are multivitamins, a growing proportion are botanical supplements. Unfortunately, there is great opportunity for fraud to take place, including selling the wrong plant species or using the wrong part of the plant (bark instead of root, leaf instead of flower, etc).

Adulteration of Botanical Supplements



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Authentication and identification of botanical supplements is difficult due to the complexity of samples, the fact that they often lack marker peaks, differences in harvesting and drying processes, and the natural variations in metabolite levels due to geographical location and season.

As discussed in Poster Note 70540 one approach that is gaining acceptance for studying botanical supplements is gradient HPLC with coulometric electrochemical array detection. Hoodia supplement is used to illustrate the capabilities of this approach.



Table 7-4. Some examples of botanical supplements and their known adulterants.

Proposed Botanical	Adulterant
Ginkgo (<i>Ginkgo biloba</i>) leaf extract standardized to flavonol glycosides and terpenes	Ginkgo (<i>Ginkgo biloba</i>) leaf extract with added flavonol glycosides or aglycones (e.g., rutin, quercetin, etc.)
Asian species of <i>Akebia</i> and <i>Clematis</i> stem	<i>Aristolochia manshuriensis</i> stem (<i>guan mu tong</i>)
Skullcap herb (<i>Scutellaria lateriflora</i>)	Germander herb (<i>Teucrium chamaedrys</i>)
Plantain leaf (<i>Plantago lanceolata</i>)	<i>Digitalis lanata</i> leaf
Black cohosh root/rhizome (<i>Actaea racemosa</i>)	Chinese cimicifuga root/rhizome (<i>Actaea</i> spp.)

Did You Know?

The leading food categories with reported cases of food fraud include:

- Olive Oil
- Milk and milk-based products
- Honey, maple syrup, and other natural sweeteners
- Fruit juice
- Coffee and tea
- Spices
- Clouding agents

Screening Botanical Supplements

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Principal Component Analysis of Hoodia

Authentic Hoodia herb, capsules, and tablets were extracted and analyzed using a gradient HPLC with coulometric electrochemical array detection

method. The resulting chromatograms, containing patterns of both known and unknown analytes were then used to form a database. This approach could distinguish between sample types and identify a sample of a common adulterant.




Figure 7-34. HPLC-ECD Array chromatograms of *Hoodia gordonii*: (A) mixed parts extract #3165 and (B) commercial product tablet extract HG6. Thanks to Drs. Ikhlas Khan and Bharathi Avula of the National Center for Natural Products Research at The University of Mississippi, Oxford, MS, for providing the plant materials and hoodigoside standards.

Screening Botanical Supplements

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)




Figure 7-35. Principal component analysis of hoodia samples showing differentiation between authentic dried plant, capsules/tablets and adulterant.

Did You Know?

The Grocery Manufacturers Association (GMA) estimates that fraud may cost the global food industry between \$10 billion and \$15 billion per year, affecting approximately 10% of all commercially sold food products. Compared to the trillions of dollars spent on food and food ingredients globally each year, however, the prevalence of food fraud is ultimately very low.



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)


[References](#)



Metabolomics

Adulteration and Authentication

TRYPTOPHAN



Chemical structure of tryptophan

Contamination of Tryptophan Supplements

Eosinophilia-Myalgia Syndrome (EMS) is an incurable and sometimes fatal flu-like neurological condition first recognized in 1989. EMS initially reported in New Mexico, was found to be related to use of L-tryptophan. Evidence suggested that the cause of the problem was improperly prepared tryptophan supplements, contaminated with potentially toxic metabolites. EMS resulted in 37 deaths and the banning of tryptophan supplements from 1991 to 2001

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Screening Tryptophan Supplements

Analysis of tryptophan supplements, commonly used for insomnia, depression, and PMS, using gradient HPLC with coulometric electrochemical array detection generated patterns of hundreds of electrochemically active analytes—both known and unknown compounds. Interrogation of these metabolomic profiles using chemometric software readily distinguished between control and contaminated supplement samples.




Figure 7-36. Structures of some possible tryptophan contaminants.




Figure 7-37. Principal component analysis of control and contaminated tryptophan supplements.



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Metabolomics

Adulteration and Authentication



Adulteration of Olive Oil

Olive oil, high in monounsaturated fatty acids, is a key component of the heart-healthy Mediterranean diet recommended by doctors. But not all olive oils are nutritionally the same. There are a number of commercial grades of olive oil differing in quality, nutrient levels, antioxidant activity, and final chemical composition.

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)


[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



Screening Olive Oil for Adulteration

“Virgin” olive oil is the juice of freshly harvested olive fruits; “extra virgin” is virgin olive oil of the highest quality based on both sensory and physical–chemical, characteristics. “Refined olive oil” is one that has been chemically purified; commercial “olive oil” is a blend of virgin and refined oils. Finally, “pomace” olive oil one that is obtained from the solid remains of olives already pressed for juice using solvent extraction.

Unfortunately, as consumer demand for extra virgin olive oil (EVOO) is set to outpace production, the inevitable increase in EVOO price makes it an ideal target for adulteration. Adulteration can take many forms. The ones that impact health the least are those that maintain the safety of

the product such as passing off inferior olive oil as superior olive oil, or blending olive oil with cheaper oils such as sunflower, soybean, sesame, rapeseed, canola, corn, palm, and hazelnut. Blending with peanut oil is of concern for those with peanut allergies.

As discussed in Poster Note 70689, to help tackle the problem of olive oil adulteration a simple HPLC-based method using charged aerosol detection was used to generate patterns of fats (triglycerides) found in oil samples. These patterns can be used to determine pure EVOO samples from olive oil blended with corn oil, hazelnut oil, or pomace.




Figure 7-38. Principal component analysis of hydrolyzed EVOO with hazelnut (HA–HF), corn (CA–CF), and pomace oil (PA–PF) adulteration at six levels (1, 5, 10, 25, 50, and 75 mass-%, respectively) using hydrolyzed oil data. Inset shows results for pomace oil analysis.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



References



Technical Collateral and Peer Reviewed Journals

Here you'll find a multitude of references using our HPLC, ion chromatography and sample preparation solutions.

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



References



HPLC and UHPLC References



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Carbohydrates

Title	Authors	Publication	Publication Date
Carbohydrate and oligosaccharide analysis with a universal HPLC detector.	Asa, D.	<i>American Laboratory</i> 38, 16.	2006
Determination of levoglucosan in atmospheric aerosols using high performance liquid chromatography with aerosol charge detection.	Dixon, R. W.; Baltzell, G.	<i>J. Chromatogr., A.</i> 1109 (2), 214–221	2006 Mar 24
Composition of structural carbohydrates in biomass: Precision of a liquid chromatography method using a neutral detergent extraction and a charged aerosol detector.	Godin, B.; Agneessens, R.; Gerin, P. A.; Delcarte, J.	<i>Talanta</i> 85 (4), 2014–2026	2011 Sep 30
Selectivity issues in targeted metabolomics: Separation of phosphorylated carbohydrate isomers by mixed-mode hydrophilic interaction/weak anion exchange chromatography.	Hinterwirth, H.; Lämmerhofer, M.; Preinerstorfer, B.; Gargano, A.; Reischl, R.; Bicker, W.; Trapp, O.; Brecker, L.; Lindner, W.	<i>J. Sep. Sci.</i> 33 (21), 3273–3282	2010 Nov
Investigation of polar organic solvents compatible with Corona charged aerosol detection and their use for the determination of sugars by hydrophilic interaction liquid chromatography.	Hutchinson, J. P.; Remenyi, T.; Nesterenko, P.; Farrell, W.; Groeber, E.; Szucs, R.; Dicinowski, G.; Haddad, P. R.	<i>Anal. Chim. Acta.</i> 750, 199–206	2012 Oct 31
Characterization of an endoglucanase belonging to a new subfamily of glycoside hydrolase family 45 of the basidiomycete <i>Phanerochaete chrysosporium</i>.	Igarashi, K.; Ishida, T.; Hori, C.; Samejima, M.	<i>Appl. Environ. Microbiol.</i> 74 (18), 5628–5634	2008 Sep
Direct detection method of oligosaccharides by high-performance liquid chromatography with charged aerosol detection.	Inagaki, S.; Min, J. Z.; Toyo'oka, T.	<i>Biomed. Chromatogr.</i> 21 (4), 338–342	2007 Apr
Differential selectivity of the <i>Escherichia coli</i> cell membrane shifts the equilibrium for the enzyme-catalyzed isomerization of galactose to tagatose.	Kim, J. H.; Lim, B. C.; Yeom, S. J.; Kim, Y. S.; Kim, H. J.; Lee, J. K.; Lee, S. H.; Kim, S. W.; Oh, D. K.	<i>Appl. Environ. Microbiol.</i> 74 (8), 2307–2313	2008 Apr
Elution strategies for reversed-phase high-performance liquid chromatography analysis of sucrose alkanolate regioisomers with charged aerosol detection.	Lie, A.; Pedersen, L. H.	<i>J. Chromatogr., A.</i> 1311, 127–133	2013 Oct 11
Design of experiments and multivariate analysis for evaluation of reversed-phase high-performance liquid chromatography with charged aerosol detection of sucrose caprate regioisomers	Lie, A.; Wimmer, R.; Pedersen, L. H.	<i>J. Chromatogr., A.</i> 1281, 67–72	2013 Mar 15
Solvent effects on the retention of oligosaccharides in porous graphitic carbon liquid chromatography	Melmer, M.; Stangler, T.; Premstaller, A.; Lindner, W.	<i>J. Chromatogr., A</i> 1217 (39) 6092–6096	2010 Sep 24
Practical preparation of lacto-N-biose I, a candidate for the bifidus factor in human milk	Nishimoto, M.; Kitaoka, M.	<i>Biosci., Biotechnol., Biochem.</i> 71 (8), 2101–2104	2007 Aug



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Carbohydrates

Title	Authors	Publication	Publication Date
Cellotriose and cellotetraose as inducers of the genes encoding cellobiohydrolases in the basidiomycete <i>Phanerochaete chrysosporium</i>	Suzuki, H.; Igarashi, K.; Samejima, M.	<i>Appl. Environ. Microbiol.</i> 76 (18), 6164–6170	2010 Sep
1,2-alpha-L-Fucosynthase: A glycosynthase derived from an inverting alpha-glycosidase with an unusual reaction mechanism	Wada, J.; Honda, Y.; Nagae, M.; Kato, R.; Wakatsuki, S.; Katayama, T.; Taniguchi, H.; Kumagai, H.; Kitaoka, M.; Yamamoto, K.	<i>FEBS Lett.</i> 582 (27), 3739–3743	2008 Nov 12
Efficient separation of oxidized cello-oligosaccharides generated by cellulose degrading lytic polysaccharide monoxygenases	Westereng, B.; Agger, J. W.; Horn, S. J.; Vaaje-Kolstad, G.; Aachmann, F. L.; Stenström, Y. H.; Eijsink, V. G.	<i>J. Chromatogr., A.</i> 1271 (1), 144–152	2013 Jan 4
Distribution of in vitro fermentation ability of lacto-N-Biose I, a major building block of human milk oligosaccharides, in bifidobacterial strains	Xiao, J. Z.; Takahashi, S.; Nishimoto, M.; Odamaki, T.; Yaeshima, T.; Iwatsuki, K.; Kitaoka, M.	<i>Appl. Environ. Microbiol.</i> 76 (1), 54–59	2010 Jan





Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Characterization of phenolic compounds in strawberry (<i>Fragaria x ananassa</i>) fruits by different HPLC detectors and contribution of individual compounds to total antioxidant capacity	Aaby, K.; Ekeberg, D.; Skrede, G.	<i>J. Agric. Food Chem.</i> 55 (11), 4395–4406	2007 May 30
Analysis of flavonoids and other phenolic compounds using high-performance liquid chromatography with coulometric array detection: relationship to antioxidant activity	Aaby, K.; Hvattum, E.; Skrede, G.	<i>J. Agric. Food Chem.</i> 52 (15), 4595–4603	2004 Jul 28
Aqueous extract of Astragali Radix induces human natriuresis through enhancement of renal response to atrial natriuretic peptide	Ai, P.; Yong, G.; Dingkun, G.; Qiuyu, Z.; Kaiyuan, Z.; Shanyan, L.	<i>J. Ethnopharmacol.</i> 116 (13), 413–421	2008 Mar 28
Antioxidant, α-amylase inhibitory and oxidative DNA damage protective property of <i>Boerhaavia diffusa</i> (Linn.) root	Akhter, F.; Hashim, A.; Khan, M. S.; Ahmad, S.; Iqbal, D.; Srivastava, A. K.; Siddiqui, M. H.	<i>S. Afr. J. Bot.</i> 88, 265–272	2013 Sep
Antioxidant activity and metabolite profile of quercetin in vitamin-E-depleted rats.	Ameho, C. K.; Chen, C. Y. O.; Smith, D.; Sánchez-Moreno, C.; Milbury, P. E.; Blumberg, J. B.	<i>J. Nutr. Biochem.</i> 19 (7), p.467–474	2008 Jul
Evaluation of tolerable levels of dietary quercetin for exerting its antioxidative effect in high cholesterol-fed rats	Azuma, K.; Ippoushi, K.; Terao, J.	<i>Food Chem. Toxicol.</i> 48 (4), 1117–1122	2010 Apr
Recent methodology in ginseng analysis	Baek, S.; Bae, O.; Park, J.	<i>J. Ginseng Res.</i> 36 (2), 119–134	2012 Apr
Sensitive determination of saponins in radix et rhizoma notoginseng by charged aerosol detector coupled with HPLC	Bai, C.; Han, S.; Chai, X.; Jiang, Y.; Li, P.; Tu, P.	<i>J. Liq. Chromatogr. Relat. Technol.</i> 32 (2), 242–260	2010 Aug 27
Comprehensive analysis of polyphenols in 55 extra virgin olive oils by HPLC-ECD and their correlation with antioxidant activities	Bayram, B.; Esatbeyoglu, T.; Schulze, N.; Ozcelik, B.; Frank, J.; Rimbach, G.	<i>Plant Foods Hum. Nutr. (N. Y., NY, U.S.)</i> 67 (4), 326–336	2012 Dec
Hydrogen sulfide mediates the vasoactivity of garlic	Benavides, G. A.; Squadrito, G. L.; Mills, R. W.; Patel, H. D.; Isbell, T. S.; Patel, R. P.; Darley-Usmar, V. M.; Doeller, J. E.; Kraus, D. W.	<i>Proc. Natl. Acad. Sci. U.S.A.</i> 104 (46), 17977–17982	2007 Nov
Analysis of selected stilbenes in <i>Polygonum cuspidatum</i> by HPLC coupled with CoulArray detection	Benová, B.; Adam, M.; Onderková, K.; Královský, J.; Krajček, M.	<i>J. Sep. Sci.</i> 31 (13), 2404–2409	2008 Jul
Rapid and complete extraction of phenols from olive oil and determination by means of a coulometric electrode array system	Brenes, M.; García, A.; García, P.; Garrido, A.	<i>J. Agric. Food Chem.</i> 48 (11), 5178–5183	2000 Nov
The real nature of the indole alkaloids in <i>Cortinarius infractus</i>: Evaluation of artifact formation through solvent extraction method development	Brondz, I.; Ekeberg, D.; Høiland, K.; Bell, D.; Annino, A.	<i>J. Chromatogr., A</i> 1148 (1), 1–7	2007 Apr 27



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Chemotaxonomic differentiation between <i>Cortinarius infractus</i> and <i>Cortinarius subtortus</i> by supercritical fluid chromatography connected to a multi-detection system	Brondz, I.; Høiland, K.	<i>Trends Chromatogr.</i> 4, 79–87	2008
Carotenoid bioavailability is higher from salads ingested with full-fat than with fat-reduced salad dressings as measured with electrochemical detection	Brown, M. J.; Ferruzzi, M. G.; Nguyen, M. L.; Cooper, D. A.; Eldridge, A. L.; Schwartz, S. J.; White, W. S.	<i>Am. J. Clin. Nutr.</i> 80 (2), 396–403	2004 Aug
Naringenin from cooked tomato paste is bioavailable in men	Bugianesi, R.; Catasta, G.; Spigno, P.; D'Uva, A.; Maiani, G.	<i>J. Nutr.</i> 132 (11), 3349–3352	2002 Nov
"Dilute-and-shoot" triple parallel mass spectrometry method for analysis of vitamin D and triacylglycerols in dietary supplements	Byrdwell, W. C.	<i>Anal. Bioanal. Chem.</i> 401 (10), 3317–3334	2011 Dec
Human skeletal muscle ascorbate is highly responsive to changes in vitamin C intake and plasma concentrations	Carr, A. C.; Bozonet, S. M.; Pullar, J. M.; Simcock, J. W.; Vissers, M. C.	<i>Am. J. Clin. Nutr.</i> 97 (4), 800–807	2013 Apr
Utilization of RP-HPLC fingerprinting analysis for the identification of diterpene glycosides from <i>Stevia rebaudiana</i>	Chaturvedula, V.; Prakash, I.	<i>Int. J. Res. Phytochem. Pharmacol.</i> 1 (2), 88–92	2011 Jun 9
Acid and alkaline hydrolysis studies of stevioside and rebaudioside A	Chaturvedula, V.; Prakash, I.	<i>J. Appl. Pharm. Sci.</i> 1 (8), 104–108	2011 Oct
Spectral analysis and chemical studies of the sweet constituent, rebaudioside A	Chaturvedula, V.; Prakash, I.	<i>Eur. J. Med. Plants</i> 2 (1), 57–65	2012 Feb
Flavonoids from almond skins are bioavailable and act synergistically with vitamins C and E to enhance hamster and human LDL resistance to oxidation	Chen, C.; Milbury, P. E.; Lapsley, K.; Blumberg, J. B.	<i>J. Nutr.</i> 135 (6), 1366–1373	2005 Jun 1
Photostability of rebaudioside A and stevioside in beverages	Clos, J. F.; Dubois, G. E.; Prakash, I.	<i>J. Agric. Food Chem.</i> 56 (18), 8507–8513	2008 Sep 24
CoulArray electrochemical evaluation of tocopherol and tocotrienol isomers in barley, oat and spelt grains	Colombo, M. L.; Marangon, K.; Bugatti, C.	<i>Nat. Prod. Commun.</i> 4 (2), 251–254	2009 Feb
Composition and stability of phytochemicals in five varieties of black soybeans (<i>Glycine max</i>)	Correa, C. R.; Li, L.; Aldini, G.; Carini, M.; Oliver Chen, C. Y.; Chun, H.; Cho, S.; Park, K.; Russell, R. M.; Blumberg, J. B.; Yeum, K.	<i>Food Chem.</i> 123 (4), 1176–1184	2010 Dec 15
Effect of UV-B light and different cutting styles on antioxidant enhancement of commercial fresh-cut carrot products	Du, W.; Avena-Bustillos, R. J.; Breksa, A. P., III.; McHugh, T. H.	<i>Food Chem.</i> 134 (4), 1862–1869	2012 Oct 15



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Phenols, lignans and antioxidant properties of legume and sweet chestnut flours	Durazzo, A.; Turfani, V.; Azzini, E.; Maiani, G.; Carcea, M.	<i>Food Chem.</i> 140 (4), 666–671	2013 Oct 15
alpha-Lipoic acid in dietary supplements: development and comparison of HPLC-CEAD and HPLC-ESI-MS methods	Durrani, A. I.; Schwartz, H.; Schmid, W.; Sontag, G.	<i>J. Pharm. Biomed. Anal.</i> 45 (4), 694–699	2007 Nov 30
Comparison between evaporative light scattering detection and charged aerosol detection for the analysis of saikosaponins	Eom, H. Y.; Park, S. Y.; Kim, M. K.; Suh, J. H.; Yeom, H.; Min, J. W.; Kim, U.; Lee, J.; Youm, J. R.; Han, S. B.	<i>J. Chromatogr., A.</i> 1217 (26), 4347–4354	2010 Jun 25
Assessment of microcystin purity using charged aerosol detection	Edwards, C.; Lawton, L. A.	<i>J. Chromatogr., A.</i> 1217 (32), 5233–5238	2010 Aug 6
Analysis of lycopene geometrical isomers in biological microsamples by liquid chromatography with coulometric array detection	Ferruzzi, M. G.; Nguyen, M. L.; Sander, L. C.; Rock, C. L.; Schwartz, S. J.	<i>J. Chromatogr., B: Biomed. Sci. Appl.</i> 760 (2), 289–299	2001 Sep 5
Charged aerosol detection to characterize components of dispersed-phase formulations	Fox, C. B.; Sivananthan, S. J.; Mikasa, T. J.; Lin, S.; Parker, S. C.	<i>Adv. Colloid Interface Sci.</i> 199–200, 59–65	2013 Nov
HPLC with charged aerosol detection for the measurement of natural products	Fukushima, K.; Kanedai, Y.; Hirose, K.; Matsumoto, T.; Hashiguchi, K.; Senda, M.; et al.	<i>Chromatography 27 (Suppl. 1)</i> , 83–86	2006
Determination of heterocyclic aromatic amines in beef extract, cooked meat and rat urine by liquid chromatography with coulometric electrode array detection	Gerbl, U.; Cichna, M.; Zsivkovits, M.; Knasmüller, S.; Sontag, G.	<i>J. Chromatogr., B: Anal. Technol. Biomed. Life Sci.</i> 802 (1), 107–113	2004 Mar 25
Determination of macrolide antibiotics in porcine and bovine urine by high-performance liquid chromatography coupled to coulometric detection	González de la Huebra, M. J.; Vincent, U.; Bordin, G.; Rodríguez, A. R.	<i>Anal. Bioanal. Chem.</i> 382 (2), 433–439	2005 May
Development and validation of HPLC-DAD-CAD-MS3 method for qualitative and quantitative standardization of polyphenols in <i>Agrimoniae eupatoriæ herba</i> (Ph. Eur)	Granica, S.; Krupa, K.; Klebowska, A.; Kiss, A. K.	<i>J. Pharm. Biomed. Anal.</i> 86, 112–122	2013 Dec
Total reducing capacity of fresh sweet peppers and five different Italian pepper recipes	Greco, L.; Riccio, R.; Bergero, S.; Del Re, A. A. M.; Trevisan, M.	<i>Food Chem.</i> 103 (4), 1127–1133	2007 Jan
Urinary 3-(3,5-dihydroxyphenyl)-1-propanoic acid, an alkylresorcinol metabolite, is a potential biomarker of whole-grain intake in a U.S. population	Guyman, L. A.; Adlercreutz, H.; Koskela, A.; Li, L.; Beresford, S. A.; Lampe, J. W.	<i>J. Nutr.</i> 138 (10), 1957–1962	2008 Oct
Multidimensional LC x LC analysis of phenolic and flavone natural antioxidants with UV-electrochemical coulometric and MS detection	Hájek, T.; Skeríková, V.; Cesla, P.; Vynuchalová, K.; Jandera, P.	<i>J. Sep. Sci.</i> 31 (19), 3309–3328	2008 Oct
Determination of the urinary aglycone metabolites of vitamin K by HPLC with redox-mode electrochemical detection	Harrington, D. J.; Soper, R.; Edwards, C.; Savidge, G. F.; Hodges, S. J.; Shearer, M. J.	<i>J. Lipid Res.</i> 46 (5), 1053–1060	2005 May



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Bioavailability and antioxidant effect of epigallocatechin gallate administered in purified form versus as green tea extract in healthy individuals	Henning, S. M.; Niu, Y.; Liu, Y.; Lee, N. H.; Hara, Y.; Thames, G. D.; Minutti, R. R.; Carpenter, C. L.; Wang, H.; Heber, D.	<i>J. Nutr. Biochem.</i> 16 (10), 610–616	2005 Oct
Procyanidin dimer B₂ [epicatechin-(4β-8)-epicatechin] in human plasma after the consumption of a flavanol-rich cocoa	Holt, R. R.; Lazarus, S. A.; Sullards, M. C.; Zhu, Q. Y.; Schramm, D. D.; Hammerstone, J. F.; Fraga, C. G.; Schmitz, H. H.; Keen, C. L.	<i>Am. J. Clin. Nutr.</i> 76 (4), 798–804	2002 Oct
Effects of natural (RRR α-tocopherol acetate) or synthetic (all-rac α-tocopherol acetate) vitamin E supplementation on reproductive efficiency in beef cows	Horn, M.; Gunn, P.; Van Emon, M.; Lemenager, R.; Burgess, J.; Pyatt, N. A.; Lake, S. L.	<i>J. Anim. Sci. (Savoy, IL, U.S.)</i> 88 (9), 3121–3127	2010 Sep
RP-HPLC analysis of phenolic compounds and flavonoids in beverages and plant extracts using a CoulArray detector	Jandera, P.; Skeifíková, V.; Rehová, L.; Hájek, T.; Baldríanová, L.; Skopová, G.; Kellner, V.; Horna, A.	<i>J. Sep. Sci.</i> 28 (9–10), 1005–1022	2005 Jun
A new application of charged aerosol detection in liquid chromatography for the simultaneous determination of polar and less polar ginsenosides in ginseng products	Jia, S.; Li, J.; Yunusova, N.; Park, J. H.; Kwon, S. W.; Lee, J.	<i>Phytochem. Anal.</i> 24 (4), 374–380	2013 Jul–Aug
A combination of aspirin and γ-tocopherol is superior to that of aspirin and α-tocopherol in anti-inflammatory action and attenuation of aspirin-induced adverse effects	Jiang, Q.; Moreland, M.; Ames, B. N.; Yin, X.	<i>J. Nutr. Biochem.</i> 20 (11), 894–900	2009 Nov
HPLC analysis of rosmarinic acid in feed enriched with aerial parts of <i>Prunella vulgaris</i> and its metabolites in pig plasma using dual-channel coulometric detection	Jirovský, D.; Kosina, P.; Myslíková, M.; Stýskála, J.; Ulrichová, J.; Simánek V.	<i>J. Agric. Food Chem.</i> 55 (19), 7631–7637	2007 Sep 19
Molar absorptivities and reducing capacity of pyranoanthocyanins and other anthocyanins	Jordheim, M.; Aaby, K.; Fossen, T.; Skrede, G.; Andersen, Ø. M.	<i>J. Agric. Food Chem.</i> 55 (26), 10591–10598	2007 Dec 26
Sensitive electrochemical detection method for alpha-acids, beta-acids and xanthohumol in hops (<i>Humulus lupulus</i> L.)	Kac, J.; Vovk, T.	<i>J. Chromatogr., B: Anal. Technol. Biomed. Life Sci.</i> 850 (1–2), 531–537	2007 May 1
Determination of phenolic compounds and hydroxymethylfurfural in meads using high performance liquid chromatography with coulometric-array and UV detection	Kahoun, D.; Rezková, S.; Veskrnová, K.; Královský, J.; Holcapek, M.	<i>J. Chromatogr., A</i> 1202 (1), 19–33	2008 Aug 15
Analysis of terpene lactones in a Ginkgo leaf extract by high-performance liquid chromatography using charged aerosol detection	Kakigi, Y.; Mochizuki, N.; Icho, T.; Hakamatsuka, T.; Goda, Y.	<i>Biosci., Biotechnol., Biochem.</i> 74 (3), 590–594	2010
Linear aglycones are the substrates for glycosyltransferase DesVII in methymycin biosynthesis: analysis and implications	Kao, C.; Borisova, S.; Kim, H.; Liu, H.	<i>J. Am. Chem. Soc.</i> 128 (17), 5606–5607	2006 May 3



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Antioxidant-rich food intakes and their association with blood total antioxidant status and vitamin C and E levels in community-dwelling seniors from the Quebec longitudinal study NuAge	Khalil, A.; Gaudreau, P.; Cherki, M.; Wagner, R.; Tessier, D. M.; Fulop, T.; Shatenstein, B.	<i>Exp. Gerontol.</i> 46 (6), 475–481	2011 Jun
Certification of a pure reference material for the ginsenoside Rg1	Kim, D.; Chang, J.; Sohn, H.; Cho, B.; Ko, S.; Nho, K.; Jang, D.; Lee, S.	<i>Accredit. Qual. Assur.</i> 15 (2), 81–87	2009 Sep
Optimization of pressurized liquid extraction for spicatoside A in <i>Liriope platyphylla</i>	Kim, S. H.; Kim, H. K.; Yang, E. S.; Lee, K. Y.; Kim, S. D.; Kim, Y. C.; Sung, S. H.	<i>Sep. Purif. Technol.</i> 71 (2), 168–172	2010
Production of surfactin and iturin by <i>Bacillus licheniformis</i> N1 responsible for plant disease control activity	Kong, H. G.; Kim, J. C.; Choi, G. J.; Lee, K. Y.; Kim, H. J.; Hwang, E. C.; Moon, B. J.; Lee, S. W.	<i>Plant Pathol. J.</i> 26 (2), 170–177	2010
Transepithelial transport of microbial metabolites of quercetin in intestinal Caco-2 cell monolayers	Konishi, Y.	<i>J. Agric. Food Chem.</i> 53 (3), 601–607	2005 Feb 9
Absorption and bioavailability of artemillin C in rats after oral administration	Konishi, Y.; Hitomi, Y.; Yoshida, M.; Yoshioka, E.	<i>J. Agric. Food Chem.</i> 53 (26), 9928–9933	2005 Dec 28
Pharmacokinetic study of caffeic and rosmarinic acids in rats after oral administration	Konishi, Y.; Hitomi, Y.; Yoshida, M.; Yoshioka, E.	<i>J. Agric. Food Chem.</i> 53 (12), 4740–4746	2005 Jun 15
Intestinal absorption of <i>p</i>-coumaric and gallic acids in rats after oral administration	Konishi, Y.; Hitomi, Y.; Yoshioka, E.	<i>J. Agric. Food Chem.</i> 52 (9), 2527–2532	2004 May 5
Microbial metabolites of ingested caffeic acid are absorbed by the monocarboxylic acid transporter (MCT) in intestinal Caco-2 cell monolayers	Konishi, Y.; Kobayashi, S.	<i>J. Agric. Food Chem.</i> 52 (21), 6418–6424	2004 Oct 20
Transepithelial transport of rosmarinic acid in intestinal Caco-2 cell monolayers	Konishi, Y.; Kobayashi, S.	<i>Biosci., Biotechnol., Biochem.</i> 69 (3), 583–591	2005 Mar
Effects of various doses of selenite on stinging nettle (<i>Urtica dioica</i> L.)	Krystofova, O.; Adam, V.; Babula, P.; Zehnalek, J.; Beklova, M.; Havel, L.; Kizek, R.	<i>Int. J. Environ. Res. Public Health</i> 7 (10), 3804–3815	2010 Oct
Biofortified cassava increases β-carotene and vitamin A concentrations in the TAG-rich plasma layer of American women	La Frano, M. R.; Woodhouse, L. R.; Burnett, D. J.; Burri, B. J.	<i>Br. J. Nutr.</i> 110 (2), 310–320	2013 Jul 28
Chlorogenic acid is absorbed in its intact form in the stomach of rats	Lafay, S.; Gil-Izquierdo, A.; Manach, C.; Morand, C.; Besson, C.; Scalbert, A.	<i>J. Nutr.</i> 136 (5), 1192–1197	2006 May
Determination of 4-ethylcatechol in wine by high-performance liquid chromatography-coulometric electrochemical array detection	Larcher, R.; Nicolini, G.; Bertoldi, D.; Nardin, T.	<i>Anal. Chim. Acta</i> 609 (2), 235–240	2008 Feb 25
Determination of volatile phenols in wine using high-performance liquid chromatography with a coulometric array detector	Larcher, R.; Nicolini, G.; Puecher, C.; Bertoldi, D.; Moser, S.; Favaro, G.	<i>Anal. Chim. Acta</i> 582 (1), 55–60	2007 Jan 16



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Acute, quercetin-induced reductions in blood pressure in hypertensive individuals are not secondary to lower plasma angiotensin-converting enzyme activity or endothelin-1: nitric oxide	Larson, A.; Witman, M. A. H.; Guo, Y.; Ives, S.; Richardson, R. S.; Bruno, R. S.; Jalili, T.; Symons, J. D.	<i>Nutr. Res. (N. Y., NY, U.S.)</i> 32 (8), 557–564	2012 Aug
High-performance liquid chromatography method for the determination of folic acid in fortified food products	Lebiedzinska, A.; Dałbrowska, M.; Szefer, P.; Marszałł M.	<i>Toxicol. Mech. Methods</i> 18 (6), 463–467	2008 Jul
Reversed-phase high-performance liquid chromatography method with coulometric electrochemical and ultraviolet detection for the quantification of vitamins B(1) (thiamine), B(6) (pyridoxamine, pyridoxal and pyridoxine) and B(12) in animal and plant foods	Lebiedzinska, A.; Marszałł, M. L.; Kuta, J.; Szefer, P.	<i>J. Chromatogr., A</i> 1173 (1–2), 71–80	2007 Nov 30
An improved method for the determination of green and black tea polyphenols in biomatrices by high-performance liquid chromatography with coulometric array detection	Lee, M. J.; Prabhu, S.; Meng, X.; Li, C.; Yang, C. S.	<i>Anal. Biochem.</i> 279 (2), 164–169	2000 Mar 15
Characterisation, extraction efficiency, stability and antioxidant activity of phytonutrients in <i>Angelica keiskei</i>	Li, L.; Aldini, G.; Carini, M.; Chen, C. Y. O.; Chun, H.; Cho, S.; Park, K.; Correa, C. R.; Russell, R. M.; Blumberg, J. B.; Yeum, K.	<i>Food Chem.</i> 115 (1), 227–232	2009 Jul
Vitamin A equivalence of the β-carotene in β-carotene-biofortified maize porridge consumed by women	Li, S.; Nugroho, A.; Rocheford, T.; White, W. S.	<i>Am. J. Clin. Nutr.</i> 92 (5), 1105–1112	2010 Nov
Phase IIa chemoprevention trial of green tea polyphenols in high-risk individuals of liver cancer: modulation of urinary excretion of green tea polyphenols and 8-hydroxydeoxyguanosine	Luo, H.; Tang, L.; Tang, M.; Billam, M.; Huang, T.; Yu, J.; Wei, Z.; Liang, Y.; Wang, K.; Zhang, Z. Q.; Zhang, L.; Wang, J. S.	<i>Carcinogenesis</i> 27 (2), 262–268	2006 Feb
Determination of four water-soluble compounds in <i>Salvia miltiorrhiza</i> Bunge by high-performance liquid chromatography with a coulometric electrode array system	Ma, L.; Zhang, X.; Guo, H.; Gan, Y.	<i>J. Chromatogr., B: Anal. Technol. Biomed. Life Sci.</i> 833 (2), 260–263	2006 Apr 3
Effect of green tea powder (<i>Camellia sinensis</i> L. cv. Benifuuki) particle size on O-methylated EGCG absorption in rats. The Kakegawa Study	Maeda-Yamamoto, M.; Ema, K.; Tokuda, Y.; Monobe, M.; Tachibana, H.; Sameshima, Y.; Kuriyama, S.	<i>Cytotechnology</i> 63 (2), 171–179	2011 Mar
Supplementation of a γ-tocopherol-rich mixture of tocopherols in healthy men protects against vascular endothelial dysfunction induced by postprandial hyperglycemia	Mah, E.; Noh, S. K.; Ballard, K. D.; Park, H. J.; Volek, J. S.; Bruno, R. S.	<i>J. Nutr. Biochem.</i> 24 (1), 196–203	2013 Jan



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Mediterranean diet reduces endothelial damage and improves the regenerative capacity of endothelium	Marin, C.; Ramirez, R.; Delgado-Lista, J.; Yubero-Serrano, E. M.; Perez-Martinez, P.; Carracedo, J.; Garcia-Rios, A.; Rodriguez, F.; Gutierrez-Mariscal, F. M.; Gomez, P.; Perez-Jimenez, F.; Lopez-Miranda, J.	<i>Am. J. Clin. Nutr.</i> 93 (2), 267–274	2011 Feb
Photodiode array (PDA) and other detection methods in HPLC of plant metabolites	Markowski, W.; Waksmundzka-Hajnos, M.	Chapter 13 in <i>High Performance Liquid Chromatography in Phytochemical Analysis</i> , Chromatographic Science Series, Markowski, W., Sherma, J., Eds.; Taylor & Francis Group, LLC: Boca Raton, FL; 331–350	2010 Nov
Determination of water-soluble vitamins in infant milk and dietary supplement using a liquid chromatography on-line coupled to a corona-charged aerosol detector	Márquez-Sillero, I.; Cárdenas, S.; Valcárcel, M.	<i>J. Chromatogr., A.</i> 1313C, 253–258	2013 Oct 25
Sensitive high-performance liquid chromatographic method using coulometric electrode array detection for measurement of phytoestrogens in dried blood spots	Melby, M. K.; Watanabe, S.; Whitten, P. L.; Worthman, C. M.	<i>J. Chromatogr., B: Anal. Technol. Biomed. Life Sci.</i> 826 (1–2), 81–90	2005 Nov 5
Phenolic acids from beer are absorbed and extensively metabolized in humans	Nardini, M.; Natella, F.; Scaccini, C.; Ghiselli, A.	<i>J. Nutr. Biochem.</i> 17 (1), 14–22	2006 Jan
High-performance liquid chromatography analysis of plant saponins: An update 2005-2010	Negi, J. S.; Singh, P.; Pant, G. J.; Rawat, M. S.	<i>Pharmacogn. Rev.</i> 5 (10), 155–158	2011 Jul
Physicochemical effect of pH and antioxidants on mono- and triglutamate forms of 5-methyltetrahydrofolate, and evaluation of vitamin stability in human gastric juice: Implications for folate bioavailability	Ng, X.; Lucock, M.; Veysey, M.	<i>Food Chem.</i> 106 (1), 200–210	2008 Jan
Practical preparation of lacto-N-biose I, a candidate for the bifidus factor in human milk	Nishimoto, M.; Kitaoka, M.	<i>Biosci., Biotechnol., Biochem.</i> 71 (8), 2101-2104	2007 Aug
Hydrophilic interaction liquid chromatography—charged aerosol detection as a straightforward solution for simultaneous analysis of ascorbic acid and dehydroascorbic acid	Nováková, L.; Solichová, D.; Solich, P.	<i>J. Chromatogr., A.</i> 1216 (21), 4574–4581	2009 May 22



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
No effect on adenoma formation in Min mice after moderate amount of flaxseed	Oikarinen, S.; Heinonen, S. M.; Nurmi, T.; Adlercreutz, H.; Mutanen, M.	<i>Eur. J. Nutr.</i> 44 (5), 273–280	2005 Aug
Measurement of isoflavones using liquid chromatography with multi-channel coulometric electrochemical detection	Ouchi, K.; Gamache, P.; Acworth, I.; Watanabe, S.	<i>BioFactors.</i> 22 (1–4), 353–356	2004
Quantitation of clovamide-type phenylpropenoic acid amides in cells and plasma using high-performance liquid chromatography with a coulometric electrochemical detector	Park, J. B.	<i>J. Agric. Food Chem.</i> 53 (21), 8135–8140	2005 Oct 19
Synthesis, HPLC measurement and bioavailability of the phenolic amide amkamide	Park, J. B.	<i>J. Chromatogr. Sci.</i> [Epub ahead of print]	2013 May 27
Synthesis of safflomide and its HPLC measurement in mouse plasma after oral administration	Park, J. B.; Chen, P.	<i>J. Chromatogr., B: Anal. Technol. Biomed. Life Sci.</i> 852 (1–2), 398–402	2007 Jun 1
Determination of lignans in human plasma by liquid chromatography with coulometric electrode array detection	Peñalvo, J. L.; Nurmi, T.; Haajanen, K.; Al-Maharik, N.; Botting, N.; Adlercreutz, H.	<i>Anal. Biochem.</i> 332 (2), 384–393	2004 Sep 15
Supercritical antisolvent fractionation of lignans from the ethanol extract of flaxseed	Perretti, G.; Virgili, C.; Troilo, A.; Marconi, O.; Regnicoli, G. F.; Fantozzi, P.	<i>J. Supercrit. Fluids</i> 75, 94–100	2013 Mar
Analysis of flavonoids in honey by HPLC coupled with coulometric electrode array detection and electrospray ionization mass spectrometry	Petrus, K.; Schwartz, H.; Sontag, G.	<i>Anal. Bioanal. Chem.</i> 400 (8), 2555–2563	2011 Jun
High-dose supplementation with natural α-tocopherol does neither alter the pharmacodynamics of atorvastatin nor its phase I metabolism in guinea pigs	Podszun, M. C.; Grebenstein, N.; Hofmann, U.; Frank, J.	<i>Toxicol. Appl. Pharmacol.</i> 266 (3), 452–458	2013 Feb 1
Application of high-performance liquid chromatography with charged aerosol detection for universal quantitation of undeclared phosphodiesterase-5 inhibitors in herbal dietary supplements	Poplawska, M.; Blazewicz, A.; Bukowska, K.; Fijalek, Z.	<i>J. Pharm. Biomed. Anal.</i> 84, 232–243	2013 Oct
Isolation and analysis of ginseng: advances and challenges	Qi, L.; Wang, C.; Yuan, C.	<i>Nat. Prod. Rep.</i> 28 (3), 467–495	2011 Mar
Folate analysis in complex food matrices: Use of a recombinant Arabidopsis γ-glutamyl hydrolase for folate deglutamylation	Ramos-Parra, P. A.; Urrea-López, R.; Diaz de la Garza, R. I.	<i>Food Res. Int.</i> 54 (1), 177–185	2013 Nov
Optimisation of gradient HPLC analysis of phenolic compounds and flavonoids in beer using a coularray detector	Rehová, L.; Skeríková, V.; Jandera, P.	<i>J. Sep. Sci.</i> 27 (15–16), 1345–1359	2004 Nov
Chiral separation of (+)/(-)-catechin from sulfated and glucuronidated metabolites in human plasma after cocoa consumption	Ritter, C.; Zimmermann, B. F.; Galensa, R.	<i>Anal. Bioanal. Chem.</i> 397 (2), 723–730	2010 May



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Analysis of alkylresorcinols in cereal grains and products using ultrahigh-pressure liquid chromatography with fluorescence, ultraviolet, and CoulArray electrochemical detection	Ross, A. B.	<i>J. Agric. Food Chem.</i> 60 (36), 8954–8962	2012 Sep 12
Rapid and sensitive analysis of alkylresorcinols from cereal grains and products using HPLC-CoulArray-based electrochemical detection	Ross, A. B.; Kochhar, S.	<i>J. Agric. Food Chem.</i> 57 (12), 5187–5193	2009 Jun 24
Analysis of soy isoflavone plasma levels using HPLC with coulometric detection in postmenopausal women	Saracino, M. A.; Raggi, M. A.	<i>J. Pharm. Biomed. Anal.</i> 53 (3), 682–687	2010 Nov 2
A biosynthetic pathway for BE-7585A, a 2-thiosugar-containing angucycline-type natural product	Sasaki, E.; Ogasawara, Y.; Liu, H. W.	<i>J. Am. Chem. Soc.</i> 132 (21), 7405–7417	2010 Jun 2
The senescence-accelerated mouse-prone 8 is not a suitable model for the investigation of cardiac inflammation and oxidative stress and their modulation by dietary phytochemicals	Schiborr, C.; Schwamm, D.; Kocher, A.; Rimbach, G.; Eckert, G. P.; Frank, J.	<i>Pharmacol. Res.</i> 74, 113–120	2013 Aug
Comprehensive impurity profiling of nutritional infusion solutions by multidimensional off-line reversed-phase liquid chromatography × hydrophilic interaction chromatography-ion trap mass-spectrometry and charged aerosol detection with universal calibration	Schiesel, S.; Lämmerhofer, M.; Lindner, W.	<i>J. Chromatogr., A.</i> 1259, 100–110	2012 Oct 12
The effect of α-tocopherol supplementation on training-induced elevation of S100B protein in sera of basketball players	Schulpis, K. H.; Moukas, M.; Parthimos, T.; Tsakiris, T.; Parthimos, N.; Tsakiris, S.	<i>Clin. Biochem.</i> 40 (12), 900–906	2007 Aug
Determination of secoisolariciresinol, lariciresinol and isolariciresinol in plant foods by high performance liquid chromatography coupled with coulometric electrode array detection	Schwartz, H.; Sontag, G.	<i>J. Chromatogr., B: Anal. Technol. Biomed. Life Sci.</i> 838 (2), 78–85	2006 Jul 11
Assessment of probiotic strains ability to reduce the bioaccessibility of aflatoxin M 1 in artificially contaminated milk using an in vitro digestive model	Serrano-Niño, J. C.; Cavazos-Garduño, A.; Hernandez-Mendoza, A.; Applegate, B.; Ferruzzi, M. G.; San Martin-González, M. F.; García, H. S.	<i>Food Control</i> 31 (1), 202–207	2013 May
Intestinal uptake of quercetin-3-glucoside in rats involves hydrolysis by lactase phlorizin hydrolase	Sesink, A. L.; Arts, I. C.; Faassen-Peters, M.; Hollman, P. C.	<i>J. Nutr.</i> 133 (3), 773–776	2003 Mar
Quercetin glucuronides but not glucosides are present in human plasma after consumption of quercetin-3-glucoside or quercetin-4'-glucoside	Sesink, A. L.; O'Leary, K. A.; Hollman, P. C.	<i>J. Nutr.</i> 131 (7), 1938–1941	2001 Jul
Co-administration of quercetin and catechin in rats alters their absorption but not their metabolism	Silberberg, M.; Morand, C.; Manach, C.; Scalbert, A.; Remesy, C.	<i>Life Sci.</i> 77 (25), 3156–3167	2005 Nov 4
Nutritional status is altered in the self-neglecting elderly	Smith, S. M.; Mathews Oliver, S. A.; Zwart, S. R.; Kala, G.; Kelly, P. A.; Goodwin, J. S.; Dyer, C. B.	<i>J. Nutr.</i> 136 (10), 2534–2541	2006 Oct



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Binding of heterocyclic aromatic amines by lactic acid bacteria: results of a comprehensive screening trial	Stidl, R.; Sontag, G.; Koller, V.; Knasmüller, S.	<i>Mol. Nutr. Food Res.</i> 52 (3), 322–329	2008 Mar
Direct separation and detection of biogenic amines by ion-pair liquid chromatography with chemiluminescent nitrogen detector	Sun, J.; Guo, H. X.; Semin, D.; Cheetham, J.	<i>J. Chromatogr., A.</i> 1218 (29), 4689–4697	2011 Jul 22
Rapid purification method for fumonisin B1 using centrifugal partition chromatography	Szekeres, A.; Lorántfy, L.; Bencsik, O.; Kecskeméti, A.; Szécsi, Á.; Mesterházy, Á.; Vágvölgyi, C.	<i>Food Addit. Contam.</i> 30 (1), 147–155	2013
Determination of coenzyme Q10 in over-the-counter dietary supplements by high-performance liquid chromatography with coulometric detection	Tang, P. H.	<i>J. AOAC Int.</i> 89 (1), 35–39	2006 Jan–Feb
α-Tocopherol supplementation restores the reduction of erythrocyte glucose-6-phosphate dehydrogenase activity induced by forced training	Tsakiris, S.; Reclus, G. J.; Parthimos, T.; Tsakiris, T.; Parthimos, N.; Schulpis, K. H.	<i>Pharmacol. Res.</i> 54 (5), 373–379	2006 Nov
Tissue distribution of isoflavones in ewes after consumption of red clover silage	Urpi-Sarda, M.; Morand, C.; Besson, C.; Kraft, G.; Viala, D.; Scalbert, A.; Besle, J. M.; Manach, C.	<i>Arch. Biochem. Biophys.</i> 476 (2), 205–210	2008 Aug 15
Performance evaluation of charged aerosol and evaporative light scattering detection for the determination of ginsenosides by LC	Wang, L.; He, W. S.; Yan, H. X.; Jiang, Y.; Bi, K. S.; Tu, P. F.	<i>Chromatographia</i> 70 (3–4), 603–608	2009 Aug
Catechins are bioavailable in men and women drinking black tea throughout the day	Warden, B. A.; Smith, L. S.; Beecher, G. R.; Balentine, D. A.; Clevidence, B. A.	<i>J. Nutr.</i> 131 (6), 1731–1737	2001 Jun
Identification and quantification of polyphenol phytoestrogens in foods and human biological fluids	Wilkinson, A. P.; Wähälä, K.; Williamson, G.	<i>J. Chromatogr., B: Anal. Technol. Biomed. Life Sci.</i> 777 (1–2), 93–109	2002 Sep 25
Bioavailability and pharmacokinetics of caffeoylquinic acids and flavonoids after oral administration of Artichoke leaf extracts in humans	Wittemer, S. M.; Ploch, M.; Windeck, T.; Müller, S. C.; Drewelow, B.; Derendorf, H.; Veit, M.	<i>Phytomedicine</i> 12 (1–2), 28–38	2005 Jan
Validated method for the determination of six metabolites derived from artichoke leaf extract in human plasma by high-performance liquid chromatography-coulometric-array detection	Wittemer, S. M.; Veit, M.	<i>J. Chromatogr., B: Anal. Technol. Biomed. Life Sci.</i> 793 (2), 367–375	2003 Aug 15
HPLC in natural product analysis: The detection issue	Wolfender, J. L.	<i>Planta Med.</i> 75 (07), 719–734	2009 Jun
Simultaneous determination of isoflavones and bisphenol A in rat serum by high-performance liquid chromatography coupled with coulometric array detection	Yasuda, S.; Wu, P. S.; Hattori, E.; Tachibana, H.; Yamada, K.	<i>Biosci., Biotechnol., Biochem.</i> 68 (1), 51–58	2004 Jan
Impurities from polypropylene microcentrifuge tubes as a potential source of interference in simultaneous analysis of multiple lipid-soluble antioxidants by HPLC with electrochemical detection	Yen, H. C.; Hsu, Y. T.	<i>Clin. Chem. Lab. Med.</i> 42 (4), 390–395	2004 Apr



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Food, Nutrition, Natural Products, and Supplements

Title	Authors	Publication	Publication Date
Simultaneous determination of triterpenoid saponins from <i>pulsatilla koreana</i> using high performance liquid chromatography coupled with a charged aerosol detector (HPLC-CAD)	Yeom, H.; Suh, J. H.; Youm, J. R.; Han, S. B.	<i>Bull. Korean Chem. Soc.</i> 31 (5), 1159–1164	2010
DPPH radical scavenging activities of 31 flavonoids and phenolic acids and 10 extracts of Chinese materia medica	Yuan, Y.; Chen, C.; Yang, B.; Kusu, F.; Kotani, A.	<i>Zhongguo Zhongyao Zazhi</i> 34 (13), 1695–1700	2009 Jul
Determination of residual clenbuterol in pork meat and liver by HPLC with electrochemical detection	Zhang, X. Z.; Gan, Y. R.; Zhao, F. N.	<i>Yaoyue Xuebao</i> 39 (4), 276–280	2004 Apr
Identification of equol producers in a Japanese population by high-performance liquid chromatography with coulometric array for determining serum isoflavones	Zhao, J. H.; Sun, S. J.; Arai, Y.; Oguma, E.; Yamada, K.; Horiguchi, H.; Kayama, F.	<i>Phytomedicine</i> 13 (5), 304–309	2006 May
Simultaneous sampling of volatile and non-volatile analytes in beer for fast fingerprinting by extractive electrospray ionization mass spectrometry	Zhu, L.; Hu, Z.; Gamez, G.; Law, W. S.; Chen, H.; Yang, S.; Chingin, K.; Balabin, R. M.; Wang, R.; Zhang, T.; Zenobi, R.	<i>Anal. Bioanal. Chem.</i> 398 (1), 405–413	2010 Sep
Comparison of various easy-to-use procedures for extraction of phenols from apricot fruits	Zitka, O.; Sochor, J.; Rop, O.; Skalickova, S.; Sobrova, P.; Zehnalek, J.; Beklova, M.; Krska, B.; Adam, V.; Kizek, R.	<i>Molecules</i> 16 (4), 2914–2936	2011 Apr 4





Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Lipids

Title	Authors	Publication	Publication Date
Development of analytical procedures to study changes in the composition of meat phospholipids caused by induced oxidation	Cascone, A.; Eerola, S.; Ritieni, A.; Rizzo, A.	<i>J. Chromatogr., A</i> 1120 (1–2), 211–220	2006 Jul 7
Evaporative light scattering and charged aerosol detector.	Chaminade, P.	Chapter 5. In <i>Hyphenated and Alternative Methods of Detection in Chromatography</i> , Chromatographic Science Series; Shalliker, A., Ed.; Taylor & Francis Group, LLC: Boca Raton, FL.; 145–160	2012
Simple and efficient profiling of phospholipids in phospholipase D-modified soy lecithin by HPLC with charged aerosol detection	Damjanovic, J.; Nakano, H.; Iwasaki, Y.	<i>J. Am. Oil Chem. Soc.</i> 90 (7), 951–957	2013 Jul
Discriminating olive and non-olive oils using HPLC-CAD and chemometrics	de la Mata-Espinosa, P.; Bosque-Sendra, J. M.; Bro, R.; Cuadros-Rodríguez, L.	<i>Anal. Bioanal. Chem.</i> 399 (6), 2083–2092	2011 Feb
Olive oil quantification of edible vegetable oil blends using triacylglycerols chromatographic fingerprints and chemometric tools	de la Mata-Espinosa, P.; Bosque-Sendra, J. M.; Bro, R.; Cuadros-Rodríguez, L.	<i>Talanta</i> 85 (1), 177–182	2011 Jul 15
Quantification of triacylglycerols in olive oils using HPLC-CAD	de la Mata-Espinosa, P.; Bosque-Sendra, J.; Cuadros-Rodríguez, L.	<i>Food Analytical Methods</i> 4 (4), 574–581	2011 Dec
Quantification of pegylated phospholipids decorating polymeric microcapsules of perfluorooctyl bromide by reverse phase HPLC with a charged aerosol detector	Díaz-López, R.; Libong, D.; Tsapis, N.; Fattal, E.; Chaminade, P.	<i>J. Pharm. Biomed. Anal.</i> 48 (3), 702–707	2008 Nov 4
Squalene emulsions for parenteral vaccine and drug delivery	Fox, C. B.	<i>Molecules</i> 14 (9), 3286–3312	2009 Sep 1
Interactions between parenteral lipid emulsions and container surfaces	Gonyon, T.; Tomaso, A.; Kotha, P.; Owen, H.; Patel, D.; Carter, P.; Cronin, J.; Green, J.	<i>PDA J. Pharm. Sci. and Tech.</i> 67 (3), 247–254	2013 May–Jun
Composition analysis of positional isomers of phosphatidylinositol by high-performance liquid chromatography	Iwasaki, Y.; Masayama, A.; Mori, A.; Ikeda, C.; Nakano, H.	<i>J. Chromatogr., A</i> 1216 (32), 6077–6080	2009 Aug 7
Determination of phospholipid and its degradation products in liposomes for injection by HPLC-charged aerosol detection (CAD)	Jiang, Q.; Yang, R.; Mei, X.	<i>Chinese Pharmaceutical Journal (Zhongguo Yaoxue Zazhi, Beijing, China)</i> 42 (23), 1794–1796	2007



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Lipids

Title	Authors	Publication	Publication Date
Rapid quantification of yeast lipid using microwave-assisted total lipid extraction and HPLC-CAD	Khoomrung, S.; Chumnanpuen, P.; Jansa-Ard, S.; Ståhlman, M.; Nookaew, I.; Borén, J.; Nielsen, J.	<i>Anal. Chem.</i> 85 (10), 4912–4919	2013 May 21
A new liquid chromatography method with charge aerosol detector (CAD) for the determination of phospholipid classes. Application to milk phospholipids	Kiełbowicz, G.; Micek, P.; Wawrzenczyk, C.	<i>Talanta</i> 105, 28–33	2013 Feb 15
An LC method for the analysis of phosphatidylcholine hydrolysis products and its application to the monitoring of the acyl migration process	Kiełbowicz, G.; Smuga, D.; Gładkowski, W.; Chojnacka, A.; Wawrzenczyk, C.	<i>Talanta</i> 94, 22–29	2012 May 30
Separation of acylglycerols, FAME and FFA in biodiesel by size exclusion chromatography	Kittirattanapiboon, K.; Krisnangkura, K.	<i>Eur. J. Lipid Sci. Technol.</i> 110 (5), 422–427	2008 Mar 17
Quantitation of triacylglycerols from plant oils using charged aerosol detection with gradient compensation	Lísa, M.; Lynen, F.; Holčápek, M.; Sandra, P.	<i>J. Chromatogr., A.</i> 1176 (1–2), 135–142	2007 Dec 28
Quantitative study of the stratum corneum lipid classes by normal phase liquid chromatography: comparison between two universal detectors	Merle, C.; Laugel, C.; Chaminade, P.; Baillet-Guffroy, A.	<i>J. Liq. Chromatogr. Relat. Technol.</i> 33, 629–644	2010 Mar
The analysis of lipids via HPLC with a charged aerosol detector	Moreau, R. A.	<i>Lipids</i> 41 (7), 727–34	2006 Jul
Lipid analysis via HPLC with a charged aerosol detector	Moreau, R. A.	<i>Lipid Technol.</i> 21 (8–9), 191–194	2009 Oct 23
Extraction and analysis of food lipids	Moreau, R. A.; Winkler-Moser, J. K.	Chapter 6 in <i>Methods of Analysis of Food Components and Additives</i> , Second Edition; Ötles, S., Ed.; Taylor & Francis Group, LLC: Boca Raton, FL.; 115–134	2011 Nov
Aerosol based detectors for the investigation of phospholipid hydrolysis in a pharmaceutical suspension formulation	Nair, L.; Werling, J.	<i>J. Pharm. Biomed. Anal.</i> 49 (1), 95–99	2009 Jan 15
Structure/function relationships of adipose phospholipase A2 containing a cys-his-his catalytic triad	Pang, X. Y.; Cao, J.; Addington, L.; Lovell, S.; Battaile, K. P.; Zhang, Rao, J. L.; Dennis, E. A.; Moise, A. R.	<i>J. Biol. Chem.</i> 287 (42), 35260–35274	2012 Oct 12
Simultaneous assessment of lipid classes and bile acids in human intestinal fluid by solid-phase extraction and HPLC methods	Persson, E.; Löfgren, L.; Hansson, G.; Abrahamsson, B.; Lennernäs, H.; Nilsson, R.	<i>J. Lipid Res.</i> 48 (1), 242–251	2007 Jan



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: HPLC and UHPLC Methods

Lipids

Title	Authors	Publication	Publication Date
The use of charged aerosol detection with HPLC for the measurement of lipids	Plante, M.; Bailey, B.; Acworth, I.	<i>Methods Mol. Biol.</i> (Totowa, NJ, U.S.) 579, 469–482	2009
Comparison between charged aerosol detection and light scattering detection for the analysis of Leishmania membrane phospholipids	Ramos, R. G.; Libong, D.; Rakotomanga, M.; Gaudin, K.; Loiseau, P. M.; Chaminade, P.	<i>J. Chromatogr., A.</i> 1209 (1–2), 88–94	2008 Oct 31
Authentication of geographical origin of palm oil by chromatographic fingerprinting of triacylglycerols and partial least square-discriminant analysis	Ruiz-Samblás, C.; Arrebola-Pascual, C.; Tres, A.; van Ruth, S.; Cuadros-Rodríguez, L.	<i>Talanta.</i> 116, 788–793	2013 Nov 15
Simple and precise detection of lipid compounds present within liposomal formulations using a charged aerosol detector	Schönherr, C.; Touchene, S.; Wilser, G.; Peschka-Süss, R.; Francese, G.	<i>J. Chromatogr., A.</i> 1216 (5), 781–786	2009 Jan 30
Determination of intraluminal individual bile acids by HPLC with charged aerosol detection	Vertzoni, M.; Archontaki, H.; Reppas, C.	<i>J. Lipid Res.</i> 49 (12), 2690–2695	2008 Dec
Neurolipids and the use of a charged aerosol detector	Waraska, J.; Acworth, I.	<i>Am. Biotechnol. Lab.</i> 26 (1), 12–13	2008





Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Technical Collateral: HPLC and UHPLC Methods

Product Number	Technique	Title
AB 119	UV	Rapid Separation of Paclitaxel and Related Compounds in Paclitaxel Injection
AB 134	MS	LC-MS Analysis of Anthocyanins in Bilberry Extract
AB 139	UV	Separation of Schizandrin, Schizandrin A, and Schizandrin B in a Tablet Sample
AB 153	UV	Save the Flavor – Robust Iso- α -Acids Assaying in Beer within Ten Minutes
AB 155	UV	Monitor the Brewing Process with LC-Transformation of Hop alpha-Acids into Beer Iso-alpha-Acids
AN 109	FLD	Determination of Glyphosate by Cation-Exchange Chromatography with Postcolumn Derivatization
AN 156	UV	The Everlasting Paradigm-Keep Beer Tradition or Prevent Beer from a Skunky Off-Flavor?
AN 196	FLD	Determination of Polycyclic Aromatic Hydrocarbons (PAHs) in Edible Oils by Donor-Acceptor Complex Chromatography (DACC)-HPLC with Fluorescent Detection
AN 207	UV	Chromatographic Fingerprinting of <i>Flos Chrysanthema indicis</i> Using HPLC
AN 213	UV/FLD	Determination of Polycyclic Aromatic Hydrocarbons (PAHs) in Tap Water Using on-Line Solid-Phase Extraction Followed by HPLC with UV and Fluorescence Detections
AN 216	UV	Determination of Water- and Fat-Soluble Vitamins in Functional Waters by HPLC with UV-PDA Detection
AN 224	UV	Determination of Melamine in Milk Powder by Reversed-Phase HPLC with UV Detection
AN 232	UV	Determination of Anthraquinones and Stilbenes in Giant Knotweed Rhizome by HPLC with UV Detection
AN 236	UV	Determination of Iodide and Iodate in Seawater and Iodized Table Salt by HPLC-UV Detection
AN 245	UV	Fast Analysis of Dyes in Foods and Beverages
AN 251	UV	Determination of Water- and Fat-Soluble Vitamins in Nutritional Supplements by HPLC with UV Detection
AN 252	UV	HPLC Assay of Water-Soluble Vitamins, Fat-Soluble Vitamins, and a Preservative in Dry Syrup Multivitamin Formulation
AN 261	UV	Sensitive Determination of Microcystins in Drinking and Environmental Waters
AN 264	UV	Fast Determination of Anthocyanins in Pomegranate Juice
AN 266	FLD	Determination of Sialic Acids Using UHPLC with Fluorescence Detection
AN 272	FLD	Faster Yet Sensitive Determination of N-Methylcarbamates in Rice, Potato, and Corn by HPLC
AN 275	UV	Sensitive Determination of Catechins in Tea by HPLC
AN 287	UV	Two-Dimensional HPLC Combined with On-Line SPE for Determination of Sudan Dyes I-IV in Chili Oil
AN 292	UV	Determination of Aniline and Nitroanilines in Environmental and Drinking Waters by On-Line SPE
AN 293	CAD and UV	Steviol Glycoside Determination by HPLC with Charged Aerosol and UV Detections Using the Acclaim Trinity P1 Column
AN 299	UV	HPLC Analysis of Six Active Components of <i>Caulis Ionicerae</i> Using a Phenyl-1 Column
AN 1008	UV	Determination of Nitidine Chloride, Toddalolactone, and Chelerythrine Chloride by HPLC



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Technical Collateral: HPLC and UHPLC Methods

Product Number	Technique	Title
AN 1020	EC, UV	Chalcinoids and Bitter Acids in Beer by HPLC with UV and ECD
AN 1023	UV	Determination of Sudan Dyes I-IV in Curry Paste
AN 1026	CAD	Fatty Acid Esters at Low Nanogram Levels
AN 1027	CAD	Ginseng
AN 1028	CAD	Ginkgo biloba
AN 1029	CAD	Black Cohosh
AN 1030	CAD	Soy Saponins
AN 1032	CAD	Unsaturated Fatty Acid: Arachidonic, Linoleic, Linolenic and Oleic Acids
AN 1033	CAD	Corn Syrup
AN 1034	CAD	Honey Sugars
AN 1035	CAD	Phenolic Acids
AN 1036	CAD	Water-Soluble Antioxidants: Ascorbic Acid, Glutathione and Uric Acid
AN 1037	CAD	Artificial Sweeteners-Global Method
AN 1039	CAD	Simultaneous Measurement of Glycerides (Mono-, Di- and Triglycerides) and Free Fatty Acids in Palm Oil
AN 1040	CAD	Analysis of Commercially Available Products Containing Stevia
AN 1041	CAD	Phytosterols
AN 1042	UV	Rapid Separation of Anthocyanins in Cranberry and Bilberry Extracts Using a Core-Shell Particle Column
AN 1045	UV	Determination of Phthalates in Drinking Water by UHPLC with UV Detection
AN 1046	UV	Determination of Phenylurea Compounds in Tap Water and Bottled Green Tea
AN 1055	CAD	Determination of Virginiamycin, Erythromycin, and Penicillin in Dried Distillers Grains with Solubles
AN 1063	ECD	Targeted Analyses of Secondary Metabolites in Herbs, Spices, and Beverages Using a Novel Spectro-Electro Array Platform
AN 1064	ECD	Product Authentication and Adulteration Determination Using a Novel Spectro-Electro Array Platform
AN 1067	UV	Determination of Carbendazim in Orange Juice
AN 1069	UV	Two-Dimensional HPLC Determination of Water-Soluble Vitamins in a Nutritional Drink
AN 1070	UV	Determination of Inositol Phosphates in Dried Distillers Grains and Solubles
AN 20583	UV	Determination of Catechins and Phenolic Acids in Red Wine by Solid Phase Extraction and HPLC
AN 20610	UV	Fast Analysis of Coffee Bean Extracts Using a Solid Core HPLC Column
AN 20663	CAD	Comparative Analysis of Cooking Oils Using a Solid Core HPLC Column
AN 20847	CAD	Analysis of a Sports Beverage for Electrolytes and Sugars Using Multi-Mode Chromatography with Charged Aerosol Detection



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Technical Collateral: HPLC and UHPLC Methods

Product Number	Technique	Title
AN 70158	CAD	Novel Universal Approach for the Measurement of Natural Products in a Variety of Botanicals and Supplements
AN 70277	CAD	Simultaneous Analysis of Glycerides and Fatty Acids in Palm Oil
AU 144	UV	Determination of Hexavalent Chromium in Drinking Water Using Ion Chromatography
AU 170	UV	Fast Determination of Vanillin and its Synthesis Precursor by HPLC
AU 182	CAD	Measuring Lactose in Milk: A Validated Method
AU 184	CAD, UV	Mogroside V Determination by HPLC with Charged Aerosol and UV Detections
CAN 106	UV	Determination of the Punicalagins Found in Pomegranate by High Performance Liquid Chromatography
CAN 111	CAD	Determination of Triterpenes in <i>Centella asiatica</i> (Gotu Kola) by HPLC-CAD
CAN 112	CAD	Determination of Ginsenosides in Panax ginseng by HPLC-CAD
CAN 115	FLD	Clean-Up and Analysis of Aflatoxins and Ochratoxin A in Herbs and Spices
LPN 2062	MS	Profiling Analysis of 15 Prominent Naturally Occurring Phenolic Acids by LC-MS
LPN 2069	FLD	Fast and Effective Determination of Aflatoxins in Grains or Food Using Accelerated Solvent Extraction followed by HPLC
LPN 2421	UV	Achieving Maximum Productivity by Combining UHPLC with Advanced Chromatographic Techniques
LPN 2818	CAD	Analysis of Fat-Soluble Vitamins and Antioxidants in Supplements by RP-HPLC
LPN 2870	FLD	Benefits of High-Speed Wavelength Switching in UHPLC Methods Using Fluorescence Detection
LPN 2930	CAD	Determination of the Composition of Natural Products by HPLC with Charged Aerosol Detection
LPN 2923	CAD	Simple and Direct Analysis of Falcarinol and Other Polyacetylenic Oxylipins in Carrots by Reversed-Phase HPLC and Charged Aerosol Detection
LPN 2931	CAD	Quantification of Underivatized Omega-3 and Omega-6 Fatty Acids in Foods by HPLC CAD
LPN 2932	ECD	A Versatile Detector for the Sensitive and Selective Measurement of Numerous Fat-Soluble Vitamins and Antioxidants in Human Plasma and Plant Extracts
LPN 2934	CAD	Sensitive Analysis of Commonly Used Artificial and Natural Sweeteners Including Stevia and Their Impurities and Degradation Products
LPN 2991	CAD	Evaluation of Methods for the Characterization and Quantification of Polysorbates and Impurities Along with Other Surfactants and Emulsifiers Used in the Food and Pharmaceutical Industries
PN 70026	CAD	Carbohydrate Analysis Using PAD, FLD, CAD and MS Detectors
PN 70037	CAD	Sensitive HPLC Method for Triterpenoid Analysis Using Charged Aerosol Detection with Improved Resolution
PN 70055	CAD	Direct Analysis of Surfactants using HPLC with Charged Aerosol Detection
PN 70138	UV	Rapid Determination of Polyphenol Antioxidants in Green Tea and Cranberry Extract Using Core Shell Columns
PN 70538	CAD	Analysis of Silicone Oils by HPLC-CAD
PN 70540	CAD, ECD	Profiling <i>Hoodia</i> Extracts by HPLC with CAD, ECD, Principal Component Analysis

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



References



Ion Chromatography References



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Technical Collateral: Ion Chromatography Methods

Product Number	Technique	Title
AB 127	IC-PAD	Determination of Carbohydrates in Fruit Juice Using Capillary High-Performance Anion-Exchange Chromatography
AB 135	IC-SC	Determination of Anions and Organic Acids in Brewed Coffee Samples Using Capillary IC
AB 137	IC-SC	Determination of Inorganic and Organic Acids in Apple and Orange Juice Samples Using Capillary IC
AN 25	IC-SC	Determination of Inorganic Ions and Organic Acids in Non-Alcoholic Carbonated Beverages
AN 37	IC-PAD	Determination of Iodide and Iodate in Soy- and Mil-Based Infant Formulas
AN 46	IC-PAD	Ion Chromatography: A Versatile Technique for the Analysis of Beer
AN 54	IC-PAD	Determination of Total and Free Sulfite in Foods and Beverages
AN 67	IC-PAD	Determination of Plant-Derived Neutral Oligo- and Polysaccharides
AN 81	IC-SC	Ion Chromatographic Determination of Oxyhalides and Bromide at Trace Level Concentrations in Drinking Water Using direct Injection
AN 82	IC-PAD	Analysis of Fruit Juice Adulterated with Medium Invert Sugar from Beets
AN 87	IC-PAD	Determination of Sugar Alcohols in Confections and Fruit Juices by High-Performance Anion-Exchange Chromatography with Pulsed Amperometric Detection
AN 101	IC-SC	Trace Level Determination of Bromate in Ozonated Drinking Water Using Ion Chromatography
AN 112	IC-UV	Determination of Nitrate and Nitrite in Meat Using High-Performance Anion-Exchange Chromatography
AN 121	IC-SC	Analysis of Low Concentrations of Perchlorate in Drinking Water and Ground Water by Ion Chromatography
AN 123	IC-SC	Determination of Inorganic Anions and Organic Acids in Fermentation Broths
AN 133	IC-SC	Determination of Inorganic Anions in Drinking Water by Ion Chromatography
AN 136	IC-SC and IC-UV	Determination of Inorganic Oxyhalide Disinfection Byproduct Anions and Bromide in Drinking Water Using Ion Chromatography with the Addition of a Postcolumn Reagent for Trace Bromate Analysis
AN 140	IC-SC	Fast Analysis of Anions in Drinking Water by Ion Chromatography
AN 143	IC-SC	Determination of Organic Acids in Fruit Juices
AN 149	IC-SC	Determination of Chlorite, Bromate, Bromide, and Chlorate in Drinking Water by Ion Chromatography with an On-Line-Generated Postcolumn Reagent for Sub- $\mu\text{g/L}$ Bromate Analysis
AN 150	IC-PAD	Determination of Amino Acids in Cell Cultures and Fermentation Broths
AN 154	IC-SC	Determination of Inorganic Anions in Environmental Waters Using a Hydroxide-Selective Column
AN 155	IC-PAD	Determination of Trans-Galactooligosaccharides in Foods by AOAC Method 2001.02



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Technical Collateral: Ion Chromatography Methods

Product Number	Technique	Title
AN 165	IC-SC	Determination of Benzoate in Liquid Food Products by Reagent-Free Ion Chromatography
AN 167	IC-SC	Determination of Trace Concentrations of Oxyhalides and Bromide in Municipal and Bottled Waters Using a Hydroxide-Selective Column with a Reagent-Free Ion Chromatography System
AN 168	IC-UV	Determination of Trace Concentrations of Disinfection By-Product Anions and Bromide in Drinking Water Using Reagent-Free Ion Chromatography Followed by Postcolumn Addition of Iol-Dianisidine for Trace Bromate Analysis
AN 169	IC-SC	Rapid Determination of Phosphate and Citrate in Carbonated Soft Drinks Using a Reagent-Free Ion Chromatography System
AN 172	IC-SC	Determination of Azide in Aqueous Samples by Ion Chromatography with Suppressed Conductivity Detection
AN 173	IC-PAD	Direct Determination of Cyanide in Drinking Water by Ion Chromatography with Pulsed Amperometric Detection (PAD)
AN 178	IC-SC	Improved Determination of Trace Concentrations of Perchlorate in Drinking Water Using Preconcentration with Two-Dimensional Ion Chromatography and Suppressed Conductivity Detection
AN 182	IC-SC and IC-PAD	Determination of Biogenic Amines in Alcoholic Beverages by Ion Chromatography with Suppressed Conductivity and Integrated Pulsed Amperometric Detections
AN 183	IC-SC and IC-PAD	Determination of Biogenic Amines in Fermented and Non-Fermented Foods Using Ion Chromatography with Suppressed Conductivity and Integrated Pulsed Amperometric Detections
AN 187	IC-SC	Determination of sub- $\mu\text{g/L}$ Bromate in Municipal Waters Using Preconcentration with Two-Dimensional Ion Chromatography and Suppressed Conductivity Detection
AN1 88	IC-PAD	Determination of Glycols and Alcohols in Fermentation Broths Using Ion-Exclusion Chromatography and Pulsed Amperometric Detection
AN 197	IC-PAD	Determination of Glucosamine in Dietary Supplements Using HPAE-PAD
AN 227	ICE-PAD	Determination of Total Cyanide in Municipal Wastewater and Drinking Water Using Ion-Exclusion Chromatography with Pulsed Amperometric Detection (ICE-PAD)
AN 248	IC-PAD	Determination of Lactose in Lactose-Free Milk Products by High-Performance Anion-Exchange Chromatography with Pulsed Amperometric Detection
AN 253	IC-PAD	HPAE-PAD Determination of Infant Formula Sialic Acids
AN 270	IC-PAD	Determination of Hydroxymethylfurfural in Honey and Biomass
AN 273	IC-SC	Determination of Organic Acids in Fruit Juices and Wines by High-Pressure IC
AN 279	IC-SC	Time Savings and Improved Reproducibility of Nitrate and Nitrite Ion Chromatography Determination in Milk Samples
AN 280	IC-PAD	Carbohydrates in Coffee: AOAC Method 995.13 vs a New Fast Ion Chromatography Method
AN 295	IC-SC	Determination of Phytic Acid in Soybeans and Black Sesame Seeds
AN 1007	IC-SC	Determination of Mono-, Di-, and Triphosphates and Citrate in Shrimp by Ion Chromatography



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Technical Collateral: Ion Chromatography Methods

Product Number	Technique	Title
AN 1044	IC-SC	Determination of Anions in Dried Distillers Grains with Solubles
AN 1068	IC-SC	Determination of Organic Acids in Fruit Juices and Wines by High-Pressure IC
AU 132	IC-UV	Determination of Nitrite and Nitrate in drinking Water by Ion Chromatography with Direct UV Detection
AU 144	IC-UV	Determination of Hexavalent Chromium in Drinking Water Using Ion Chromatography
AU 148	IC-SC	Determination of Perchlorate in Drinking Water Using Reagent-Free Ion Chromatography
AU 150	IC-PAD	Determination of Plant-Derived Neutral Oligo- and Polysaccharides Using the CarboPac PA200
AU 151	IC-PAD	Determination of Sucralose in Reduced- Carbohydrate Colas using High-Performance Anion-Exchange Chromatography with Pulsed Amperometric Detection
AU 189	IC-SC	Determination of Choline in Infant Formula and Other Food Samples by IC
LPN 2982	IC-SC	Determination of Inorganic Anions and Organic Acids in Beverages Using a Capillary IC on a Monolith Anion-Exchange Column
PN 70743	IC-SC	Determination of Perchlorate Levels in Food and Soil Samples Using Accelerated Solvent Extraction and Ion Chromatography
TN 20	IC-PAD	Analysis of Carbohydrates by High-Performance Anion-Exchange Chromatography with Pulsed Amperometric Detection (HPAE-PAD)
TN 126	IC-SC	Determination of Organic Acids in Beer Samples Using a High-Pressure Ion Chromatography System
TN 135	IC-PAD	Determinations of Monosaccharides and Disaccharides in Beverages by Capillary HPAE-PAD

Table of Contents

Introduction

Analytical Technologies

Targeted Methods

Coffee

Controlled Substances

Fruit Juice

Melamine

Mixtures of Complex Materials

Non-targeted Methods

Herb Blending

Beer

Wine

Fruit Juice

Tea

Botanical Supplements

Tryptophan Supplements

Olive Oil

References



References



Sample Preparation References



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: Sample Preparation Methods

Title	Authors	Publication	Publication Date
Accelerated, microwave-assisted, and conventional solvent extraction methods affect anthocyanin composition from colored grains	Abdel-Aal el-SM; Akhtar, H.; Rabalski, I.; Bryan, M.	<i>J. Food Sci.</i> 79 (2), C138–46	2014 Feb
Multiresidue method for the analysis of pesticide residues in fruits and vegetables by accelerated solvent extraction and capillary gas chromatography	Adou, K.; Bontoyan, W. R.; Sweeney, P. J.	<i>J. Agric. Food Chem.</i> 49 (9), 4153–4160	2001 Sep
The development of an optimized sample preparation for trace level detection of 17α-ethinylestradiol and estrone in whole fish tissue	Al-Ansari, A. M.; Saleem, A.; Kimpe, L. E.; Trudeau, V. L.; Blais, J. M.	<i>J. Chromatogr. B Analyt. Technol. Biomed. Life Sci.</i> 879 (30), 3649–52	2011 Nov
Determination of polyphenolic profiles of basque cider apple varieties using accelerated solvent extraction	Alonso-Salces, R. M.; Korta, E.; Barranco, A.; Berrueta, L.A.; Gallo, B.; Vicent, F.	<i>J. Agric. Food Chem.</i> 49 (8), 3761–376	2001
Pressurized liquid extraction for the determination of polyphenols in apple	Alonso-Salces, R. M.; Korta, E.; Barranco, A.; Berrueta, L. A.; Gallo, B.; Vicente, F.;	<i>J. Chromatogr., A.</i> 933 (1–2), 37–43	2001 Nov
Methods for extraction and determination of phenolic acids in medicinal plants: a review	Arceusz, A.; Wesolowski, M.; Konieczynski, P.	<i>Nat. Prod. Commun.</i> 8 (12), 1821–9	2013 Dec
Study of an accelerated solvent extraction procedure for the determination of acaricide residues in honey by high-performance liquid chromatography-diode array detector	Bakkali, A.; Korta, E.; Berrueta, L. A.	<i>J. Food Protection</i> 65 (1), 161–166	2002
Pressurized liquid extraction of medicinal plants	Benthin, B.; Danz, H.; Hamburger, M.	<i>J. Chromatogr., A.</i> 837 (1-2), 211–9	1999 Apr
Comparison of the chemical composition of extracts from <i>Scutellaria lateriflora</i> using accelerated solvent extraction and supercritical fluid extraction versus standard hot water or 70% ethanol extraction	Bergeron, C.; Gafner, S.; Clausen, E.; Carrier, D. J.	<i>J. Agric. Food Chem.</i> 53 (8), 3076–80	2005 Apr
Polybrominated diphenyl ethers (PBDEs) in Mediterranean mussels (<i>Mytilus gallo-provincialis</i>) from selected Apulia coastal sites evaluated by GC-HRMS	Bianco, G.; Novario, G.; Anzilotta, G.; Palma, A.; Mangone, A.; Cataldi, T. R.	<i>J. Mass Spectrom.</i> 45 (9), 1046–55	2010 Sep
Free and bound phenolic compounds in barley (<i>Hordeum vulgare</i> L.) flours. evaluation of the extraction capability of different solvent mixtures and pressurized liquid methods by micellar electrokinetic chromatography and spectrophotometry	Bonoli, M.; Marconi, E.; Caboni, M. F.	<i>J. Chromatogr., A.</i> 19; 1057 (1-2), 1–12	2004 Nov
Pressurized liquid extraction of lipids for the determination of oxysterols in egg-containing food	Boselli, E.; Velazco, V.; Caboni, M. F.; Lercker, G.	<i>J. Chromatogr., A.</i> 11; 917 (1-2), 239–44	2001 May
Optimisation of accelerated solvent extraction of cocaine and benzoylecgonine from coca leaves	Brachet, A.; Rudaz, S.; Mateus, L.; Christen, P.; Veuthey, J-L.	<i>J. Sep. Sci.</i> 24 (10-11), 865–873	2001 Nov



Peer Reviewed Journals: Sample Preparation Methods

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Title	Authors	Publication	Publication Date
Multi-residue determination of 130 multiclass pesticides in fruits and vegetables by gas chromatography coupled to triple quadrupole tandem mass spectrometry	Cervera, M.I.; Medina, C.; Portolés, T.; Pitarch, E.; Beltrán, J.; Serrahima, E.; Pineda, L.; Muñoz, G.; Centrich, F.; Hernández, F.	<i>Anal. Bioanal. Chem.</i> 397 (7), 2873–91	2010 Aug
Influence of extraction methodologies on the analysis of five major volatile aromatic compounds of citronella grass (<i>Cymbopogon nardus</i>) and lemongrass (<i>Cymbopogon citratus</i>) grown in Thailand	Chanthai, S.; Prachakoll, S.; Ruangviriyachai, C.; Luthria, D. L.	<i>J. AOAC Int.</i> 95 (3), 763–72	2012 May-Jun
Accelerated solvent extraction of vitamin K₁ in medical foods in conjunction with matrix solid-phase dispersion	Chase, G. W.; Thompson, B.	<i>J. AOAC Int.</i> 83 (2), 407–10	2000
Development of a liquid chromatography-tandem mass spectrometry with pressurized liquid extraction method for the determination of benzimidazole residues in edible tissues	Chen, D.; Tao, Y.; Zhang, H.; Pan, Y.; Liu, Z.; Huang, L.; Wang, Y.; Peng, D.; Wang, X.; Dai, M.; Yuan, Z.	<i>J. Chromatogr. B Analyt. Technol. Biomed. Life Sci.</i> 879 (19), 1659–67	2011 Jun
Determination of 88 pesticide residues in tea using gas chromatography-tandem mass spectrometry	Chen, H.; Liu, X.; Wang, Q.; Jiang, Y.	<i>Se Pu.</i> 29 (5), 409–16	2011 May
Optimization of accelerated solvent extraction for the determination of chlorinated pesticides from animal feed	Chen, S.; Gfrerer, M.; Lankmayr, E.; Quan, X.; Yang, F.	<i>Chromatographia</i> 58, 631–636	2003
Uptake of oxytetracycline, sulfamethoxazole and ketoconazole from fertilised soils by plants	Chitescu, C. L.; Nicolau, A. I.; Stolker, A. A.	<i>Food Addit. Contam. Part A Chem. Anal. Control Expo. Risk Assess.</i> 30 (6), 1138–46	2013
Ultrasonic or accelerated solvent extraction followed by U-HPLC-high mass accuracy MS for screening of pharmaceuticals and fungicides in soil and plant samples	Chitescu, C. L.; Oosterink, E.; de Jong, J.; Stolker, A. A.	<i>Talanta</i> 2012; 88, 653–62	2011 Jan
Evaluation of analytical methods for determining pesticides in baby foods and adult duplicate-diet samples	Chuang, J. C.; Hart, K.; Chang, J. S.; Boman, L. E.; Van Emon, J. M.; Reed, A. W.	<i>Anal. Chim. Acta.</i> 444 (1), 87–95	2001 Oct
Comparison of extraction techniques and modeling of accelerated solvent extraction for the authentication of natural vanilla flavors	Cicchetti, E.; Chaintreau, A.	<i>J. Sep. Sci.</i> 32 (11), 1957–64	2009 Jun
Development of a fast and convenient method for the isolation of triterpene saponins from <i>Actaea racemosa</i> by high-speed countercurrent chromatography coupled with evaporative light scattering detection	Cicek, S. S.; Schwaiger, S.; Ellmerer, E. P.; Stuppner, H.	<i>Planta. Med.</i> 76 (5), 467–73	2010 Mar
Extraction of bitter acids from hops and hop products using pressurized solvent extraction (PSE)	Culík, J.; Jurková, M.; Horák, T.; Cejka, P.; Kellner, V.; Dvorák, J.; Karásek, P.; Roth, M.	<i>J. Inst. Brew.</i> 115 (3), 220–225	2009
Comparison of methods for extraction of flavanones and xanthenes from the root bark of the osage orange tree using liquid chromatography	da Costa, C. T.; Margolis, S. A.; Benner, Jr. B.A.; Horton, D.	<i>J. Chromatogr., A.</i> 831 (2), 167–178	1999 Jan
Pressurized liquid extraction prior to liquid chromatography with electrochemical detection for the analysis of vitamin E isomers in seeds and nuts	Delgado-Zamarreño, M. M.; Bustamante-Rangel, M.; Sánchez-Pérez, A.; Carabias-Martínez, R.	<i>J. Chromatogr., A.</i> 12; 1056 (1-2), 249–52	2004 Nov



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: Sample Preparation Methods

Title	Authors	Publication	Publication Date
Pressurized fluid extraction of carotenoids from <i>Haematococcus pluvialis</i> and <i>Dunaliella salina</i> and kavalactones from <i>Piper methysticum</i>	Denery, J. R.; Dragull, K.; Tang, C. S.; Li, Q. X.	<i>Anal. Chim. Acta.</i> 501 (2), 175–181	2004 Jan
Development and comparison of two multiresidue methods for the analysis of 17 mycotoxins in cereals by liquid chromatography electrospray ionization tandem mass spectrometry	Desmarchelier, A.; Oberson, J. M.; Tella, P.; Gremaud, E.; Seefelder, W.; Mottier, P.	<i>J. Agric. Food Chem.</i> 58 (13), 7510–9	2010 Jul
Identification, extraction and quantification of the synthetic cannabinoid JWH-018 from commercially available herbal marijuana alternatives	Dunham, S. J.; Hooker, P. D.; Hyde, R. M.	<i>Forensic Sci. Int.</i> 223 (1-3), 241–4	2012 Nov
Evaluation of polyphenol contents in differently processed apricots using accelerated solvent extraction followed by high-performance liquid chromatography-diode array detector	Erdogan, S.; Erdemoglu, S.	<i>Int. J. Food Sci. Nutr.</i> 62 (7), 729–39	2011 Nov
Determination of 2,4,6-trichloroanisole and guaiacol in cork stoppers by pressurised fluid extraction and gas chromatography–mass spectrometry	Ezquerro, Ó.; Garrido-López, Á.; Tena, M. T.	<i>J. Chromatogr., A.</i> 1102 (12), 18–24	2006 Jan
Multiwalled carbon nanotubes as matrix solid-phase dispersion extraction absorbents to determine 31 pesticides in agriculture samples by gas chromatography–mass spectrometry	Fang, G.; Min, G.; He, J.; Zhang, C.; Qian, K.; Wang, S.	<i>J. Agric. Food Chem.</i> 57 (8), 3040–5	2009 Apr
High-anthocyanin strawberries through cultivar selection	Fredericks, C. H.; Fanning, K. J.; Gidley, M. J.; Netzel, G.; Zabar, D.; Herrington, M.; Netzel, M.	<i>J. Sci. Food Agric.</i> 93 (4), 846–52	2013 Mar
Optimal extraction and fingerprint analysis of <i>Cnidii fructus</i> by accelerated solvent extraction and high performance liquid chromatographic analysis with photodiode array and mass spectrometry detections	Gao, F.; Hu, Y.; Ye, X.; Li, J.; Chen, Z.; Fan, G.	<i>Food Chem.</i> 1 141 (3), 1962–71	2013 Dec
Simultaneous analysis of seven alkaloids in <i>Coptis-evodia</i> herb couple and Zuojin pill by UPLC with accelerated solvent extraction	Gao, X.; Yang, X. W.; Marriott, P. J.	<i>J. Sep. Sci.</i> 33 (17-18), 2714–22	2010 Sep
Determination of chromones in <i>Dysophylla stellata</i> by HPLC: method development, validation and comparison of different extraction methods	Gautam, R.; Srivastava, A.; Jachak, S. M.	<i>Nat. Prod. Commun.</i> 5 (4), 555–8	2010 Apr
Comparison of different extraction techniques for the determination of chlorinated pesticides in animal feed	Gfrerer, M.; Chen, S.; Lankmayr, E.; Xie, Q.; Yang, F.	<i>Anal. Bioanal. Chem.</i> 378 (7), 1861–1867	2004
Speciation analysis of selenium compounds in yeasts using pressurised liquid extraction and liquid chromatography–microwave-assisted digestion–hydride generation–atomic fluorescence spectrometry	Gómez-Ariza, J. L.; Caro de la Torre, M. A.; Giráldez, I.; Morales, E.	<i>Anal. Chim. Acta.</i> 524, (1–2), 305–314	2004 Oct



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: Sample Preparation Methods

Title	Authors	Publication	Publication Date
Multianalysis of 35 mycotoxins in traditional Chinese medicines by ultra-high-performance liquid chromatography-tandem mass spectrometry coupled with accelerated solvent extraction	Han, Z.; Ren Y.; Zhu, J.; Cai, Z.; Chen, Y.; Luan, L.; Wu, Y.	<i>J. Agric. Food Chem.</i> 60 (33), 8233-47.	2012 Aug
Pressurized liquid extraction-capillary electrophoresis-mass spectrometry for the analysis of polar antioxidants in rosemary extracts	Herrero, M.; Arráez-Román, D.; Segura A.; Kenndler, E.; Gius, B.; Raggid, M. A.; Ibáñez, E.; Cifuentes, A.	<i>J. Chromatogr., A</i> 1084 (1-2), 54-62.	2005 Aug
Accelerated solvent extraction of alkylresorcinols in food products containing uncooked and cooked wheat	Holt, M D.; Moreau, R A.; DerMarderosian, A.; McKeown, N.; Jacques, P. F.	<i>J. Agric. Food Chem.</i> 60 (19), 4799-802	2012 May
Application of response surface methodology to optimize pressurized liquid extraction of antioxidant compounds from sage (<i>Salvia officinalis</i> L.), basil (<i>Ocimum basilicum</i> L.) and thyme (<i>Thymus vulgaris</i> L.)	Hossain, M. B.; Brunton, N. P.; Martin-Diana, A. B.; Barry-Ryan, C.	<i>Food Funct.</i> 1(3), 269-77	2010 Dec
A review of modern sample-preparation techniques for the extraction and analysis of medicinal plants	Huie, C. W.	<i>Anal. Bioanal. Chem.</i> 373 (1-2), 23-30.	2002 May
Polychlorinated dioxins, furans, and biphenyls, and polybrominated diphenyl ethers in a U.S. meat market basket and estimates of dietary intake	Huwe, J. K.; Larsen, G. L.	<i>Environ. Sci. Technol.</i> 39 (15), 5606-5611	2005
Study of the effect of sample preparation and cooking on the selenium speciation of selenized potatoes by HPLC with ICP-MS and electrospray ionization MS/MS	Infante, H. G.; Borrego, A. A.; Peachey, E.; Hearn, R.; O'Connor, G.; Barrera, T G.; Ariza, J. L.	<i>J. Agric. Food Chem.</i> 57(1), 38-45.	2009 Jan
Pentacyclic triterpene distribution in various plants – rich sources for a new group of multi-potent plant extracts	Jäger, S.; Trojan, H.; Kopp, T.; Laszczyk, M. N.; Scheffler, A.	<i>Molecules.</i> 14 (6), 2016-31.	2009 Jun
Comprehensive multiresidue method for the simultaneous determination of 74 pesticides and metabolites in traditional Chinese herbal medicines by accelerated solvent extraction with high-performance liquid chromatography/tandem mass spectrometry	Jia, Z.; Mao, X.; Chen, K.; Wang, K.; Ji S.	<i>J. AOAC Int.</i> ; 93(5), 1570-88.	2010 Sep-Oct
Gas chromatography-mass spectrometry (GC-MS) method for the determination of 16 European priority polycyclic aromatic hydrocarbons in smoked meat products and edible oils	Jira, W.; Ziegenhals, K.; Speer, K.	<i>Food Addit. Contam. Part A Chem. Anal. Control Expo. Risk Assess.</i> 25 (6), 704-13.	2008 Jun
Assessing pressurized liquid extraction for the high-throughput extraction of marine-sponge-derived natural products	Johnson, T. A.; Morgan, M. V.; Aratow, N. A.; Estee, S. A.; Sashidhara, K. V.; Loveridge, S. T.; Segraves, N L.; Crews, P.	<i>J. Nat. Prod.</i> 73 (3), 359-64.	2010 Mar
Lipophilic stinging nettle extracts possess potent anti-inflammatory activity, are not cytotoxic and may be superior to traditional tinctures for treating inflammatory disorders	Johnson, T. A.; Sohn, J.; Inman, W. D.; Bjeldanes, L. F.; Rayburn, K.	<i>Phytomedicine</i> 20(2), 143-7.	2013 Jan
Effects of solvent and temperature on pressurized liquid extraction of anthocyanins and total phenolics from dried red grape skin	Ju Z. Y.; Howard, L. R.	<i>J. Agric. Food Chem.</i> 51 (18), 5207-13.	2003 Aug



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: Sample Preparation Methods

Title	Authors	Publication	Publication Date
Accelerated solvent extraction of ochratoxin A from rice samples	Juan, C.; González, L.; Soriano, J. M.; Moltó, J. C.; Mañes, J.	<i>J. Agric. Food Chem.</i> 53 (24), 9348–9351	2005
Accelerated solvent extraction of paclitaxel and related compounds from the bark of <i>Taxus cuspidate</i>	Kawamura, F.; Kikuchi, Y.; Ohira, T.; Yatagai, M.	<i>J. Nat. Prod.</i> 62 (2), 244–7.	1999 Feb
Determination of polybromodiphenyl ethers (PBDEs) in milk cream by gas chromatography-mass spectrometry	Kinani, S.; Bouchonnet, S.; Abjean, J.; Campargue, C.	<i>Food Addit. Contam. Part A Chem. Anal. Control Expo. Risk Assess.</i> 25 (8), 1007–14	2008 Aug
Determination of isoflavones in soy bits by fast column high-performance liquid chromatography coupled with UV-visible diode-array detection	Klejduš, B.; Mikelová, R.; Petřlová, J.; Potešil, D.; Adam, V.; Stiborová, J.; Hodek, P.; Vacek, J.; Kizek, R.; Kubán, V.	<i>J. Chromatogr., A.</i> 1084 (1–2), 19, 71–79	2005 Aug
Accelerated solvent extraction of lignin from <i>Aleurites moluccana</i> (candlenut) nutshells	Klein, A. P.; Beach, E. S.; Emerson, J. W.; Zimmerman, J. B.	<i>J. Agric. Food Chem.</i> 58 (18), 10045–8	2010 Sep
Application of TLC method with video scanning in estimation of daily dietary intake of specific flavonoids – preliminary studies	Koch, W.; Kukuła-Koch, W.; Marzec, Z.; Marc, D.	<i>Acta Pol. Pharm.</i> 70 (4), 611–20	2013 Jul-Aug
Evaluation of a fibrous cellulose drying agent in supercritical fluid extraction and pressurized liquid extraction of diverse pesticides	Lehotay, S. J.; Lee, C. H.	<i>J. Chromatogr., A.</i> 785 (1–2), 313–27	1997 Oct
Application of accelerated solvent extraction to the investigation of saikosaponins from the roots of <i>Bupleurum falcatum</i>	Li, W.; Liu, Z.; Wang, Z.; Chen, L.; Sun, Y.; Hou, J.; Zheng, Y.	<i>J. Sep. Sci.</i> 33 (12), 1870–6	2010 Jun
Applicability of accelerated solvent extraction for synthetic colorants analysis in meat products with ultrahigh performance liquid chromatography-photodiode array detection	Liao, Q. G.; Li, W. H.; Luo, L. G.	<i>Anal. Chim. Acta.</i> 716, 128–32	2012 Feb
Extraction, isolation, and purification of analytes from samples of marine origin – a multivariate task	Liguori, L.; Bjørsvik, H. R.	<i>J. Chromatogr. B Analyt. Technol. Biomed. Life Sci.</i> 910, 46–53	2012 Dec
Investigation on levels of polybrominated diphenyl ethers in retail fish and egg products in Shenzhen	Liu, B.; Zhang, L. S.; Zhang, J. Q.; Jiang, Y. S.; Zhou, J.; Huang, H. Y.	<i>Zhonghua Yu Fang Yi Xue Za Zhi.</i> 45 (12), 1068–72	2011 Dec
Characterization of secondary volatile profiles in <i>Nigella sativa</i> seeds from two different origins using accelerated solvent extraction and gas chromatography-mass spectrometry	Liu, X.; Abd El-Aty, A. M.; Cho, S. K.; Yang, A.; Park, J. H.; Shim, J. H.	<i>Biomed. Chromatogr.</i> 26 (10), 1157–62	2012 Oct
Accelerated solvent extraction of monacolin K from red yeast rice and purification by high-speed counter-current chromatography	Liu, Y.; Guo, X.; Duan, W.; Wang, X.; Du, J.	<i>J. Chromatogr. B Analyt. Technol. Biomed. Life Sci.</i> 878 (28), 2881–5	2010 Oct
Multiresidue determination of organophosphorus pesticides in ginkgo leaves by accelerated solvent extraction and gas chromatography with flame photometric detection	Lu, Y.; Yi, X.	<i>J. AOAC Int.</i> 88 (3), 729–735	2005



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: Sample Preparation Methods

Title	Authors	Publication	Publication Date
Influence of sample preparation on assay of phenolic acids from eggplant	Luthria, DL.; Mukhopadhyay, S.	<i>J. Agric. Food Chem.</i> 54 (1), 41–47	2006
Pressurised solvent extraction for organotin speciation in vegetable matrices	Marcic, C.; Lespe, S G.; Potin-Gautier, M.	<i>Anal. Bioanal. Chem.</i> 382 (7), 1574–83	2005 Aug
Comparison of different methods for the determination of the oil content in oilseeds	Matthäus, B.; Brühl, L.	<i>J. AOCS</i> 78 95–102.	2001 Jan
A comparison of automated and traditional methods for the extraction of arsenicals from fish	McKiernan, J. W.; Creed, J. T.; Brockhoff, C. A.; Caruso, J. A.; Lorenzana, R. M.	<i>J. Anal. At. Spectrom.</i> 14, 607–613	1999
Subcritical solvent extraction of anthocyanins from dried red grape pomace	Monrad, J. K.; Howard, L. R.; King, J.; Srinivas, K.; Mauromoustakos, A.	<i>J. Agric. Food Chem.</i> 58 (5), 2862–8	2010 Mar
Subcritical solvent extraction of procyanidins from dried red grape pomace	Monrad, J. K.; Howard, L. R.; King, J. W.; Srinivas, K.; Mauromoustakos, A.	<i>J. Agric. Food Chem.</i> 58 (7), 4014–21	2010 Apr
Pressurized liquid extraction of polar and nonpolar lipids in corn and oats with hexane, methylene chloride, isopropanol, and ethanol	Moreau, R. A.; Powell, M. J.; Singh, V.	<i>J. Oil Fat Industr.</i> 80 (11), 1063–1067	2003 Jan
Accelerated solvent extraction for natural products isolation	Mottaleb, M. A.; Sarker, S. D.	<i>Methods Mol. Biol.</i> 864, 75–87	2012
Optimization of extraction process for phenolic acids from black cohosh (<i>Cimicifuga racemosa</i>) by pressurized liquid extraction	Mukhopadhyay, S.; Luthria, D. L.; Robbins, R. J.	<i>J. Sci. Food Agric.</i> 86 (1), 156–162, 15	2006 Jan
Anxiolytic activity of a supercritical carbon dioxide extract of <i>Souroubea sympetala</i> (Marcgraviaceae)	Mullally, M.; Kramp, K.; Cayer, C.; Saleem, A.; Ahmed, F; McRae, C.; Baker, J.; Goulah, A.; Otorola, M.; Sanchez, P.; Garcia, M.; Poveda, L.; Merali, Z.; Durst, T.; Trudeau, V. L.; Arnason, J. T.	<i>Phytother. Res.</i> 25 (2), 264–70	2011 Feb
On-line clean-up of pressurized liquid extracts for the determination of polychlorinated biphenyls in feedingstuffs and food matrices using gas chromatography–mass spectrometry	Müller, A.; Björklund, E.; von Holst, C.	<i>J. Chromatogr., A.</i> 925 (1–2), 197–205	2001 Aug
Analysis of multiple herbicides in soybeans using pressurized liquid extraction and capillary electrophoresis	Nemoto, S.; Lehotay, S. J.	<i>J. Agric. Food Chem.</i> ; 46 (6), 2190–2199	1998
Comparison of sample preparation methods, validation of an UPLC-MS/MS procedure for the quantification of tetrodotoxin present in marine gastropods and analysis of pufferfish	Nzoughet, J. K.; Campbell, K.; Barnes, P.; Cooper, K. M.; Chevallier, O. P; Elliott, C. T.	<i>Food Chem.</i> 15; 136 (3-4), 1584–9	2013 Feb
Multiresidue analysis of pesticides in vegetables and fruits using two-layered column with graphitized carbon and water absorbent polymer	Obana, H.; Akutsu, K.; Okihashi, M.; Hori, S.	<i>The Analyst</i> 123, 711–714	1998
Analysis of 2-alkylcyclobutanones with accelerated solvent extraction to detect irradiated meat and fish	Obana, H.; Furuta, M.; Tanaka, Y.	<i>J. Agric. Food Chem.</i> 53 (17), 6603–8	2005 Aug



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: Sample Preparation Methods

Title	Authors	Publication	Publication Date
Determination of organophosphorus pesticides in foods using an accelerated solvent extraction system	Obana, H.; Kikuchi, K.; Okihashi, M.; Hori, S.	<i>Analyst</i> 122 (3), 217–20	1997 Mar
Pressurized hot water extraction of berberine, baicalin and glycyrrhizin in medicinal plants	Ong, E. S.; Shea Mei, L.	<i>Anal. Chim. Acta.</i> 482 (1), 81–89	2003 Apr
Pressurized liquid extraction of berberine and aristolochic acids in medicinal plants	Ong E. S.; Woo S. O.; Yong, Y. K.	<i>J. Chromatogr., A.</i> 904 (1), 57–6422	2000 Dec
Rapid determination of pesticide multiresidues in vegetables and fruits by accelerated solvent extraction coupled with online gel permeation chromatography-gas chromatography-mass spectrometry	Ouyang, Y.; Tang, H.; Wu, Y.; Li, G.	<i>Se Pu.</i> 30(7), 654–9	2012 Jul
Determination of zearalenone from wheat and corn by pressurized liquid extraction and liquid chromatography-electrospray mass spectrometry	Pallaroni, L.; von Holst, C.	<i>J. Chromatogr., A.</i> 993, 39–45	2003
Development of an extraction method for the determination of zearalenone in corn using less organic solvents	Pallaroni, L.; von Holst, C.	<i>J. Chromatogr., A.</i> 5 1055 (1-2), 247–9	2004 Nov
Stability of phenolic compounds during extraction with superheated solvents	Palma, M.; Piñeiro, Z.; Barroso, C. G.	<i>J. Chromatogr., A.</i> 6 921 (2), 169–74	2001 Jul
Extraction and analysis of trace amounts of cyclonite (RDX) and its nitroso-metabolites in animal liver tissue using gas chromatography with electron capture detection (GC-ECD)	Pan, X.; Zhang, B.; Cobb, G. P.	<i>Talanta</i> 67 (4), 816–23	2005 Oct
Simultaneous determination of 405 pesticide residues in grain by accelerated solvent extraction then gas chromatography-mass spectrometry or liquid chromatography-tandem mass spectrometry	Pang, G.; Liu, Y.; Fan, C.; Zhang, J.; Cao, Y.; Li, X.; Li, Z.; Wu, Y.; Guo, T.	<i>Anal. Bioanal. Chem.</i> 384, 1366–1408	2006 Mar
Automated sample preparation by pressurized liquid extraction-solid-phase extraction for the liquid chromatographic-mass spectrometric investigation of polyphenols in the brewing process	Papagiannopoulos, M.; Mellenthin, A.	<i>J. Chromatogr., A.</i> 8 976 (1-2), 345–8	2002 Nov
Online coupling of pressurized liquid extraction, solid-phase extraction and high-performance liquid chromatography for automated analysis of proanthocyanidins in malt	Papagiannopoulos, M.; Zimmermann, B.; Mellenthin, A.; Krappe, M.; Maio, G.; Galensa, R.	<i>J. Chromatogr., A.</i> 7 958 (1-2), 9–16	2002 Jun
Simultaneous determination of 13 quinolones from feeds using accelerated solvent extraction and liquid chromatography	Pecorelli, I.; Galarini, R.; Bibi, R.; Floridi, A. I.; Casciarri, E.; Floridi, A.	<i>Anal. Chim. Acta.</i> 483 (1-2), 81–89	2003 April
Comparison of soxhlet, ultrasound-assisted and pressurized liquid extraction of terpenes, fatty acids and Vitamin E from <i>Piper gaudichaudianum</i> Kunth	Péres, V. F.; Saffi, J.; Melecchi, M. I.; Abad, F. C.; de Assis Jacques, R.; Martinez, M. M.; Oliveira, E. C.; Caramão, E. B.	<i>J. Chromatogr., A.</i> 1105 (1-2), 115–8	2006 Feb
Pressurised fluid extraction (PFE) as an alternative general method for the determination of pesticide residues in rape seed	Pihlström, T.; Isaac, G.; Waldebäck, M.; Osterdahl, B. G.; Markides, K. E.	<i>Analyst</i> 127 (4), 554–9	2002 Apr
Determination of catechins by means of extraction with pressurized liquids	Piñeiro, Z.; Palma, M.; Barroso C. G.	<i>J. Chromatogr., A.</i> 13 1026 (1-2), 19–23.	2004 Feb



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: Sample Preparation Methods

Title	Authors	Publication	Publication Date
An improved clean-up strategy for simultaneous analysis of polychlorinated dibenzo-p-dioxins (PCDD), polychlorinated dibenzofurans (PCDF), and polychlorinated biphenyls (PCB) in fatty food samples	Pirard, C.; Focant, J. F.; De, P. E.	<i>Anal. Bioanal. Chem.</i> 372 (2), 373–81.	2002 Jan
Extraction of polar and hydrophobic pollutants using accelerated solvent extraction (ASE)	Pörschmann, J., Plugge, J.	<i>Fresen. J. Anal. Chem.</i> 364 (7), 643–645	1999
Quantification of the total amount of artemisinin in leaf samples by thin layer chromatography	Quennoz, M.; Bastian, C.; Simonnet, X.; Grogg, A. F.	<i>Chimia (Aarau)</i> 64 (10), 755–7.	2010
Determination of fat in dairy products using pressurized solvent extraction	Richardson, R. K.	<i>J. AOAC Int.</i> 84 (5), 1522–1533	2001
Influence of altitudinal variation on the content of phenolic compounds in wild populations of <i>Calluna vulgaris</i>, <i>Sambucus nigra</i>, and <i>Vaccinium myrtillus</i>	Rieger, G.; Müller, M.; Guttenberger, H.; Bucar, F.	<i>J. Agric. Food Chem.</i> 56 (19), 9080–6.	2008 Oct
Pressurized liquid extraction of isoflavones from soybeans	Rostagno, M. A.; Palma, M.; Barroso, C. G.	<i>Anal. Chim. Acta.</i> 522 (2), 169–177.	2004 Sep
A multi-residue method for the analysis of organophosphorus residues in cooked and polished rice using accelerated solvent extraction and dispersive-solid phase extraction (D-SPE) technique and uncertainty measurement	Sanyal, D.; Rani, A.; Alam, S.	<i>J. Environ. Sci. Health, B</i> 44 (7), 706–16.	2009 Sep
Accelerated solvent extraction of lipids for determining the fatty acid composition of biological material	Schäfer, K.	<i>Anal. Chim. Acta.</i> 358 (1), 69–77	1998 Jan
HPLC analysis of kaempferol and quercetin derivatives isolated by different extraction techniques from plant matrix	Skalicka-Wozniak, K.; Szypowski, J.; Glowniak, K.	<i>J. AOAC Int.</i> 94 (1), 17–21.	Jan-Feb 2011
Statistical evaluation of fatty acid profile and cholesterol content in fish (common carp) lipids obtained by different sample preparation procedures	Spiric, A.; Trbovic, D.; Vranic, D.; Djinic, J.; Petronijevic, R.; Matekalo-Sverak, V.	<i>Anal. Chim. Acta.</i> 672 (1-2), 66–71.	2010 Jul
Application of accelerated solvent extraction in the analysis of organic contaminants, bioactive and nutritional compounds in food and feed	Sun, H.; Ge, X.; Lv, Y.; Wang, A.	<i>J. Chromatogr., A.</i> 1237, 1–23.	2012 May
Development of an accelerated solvent extraction, ultrasonic derivatisation LC-MS/MS method for the determination of the marker residues of nitrofurans in freshwater fish	Tao, Y.; Chen, D.; Wei, H.; Yuanhu, P.; Liu, Z.; Huang, L.; Wang, Y.; Xie, S.; Yuan, Z.	<i>Food Addit. Contam. Part A Chem. Anal. Control Expo. Risk Assess.</i> 29 (5), 736–45.	2012
Simultaneous determination of lincomycin and spectinomycin residues in animal tissues by gas chromatography-nitrogen phosphorus detection and gas chromatography-mass spectrometry with accelerated solvent extraction	Tao, Y.; Chen, D.; Yu, G.; Yu, H.; Pan, Y.; Wang, Y.; Huang, L.; Yuan, Z.	<i>Food Addit. Contam. Part A Chem. Anal. Control Expo. Risk Assess.</i> 28 (2), 145–54.	2011 Feb
Determination of 17 macrolide antibiotics and avermectins residues in meat with accelerated solvent extraction by liquid chromatography-tandem mass spectrometry	Tao, Y.; Yu, G.; Chen, D.; Pan, Y.; Liu, Z.; Wei, H.; Peng, D.; Huang, L.; Wang, Y.; Yuan, Z.	<i>J. Chromatogr. B Analyt. Technol. Biomed. Life Sci.</i> 897, 64–71.	2012 May
Determination of seven toxaphene congeners in ginseng and milkvetch root by gas chromatography tandem mass spectrometry	Tian, S.; Mao, X.; Miao, S.; Jia, Z.; Wang, K.; Ji, S.	<i>Se Pu.</i> 30 (1), 14–20.	2012 Jan



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: Sample Preparation Methods

Title	Authors	Publication	Publication Date
A consecutive preparation method based upon accelerated solvent extraction and high-speed counter-current chromatography for isolation of aesculin from <i>Cortex fraxinus</i>	Tong, X.; Zhou, T, Xiao, X.; Li, G.	<i>J. Sep. Sci.</i> 35 (24), 3609–14	2012 Dec
Characterization of anthocyanins and anthocyanidins in purple-fleshed sweetpotatoes by HPLC-DAD/ESI-MS/MS	Truong, V. D.; Deighton, N.; Thompson, R. T.; McFeeters, R. F.; Dean, L. O.; Pecota, K. V.; Yencho, G. C.	<i>J. Agric. Food Chem.</i> 58 (1), 404–10	2010 Jan
Fat extraction from acid- and base-hydrolyzed food samples using accelerated solvent extraction	Ullah, S. M.; Murphy, B.; Dorich, B.; Richter, B.; Srinivasan, K.	<i>J. Agric. Food Chem.</i> 59 (6), 2169–74.	2011 Mar
Analysis of zearalenone in cereal and swine feed samples using an automated flow-through immunosensor	Urraca, J. L.; Benito-Peña, E.; Pérez-Conde, C.; Moreno-Bondi, M. C.; Pestka, J. J.	<i>J. Agric. Food Chem.</i> 53 (9), 3338–3344	2005
Accelerated solvent extraction and gas chromatography/mass spectrometry for determination of polycyclic aromatic hydrocarbons in smoked food samples	Wang, G.; Lee, A. S.; Lewis, M.; Kamath, B.; Archer, R. K.	<i>J. Agric. Food Chem.</i> 47 (3), 1062–6.	1999 Mar
Subcritical water extraction of alkaloids in <i>Sophora flavescens</i> Ait. and determination by capillary electrophoresis with field-amplified sample stacking	Wang, H.; Lu, Y.; Chen, J.; Li, J.; Liu, S.	<i>J. Pharm. Biomed. Anal.</i> 58, 146–51.	2012 Jan
Evaluation of Soxhlet extraction, accelerated solvent extraction and microwave-assisted extraction for the determination of polychlorinated biphenyls and polybrominated diphenyl ethers in soil and fish samples	Wang, P.; Zhang, Q.; Wang, Y.; Wang, T.; Li X.; Ding, L.; Jiang, G.	<i>Anal. Chim. Acta.</i> 663 (1), 43–8.	2010 Mar
Determination of ten pesticides of pyrazoles and pyrroles in tea by accelerated solvent extraction coupled with gas chromatography-tandem mass spectrometry	Xu, D.; Lu, S.; Chen, D.; Lan, J.; Zhang, Z.; Yang, F.; Zhou, Y.	<i>Se Pu.</i> ; 31 (3), 218–22.	2013 Mar
Online cleanup of accelerated solvent extractions for determination of adenosine 5'-triphosphate (ATP), adenosine 5'-diphosphate (ADP), and adenosine 5'-monophosphate (AMP) in royal jelly using high-performance liquid chromatography	Xue, X.; Wang, F.; Zhou, J.; Chen, F.; Li, Y.; Zhao, J.	<i>J. Agric. Food Chem.</i> 57 (11), 4500–5.	2009 Jun
Identification and quantitation of eleven sesquiterpenes in three species of <i>Curcuma</i> rhizomes by pressurized liquid extraction and gas chromatography–mass spectrometry	Yang, F. Q.; Li, S.; Chen, Y.; Lao, S. C.; Wang, Y.T.; Dong, T. T. X.; Tsim, K. W. K.	<i>J. Pharm. Biomed. Anal.</i> 39 (3/4), 552–558	2005 Sep
Dispersive solid-phase extraction cleanup combined with accelerated solvent extraction for the determination of carbamate pesticide residues in <i>Radix glycyrrhizae</i> samples by UPLC-MS-MS	Yang, R. Z.; Wang, J. H.; Wang, M. L.; Zhang, R.; Lu, X. Y.; Liu, W. H.	<i>J. Chromatogr. Sci.</i> 49 (9), 702–8.	2011 Oct
Simultaneous determination of amitraz and its metabolite residue in food animal tissues by gas chromatography-electron capture detector and gas chromatography-mass spectrometry with accelerated solvent extraction	Yu, H.; Tao, Y.; Le, T.; Chen, D.; Ishsan, A.; Liu, Y.; Wang, Y.; Yuan, Z.	<i>J. Chromatogr. B Analyt. Technol. Biomed. Life Sci.</i> 878 (21), 1746–52.	2010 Jul
Simultaneous determination of fluoroquinolones in foods of animal origin by a high performance liquid chromatography and a liquid chromatography tandem mass spectrometry with accelerated solvent extraction	Yu, H.; Tao, Y.; Chen, D.; Pan, Y.; Liu, Z.; Wang, Y.; Huang, L.; Dai, M.; Peng, D.; Wang, X.; Yuan, Z.	<i>J. Chromatogr. B Analyt. Technol. Biomed. Life Sci.</i> 885-886, 150–9.	2012 Feb



Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)

Peer Reviewed Journals: Sample Preparation Methods

Title	Authors	Publication	Publication Date
Determination of pentachlorophenol residue in meat and fish by gas chromatography-electron capture detection and gas chromatography-mass spectrometry with accelerated solvent extraction	Zhao, D.	<i>J. Chromatogr. Sci.</i>	2013 May
Response surface modeling and optimization of accelerated solvent extraction of four lignans from <i>fructus schisandrae</i>	Zhao, L. C.; He, Y. Deng.; X, Yang, G. L.; Li, W.; Liang, J.; Tang, Q. L.	<i>Molecules</i> . 17 (4), 3618–29	2012 Mar
Determination of acetanilide herbicides in cereal crops using accelerated solvent extraction, solid-phase extraction and gas chromatography-electron capture detector	Zhang, Y.; Yang, J.; Shi, R.; Su, Q.; Yao, L.; Li, P.	<i>J. Sep. Sci.</i> 34 (14), 1675–82	2011 Jul
Application of accelerated solvent extraction coupled with high-performance counter-current chromatography to extraction and online isolation of chemical constituents from <i>Hypericum perforatum</i> L	Zhang, Y.; Liu, C.; Yu, M.; Zhang, Z.; Qi, Y.; Wang, J.; Wu, G.; Li, S.; Yu, J.; Hu, Y.	<i>J. Chromatogr., A</i> . 1218 (20), 2827–34	2011 May
Analysis of volatile components in Qingshanlvshui tea using solid-phase microextraction/accelerated solvent extraction-gas chromatography-mass spectrometry	Zhan, J.; Lu, S.; Meng, Z.; Xiang, N.; Cao, Q.; Miao, M.	<i>Se Pu</i> . 26 (3), 301–5.	2008 May





Technical Collateral: Sample Preparation Methods

Product Number	Technique	Title
AN 326	HPLC-UV	Extraction of Drugs from Animal Feeds Using Accelerated Solvent Extraction (ASE)
AN 335	HPLC-UV	Accelerated Solvent Extraction (ASE) of Active Ingredients from Natural Products
AN 356	IC-conductivity	Determination of Perchlorate in Vegetation Samples Using Accelerated Solvent Extraction and Ion Chromatography
AN 357	HPLC	Extraction of Phenolic Acids from Plant Tissue Using Accelerated Solvent Extraction (ASE)
AN 363	HPLC	Extraction of Herbal Marker Compounds Using Accelerated Solvent Extraction Compared to Traditional Pharmacopoeia Protocols

Table of Contents

[Introduction](#)

[Analytical Technologies](#)

[Targeted Methods](#)

[Coffee](#)

[Controlled Substances](#)

[Fruit Juice](#)

[Melamine](#)

[Mixtures of Complex Materials](#)

[Non-targeted Methods](#)

[Herb Blending](#)

[Beer](#)

[Wine](#)

[Fruit Juice](#)

[Tea](#)

[Botanical Supplements](#)

[Tryptophan Supplements](#)

[Olive Oil](#)

[References](#)



www.thermoscientific.com/liquidchromatography

©2015 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific,
Sunnyvale, CA USA is
ISO 9001:2008 Certified.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
S C I E N T I F I C

A Thermo Fisher Scientific Brand