

# What is Eluent Generation?

Monika Verma  
Thermo Fisher Scientific, Sunnyvale, CA, USA

## Executive Summary

Eluent generation is an electrolytic process that allows the production of high purity eluents for ion chromatography, and is a part of Reagent-Free Ion Chromatography™ (RFIC™). The entire process only requires deionized water, eliminating the need to handle any acids or bases.

## Keywords

RFIC-EG, Electrolytic Process,  
Eluent Generation, Deionized Water

## What exactly does eluent generation mean?

Eluent generation allows the automatic production of high purity ion chromatography (IC) eluents. This is made possible through precise control of the electric current applied to the electrolysis of water to generate hydroxide and hydronium ions. Eluent generation eliminates the need to manually prepare eluents from concentrated acids and bases. The only routine reagent needed is deionized water. Furthermore, since the instrument pump seals and pistons only come in contact with deionized water instead of acids and bases which can precipitate, overall pump maintenance is significantly reduced.

## How does the process of eluent generation work?

Figure 1 illustrates how eluent generation is an electrolytic process by which eluent is generated. A pair of electrodes is positioned with an ion exchange membrane separating them; when a current is applied to the electrodes, electrolysis of water generates hydroxide at the cathode and hydronium at the anode. The ion exchange membrane prevents the species from recombining into water, and allows a counter-ion from the Eluent Generation Cartridge to migrate across the membrane to form the eluent. The eluent concentration is varied by changing the applied current to the electrodes; thus the eluent concentration is infinitely variable throughout a single run. This entire process can be done without the use of extra pumps, fittings, valves or any moving parts.

It is important to understand that eluent generation is NOT the dilution of an eluent concentrate with water by the use of a burette type device. A dilution procedure does not allow the ability to change concentration throughout a run without the use of additional mechanical devices such as valves nor can it further purify the eluent. In addition, dilution is not as accurate and reproducible as Eluent Generation.




Figure 1. The Thermo Scientific™ Dionex™ EluGen™ EGC III KOH eluent generator cartridge consists of a KOH generation chamber and a  $K^+$  electrolyte reservoir, connected by a cation exchange connector. A high-pressure connector permits the passage of  $K^+$  ions from the  $K^+$  electrolyte reservoir into the electrolytic chamber.

### What applicable solutions exist on the market?

The Dionex Eluent Generator Cartridge (EGC), shown in Figure 2, is at the core of the patented eluent generation technology used in RFIC-EG systems. The EGC KOH, NaOH, or LiOH cartridges produce potassium, sodium, or lithium hydroxide eluents respectively for anion exchange applications. With hydroxide-selective columns, such as the Thermo Scientific™ Dionex™ IonPac™ AS18, common anions can be rapidly separated using the simplicity of an RFIC-EG system. The EGC MSA cartridge produces methanesulfonic acid eluent for cation-exchange chromatography applications. When the MSA eluent is used with the Dionex IonPac CS12A and CS16 columns, it allows for isocratic or gradient separation of Group I and Group II cations and ammonium. The EGC  $K_2CO_3$  cartridge produces carbonate-only eluent. An Electrolytic pH Modifier (EPM) can be added for generating these carbonate/bicarbonate eluents. It is also important to note that these systems eliminate variability and errors associated with manual eluent preparation which directly improves method reproducibility.



Figure 2. Thermo Scientific Dionex EGC-500 cartridge.

### Which laboratories will benefit from these solutions?

RFIC-EG systems facilitate drinking, waste, and groundwater analyses for regulatory compliance. Furthermore, they provide the accuracy and reproducibility needed for the analysis of high-purity water in multiple industries.

For example, the EPA's Statistics and Analytical Support Branch and the Office of Ground Water and Drinking Water's Technical Support Center have concluded that the use of hydroxide eluents falls within the method flexibility allowed in EPA Methods 300.0 Part A and 300.1 Part A for determining inorganic anions. These methods are critical for compliance monitoring under the Clean Water Act and Safe Drinking Water Act.

RFIC-EG systems have redefined ion chromatography by making it possible to just add water to completely operate an IC. This system allows for a simpler and more reliable way to help deliver superior results while simultaneously saving time and labor.

### References

- ICS-2100 Bid Specifications
- RFIC System Capabilities Brochure
- <http://www.dionex.com/en-us/products/ion-chromatography/ic-rfic-solutions/rfic-rfic-eg/lp-72708.html>

### [www.thermoscientific.com/dionex](http://www.thermoscientific.com/dionex)

©2013 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Australia** +61 3 9757 4486  
**Austria** +43 1 333 50 34 0  
**Belgium** +32 53 73 42 41  
**Brazil** +55 11 3731 5140  
**China** +852 2428 3282

**Denmark** +45 70 23 62 60  
**France** +33 1 60 92 48 00  
**Germany** +49 6126 991 0  
**India** +91 22 6742 9494  
**Italy** +39 02 51 62 1267

**Japan** +81 6 6885 1213  
**Korea** +82 2 3420 8600  
**Netherlands** +31 76 579 55 55  
**Singapore** +65 6289 1190  
**Sweden** +46 8 473 3380



Thermo Fisher Scientific, Sunnyvale, CA  
 USA is ISO 9001:2008 Certified.

**Switzerland** +41 62 205 9966  
**Taiwan** +886 2 8751 6655  
**UK/Ireland** +44 1442 233555  
**USA and Canada** +847 295 7500

**Thermo**  
 SCIENTIFIC

Part of Thermo Fisher Scientific