



Thermo Scientific
Pesticide Explorer Collection


Start-to-finish

workflows for pesticide analysis

Comprehensive Pesticide Analysis Solutions

Pesticide Explorer Collection

Selection table	Triple Quadrupole Solutions		Orbitrap Solutions	
	1-Standard Quan	2-Premium Quan	3-HRAM Quan	4-HRAM Quan/Qual
Lab Profile				
<i>Routine - Standard Quantitation</i>	●			
<i>Routine - High Sensitivity Quantitation</i>		●		
<i>Routine - Targeted HRAM Quantitation</i>			●	
<i>Routine Plus - Targeted and Non-targeted Screening and Quantitation</i>				●
Workflow and Components	Targeted	Targeted	Targeted	Targeted + Non-targeted
<i>Thermo Scientific™ QuEChERS kit</i>	●	●	●	●
<i>Pre-set Methods (LC-MS)</i>	●	●	●	●
<i>Thermo Scientific™ Accucore™ aQ C18 LC columns (100 x 2.1 mm, 2.6µm)</i>	●	●	●	●
<i>Thermo Scientific™ UltiMate™ 3000 LC system (HPG-3400RS)</i>	●	●	●	●
MS System				
<i>Thermo Scientific™ TSQ Endura MS</i>	●			
<i>Thermo Scientific™ TSQ Quantiva MS</i>		●		
<i>Thermo Scientific™ Q Exactive™ Focus MS</i>			●	●
Software				
<i>Thermo Scientific™ TraceFinder™ software</i>	●	●	●	●
<i>Thermo Scientific™ Compound Discoverer™ software</i>				●
<i>Thermo Scientific™ SIEVE™ software</i>				●
<i>Thermo Scientific™ HRAM MS/MS Spectral Library</i>			●	●
<i>Quick-Start Manual</i>	●	●	●	●
<i>3-day training</i>	●	●	●	●



The Thermo Scientific™ Pesticides Explorer Collection is a comprehensive set of liquid chromatography-mass spectrometry (LC-MS) solutions designed for laboratories performing routine quantitation, targeted and non-targeted screening of pesticide residues in food matrices. Each configuration includes all the workflow components needed—consumables, hardware, software and built-in instrument and data processing methods—pre-configured and tested, from your single trusted supplier, Thermo Fisher Scientific.

Choice of configurations simplify complex pesticide methods

To ensure your immediate success, the Pesticide Explorer Collection is available in four pre-configured, pre-tested solutions.

For laboratories performing routine targeted quantitation, the Pesticide Explorer Collection offers two triple quadrupole mass spectrometer workflow-based configurations for Standard and Premier Quantitation. The Standard Quantitation configuration includes the TSQ Endura triple quadrupole mass spectrometer and ensures compliance against regulated levels of detection. For quantitative applications demanding the highest possible sensitivity, the Premier Quantitation configuration features the TSQ Quantiva triple quadrupole mass spectrometer to meet or exceed regulatory limits.

For laboratories offering specialized analytical services, there is a choice of two Thermo Scientific™ HRAM Orbitrap configurations for quantitation and screening employing the Thermo Scientific™ Q Exactive™ Focus mass spectrometer. HRAM capability significantly enhances quantitative accuracy even when analyzing complex sample matrices with limited prior clean up as well as the ability to quantify and confirm in one single analysis.

Configurations available in the Pesticide Explorer Collection

- Pesticide Explorer Collection – Standard Quantitation
- Pesticide Explorer Collection – Premier Quantitation
- Pesticide Explorer Collection – Orbitrap HRAM Quantitation
- Pesticide Explorer Collection – Orbitrap HRAM Screening and Quantitation

Robust, Routine Workflows for Triple Quadrupo

The Standard and Premier Quantitation configurations provide everything needed to perform robust, routine workflows for SRM-based quantitation of pesticides, from the QuEChERS sample extraction kit to proven multi-class pesticide residue analysis methods. Pre-configured methods for each solution are easy to access on the USB drive included with the Pesticide Explorer Collection installation guide.




Select compounds from the compound database to automatically create the instrument and processing method.

Users can quickly and easily set up a new method using the TraceFinder software compound database, or by uploading and modifying pre-configured methods.



UltiMate 3000 LC System with TSQ Endura MS or TSQ Quantiva MS

SELECT




Flexibility in the instrument method allows you to create or edit pre-configured methods with SRM transitions and retention times with ease.



Select a pre-configured TSQ-tested method with column and conditions. Load and go!


File-Based Quantitation




The color flagging features in TraceFinder software enable the analyst to quickly review and confirm results.

LOAD METHOD AND ACQUIRE DATA

REVIEW AND REPORT



Data is acquired under optimal conditions for large multi-component pesticide residue analysis.




After data review, generate high-quality standard and custom reports to rapidly turn your samples into results.

Routine Workflows for Orbitrap HRAM Targeted

The Orbitrap HRAM Quantitation configuration provides high resolution accurate mass analysis, a unique capability that enables quantitation without compromise in sensitivity, accuracy, precision, or linear dynamic range. High resolving power is particularly useful for the analysis of contaminants in complex matrices and can overcome the masking effects of isobaric interferences, allowing detection of analytes at very low concentrations.

Reduce method development time by using the HRAM spectral fragmentation library and compound database to identify compounds with speed and confidence.




SELECT

Pre-configured instrument methods for targeted quantitation enable the user to quickly start acquiring data.

LOAD METHOD AND ACQUIRE DATA




Q Exactive Focus MS and UltiMate 3000 LC System



Master Method templates available with the Pesticide Explorer Collection instantly allows users to acquire and process data for quantitation.

Quantitation

Review results rapidly by flagging method criteria and/or regulatory levels. Choose from a large number of pre-defined and customizable report templates meeting a variety of regulatory requirements. View the result of changes to the data processing method immediately without manual intervention.



In data review, the analyst is able to quickly glance through all compounds and samples in a single view.

REVIEW AND REPORT


Quantitation Report - 2						
#	Compound Name	RT	Scan Peak	Response	Conc. Units	Det. (ppm)
14	1 Ergosterol	8.72	189 1313	336915070	1 ppm	0.00
16	12 Abiraterone-Cl	7.16	221 1528	2016276251	1 ppm	0.00
17	10 Metoprolol-Cl	7.54	213 1885	306262660	1 ppm	0.02
Surrogate						
#	Compound Name	RT	Scan Peak	Response	Conc. Units	Det. (ppm)
22	16 Surrogate	8.57	189 1313	336915070	1 ppm	0.00
Target Compounds						
#	Compound Name	RT	Scan Peak	Response	Conc. Units	Det. (ppm)
28	2 Acetophenone	2.50	142 8626	54803331	10 416 ppm	0.00
31	7 Carbamazepine	3.36	182 8788	238621382	0 647 ppm	0.01
32	4 Acetaminophen	3.84	233 8146	225653736	11 786 ppm	0.00
33	5 Aspirin	4.58	116 8528	49688223	10 287 ppm	0.00
34	6 Benzocaine	5.65	224 8817	78140018	12 336 ppm	0.00
35	7 Paracetamol-sulfate	5.86	320 108	476588873	10 887 ppm	-0.21
36	8 Dext	5.97	213 8849	13622426	11 219 ppm	0.00
37	9 Paracetamol-sulfate	6.21	336 1629	168820360	10 288 ppm	0.00
38	10 Theophylline	6.86	368 8863	8468204	10 651 ppm	0.21
39	11 Phenylephrine	6.96	277 816	948188177	6 825 ppm	-0.01
40	12 Abacavir	7.22	276 129	308818422	10 134 ppm	0.00
41	14 Phorbol-sulfate	7.51	295 8088	92228147	10 28 ppm	0.00
42	16 Ibuprofen	7.68	287 1482	362514270	10 763 ppm	0.01
43	17 Benzocaine	8.50	360 824	60061642	9 887 ppm	0.01

Customizable reports can easily be exported in MS Excel format for possible LIMS integration.


New Analytical Capability Workflows For

The Q Exactive Focus mass spectrometer produces data that can be used for highly sensitive and selective quantitation as well for in-depth screening. Built-in databases designed for food safety and environmental analyses make quantitation, and targeted as well as non-targeted screening from a single data set seamless, with the option to reanalyze data retrospectively at a future date without the need for sample reinjection. Even without using reference standards, TraceFinder software allows the user to perform relative quantitation between multiple samples in their workflow.

Targeted Screening




Q Exactive Focus MS and UltiMate 3000 LC System



With full-scan, targeted data-dependant MS/MS, the analyst enjoys greater flexibility to detect large sets of compounds in a single run.

METHOD SELECTION

Non-targeted Screening



With variable data independent acquisition (vDIA), full-scan MS and MS/MS HRAM analysis, no sample-specific method optimization is necessary, and the risk of missing important non-targeted compounds is greatly reduced.

Orbitrap HRAM Screening and Quantitation

Built-in databases in TraceFinder software allow the user to perform quantitation.




Target Screening Summary Report														
1	2	3	4	5	6	7	8	9	10	11	12			
Batch	Method	Inst. Method	Sample ID	Raw File	Sample Folder	Viol. Position	Acquisition Date	Lab Name	Agency	Default Laboratory	Study Name			
1	MPV Screening Focus	MPV Screening	MPV_Screening_V04_A030416_011114.rmg		Unknown_2	181	2012114 6 14 15.348							
C:\Thermo\TraceFinder\3.1\FF\Orbitrap\Completed\FFS-New DB.cdf														
13	14	15	16	17	18	19	20	21	22	23	24			
Found	Confirmed	Target Name	M. Area	RT (Min)	Formula	Expected m/z	Measured m/z	Coll. m/z	Isotopic Pattern (%)	Num. Peaks Matched	Library Name (ID)	Library Score (ID)	Fragments Found	Active
1	1	Cytosine	1.19E+07	0.55	C1H2N2O2	108.0546	108.0546	108.0546	100	2 of 2	Cytosine (10)	5	0	Hydrogen
1	1	Methylammonium	8.48E+06	0.65	C14H15N2O2S	262.0916	262.0916	262.0916	100	2 of 2	Methylammonium (10)	5	0	Hydrogen
1	1	Deoxyfluoramide	3.18E+06	0.65	C8H10F2NO2S	214.0646	214.0646	214.0646	100	2 of 2	Deoxyfluoramide (10)	5	0	Hydrogen
1	1	Progesterone, M201H	5.54E+07	0.70	C21H30O2	314.2246	314.2246	314.2246	100	2 of 2	Progesterone (10)	5	0	Hydrogen
1	1	Progesterone	1.43E+07	0.70	C21H30O2	314.2246	314.2246	314.2246	100	2 of 2	Progesterone (10)	5	0	Hydrogen
1	1	Debutolan	5.48E+07	0.75	C12H15N2O2S	232.1104	232.1104	232.1104	100	2 of 2	Debutolan (10)	5	0	Hydrogen
1	1	Expiprin H2	3.43E+06	0.75	C9H12O2N2S	185.1016	185.1016	185.1016	100	2 of 2	Expiprin H2 (10)	5	0	Hydrogen
1	1	Progiprin	1.18E+07	0.80	C11H15N2O2S	214.1176	214.1176	214.1176	100	1 of 2	Progiprin (10)	5	0	Hydrogen
1	1	Progiprin, M201H	3.10E+06	0.80	C11H15N2O2S	214.1176	214.1176	214.1176	100	2 of 2	Progiprin (10)	5	0	Methylammonium
1	1	Progiprin, M201H	2.78E+07	0.80	C11H15N2O2S	214.1176	214.1176	214.1176	100	2 of 2	Progiprin (10)	5	0	Methylammonium
1	1	Tamoxifen, 2,3,5-	1.38E+07	0.80	C18H19NO2	273.1416	273.1416	273.1416	100	1 of 2	Tamoxifen (10)	5	0	Hydrogen
1	1	2,3,5-Tamoxifen	1.38E+07	0.80	C18H19NO2	273.1416	273.1416	273.1416	100	1 of 2	2,3,5-Tamoxifen (10)	5	0	Hydrogen
1	1	Deoxyfluoramide	9.95E+05	0.65	C8H10F2NO2S	214.0646	214.0646	214.0646	100	1 of 1	Deoxyfluoramide (10)	5	0	Hydrogen

View positive hits during data review; confirmation is easily done using fragmentation, spectra, retention time, and isotope pattern information.

Detailed reports are generated and exported for printing and also customizable for individual customers.

CONFIRM RESULTS

REVIEW AND REPORT



When integrated into the TraceFinder workflow, SIEVE software extracts compound spectra from complex data and performs differential and statistical analysis to determine the compounds that vary significantly within the sample set.

Compound Discoverer software can be used to identify the degradants and metabolites of known target compounds.

Power of HRAM MS/MS Spectral Libraries and Compound Databases

The Thermo Scientific Pesticide Explorer Collection includes the High-Resolution Accurate-Mass MS/MS Spectral Libraries and TraceFinder compound databases. These were designed specifically for targeted screening, quantitation and analysis of non-targeted contaminants in food, environmental, clinical research and forensic toxicology sample matrices. Fully integrated and searchable using TraceFinder software with over 2,600 compounds and more than 15,000 spectra, the HRAM MS/MS Spectral Libraries can be used to screen and identify a variety of known and unknown compounds with speed and confidence. To reduce method development time and increase confidence in compounds detected, the HRAM MS/MS compound database includes retention times that were determined using the same internal standards.

Laboratories using The Pesticide Explorer Collection comprising of either the Orbitrap HRAM Quantitation or Orbitrap HRAM Screening and Quantitation configuration will enjoy the time-savings provided by an “off-the-shelf” database of MS/MS spectra – no need to build your own libraries! To meet specific laboratory requirements, the libraries can be expanded and customized by adding new compounds and spectra.

COMPOUND CLASS	
Food Safety and Environmental	Forensic Toxicology
Emerging Environmental Contaminants	Drugs of Abuse
Pesticides	Natural and Industrial Toxins
Veterinary Drugs	Prescription Drugs
Mycotoxins	Performance Enhancing Drugs
Perfluorinated Compounds (PFCs)	Other Drug Monitoring Research

Table 1. Compound classes provided in the HRAM MS/MS spectral libraries.


Compound Groupings	Unique Entries	Total Spectra
Environmental and Food Safety	1,634	8,906
Clinical Research and Forensic Toxicology	926	4,630

Table 2. Unique entries and total spectra in the HRAM MS/MS spectral libraries.




Thermo Scientific m/z Cloud Library assists with unknown identification

The m/z Cloud library is a continuously expanding, highly curated mass spectral fragmentation database that assists in the identification of compounds when they are not in local spectral fragmentation libraries. The library features a freely searchable collection of HRAM spectra that can be accessed using an advanced spectral correlation algorithm.



Q EXACTIVE FOCUS

Comprehensive, High Productivity Reporting



As the identification of pesticide residues becomes increasingly important, so does that ability to review and report results quickly and accurately. In particular, high productivity multi-residue analyses require accurate and robust reporting across multiple regulatory agencies. The TraceFinder software reporting system allows creation and automated generation of quantitative and qualitative reports covering a broad array of pesticide residue analysis methods. The reports can be tailored to meet specific reporting needs.

The heat map reporting templates allow for the quick evaluation of the presence of components and their relative abundance across the sample set.

Powerful Software Empowers Your Pesticide Analyses



TraceFinder software provides an extensive suite of targeted quantitation and screening LC-MS workflows, experiments and reports for environmental and food safety applications as well as protects the integrity of the analytical data.

Further updates on the TraceFinder compound database are offered to registered users at <https://thermo.flexnetoperations.com>



SIEVE software enables label-free, semi-quantitative differential analysis of complex LC-MS data sets. Whether your analytes are large or small molecules or whether you are comparing two or hundreds of samples, SIEVE software confidently and reproducibly identifies components with statistically significant inter-sample differences in abundance.



Designed to save time, Compound Discoverer software includes an extensive set of tools to ensure confident compound identification and structural elucidation in pharmaceutical metabolism, impurity analysis, forensic toxicology, and food and environmental applications.

Single Provider Solution Ensures Your Success


The Pesticide Explorer Collection provides start-to-finish workflows tailored to help food monitoring and testing laboratories reduce startup time and cost. The collection provides compelling productivity and efficiency enhancements for both startup laboratories and laboratories adding new analytical capabilities to address evolving customer and industry demands. Regardless of staff or laboratory expertise, the Pesticide Explorer Collection is designed to make it easier for both new and experienced users to obtain reliable, unambiguous, high-quality LC-MS/MS results. At Thermo Fisher Scientific, we'll match you to the right products - through experienced commercial support combined with the most comprehensive line of instruments, equipment, consumables, chemicals, reagents and software available.



QuEChERs products are a convenient and effective approach for determining pesticide residues in fruit, vegetables and other foods. The extraction and clean-up products are easy to use and enable researchers to determine greater numbers of pesticides than with standard SPE.



Perform the full range of HPLC or nano HPLC applications on one system platform using UltiMate 3000 RSLC systems. For added mobile phase flexibility, the included Thermo Scientific™ Dionex™ Ultimate™ HPG-3400RS Rapid Separation Binary Pump is included, provides true ballistic gradients with excellent retention time precision.



Simplify and improve your analytical results with Thermo Scientific™ HPLC columns. Available in particle sizes and column designs to meet all separation needs, they improve resolution, enhance sensitivity, and deliver faster analysis and consistent performance.



Meet today's challenges with us to safeguard the global food supply. Triple quadrupole MS delivers SRM sensitivity and speed to detect targeted compounds more quickly and HRAM solutions using Thermo Scientific™ Orbitrap™ MS enables screening with accurate quantitation of hundreds of contaminants.



Screening, identification, and quantification of food and environmental contaminants are faster and more confident with High-Resolution Accurate-Mass MS/MS Spectral Libraries. The libraries provide high-resolution, accurate-mass (HRAM) MS/MS spectra for identification, confirmation, and quantification of thousands of compounds.

www.thermofisher.com/Pesticides-LCMS

©2015-16 Thermo Fisher Scientific Inc. All rights reserved. mzCloud and the mzCloud logo are trademarks of HighChem LLC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. *Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.*

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic) 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 10 3292 200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand