



**ThermoFisher**  
S C I E N T I F I C


# Modular automation solution for targeted genotyping by sequencing for animal breeding

Christina Buchanan-Wright, PhD  
Sr Product Manager

Chris Willis  
Staff Scientist

The world leader in serving science

# AgriSeq™ targeted GBS is a flexible, powerful genotyping system



Select targets



Construct library



Run sequence



Analyze data



Custom BFX design targeting 50-5000 SNPs, MNPs, or In/Dels

Single tube amplicon based multiplexing

Barcode samples & combine different panels or species in same run

Automated genotype calling & flanking 200 bp SNP identification

**Low sample cost \$5-15/sample**

# AgriSeq targeted GBS generates consistent, complete results

## Metrics

- Sample call rates: ave number of markers called across diverse sample set
- Uniformity: % of bases covered at least 0.2X of the average coverage
- % On Target Read: Specificity of panel to targeted amplicons

Species	Markers	Sample Call Rate	% Uniformity	% on Target reads
Bovine	190	99.7%	98%	98%
Canine	229	99.2%	99%	99%
Feline	62	99.8%	99%	97%
Porcine	3000	96.2%	98%	99%
Salmon	3152	93.9%	94%	99%
Equine	204	99.2%	99%	99%
Cucumber	2804	91.4%	97%	100%
Maize	1079	87.5%	87%	98%
Soybean	1134	98.3%	97%	99%

# AgriSeq can generate 2.6M genotypes/day

## Ion™ 550 Chip & GeneStudio™ Systems

- Reads: 100 -130M
- Read-length: Up to 200 bp
- Run time: 2.5 hours




100X average coverage to ensure maximal genotyping

Flexible to scale different marker vs. sample combinations

How to calculate average coverage:


$$\frac{\text{Reads/chip}}{\text{Samples} \times \text{SNPs}}$$





Recommended coverage:  
 - Average = >100  
 - Min = >20\*

		Number of Samples		
		per Chip	per Day	per Week
Number of markers	1500	768	1536	7680
	2500	460	920	4600
	5000	230	460	2300

# 2-3 day AgriSeq™ workflow with minimal hands on time




10 ng  
of DNA →



<b>Total time</b>	<b>6-7 hours</b>	<b>overnight</b>	<b>2.5 hours</b>	<b>6-24 hours</b>
<b>Hands on time</b>	<b>&lt;3 hours</b>	<b>&lt;15 minutes</b>	<b>&lt;15 minutes</b>	<b>&lt;15 minutes</b>

Workflow is achievable with 1 FTE, but opportunity for automation to reduce hands-on-time, operator fatigue and potential for error

# AgriSeq Library Prep Workflow



- 96 and 384-well protocols
- Best performance with normalized DNA
- Requires thermocycler, plate magnet and AMP Pure beads
- Pool samples post ligation to simplify sample handling, reduce tip usage and clean up reagent costs (<75% AmpPure Bead usage)




# Automation Platform: Gilson PipetMax



- Small footprint: Easily fits on most benchtops
- Low cost: less than \$30,000
- Gilson PipetMan technology: precise, accurate pipetting
- Throughput: Moderate
- Normalization assistant available

# AgriSeq PipetMax Workflow

AgriSeq PipetMAX Script Modules
Amplification Setup
Sample Addition
Pre-Ligation Enzyme Addition
Barcode Addition
Barcode Enzyme Mix Addition
Post Ligation Pooling
Library Purification
Normalization Cleanup





- Modular workflow allows for the separation of processes, flexibility and ease of use
- Step by step setup guide for easy tracking
- Scripts developed for processing either 2 -96 well plates or 1-384 plate
- Hands-on time reduced from 3 hours to 1 hour



- Each AgriSeq automated workflow was validated with two different custom GBS panels:
  1. Bovine ISAG SNP Parentage Panel (2013)
    - 200 amplicon panel
    - Marker chosen by ISAG for bovine parentage verification
    - 96 diverse samples
  2. Porcine Imputation Panel
    - 1500 amplicon SNP panel
    - Markers chosen by a collaborator for porcine genotyping
    - 96 diverse samples

# High Performance for Manual Workflows


- 96 diverse samples were tested with the AgriSeq 384-well and AgriSeq 96-well protocols in replicates (n=4).




- 99.9% reproducible genotype calls between runs
- 99.4% concordance with Illumina Array Data (n=8800 datapoints)

# Similar performance between manual and automated workflows

## Bovine 200 Marker



## Porcine 1500 Marker




Run	Mean Genotyping Call Concordance	Stdev
Manual run	99.7%	0.7%
Gilson Agriseq-96 Run1	99.7%	0.7%
Gilson AgriSeq-96 Run2	99.1%	1.7%
Gilson AgriSeq-96 Run3	99.4%	1.2%

Run	Mean Genotyping Call Concordance	Std Dev
Manual run	99.5%	1.0%
Gilson Agriseq-384 Run1	99.1%	2.0%
Gilson AgriSeq-384 Run2	98.3%	2.5%
Gilson AgriSeq-384 Run3	98.5%	2.1%

\*No significant difference was observed in other sequencing metrics included % on target reads and uniformity between manual and automated processing

# Identification of Flanking SNPs by AgriSeq Workflow

AgriSeq™ GBS generates ~200 bp of flanking sequence



	Bovine	Porcine
Number of Targeted Markers	200	1500
Number of Samples	96	96
Number of Non-Hot Spot Variants	38,888	468,864
Additional Unique Markers	495	1221

**If desired, additional variants can be utilized for additional analysis**

- Discovery of novel variants
- Utility for additional linkage analysis
- Generation of micro-haplotypes for improved marker specificity or enhanced discrimination in parentage/traceability analysis


# Conclusions

- The AgriSeq library prep workflow along with custom GBS panels provide a high throughput, cost-effective method for animal genotyping.
- AgriSeq library prep can be easily automated on the Gilson PipetMax (or other platforms) to improve lab efficiency while reducing risk of technical errors while maintaining high call rates & genotyping concordance.
- Thermo Fisher Scientific offers bioinformatics services, wet lab validation and custom automation support for AgriSeq GBS.
- Visit our booth to learn more about our portfolio of genotyping tools to advance your research.





# Acknowledgements



**Thermo Fisher Scientific Agrigenomics Team (Austin, Texas)**

<b>Chris Willis</b>	<b>Haktan Suren</b>	<b>Prasad Siddavatam</b>	<b>Rick Conrad</b>	<b>Michelle Swimley</b>	<b>Jie Lu</b>	<b>Angela Burrell</b>
---------------------	---------------------	--------------------------	--------------------	-------------------------	---------------	-----------------------

Ryan Ferretti (Neogen Corporation) - Concordance Data

© 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.  
For research use only.