



Molecular biology

Prepare for discovery

Molecular biology workflow solutions





Molecular biology solutions fit for discovery

In your pursuit to advance science, every experiment matters. There's no time to start again.

This handbook is intended to guide you by providing technical information and clear choices across the molecular biology workflow. Applied Biosystems™ and Invitrogen™ products incorporate the latest innovations to enable faster results, more assurance, and less optimization in your lab studies.

Set the foundation and explore solutions fit for discovery, from sample preparation to reverse transcription, PCR, and cloning. With this guide, there's no need to wonder if the products you choose will set you back or propel you forward.

Find additional information at [thermofisher.com/pcrandcloning](https://www.thermofisher.com/pcrandcloning).

Contents



Sample preparation

Tools and methods	7
Automation platform	10
DNA isolation	11
DNA and RNA isolation	14



Reverse transcription

Considerations	17
Reagent selection	18
Genomic DNA removal	19
Primers	20



PCR

Thermal cycler considerations	23
Plastics essentials	26
PCR enzymes	28
Oligo design and selection	32



Electrophoresis

Workflow	35
E-Gel product selection	36
Electrophoresis reagents	38



Cloning

Technologies overview	41
Restriction enzyme cloning	42
PCR cloning	44
Cloning with synthetic DNA	46
Transformation	47



Isothermal amplification

LAMP overview	51
Product selection	51



Resources

Educational resources	55
Mobile apps	56
Custom Commercial Supply	56
FAQs	57
Ordering information	59

Sample preparation






Nucleic acid isolation is a crucial first step in the molecular biology workflow, whether you are isolating genomic DNA (gDNA) or RNA. Selecting nucleic acid purification products that are optimized to provide maximum yield, purity, and integrity from virtually any sample type and application is important for your research success.

Advance your research at thermofisher.com/kingfisher




Portfolio of solutions for your nucleic acid isolation



DNA type

For gDNA extraction, plasmid isolation, and DNA cleanup



RNA type

For purification of total RNA, transcriptome RNA, messenger RNA (mRNA), microRNA (miRNA) and other small RNA, and sequence-specific RNA capture

Learn more at [thermofisher.com/kingfisherkits](https://www.thermofisher.com/kingfisherkits)



Common nucleic acid isolation methods




Automated purification instruments: automated processing of magnetic particles in a microplate format (e.g., Thermo Scientific™ KingFisher™ purification systems)

Samples are processed by moving magnetic beads (not liquid). The system utilizes magnetic rods covered with a disposable, specially designed tip comb and plates. The instrument functions without any dispensing or aspiration parts or devices. Before the run, samples and reagents, including magnetic particles, are dispensed into plates according to default protocols that are installed on the instrument.

Benefits:

- Process 6–96 samples/run
- 24- or 96-well plates for different input volumes
- Easily edit, modify, or create new protocols
- All benefits of magnetic beads (below)



Magnetic beads: 0.5–1.0 μm particles with a paramagnetic core and modified shell (e.g., Applied Biosystems™ MagMAX™ kits and Invitrogen™ Dynabeads™ magnetic beads)

Samples are lysed in solution and allowed to bind nucleic acid to magnetic particles based on specific surface modifications. Application of an external magnetic field rapidly collects the particles. Rounds of release, wash, and recapture enable purification of the desired nucleic acid.

Benefits:

- No risk of clogging
- Increased target capture efficiency
- Rapid collection and concentration of sample
- Specialized equipment not required
- Scalability

Learn more at thermofisher.com/sampleprep





Spin columns: Glass fiber, derivatized silica, or ion exchange membrane in column (e.g., Thermo Scientific™ GeneJET™ and Invitrogen™ PureLink™ kits)

Samples are lysed and passed through the membrane using centrifugal or vacuum force. Wash and elution solutions are subsequently passed through the membrane, and the sample is collected into a tube by centrifugation.

Benefits:

- Convenience
- Ease of use
- Throughput flexibility
- Specialized equipment not required



Organic extraction: Phenol-chloroform solution (e.g., Invitrogen™ DNAzol™ and TRIzol™ Reagents)

After homogenizing the sample with TRIzol Reagent, chloroform is added, and the mixture separates into a clear upper aqueous layer containing RNA, an interphase layer, and a pink lower organic layer containing the DNA and protein. RNA is precipitated from the upper aqueous layer with isopropanol. DNA is precipitated from the interphase and organic layers with ethanol. Protein is precipitated from the phenol-ethanol supernatant with isopropanol.

Benefits:

- Efficient lysis of cells and tissue
- Rapid denaturation of nucleases
- Stabilization of nucleic acids
- Great for fatty and cartilaginous samples

Automation platform: Find a model that meets your needs

Optimize and automate your DNA and RNA purification workflow with KingFisher systems. When used with compatible bead-based reagents, such as MagMAX and Dynabeads products, these versatile instruments enable the automation of DNA, RNA, protein, and cell isolation procedures. Learn more and request a demo at thermofisher.com/kingfisher.



KingFisher instrument:	Duo Prime	Flex	Apex	Presto
Instrument size	Compact benchtop	Benchtop	Benchtop	Benchtop—integrates with robotic liquid handler
Throughput level	Low to medium	High	High	Ultrahigh
Processing volume range	<ul style="list-style-type: none"> • 50–1,000 µL: 12-pin magnet head • 200–5,000 µL: 6-pin magnet head 	<ul style="list-style-type: none"> • PCR plate (20–100 µL*), skirted • 20–200 µL: 96-well plate • 50–1,000 µL: 96 deep-well plate • 200–5,000 µL: 24 deep-well plate 	<ul style="list-style-type: none"> • 15–1,000 µL: 96 deep-well plate • 15–200 µL: 96-well KingFisher standard plate • 10–80 µL: 96-well PCR plate • 30–5,000 µL: 24 deep-well plate • 30–200 µL: 96 storage tubes • 200–1,000 µL: 24 storage tubes 	<ul style="list-style-type: none"> • 50–1,000 µL: 96 deep-well plate • 200–5,000 µL: 24 deep-well plate • KingFisher 96 plate: 50–150 µL
Samples per run	6 or 12	24 or 96	24 or 96	24 or 96
Customizable protocols	Yes	Yes	Yes, with touchscreen or PC software	Yes
Heating/cooling	<ul style="list-style-type: none"> • 10°C to 75°C (plate row block A) • 4°C to 75°C (elution strip block) 	<ul style="list-style-type: none"> • From 5°C above ambient temperature to 115°C 	<ul style="list-style-type: none"> • From 4°C above ambient temperature to 100°C • Cooling down to 4°C 	<ul style="list-style-type: none"> • From 5°C above ambient temperature to 115°C
Ultraviolet lamp	8 watts (up to 16 hr)	No	2 UV lamps, max 23 h 59 min	No
Additional details	For Research Use Only	For Laboratory Use	For Laboratory Use	For Laboratory Use

* Or similar skirted PCR plate.

Typical run with a KingFisher instrument



Prep plates



Select program



Load plates



Press start



25–120 min**

Run time



Resource

Use our selection tool to find the right magnetic bead-based kit for your automated sample preparation.

Find out more at

thermofisher.com/kingfisherkits

** Can vary depending on application and instrument.

The above graphic shows the expected run time for sample preparation when using a KingFisher instrument and prefill plates. It takes less than 15 minutes to prep plates, 1 minute to select the program or protocol, 1 minute to load plates, press start, and walk away for 25–120 minutes, depending on the sample and analyte type.

Learn more at thermofisher.com/kingfisher



Selecting the right DNA isolation kits for your downstream research

Cancer research: Applied Biosystems™ MagMAX™ cell-free nucleic acid kits on KingFisher instruments are ideal for liquid biopsy research. They are optimized specifically for enrichment of cfDNA and total nucleic acid (cfTNA), not gDNA, which means increased recovery and lower starting volumes.

Application	Automation-ready extraction kit and reagents	Cat. No.
Cancer research	MagMAX Cell-Free DNA Isolation Kit	A29319
	MagMAX <i>mirVana</i> Total RNA Isolation Kit	A27828
	MagMAX FFPE DNA/RNA Ultra Kit	A31881
	MagMAX DNA Multi-Sample Ultra 2.0 Kit	A36570
	Dynabeads FlowComp Human CD3 Kit	11365D
	Dynabeads FlowComp Human CD4 Kit	11361D
	Dynabeads FlowComp Human CD8 Kit	11362D
	Dynabeads Untouched Human T Cells Kit	11344D
	Dynabeads Untouched Human CD4 T Cells Kit	11346D
	Dynabeads CD15	11137D
Dynabeads CD14	11145D	

DNA research: Using Applied Biosystems™ MagMAX™ DNA Multi-Sample Ultra 2.0 Kits and KingFisher instruments, you can isolate gDNA from 50 µL to 2 mL of whole blood, saliva, buffy coat, buccal swabs, or other biological samples. The resulting purified gDNA is ideal for many downstream molecular biology applications such as real-time PCR (qPCR), next-generation sequencing (NGS), and microarray analysis.

Application	Automation-ready extraction kit and reagents	Cat. No.
Genomics	MagMAX DNA Multi-Sample Ultra 2.0 Kit	A36570
	MagMAX Viral/Pathogen II Nucleic Acid Isolation Kit	A48383R
	MagMAX <i>mirVana</i> Total RNA Isolation Kit	A27828
NGS	Dynabeads Streptavidin for Target Enrichment	65606D



Go directly to our kit selection tool at thermofisher.com/kingfisherkits



Selecting the right DNA isolation kits for your downstream research

Infectious disease research: Applied Biosystems™ MagMAX™ viral/pathogen kits on KingFisher instruments provide a sensitive and simple method for nucleic acid extraction from samples containing viruses or other pathogens.

Application	Automation-ready extraction kit and reagents	Cat. No.
Infectious disease research	MagMAX Viral/Pathogen Nucleic Acid Isolation Kit for HIV-1 dried blood spots	A53770
	MagMAX Viral/Pathogen II Nucleic Acid Isolation Kit	A48383R
	MagMAX Microbiome Ultra Nucleic Acid Isolation Kit, with bead tubes	A42358
	MagMAX DNA Multi-Sample Ultra 2.0 Kit	A36570
	Dynabeads Intact Virus Enrichment	10700D

Plasmid purification: Choose from our wide range of high-performing, cost-effective Thermo Scientific™ and Invitrogen™ kits for plasmid DNA isolation designed to isolate plasmid DNA at the purity and scale you need.

Application	Technologies for plasmid DNA isolation	Cat. No.
Plasmid purification	GeneJET Plasmid Miniprep Kit	K0502
	PureLink HiPure Plasmid Filter Maxiprep Kit	K210016
	PureLink Fast Low-Endotoxin Midi Plasmid Purification Kit	A35892
	GeneJET Endo-Free Plasmid Maxiprep Kit	K0861
	PureLink Expi Endotoxin-Free Maxi Plasmid Purification Kit	A33073



Go directly to our kit selection tool at thermofisher.com/sampleprep



Selecting the right DNA isolation kits for your downstream research

Wastewater DNA and RNA for disease surveillance: Applied Biosystems™ MagMAX™ wastewater kits offer an efficient and simple method for extracting high-quality nucleic acids from wastewater, sewage, or sludge samples for disease surveillance workflows. The purified DNA and RNA is ideal for use in a variety of downstream applications such as quantitative qPCR, digital PCR, or NGS.

Application	Automation-ready extraction kit and reagents	Cat. No.
Surveillance	MagMAX Wastewater Ultra Nucleic Acid Isolation Kit with Virus Enrichment	A52610
	MagMAX Viral/Pathogen II Nucleic Acid Isolation Kit	A48383R
	MagMAX <i>mirVana</i> Total RNA Isolation Kit	A27828
	MagMAX Microbiome Ultra Nucleic Acid Isolation Kit, with bead tubes	A42358

Streptavidin: Utilization of Invitrogen™ Dynabeads™ streptavidin magnetic beads offers the best balance of capacity and yield, reproducibility, purity, and cost for smaller-scale isolation of specific proteins (e.g., immunoprecipitation (IP)) and protein complexes (co-immunoprecipitation (co-IP)).

Application	Automation-ready extraction kit and reagents	Cat. No.
IP / co-IP	Dynabeads Streptavidin for Target Enrichment	65606D
	Dynabeads M-270 Streptavidin	65305
	Dynabeads M-280 Streptavidin	11205D
	Dynabeads MyOne Streptavidin T1	65601
	Dynabeads MyOne Streptavidin C1	65001

Go directly to our kit selection tool at thermofisher.com/kingfisherkits



Selecting the right DNA and RNA isolation kits

For the quality and performance you need, a full suite of products for DNA and RNA isolation is available for a wide range of sample types, throughputs, and input quantities. To use our online kit selection guide, go to thermofisher.com/rnaselection.



Applied Biosystems™ and Invitrogen™ technologies for DNA and total RNA isolation

Capabilities	Process a large amount of tissue	Fast isolation of RNA from a variety of samples	High-throughput purification of RNA and DNA	Process cells for gene expression
Kits	TRizol reagents	PureLink kits	MagMAX kits	Cells-to-C_T kits
Prep time	30–60 min	<20 min	45 min	≤10 min
Sample types	Most samples, particularly those more difficult to lyse	Bacteria, liquid, blood, cells, yeast, plants, tissue	Blood, plants, saliva, urine, stool, soil, plasma, serum*	Cultured cells
Starting material	100 mg of tissue or 10 ⁷ cells	10 ⁸ cells, 200 mg of tissue, 250 mg of plant tissue, 0.2 mL of blood, 5 x 10 ⁸ yeast, 10 ⁹ bacteria	Variable depending on sample	1–100,000 cells
Yield	10 ⁶ epithelial cells: 8–15 µg 100 mg tobacco leaf: 73 µg (variable depending on sample)	Up to 350 µg	Variable depending on sample	NA
High throughput-compatible	No	Yes	Yes	Yes
Technology	Organic extraction	Silica membrane spin column/filter plate	Magnetic beads	Crude lysate

* Specialty kits with optimized chemistry are available for extraction of cell-free DNA/RNA, total RNA, gDNA, and total nucleic acid.



Helpful tip

If you are not ready to process your RNA sample, simply store it in Invitrogen™ RNA_{later}™ Stabilization Solution for use at a later time. Visit thermofisher.com/stabilizerna.



Learn more at thermofisher.com/rnapreps





Reverse transcription





Reverse transcription is the synthesis of complementary DNA (cDNA) from an RNA template using reverse transcriptase enzymes. The cDNA can be used as a template for PCR amplification, cDNA library construction, RNA sequencing, and more. Selecting the right reverse transcriptase is critical to detecting low-abundance RNAs in a sample and obtaining high yields of full-length cDNA.

Find technical resources on reverse transcription at
[thermofisher.com/rteeducation](https://www.thermofisher.com/rteeducation)



Considerations for selecting the right reverse transcriptase

Sensitivity, thermostability, processivity, and inhibitor tolerance of reverse transcriptases all affect the quantity and length of cDNA synthesized.

Sensitivity

The ability of a reverse transcriptase to generate cDNA from the least amount of input RNA is an important attribute when working with low-copy genes or difficult sample sources where RNA may have already degraded.

Thermostability


Thermostable reverse transcriptases allow reactions to occur at higher temperatures, which help denature RNA with complex secondary structures or high GC content, for generation of longer cDNA, higher cDNA yields, and better coverage of RNA populations in the cDNA.

Processivity

Processivity is the enzyme's ability to add consecutive nucleotides without releasing the template. Highly processive reverse transcriptases allow synthesis of longer cDNA strands in a shorter reaction time, and overall better efficiency in making full-length cDNA.

Inhibitor tolerance

Compounds that have inhibitory effects on reverse transcriptases are common in RNA samples even after purification. Their sources include reagents used for RNA isolation and contaminants carried over from biological samples. Reverse transcriptases resistant to common inhibitors help minimize inconsistent or suboptimal results in cDNA-based assays.



(A) Sensitivity, (B) thermostability, (C) processivity, and (D) inhibitor tolerance of reverse transcriptases can affect the quantity and length of cDNA.

Learn more at thermofisher.com/reverse-transcription



Reverse transcription reagent selection guide

	Ability to optimize reaction components and conditions	Complete cDNA synthesis kit with all reaction components	Most convenient and fewest pipetting steps for RT-qPCR applications	Most convenient and fewest pipetting steps for RT-PCR applications	Go from mammalian cell lysate to cDNA synthesis without isolating RNA	cDNA synthesis and amplification directly from intact single cells or low amounts of total RNA
Product format	Stand-alone enzyme	First-strand cDNA synthesis kit	First-strand cDNA synthesis master mix for RT-qPCR	One-step RT-PCR kit	Direct RT kit	cDNA PreAmp kit
Recommended product	Invitrogen™ SuperScript™ IV Reverse Transcriptase	Invitrogen™ SuperScript™ IV First-Strand Synthesis System	Invitrogen™ SuperScript™ IV VIL0™ Master Mix	Invitrogen™ SuperScript™ IV UniPrime™ One-Step RT-PCR System	Invitrogen™ SuperScript™ IV CellsDirect™ cDNA Synthesis Kit	Invitrogen™ SuperScript™ IV Single Cell/Low Input cDNA PreAmp Kit
Applications	RT-PCR, RT-qPCR, cloning, cDNA library construction, RACE, RNA-Seq	RT-PCR, RT-qPCR, cloning, cDNA library construction, RACE, RNA-Seq	RT-qPCR	RT-PCR	RT-PCR, RT-qPCR	RT-qPCR, RNA-Seq
Total RNA input	1 pg–5 µg	1 pg–5 µg	0.01 pg–2.5 µg	0.01 pg–1 µg	1–10,000 cells	1–1,000 cells or 2 pg–10 ng
Optimal reaction temperature	50–55°C	50–55°C	50–55°C	50–55°C	50–55°C	50°C
Reverse transcription time	10 min	10 min	10 min	10 min	10 min	30 min
High cDNA yield with challenging or degraded RNA	•	•	•	•	•	
Includes PCR step				•		•
Available formats	<ul style="list-style-type: none"> • Stand-alone enzyme 	<ul style="list-style-type: none"> • cDNA synthesis kit • cDNA synthesis kit with ezDNase enzyme 	<ul style="list-style-type: none"> • SuperScript IV VIL0 Master Mix • SuperScript IV VIL0 Master Mix with ezDNase enzyme 	<ul style="list-style-type: none"> • SuperScript IV UniPrime One-Step RT-PCR System (colored) • SuperScript IV UniPrime One-Step RT-PCR System (dye-free) 	<ul style="list-style-type: none"> • SuperScript IV CellsDirect cDNA Synthesis Kit 	<ul style="list-style-type: none"> • SuperScript IV Single Cell/Low Input cDNA PreAmp Kit



Did you know?

The standard enzyme format is incompatible with lyophilization because of the glycerol in the storage buffer. The lyo-ready (lyophilization-ready) format of SuperScript reverse transcriptases has a glycerol content below 0.1% and offers greater stability for lyophilized molecular assay kits. Learn more at thermofisher.com/lyoreadyenzymes.

Learn more at thermofisher.com/superscript




Genomic DNA removal

RNA purification methods, including protocols with DNase digestion on column, often fail to completely remove gDNA. Amplification of contaminating gDNA can cause inaccurate results. Traditional gDNA decontamination protocols with DNase I include time-consuming DNase inactivation or gDNA removal steps under conditions that can damage RNA and affect results.

SuperScript IV VILO Master Mix is available in a format with the novel dsDNA-specific Invitrogen™ ezDNase™ Enzyme, which enables efficient, fast, and gentle gDNA removal from RNA samples to help ensure high confidence in RT-PCR and RT-qPCR results.


SuperScript IV VILO Master Mix cDNA synthesis workflow with ezDNase enzyme



~27 minutes



Traditional cDNA synthesis workflow with DNase I



~105 minutes

Learn more at thermofisher.com/4vilo



Reverse transcription primers


The priming strategy you choose for reverse transcription is important for cDNA synthesis efficiency, consistency, and yield. Each primer type has its benefits and drawbacks, depending on the individual target RNA.

For full-length first-strand cDNA synthesis, oligo(dT) primers are recommended because of their specificity for eukaryotic mRNA, and they allow many different targets to be studied from the same cDNA pool. Typically, oligo(dT) primers are strings of 12–20 deoxythymidines. We offer oligo(dT) in different lengths and formats for flexibility in your reverse transcription experiments.

For target mRNA containing strong transcriptional pauses, random primers are better suited because they anneal throughout the target molecules. They are also ideal for nonpolyadenylated RNA, such as bacterial RNA.



Two most common primers used in reverse transcription



Helpful tip

To avoid poly(A) slippage during priming, anchored oligo(dT) primers can be used to anneal to the 5' end of the poly(A) tail of mRNA and prevent priming within the poly(A) tail. Learn more about selection of primers for reverse transcription at thermofisher.com/rteducation.

Learn more at thermofisher.com/rtparameters





PCR



The polymerase chain reaction (PCR) is a scientific method used to make many copies of a specific piece of DNA. The process involves the following steps:

1) Denaturation: The double-stranded DNA is heated to separate the strands into single strands.

2) Annealing: Primers, which are short DNA molecules, bind to specific regions of the target DNA.

3) Extension: DNA polymerase is used to extend the primers along the single strands, creating new copies of the target DNA. This step is performed in a thermal cycler, alternating between high and low temperatures, making millions of copies of the target DNA.

PCR is a powerful tool for various applications in molecular biology.

Find technical and educational resources about PCR at [thermofisher.com/pcreducation](https://www.thermofisher.com/pcreducation)

Thermal cyclers

Thermal cyclers, which automate the heating and cooling cycles required to amplify DNA, play a critical role in the success of PCR. The following are things to consider when selecting a thermal cycler.

Precise temperature control

Thermal cyclers with precise temperature control enable you to quickly and accurately determine optimal annealing temperatures. Several block technologies, including gradient and Applied Biosystems™ VeriFlex™ Blocks temperature control, are available. A VeriFlex Block employs a separate heating and cooling element in each temperature zone, allowing better control and precision of temperatures. Learn more about the technology at thermofisher.com/veriflextechnology.

Reliability

Thermal cyclers should be able to withstand repeated use, environmental stress, and shipping conditions. Component reliability can be tested using robotic assemblies in repeated testing of frequently used instrument components such as the heated lid, touchscreens, and temperature cycling modules. Applied Biosystems™ thermal cyclers adhere to stringent reliability criteria, which are reported at thermofisher.com/thermalcyclerreliability.

Temperature accuracy

Thermal cycler temperature accuracy is a key factor in the success or failure of a PCR reaction. It is particularly important during annealing temperature optimization,

which requires both accuracy and consistency in the thermal cycler block. If the temperature set point of the instrument does not correspond to the actual temperature of the block, further temperature optimization could be required. Review a study of temperature accuracy in a number of models, available at thermofisher.com/thermalcycleraccuracy.

Features

A variety of Applied Biosystems thermal cyclers are available to fit your applications and budget. Certain features may be important to you, depending on your needs. If you perform PCR optimization frequently, you will likely benefit from an instrument with a VeriFlex Block. If you would like to run optimized assays on a new or different thermal cycler, you can save re-optimization time by using a simulation mode.

If you want remote access to your instrument, you will appreciate the convenience of cloud-enabled thermal cyclers. They allow you to design and share protocols, schedule an instrument, start or stop a run, and check run status from anywhere, on any mobile device or desktop computer.

Fleet control



If you manage multiple thermal cyclers and users, you may benefit from a single interface for viewing all instruments at a glance and setting custom permissions by instrument, user, and method. Learn more at thermofisher.com/fleetcontrol.



Helpful tip

Using the right PCR plastics for your application and instrument can improve the reliability of your PCR results. Go to thermofisher.com/findplastics to determine the right PCR plastics for you.

VeriFlex Block temperature zones



Learn more at thermofisher.com/thermalcyclers



Select the Applied Biosystems™ thermal cycler that's right for you



Cloud-enabled



Cloud-enabled

ProFlex™ PCR System

- > Do you share the device with colleagues?
- > Do you expect your throughput needs to change?
- > Do you want to access your instrument remotely?

VeritiPro™ Thermal Cycler

- > Do you perform a lot of optimizations?
- > Do you want to access your instrument remotely?

Key benefits

Ultimate flexibility and throughput

Ultimate performance

Max sample throughput

480,000 reactions

384 reactions

Max block ramp rate

6.0°C/sec

6.0°C/sec

Temperature optimization

6-zone VeriFlex Block on 96-well system
2-zone VeriFlex Block on 3 x 32-well system

6-zone VeriFlex Block on 96-well system

Compatible with Fleet Control Software

Yes

Yes

Looking for a CE-IVD-labeled thermal cycler?
Learn more at [thermofisher.com/veritiprodx](https://www.thermofisher.com/veritiprodx)





SimpliAmp™ Thermal Cycler

- > Do you need an intuitive interface?
- > Do you train new technicians often?
- > Do you want to access your instrument remotely?

Elegantly simple and precise

96 reactions

4.0°C/sec

3-zone VeriFlex Block on 96-well system

Yes

MiniAmp™ Thermal Cycler

- > Do you want an instrument with just the features needed for routine PCR?
- > Do you want to access your instrument remotely?

Routine PCR, elevated

96 reactions

3.5°C/sec

3-zone VeriFlex Block on MiniAmp™ Plus model

Yes

Automated Thermal Cycler

- > Do you want to place your instrument on a robotic platform now or in the future?

Designed for easy robotic integration

384 reactions

3.5°C/sec

None

Yes

Find out more at [thermofisher.com/thermalcyclers](https://www.thermofisher.com/thermalcyclers)



PCR and qPCR plastics, seals, and accessories

Since PCR is a sensitive detection method, PCR plastics must be of high quality and free of contaminants and inhibitors, to help enable optimal performance. Regardless of the plastics format you select, proper fit and uniform heat transfer during thermal cycling are essential.

Manufacturing quality control

Applied Biosystems™ PCR and qPCR plastic consumables are manufactured in world-class facilities dedicated to the production of high-quality molecular biology-grade plastics. After manufacturing, all plastics undergo stringent quality control.

Integrity testing: Every well of every plate is visually inspected and leak-tested. This thorough screening verifies every well is intact to protect all reactions.

Evaporation testing: Samples are run through PCR to test sealing performance. Well liquid volumes are analyzed post-PCR to verify seal integrity. This helps ensure that every production lot conforms to strict tolerances.

Biological testing: Our plastics are biologically tested to certify them as free of DNA, RNase, and PCR inhibitors. We offer plastics for laboratory use that are provided with a PCR certificate for your convenience and documentation.


Construction materials

Applied Biosystems™ MicroAmp™ optical microplates are made of polypropylene for optimal transfer of thermal energy for efficient PCR. A select medical-grade polypropylene is chosen for its exceptional biocompatibility and inert properties.

Applied Biosystems™ MicroAmp™ EnduraPlate™ microplates are constructed with a stronger polycarbonate frame to resist distortion caused by robotic grippers and to better tolerate rapid heating and cooling, while retaining thin-walled polypropylene wells for efficient heat transfer to the reaction mixture. The polycarbonate frames of the plates are available in multiple colors to help with organization and visual monitoring of assays in a high-throughput setting.



MicroAmp optical microplates



Find out more at [thermofisher.com/pcrplastics](https://www.thermofisher.com/pcrplastics)



Applied Biosystems PCR and qPCR plastics are validated and tested for reliability and optimal performance. They are “Engineer Approved” for use with all Applied Biosystems thermal cyclers and qPCR instruments, and are available in a variety of 32-, 48-, 96-, and 384-well plates; tube strips; single tubes; caps; and seals. The table below provides a detailed comparison of each product. Easily find the PCR and qPCR plastics compatible with your instrument using the online selection tool at [thermofisher.com/findplastics](https://www.thermofisher.com/findplastics).

	Small-scale experiments with a few samples	Daily experiments	Complete-workflow experiments—ideal for automation	Automation-compatible
	Single tubes, strips, caps, adhesive film, and accessories	MicroAmp optical microplates	MicroAmp EnduraPlate optical microplates	MicroAmp EnduraPlate optical reaction plates
Formats	<ul style="list-style-type: none"> • Single tubes • Single tubes with caps • 8-strip tubes with caps • 12-strip caps 	<ul style="list-style-type: none"> • 32-well • 48-well Fast • 96-well* • 96-well Fast* • 384-well* 	<ul style="list-style-type: none"> • 96-well • 96-well Fast • 384-well • 96-well full skirted 	<ul style="list-style-type: none"> • 96-well • 96-well Fast • 384-well
DNA-, RNase-, PCR inhibitor-free	Yes	Yes	Yes	Yes
Colors available	Clear, or mixed packs containing red, orange, blue, and green	Clear	Single-color packs (red, blue, green, yellow, or clear) and 5-plate sampler (one of each color)	Clear
Barcode available	No	Yes (1 or 2 sides)	Yes (3 sides)	Yes (3 sides)
Automation-compatible	No	Yes (for those with * above)	Yes	Yes



Did you know?

Low-profile plastics, also referred to as “Fast” tubes or plates, are generally required for fast (0.1 mL) thermal blocks. Fast plastics utilize lower volumes (0.1 mL) than the standard (0.2 mL) tubes or plates. The low profile minimizes the air space above the reaction, helping reduce the effects of evaporation and enhancing thermal conductivity. Learn more about PCR and qPCR plastics at [thermofisher.com/pcrplastics-education](https://www.thermofisher.com/pcrplastics-education).

PCR reagents

DNA polymerase is an essential component for PCR because of its key role in synthesizing new DNA strands. Because of the sensitive and specific nature of PCR, it is important to choose high-quality enzymes and reagents to produce optimal results. The following are things to consider when choosing PCR enzymes.

Specificity

Nonspecific amplification is one of the major hurdles in PCR, since it can drastically impact yield and sensitivity of target amplification. One way to help reduce nonspecific amplification is through the use of a hot-start DNA polymerase, which utilizes an antibody or chemical modification so that the polymerase becomes active only at the high temperature of the denaturation step. In addition to improving specificity, a hot-start DNA polymerase increases yield and allows convenient room-temperature setup for high-throughput applications.

Thermostability

Since thermal cycling is a key feature of the conditions that enable the repetitive chain reaction of amplifying DNA, thermostability of the DNA polymerase to be used is also an important feature. Highly thermostable DNA polymerases are recommended for amplifying GC-rich or long templates that often require prolonged high-temperature reactions.

Fidelity

The fidelity, or proofreading capability, of a DNA polymerase is based on its 3' to 5' exonuclease activity, which corrects misincorporated nucleotides. This function is critical in applications such as cloning, sequencing, and site-directed mutagenesis, for accurate replication of DNA sequences.

Processivity

A DNA polymerase's processivity is defined as the number of nucleotides being incorporated in a single binding event. This property often reflects synthesis rate and speed, as well as affinity for its substrates. Therefore, highly processive DNA polymerases are beneficial to amplify challenging templates such as long, GC-rich, or inhibitor-containing DNA.

Primer annealing temperature

The primer annealing temperature of each DNA fragment to be amplified often needs optimization when designing a PCR protocol. To help simplify annealing and enable co-cycling of PCR assays, consider a DNA polymerase with a reaction buffer that allows a universal annealing temperature of 60°C for primers.



Did you know?

The residual bacterial DNA in recombinant PCR enzymes poses challenges in microbial genome analysis, such as accurately detecting bacterial strains by 16S rRNA gene sequences. To enable confidence and success in microbial PCR assays, choose PCR enzymes with controlled low levels of residual bacterial and human genomic DNA.

Find out more at

thermofisher.com/broad-range-pcr



Helpful tip

Direct PCR is a way to help simplify PCR experiments, save time, and prevent sample loss in the workflow. Direct PCR allows you to amplify target sequences directly from the samples without the need to first isolate and purify the DNA.




Find out more at

thermofisher.com/direct-pcr

Find out more at thermofisher.com/pcrenzymes



Antibody-based hot-start DNA polymerase and its activation in PCR for enhanced specificity



Helpful tip

One of the most common PCR troubleshooting issues is the presence of unwanted bands, or nonspecific amplification. To reduce nonspecific amplification:

- Use hot-start PCR
- Optimize annealing temperature
- Check primer design
- Prevent DNA cross-contamination
- Decrease template and/or primer concentration
- Optimize Mg^{2+} concentration

Choose the right PCR reagent for your research needs

A comprehensive portfolio of PCR enzymes and master mixes is available with the high performance and consistency you need. Start with the selection guide below to find the best enzyme for common PCR applications.

DNA polymerase	Invitrogen™ Platinum™ SuperFi™ II DNA Polymerase	Invitrogen™ Platinum™ II Taq Hot-Start DNA Polymerase	Applied Biosystems™ AmpliTaq Gold™ 360 DNA Polymerase	Invitrogen™ Platinum™ Direct PCR Universal Master Mix
PCR type	High-fidelity PCR	Hot-start PCR	Hot-start PCR	Direct PCR
Capabilities	Highly accurate amplicon sequences, universal primer annealing, robust amplification of difficult targets	Universal primer annealing, fast DNA synthesis, detection of low-abundance targets	Chemical hot start	Detection of target DNA without genomic DNA purification
Technical specifications				
Fidelity compared to Taq polymerase	>300x	1x	1x	1x
Target length	Up to 20 kb*	Up to 5 kb	Up to 5 kb	Up to 8 kb
Hot-start modification	Antibody-mediated	Antibody-mediated	Chemical modification	Antibody-mediated
Speed	15–30 sec/kb	15 sec/kb	60 sec/kb	20 sec/kb
Universal primer annealing	Yes	Yes	No	Yes
Inhibitor tolerance	Yes	Yes	No	Yes
Blunt or 3'-A end	Blunt	3'-A	3'-A	3'-A
Compatible with Applied Biosystems™ TaqMan™ probes	No	Yes	Yes	No
Certified low level of bacterial gDNA	Yes	Yes	Yes	No
Applications				
Cloning and subcloning	•			
Site-directed mutagenesis	•			
GC-rich amplification	•	•	•	•
Template generation for sequencing	•	•	•	•
High-throughput PCR	•	•		•
Long PCR (up to 20 kb)	•			
Genotyping	•	•	•	•
Amplification of samples with suboptimal purity	•	•		•
Colony PCR	•	•	•	•
Multiplex PCR	•	•	•	•
Fast PCR	•	•		•

* Amplification of up to 40 kb fragment sizes is possible, but may require additional optimization of reaction conditions and primer design.



Innovations for superior PCR

PCR enzymes and reagents are continually being improved to help you get to your research destination faster. For example, the latest Platinum DNA polymerases are designed with the following key innovative features.

More robust and versatile


Advanced enzymatic engineering and methodology provide DNA polymerase with fast cycling, high tolerance of PCR inhibitors, and efficient amplification of challenging DNA like GC-rich sequences. These features help you amplify DNA targets confidently with speed and simplicity.

Find out more at thermofisher.com/platinumenzymes

Universal primer annealing

The innovative Platinum PCR buffers enable universal primer annealing at 60°C. This design allows you to co-cycle different PCR assays (instead of running them sequentially), drastically reducing tedious optimization steps and saving time.

Find out more at thermofisher.com/universalannealing




PCR

Direct gel loading

The latest Platinum DNA polymerases are available in a green buffer format that allows direct gel loading and eliminates tedious steps of dye addition, helping reduce pipetting errors. DNA migration is easily tracked with two dyes (blue and yellow) that are readily visible during electrophoresis (the lanes for 5 and 15 min in the figure to the right).



+



Custom DNA oligos


With a complete range of custom-synthesized oligonucleotide primers, probes, and genes, we recognize the need for quality, reliability, and convenience.

Invitrogen™ custom DNA oligos are synthesized on state-of-the-art, automated systems to increase performance, speed, and capacity. Available in a range of synthesis scales, purification options, and modifications, oligos are analyzed by mass spectrometry or capillary electrophoresis to help ensure the quality of the process and end products. This means you will receive high-quality custom DNA oligos quickly and efficiently.

Choose the right oligos and purification methods for your applications at thermofisher.com/oligos.

Best practices for primer design

Good primer design is essential for a successful PCR assay. For design tips, review the infographic below or go to thermofisher.com/primerdesign.



Primer design made easy

Whether you're performing PCR, cloning, or capillary electrophoresis (CE) sequencing, take advantage of the benefits offered by our robust and easy-to-use Primer3-based Invitrogen™ OligoPerfect™ Designer.

- **Speed up**—design primers for up to 50 genes at the same time
- **Store your data**—ability to save your projects
- **Work smarter**—recognizes .txt and .fasta file types
- **Order with ease**—seamlessly integrates with the Invitrogen™ ordering portal

Try the OligoPerfect Designer at thermofisher.com/oligoperfect or visit the oligo utility hub for our full suite of tools and calculators at thermofisher.com/oligotools.

Purification method	Description	Benefit	Application(s)	25 nmol	50 nmol	200 nmol	1 μmol	10 μmol
Desalt (25 nmol: 10–100 nt; 50 nmol: 5–100 nt)	Oligos are processed through a normal-phase chromatography column, which removes salts but not failure sequences	A salt-free DNA solution, ready to use; suitable for many PCR and sequencing applications without further purification	<ul style="list-style-type: none"> • Endpoint PCR • Isothermal sequencing • Fluorescent sequencing 	•	•	•	•	•
Cartridge (50 nmol–1 μmol, 7–55 nt)	Based on reverse-phase chromatography; removes failure sequences from the completed synthesis	Provides full-length sequences needed in some applications	<ul style="list-style-type: none"> • Antisense oligos (ASO) 	NA	•	•	•	NA
HPLC (≥50 nmol, 10–55 nt; long oligo HPLC available)	Reverse-phase high-performance liquid chromatography (HPLC) removes failure sequences or unincorporated labels the same way as cartridge purification	Guarantees highly purified primer required in some applications (≥85% full length)	<ul style="list-style-type: none"> • First-strand cDNA synthesis for generation of libraries • Fluorescent sequencing 	NA	•	•	•	•
PAGE (≥200 nmol, 7–100 nt)	Polyacrylamide gel electrophoresis (PAGE) is a method used to differentiate full-length product from failure sequences based on size and conformation	Provides the highest percentage of full-length oligos (≥85%) required for certain demanding applications such as mutagenesis or adapter production	<ul style="list-style-type: none"> • Gel shift assays • GeneTrapper screening 	NA	NA	•	•	•

Did you know?


Thermo Fisher Scientific offers scales beyond 10 μmol up to kilograms, and employs a team of manufacturing scientists dedicated to custom method development for unique modifications. For information on large-scale and complex project capabilities, visit thermofisher.com/largescaleoligos.

Find out more at thermofisher.com/oligos



Electrophoresis





Nucleic acid electrophoresis is a common technique in molecular biology to separate, identify, quantify, and/or purify nucleic acids. Setting up electrophoresis involves a number of steps to achieve optimal separation and analysis of nucleic acid samples, such as gel preparation, ladder selection, sample visualization, and gel documentation.

Find technical resources on nucleic acid electrophoresis at thermofisher.com/na-electrophoresis-education



Nucleic acid electrophoresis


Choosing the right tools for nucleic acid electrophoresis can significantly improve and accelerate results, enabling you to address downstream applications sooner.

Determining the appropriate gel type and gel concentration is an essential step that will help streamline the separation of nucleic acids. Learn more about convenient reagents for agarose gel electrophoresis, including hassle-free precast Invitrogen™ E-Gel™ agarose gels and pour-your-own Invitrogen™ UltraPure™ agarose reagents, in this section.


If you need ...	Rapid results, quality control, and a safer workflow	High-quality reagents, a versatile workflow, and cost savings
Product	E-Gel precast agarose gels	UltraPure Agarose
Product format	Precast agarose cassettes	Powder
Buffer	Dry—none*	TBE or TAE
Protocol time (approx.)	18 min	120 min
Ready to use	Yes	No
Get more information at	thermofisher.com/egel	thermofisher.com/ultrapure

* Note: This is a dry precast electrophoresis system.

E-Gel electrophoresis system workflow



Traditional DNA electrophoresis workflow



Find out more at thermofisher.com/electrophoresis



Simplify electrophoresis with E-Gel precast agarose cassettes

E-Gel precast gels

Using precast agarose gels can simplify the nucleic acid electrophoresis workflow. E-Gel precast gels are self-contained and ready for use with the agarose, electrodes, and the DNA stain packaged inside a disposable cassette. There are no gels to pour, buffers to make, staining or destaining steps to perform, or gel boxes to assemble. Just load your samples and run.

E-Gel precast gels offer excellent resolution and clarity in ≤ 18 minutes and are ideal for analyzing PCR products, restriction digests, plasmid preparations, and genotyping products. To help simplify cloning workflows, Invitrogen™ E-Gel™ CloneWell™ II gels use a double-comb design to enable recovery of purified DNA for downstream applications, without the need for additional purification kits or steps.

Find out more at thermofisher.com/egel



E-Gel gels—faster and safer workflow


Pathogen detection:

- Detection of potato spindle tuber viroid (PSTVd) RNA using Invitrogen™ SuperScript™ IV RT-LAMP Master Mix and E-Gel gels



RNA analysis:

- Optimizing loading conditions of RNA markers on an Invitrogen™ E-Gel™ system



Mouse genotyping:

- E-Gel electrophoresis systems can aid in genotyping by providing a convenient and accurate platform for visualizing DNA fragments, such as universal target fragments



Helpful tips

E-Gel precast gels are available in a variety of formats for routine and high-throughput applications, with different stains (see page 38) and agarose percentages (0.8%, 1%, 2%, and 4%). To find the right gel for your needs, see the selection guide at thermofisher.com/egelselection.

Choose Invitrogen™ E-Gel™ DNA ladders for precise electrophoresis band analysis with exceptional DNA fragment purity and quality, reduced dye masking, and improved ladder migration on E-Gel precast agarose gels. Find out more at thermofisher.com/egel-ladders.

Find out more at thermofisher.com/egels



Select your E-Gel device

Integrate your electrophoresis running and imaging into a single small device.

Enjoy speed and convenience.




Cloud-enabled
Ethernet connectivity



	Invitrogen™ E-Gel™ Power Snap Plus Electrophoresis System	Invitrogen™ E-Gel™ Power Snap Electrophoresis System
	<ul style="list-style-type: none"> > Do you use electrophoresis often? > Would you like to store your gel images on internal servers or the cloud? > Do you want to perform quantitative analysis of your gels? > Do you value speed? > Would you like to store, share, and analyze gel images online? 	<ul style="list-style-type: none"> > Do you run less than 25 samples at a time? > Do you want a simple and fast solution for your electrophoresis? > Do you need to save bench space?
Key difference	Low- to high-throughput analysis	Low-throughput analysis
Memory	64 GB	32 GB
Connectivity	USB drive, ethernet, Wi-Fi, printer	USB drive
Sample throughput	Up to 96 samples	Up to 22 samples
Applications	Genotyping, fast PCR analysis, routine electrophoresis, and cloning	Routine electrophoresis, cloning
Software	Invitrogen™ iBright™ analysis software	NA
Find out more	thermofisher.com/powersnapplus	thermofisher.com/powersnap

Electrophoresis



Electrophoresis reagents

For pouring your own agarose gels, choosing high-quality agarose, optimized DNA ladders, and improved DNA stains can help you achieve optimal electrophoresis results.

DNA stains

Detection of nucleic acid samples in gels can be improved using fluorescent dyes that are safer and/or more sensitive than ethidium bromide. The Invitrogen™ SYBR™ Safe and SYBR™ Gold stains provide greater safety and/or sensitivity with lower background fluorescence than the conventional ethidium bromide stain.

SYBR Safe stain is specifically formulated to be less hazardous than ethidium bromide and reduces your exposure to UV light. SYBR Safe stain also carries the ACT™ label.

Find out more at thermofisher.com/stains and thermofisher.com/sybrsafe

UltraPure reagents for electrophoresis

Invitrogen™ UltraPure™ reagents are specifically formulated to meet your nucleic acid analysis and purification needs. UltraPure agarose and reagents are made from highly pure biochemicals for maximum reliability and superior performance.

Find out more at thermofisher.com/ultrapure

DNA ladders

Invitrogen™ DNA ladders are available in a wide variety of size ranges (10 to 48,502 bp) and formats for different applications. To create DNA ladders of superior quality, each fragment is purified individually using proprietary chromatography-based technology. Our DNA ladders are stable during prolonged storage at room temperature and after multiple freeze-thaw cycles.

Find out more at thermofisher.com/ladders



Did you know?

Chromatographically purified nucleic acid fragments are considered the gold standard for ladders, since the technology provides higher control over quality, banding pattern, intensity, and quantity for ladder composition.

Learn more at thermofisher.com/na-electrophoresis-education


Fluorescent nucleic acid gel stains

	Safer detection	Ultimate detection
	SYBR Safe stain	SYBR Gold stain
Sensitivity (dsDNA)	Sensitive (>3 ng)	Ultrasensitive (>0.1 ng)
Less hazardous and more environmentally friendly	•	
Improved cloning efficiency	•	•



Cloning





Molecular cloning involves recombinant DNA technologies that insert a DNA sequence of interest into a vector to generate a large number of copies. Traditionally, cloning has been carried out with restriction enzymes and a DNA ligase to form a new vector capable of expressing the gene of interest. In the case of gene synthesis, researchers can obtain their desired DNA directly in a specified vector with just sequence information. Other cloning methods, such as PCR cloning, Invitrogen™ TOPO™ cloning, and gene assembly, are commonplace, exploiting unique characteristics of other DNA-modifying enzymes.

Find technical resources on molecular cloning at
[thermofisher.com/cloningeducation](https://www.thermofisher.com/cloningeducation)



Cloning and gene synthesis

From restriction enzymes to gene synthesis, a large portfolio of tools and resources is available to help you obtain high-quality cloned DNA for your next discovery.

Method	Thermo Scientific™ FastDigest™ restriction enzymes	TOPO cloning	Invitrogen™ Gateway™ cloning	Invitrogen™ GeneArt™ seamless cloning and GeneArt™ Gibson Assembly® cloning kits	Invitrogen™ GeneArt™ Type IIS assembly	Invitrogen™ GeneArt™ Strings™ DNA Fragments	Invitrogen™ GeneArt™ Gene Synthesis
Key benefits/ description	<ul style="list-style-type: none"> Familiarity, flexibility, convenience, time savings Universal protocol and complete digestion in 5–15 minutes in one buffer 100% buffer compatibility with downstream applications Direct loading on gels 	<ul style="list-style-type: none"> >95% efficiency, 5-minute PCR cloning Compatible with many other cloning systems 	<ul style="list-style-type: none"> High-throughput and high-efficiency shuttling among multiple expression vectors 	<ul style="list-style-type: none"> Seamless multifragment assembly by homologous recombination Directional cloning of up to 15 fragments Up to 95% efficiency and 15-minute cloning 	<ul style="list-style-type: none"> One-tube seamless multifragment assembly by simultaneous restriction digestion and ligation Directional cloning of up to 8 fragments, for up to 20 kb total Efficient for repetitive and very small sequences 	<ul style="list-style-type: none"> Synthesized DNA fragments ready to clone via the method of your choice No starting DNA required Pool sequence-verified 	<ul style="list-style-type: none"> Custom-cloned genes in your choice of vector Sequence-verified Can be optimized for a specific host for maximal protein expression
Technology basics	<ul style="list-style-type: none"> Restriction digestion and ligation 	<ul style="list-style-type: none"> Topoisomerase-based, ligase-free cloning 	<ul style="list-style-type: none"> Single-step, directional, and site-specific DNA recombination Restriction enzyme- and ligase-free 	<ul style="list-style-type: none"> End-terminal homology recombination using overlapping sequences Transformation-associated recombination (TAR) in <i>Saccharomyces cerevisiae</i> 	<ul style="list-style-type: none"> Type IIS restriction and ligation in a single reaction 	<ul style="list-style-type: none"> Linear dsDNA assembled from pooled synthetic oligonucleotides 200–3,000 bp, also available in library format with randomized bases 	<ul style="list-style-type: none"> DNA of interest cloned in vector 100% sequence-verified with quality assurance documentation
Needs DNA source material (gene in plasmid, library, etc.)	•	•	•	•	•		
Use your own vector	•		*	•	•	•	•

* Vector needs to be converted with Invitrogen™ Gateway™ Vector Conversion System with One Shot™ *ccdB* Survival™ 2 T1[®] Competent Cells.

Discover more at thermofisher.com/cloning



Restriction enzyme cloning

Found naturally in bacteria, restriction enzymes recognize and cleave specific DNA sequences, resulting in sticky ends (5' or 3' protruding ends) or blunt ends, enabling DNA inserts to be cloned into vectors with compatible ends. Star activity, buffer compatibility, and varying protocols for complete digestion are some common hurdles in restriction digestion.


FastDigest restriction enzymes

To simplify cloning, we offer FastDigest enzymes—an advanced line of restriction enzymes that share buffer compatibility with downstream modifying enzymes.

Benefits include:

- Complete digestion in 5–15 min
- Universal buffer allows for multiple digestions for any combination of enzymes
- No sequential digestions and buffer changes
- 176 unique specificities
- Direct loading of reaction mixture on gels

Find out more at thermofisher.com/fastdigest



Type IIS restriction enzymes

A specific group of restriction enzymes called Type IIS endonucleases cleave DNA outside of their recognition sequences. In combination with DNA ligase, Type IIS restriction enzymes are utilized to drive the insertion of one or several DNA fragments into a recipient vector without the inclusion of residual restriction enzyme sites and other unwanted DNA sequences at fragment junctions (scarless cloning).

Find FastDigest Type IIS enzymes at [thermofisher.com/fastdigesttypeiis](https://www.thermofisher.com/fastdigesttypeiis)

For GeneArt Type IIS Assembly Kits, go to [thermofisher.com/typeiis](https://www.thermofisher.com/typeiis)




PCR cloning

PCR cloning is a method in which double-stranded DNA fragments amplified by PCR are ligated into a vector. With PCR amplification, this cloning technique requires much less starting material for the insert sequence and allows introduction of new restriction and/or recombination sites to the 5' end of the inserts.

TOPO cloning

TOPO PCR cloning technology was developed to help improve cloning efficiency, simplify protocol setup, and accommodate a wide range of PCR insert sizes. TOPO cloning vectors are linearized by the activity of topoisomerase I (which also has a ligase function) that is covalently bound to the 3' phosphate on each end (see figure below). This system enables the vectors to be joined to PCR inserts with compatible ends (with up to 95% efficiency), without the need for additional ligation steps, in 5 minutes.

Find out more at thermofisher.com/topo



Quickly find your TOPO cloning kit with our interactive selection tool. Search by application, vector, or desired competent cells at thermofisher.com/topoguide.




Did you know?

The Invitrogen™ TOPO™ XL-2 Complete PCR Cloning Kit provides all the necessary elements for highly efficient cloning of extra-long PCR products from 1–13 kb. thermofisher.com/topoxl2

Gateway cloning

To shuttle a PCR insert among vectors, the Gateway cloning system offers site-specific, recombinase-based cloning. It maintains the insert's proper orientation and reading frame during shuttling using the Gateway vectors. Once a gene is cloned into an entry clone, you can then move the DNA fragment into one or more destination vectors simultaneously.

Find out more at thermofisher.com/gateway



Gateway cloning system reactions. The scheme shows the four types of plasmids and enzyme mixes involved in Gateway cloning reactions. Red arrows represent the fragment of interest. Adapted from Katzen F (2007) *Expert Opin Drug Discov* 2(4):571–589.



Seamless cloning and GeneArt Gibson Assembly cloning kits


To assemble multiple PCR fragments by end-terminal homologous recombination, several seamless cloning technologies are available for scarless and directional cloning into any vector. GeneArt seamless cloning kits offer the option of building constructs using *E. coli* and *Saccharomyces cerevisiae*.

Invitrogen™ GeneArt™ Gibson Assembly® kits allow for the simultaneous assembly of up to 15 very large DNA fragments to create precise constructs with no additional sequences, in highly efficient reactions. This cloning method circumvents the need for multiple rounds of restriction enzyme analysis and digestion, DNA end repair, dephosphorylation, ligation, enzyme inactivation, and cleanup, and is a powerful tool in synthetic biology.

GeneArt Gibson Assembly kits offer these benefits:

- Assembly of up to 15 fragments to build seamless clones
- Cloning efficiencies up to >95%
- Choice of complete kits with competent cells or master mixes

1. Enzymatic assembly of DNA molecules up to several hundred kilobases. Gibson DG et al. (2009) *Nat Methods* 6(5):343-5.



Did you know?

The Gibson Assembly method has been referenced in thousands of peer-reviewed publications and is a powerful method that can be used to seamlessly construct synthetic and natural genes, genetic pathways, and entire genomes [1].

Discover more at thermofisher.com/seamless



Table of contents

Cloning with synthetic DNA

If you lack the time to generate and clone insert DNA, including optimization and troubleshooting, our synthetic DNA fragments and cloning service might be right for you. GeneArt Strings DNA Fragments and GeneArt Gene Synthesis offer genes analogous to optimized, error-free PCR products.

GeneArt Strings DNA Fragments

A time-saving alternative to PCR, GeneArt Strings DNA Fragments are available in lengths up to 3 kb and are compatible with any downstream cloning method of choice, providing:


- Synthetic, ready-to-use DNA fragments
- DNA with your specified ends to facilitate the cloning method of choice
- No starting DNA required
- Free optimization of gene with Invitrogen™ GeneArt™ GeneOptimizer™ software for maximum protein expression
- Option of Strings DNA Libraries with mixed, randomized nucleotides using full IUPAC code

Find out more at thermofisher.com/strings

GeneArt Gene Synthesis

A reliable and cost-effective method for obtaining customized DNA constructs with 100% sequence accuracy, GeneArt Gene Synthesis offers:

- Synthetic, ready-to-transfect genes
- Cloning into several available vectors (custom options available)
- 100% sequence-verified and ready for downstream applications
- No starting DNA required
- Free optimization of gene with GeneOptimizer software for maximum protein expression



Find out more at thermofisher.com/genesyntesis




Transformation

Once the DNA fragment is cloned into a vector, transformation into bacteria is performed to enable propagation of sufficient quantities of the cloned DNA for downstream experiments. Selection of competent cells for transformation depends upon the transformation methods, strain genotypes, plasmid characteristics, and desired applications.

Visit thermofisher.com/compcells-education for technical resources on competent cells.

Choosing Invitrogen™ competent cells based on the application



Find out more at thermofisher.com/compcells

Transformation (cont.)

Medium- and high-throughput transformation

Performing bacterial transformations one by one can be very time-consuming and create a bottleneck in your experimental workflow. There are times when medium- and high-throughput transformation options are desired. Invitrogen™ MultiShot™ chemically competent cells provide three flexible product formats to meet your throughput needs.

Find out more at thermofisher.com/multishot



StripWell format

- Medium-throughput option
- Twelve 8-tube strips
- Suitable for 1–96 transformations
- Five *E. coli* strains available

FlexPlate format

- High-throughput option
- 96-well plate separates into 12 x 8-well segments
- Manual and automated platform transformations
- Six *E. coli* strains available



96-well plate

- Highest-throughput option
- Five 96-well plates
- Available with the TOP10 strain
- Stable replication of high copy number plasmids



Did you know?

Invitrogen competent cells can be provided in custom configurations per your request. Large and custom volumes as well as multiple formats are at your fingertips. Simply email us at customorders@thermofisher.com.





Isothermal amplification



Isothermal amplification is a technique that utilizes enzymes, typically strand-displacing polymerases, to amplify nucleic acid sequences at a constant temperature. This results in continuous and exponential amplification that is not limited by the constraint of thermal cycling.


Find technical and educational resources about isothermal amplification at [thermofisher.com/isothermal](https://www.thermofisher.com/isothermal)



Loop-mediated isothermal amplification

LAMP overview

One commonly used isothermal amplification technique is loop-mediated isothermal amplification (LAMP), which utilizes a set of four to six primers and strand-displacing polymerase, such as Bst DNA polymerase, to amplify target DNA at a constant temperature. The LAMP reaction produces a large amount of DNA amplicons with a characteristic ladder-like pattern that can be visualized by gel electrophoresis or detected by turbidity, fluorescence, or colorimetry.



The amplification process begins with the invasion of an inner primer into the target nucleic acid sequence, followed by extension via a strand-displacing DNA polymerase. As the extension proceeds, the first product is displaced, and an outer primer anneals to the newly synthesized strand, forming a self-hybridizing loop structure. This structure contains multiple sites for amplification initiation and serves as a seed for exponential LAMP reactions.

Find out more at thermofisher.com/lamp

SuperScript IV RT-LAMP Master Mix

SuperScript IV RT-LAMP Master Mix is a reverse transcription loop-mediated isothermal amplification (RT-LAMP)-based solution for faster and simpler detection of various pathogens, including influenza virus, measles virus, *S. enterica*, *S. aureus*, SARS-CoV-2, and other pathogens. Our master mix reagents provide maximum flexibility to optimize and accelerate your pathogen research and surveillance.

Product highlights

- **Fast**—pathogen detection in as little as 5 minutes with evolved Bst DNA polymerase
- **Efficient**—one-step reaction for reverse transcription of RNA to cDNA with SuperScript IV Reverse Transcriptase
- **Sensitive**—greater sensitivity and specificity utilizing Invitrogen™ RNaseOUT™ Recombinant Ribonuclease Inhibitor and an optimized buffer
- **Simple**—streamlined workflow: single tube format, only requires a 65°C heating block
- **Flexible**—several options for evaluating results, including real-time and endpoint detection methods



Find out more [here](#)



Lyo-ready Bst DNA Polymerase

Invitrogen™ Lyo-ready Bst DNA Polymerase is an engineered version of Bst DNA polymerase, large fragment, which shows a significantly faster reaction speed, increased sensitivity, and tolerance to inhibitors.

Lyo-ready Bst DNA Polymerase provides maximum flexibility to optimize your LAMP reaction and works with various types of pathogens, including human adenovirus, measles virus, SARS-CoV-2, and other pathogens.

Product highlights

- **Fast**—amplifies targets in as little as 10 minutes
- **Sensitive**—achieves sensitivity down to 50 copies
- **Robust**—amplifies even from inhibitor-containing RNA/DNA samples
- **Flexible**—provides the ability to optimize your LAMP or RT-LAMP reaction



Tips:

1. Did you know? The term “lyo-ready” refers to an enzyme that is provided in a liquid formulation without glycerol, making it compatible with microfluidics-based systems and various downstream applications such as lyophilization. Furthermore, it maintains the necessary stability and activity levels for direct enzymatic reactions.
2. To minimize nonspecific amplification in LAMP, follow these steps:
 - a. **Prevent cross-contamination:** Use uncontaminated reagents and maintain a clean work environment.
 - b. **Enhance primer design:** Optimize primer sequences for improved specificity.
 - c. **Optimize reaction conditions:** Adjust Lyo-Ready Bst DNA Polymerase amount and reaction time.



Find out more [here](#)





Resources




We have compiled educational resources, mobile apps, frequently asked questions, and other information to help you achieve success in your research.



Educational resources

Suitable for new and experienced molecular biologists alike, our free online educational resources are designed to help you review the basics, build your expertise, or discover our latest innovative technologies. Explore our educational resources in the following areas of molecular biology.

Reverse transcription




thermofisher.com/rteducation

PCR enzymes, plastics, and thermal cyclers




thermofisher.com/pcreducation

Nucleic acid electrophoresis



thermofisher.com/na-electrophoresis-education

Restriction enzymes, molecular cloning, and competent cells



thermofisher.com/cloningeducation



Resources

Webinars: Watch live and recorded webinars for in-depth understanding of molecular and synthetic biology techniques and tools to help elevate your research.

thermofisher.com/mbwebinars

Videos: Experience entertaining and visual learning with our educational videos on molecular biology techniques, how-to guides, tips and tricks, and more.

thermofisher.com/mbvideos

Application notes: Read white papers and application notes from our R&D scientists on our product innovations.

thermofisher.com/mbliterature

Online tools: Use our interactive online tools for PCR annealing temperature, restriction enzyme information, product selection, and more.

thermofisher.com/mbtools

Find out more at thermofisher.com/mbschool



Table of contents

Podcast: Speaking of Mol Bio



Listen now to Speaking of Mol Bio. This podcast highlights trending applications in science and the molecular biology aspects of those applications. Our hosts delve into deep discussions with CEOs, R&D scientists, researchers, and key opinion leaders across the globe.

This podcast helps scientifically curious people—from all scientific and nonscientific backgrounds—understand how modern molecular biology applications can help push the boundaries in medicine, science, drug discovery, and in the cure and treatment of diseases.

Listen now at thermofisher.com/molbiopodcast

Mobile apps



DailyCalcs—science calculator

The DailyCalcs app turns your phone into a science calculator to help simplify everyday tasks in the lab. The app features eight calculators: molarity, dilution, formula weight, transfection, unit conversions, culture vessel data, media conversions, and specific productivity. Find it in your preferred mobile application store.



Instrument Connect—remote monitoring

Instrument Connect allows you to view instrument status, monitor or schedule a run, and more on any cloud-enabled instrument, including the ProFlex, SimpliAmp, and MiniAmp PCR instruments. Access at thermofisher.com/connect.



PCR Quest—match-3 lab game

Test your PCR knowledge with our lab game—PCR Quest—where you travel from lab to lab crushing the world's toughest diseases. Download at thermofisher.com/pcrquest.

Custom Commercial Supply

As a leading supplier of molecular biology products, we offer tailored solutions for companies developing new molecular assays. Whether you're just starting or need a specific solution, we have what you need. Work with an experienced supplier that knows both raw materials and new technologies. Our dedicated business team is here to provide value beyond just our products.

What do our OEM solutions mean to you?

- Customization of products and services
- Consultation, partnership, and expertise
- Negotiated business terms
- Warranties and indemnification
- Commercial-use rights and obligations
- Risk and liability management



Find out more at thermofisher.com/mdx



Frequently asked questions

Below are some common questions and answers to help you start or troubleshoot molecular biology experiments.

Sample preparation

Which kit should I use to isolate nucleic acids from my sample?

Choosing the right product is fundamental to ensuring proper lysis of cells and tissue, as well as sufficient yield and quality of isolated nucleic acids. Look to our selection guides ([see pages 11–13](#)) to help you decide according to nucleic acid type, sample source, experimental throughput, and format as well as downstream applications.

What are the key steps to preventing RNA degradation?

The basic lab precautions listed below can help minimize RNA degradation and avoid experimental inconsistency and failure.

- Use nuclease-free pipette tips and tubes
- Use nuclease-free water and reagents
- Regularly decontaminate work surfaces
- Properly stabilize RNA sources before storage

For more tips and troubleshooting advice on sample prep, visit [thermofisher.com/rnabasics](https://www.thermofisher.com/rnabasics) and [thermofisher.com/napsupport](https://www.thermofisher.com/napsupport).

Reverse transcription

How do I improve the efficiency of cDNA synthesis when working with challenging samples (e.g., low-abundance, degraded, inhibitor-containing, or GC-rich RNA)?

When working with challenging RNA samples, select a reverse transcriptase that is highly sensitive, processive, thermostable, and resistant to common inhibitors, to help you obtain the highest cDNA yield ([see page 17](#)).

What are the benefits of using random primers, oligo(dT) primers, gene-specific primers, or oligo(dT)/random mixed primers in reverse transcription?

- Random primers are good to use with degraded RNA, RNA with high secondary structure, nonpolyadenylated RNA, or prokaryotic RNA.
- Oligo(dT) primers are an optimal choice for synthesis of full-length cDNA from eukaryotic mRNA. Applications include cDNA cloning, cDNA library construction, and 3' rapid amplification of cDNA ends (3' RACE).
- Gene-specific primers are designed based on known sequences of the target RNA. These primers offer the most specific priming and are commonly used in one-step RT-PCR.
- A mixture of oligo(dT) and random primers is often used in two-step RT-PCR to achieve the benefits of each primer type ([see page 20](#)).

For more tips and troubleshooting advice on reverse transcription, visit [thermofisher.com/rteducation](https://www.thermofisher.com/rteducation) and [thermofisher.com/rtsupport](https://www.thermofisher.com/rtsupport).

PCR amplification

How can I optimize primer annealing for PCR?

Traditionally, gradient thermal cyclers have been used to simultaneously assess a number of temperatures around the theoretical annealing point. Compared to gradient thermal cyclers, instruments with the VeriFlex technology allow more precise temperature control for faster optimization of primer annealing ([see page 23](#)).

Tedious optimization steps may be circumvented using the novel Platinum DNA polymerases. Their innovative buffers enable specific annealing at 60°C for most primers when they are designed following general primer design rules ([see pages 30–31](#)).



Frequently asked questions (cont.)

What do I need to run fast PCR?

PCR amplicons shorter than 1 kb can be amplified in as little as 40 minutes using “fast” enzymes (high processivity; **see page 30**), “fast” plastics (low profile and ultra-thin walls; **see page 27**), and “fast” thermal cyclers (fast ramp rate; **see pages 24–25**).

How can I prevent sample evaporation during PCR?

Proper sealing of your reactions will help prevent evaporation during PCR.

- When using adhesive film to seal a plate, be sure to properly align the seal to cover all wells and press firmly along all edges of the plate using an applicator tool.
- When sealing a plate using cap strips, ensure that the cap strips are compatible with the plate and thermal cycler being used. Be sure to align cap strips with each well of the plate and place firmly across the plate for a secure fit.
- Use the applicator tool (Cat. No. 4333183 or 4330015) or other comparable sealing tools as needed.

For more tips and troubleshooting advice on PCR, visit [thermofisher.com/pcreducation](https://www.thermofisher.com/pcreducation) and [thermofisher.com/pcrsupport](https://www.thermofisher.com/pcrsupport).

Nucleic acid electrophoresis

Why is it important to choose the right ladder when using E-Gel precast agarose gels?

Accurate analysis of electrophoresis bands often depends on the DNA ladder chosen for your gel run. E-Gel DNA ladders are formulated with ready-to-use buffers unique for E-Gel precast agarose gels, and DNA standards designed for optimal separation (**see page 36**).

Are there safer alternatives to ethidium bromide for staining nucleic acids in gel electrophoresis?

SYBR Safe DNA Gel Stain is a safer alternative to ethidium bromide and is commonly used in gel electrophoresis. SYBR Safe DNA stain is not classified as hazardous waste or as a pollutant under US federal regulations (**see page 38**).

Do I need a buffer to run the E-Gel system?

No. The E-Gel electrophoresis system does not require electrophoresis buffers like TBE or TAE. E-Gel cassettes already contain everything you will need and are classified as dry electrophoresis.

For more tips and troubleshooting advice on nucleic acid electrophoresis, visit [thermofisher.com/na-electrophoresis-education](https://www.thermofisher.com/na-electrophoresis-education) and [thermofisher.com/na-electrophoresis-support](https://www.thermofisher.com/na-electrophoresis-support).

Cloning

Do you have a buffer compatibility chart for restriction enzymes?

All FastDigest restriction enzymes are 100% active in one universal FastDigest buffer (**see page 42**). Hence, there is no buffer compatibility chart for FastDigest restriction enzymes.

What is the main difference between GeneArt Strings DNA Fragments and GeneArt Gene Synthesis?

GeneArt Strings DNA Fragments are custom-made, uncloned, double-stranded linear DNA fragments. GeneArt Gene Synthesis is a service offered for chemical synthesis, cloning, and sequence verification of genetic sequences (**see page 46**).

What are some key considerations for choosing competent cells for my cloning applications?

Genotype, transformation efficiency, growth rate, and throughput format are important factors in choosing competent cells for cloning. The genotype of a cell strain may determine growth conditions and suitability for transformation with specific DNA types (**see page 47**).

For more tips and troubleshooting advice on cloning, visit [thermofisher.com/cloningeducation](https://www.thermofisher.com/cloningeducation) and [thermofisher.com/cloningsupport](https://www.thermofisher.com/cloningsupport).



Ordering information

	Quantity	Cat. No.
Nucleic acid isolation		
PureLink Quick Plasmid Miniprep Kit	50 preps	K210010
PureLink HiPure Plasmid Filter Midiprep Kit	25 preps	K210014
PureLink HiPure Plasmid Filter Maxiprep Kit	10 preps	K210006
PureLink <i>Pro</i> Quick96 Plasmid Purification Kit	4 x 96 preps	K211004A
PureLink Quick Gel Extraction Kit	50 preps	K210012
TRIzol Plus RNA Purification Kit	50 preps	12183555
PureLink RNA Mini Kit	10 preps	12183020
PureLink Genomic DNA Mini Kit	10 preps	K182000
PureLink <i>Pro</i> 96 Genomic DNA Kit	4 x 96 preps	K182104A
PureLink <i>Pro</i> 96 Viral RNA/DNA Purification Kit	4 plates	12280096A
PureLink Viral RNA/DNA Mini Kit	50 preps	12280050
PureLink Genomic Plant DNA Purification Kit	50 preps	K183001
MagMAX DNA Multi-Sample Ultra Kit	500 preps	A25597
MagMAX Microbiome Ultra Nucleic Acid Isolation Kit, with bead tubes	100 preps	A42358
MagMAX Microbiome Ultra Nucleic Acid Isolation Kit, with bead plate	100 preps	A42357
KingFisher Flex Purification System with 96 Deep-Well Head	1 system	5400630
KingFisher Duo Prime Purification System	1 system	5400110
PureLink PCR Purification Kit	50 preps	K310001
PureLink Quick Gel Extraction and PCR Purification Combo Kit	50 preps	K220001
PureLink Quick Gel Extraction Kit	50 preps	K210012
Dynabeads M-270 Streptavidin	2 mL	65305
Dynabeads MyOne Streptavidin C1	2 mL	65001
KingFisher instruments		
KingFisher Apex Purification System with 96 PCR Head**	1 system	5400910
KingFisher Apex Purification System with 96 Combi Head**	1 system	5400920
KingFisher Apex Purification System with 96 Deep-Well Head**	1 system	5400930
KingFisher Apex Purification System with 24 Combi Head**	1 system	5400940
KingFisher Flex Purification System with 96 PCR Head**	1 system	5400610
KingFisher Flex Purification System with 96 KingFisher Head**	1 system	5400620
KingFisher Flex Purification System with 24 Deep-Well Head**	1 system	5400640
KingFisher Flex Purification System with 96 Deep-Well Head**	1 system	5400630
KingFisher Duo Prime Purification System	1 system	5400110
KingFisher Presto Purification System with 24 Deep-Well Head**	1 system	5400840
KingFisher Presto Purification System with 96 Deep-Well Head**	1 system	5400830

** For Laboratory Use.

	Quantity	Cat. No.
Reverse transcription		
SuperScript IV Reverse Transcriptase	2,000 units	18090010
	10,000 units	18090050
SuperScript IV First-Strand Synthesis System	50 reactions	18091050
	200 reactions	18091200
SuperScript IV VIL0 Master Mix	50 reactions	11756050
	500 reactions	11756500
SuperScript IV VIL0 Master Mix with ezDNase Enzyme	50 reactions	11766050
	500 reactions	11766500
SuperScript IV One-Step RT-PCR System	25 reactions	12594025
	100 reactions	12594100
SuperScript IV UniPrime One-Step RT-PCR System	100 reactions	12597100
	50 reactions	11750150
SuperScript IV CellsDirect cDNA Synthesis Kit	500 reactions	11750350
	500 reactions	11750550
SuperScript IV CellsDirect Lysis Reagents	500 reactions	11750550
RNaseOUT Recombinant Ribonuclease Inhibitor	5,000 units	10777019
Ribonuclease H	30 units	18021014
Random Hexamers (50 µM)	5 nmol	N8080127
Random Primers	9 A ₂₆₀ units	48190011
Oligo(dT) ₁₂₋₁₈ Primer	25 µg	18418012
Oligo(dT) ₂₀ Primer	15 µg	18418020
DNase I, Amplification Grade	100 units	18068015
Isothermal amplification		
Lyo-ready Bst DNA Polymerase	1,200 units (6 U/µL)	A56655
	6,000 units (6 U/µL)	A56656
SuperScript IV RT-LAMP Master Mix	1,200 units (40 U/µL)	A56657
	100 reactions	A51801
	400 reactions	A51802
	1,000 reactions	A51803



Ordering information (cont.)

	Quantity	Cat. No.
PCR		
Thermal cyclers		
ProFlex 3 x 32-Well PCR System	1 instrument	4484073
ProFlex 96-Well PCR System	1 instrument	4484075
SimpliAmp Thermal Cycler	1 instrument	A24811
VeritiPro Thermal Cycler, 96 well	1 instrument	A48141
MiniAmp Plus Thermal Cycler	1 instrument	A37835
MiniAmp Thermal Cycler	1 instrument	A37834
Automated Thermal Cycler, 96 well	1 instrument	A31486
Plastics		
MicroAmp EnduraPlate Optical 96-Well Fast Multicolor Reaction Plates with Barcode	5 plates	4483493
MicroAmp Optical Adhesive Film	100 covers	4311971
MicroAmp Optical 96-Well Reaction Plate	10 plates	N8010560
MicroAmp Optical 8-Cap Strips	300 strips	4323032
MicroAmp Fast Optical 96-Well Reaction Plate, 0.1 mL	10 plates	4346907
MicroAmp Fast Reaction Tube with Cap, 0.1 mL	1,000 tubes	4358297
MicroAmp EnduraPlate Optical 384-Well Multicolor Reaction Plates with Barcode	5 plates	4483316
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plates with Barcode	20 plates	4483354
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	20 plates	A32811
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	125 strips	A30589
MicroAmp EnduraPlate Optical 96-Well Full-Skirted Plates with Barcode, clear	50 plates	A31728
PCR enzymes		
Platinum II <i>Taq</i> Hot-Start DNA Polymerase	100 reactions	14966001
	500 reactions	14966005
Platinum II Hot-Start PCR Master Mix (2X)	50 reactions	14000012
	200 reactions	14000013
Platinum II Hot-Start Green PCR Master Mix (2X)	50 reactions	14001012
	200 reactions	14001013
AmpliAmp Gold 360 DNA Polymerase	100 units	4398813
	250 units	4398823
AmpliAmp Gold 360 Master Mix	1 mL	4398876
	5 mL	4398881
Platinum SuperFi II DNA Polymerase	100 units	12361010
	500 units	12361050

* North America, Europe, Middle East, and Africa.

** Asia Pacific, Japan, Latin America, and greater China.

	Quantity	Cat. No.
Platinum SuperFi II PCR Master Mix	100 reactions	12368010
	500 reactions	12368050
Platinum SuperFi II Green PCR Master Mix	100 reactions	12369010
	500 reactions	12369050
Platinum Direct PCR Universal Master Mix	100 reactions	A44647100
	500 reactions	A44647500
dNTP Set (100 mM)	4 x 250 µL	10297018
	8 x 1.25 mL	10297117
Oligos		
DNA Oligo, Desalted, Dry	25 nmol	A15612
DNA Oligo, Desalted, Dry, next-day (ordered before 1 PM Eastern Time)	25 nmol	A15613
DNA Oligo, Desalted, Liquid	25 nmol	A15611
DNA Oligo, Desalted, Dry	50 nmol	A15610
DNA Oligo, Desalted, Liquid	50 nmol	A15609
DNA Oligo, Cartridge, Dry	50 nmol	A15614
DNA Oligo, Cartridge, Liquid	50 nmol	A15608
DNA Oligo, HPLC, Dry	50 nmol	A15607
DNA Oligo, HPLC, Liquid	50 nmol	A15606
DNA Oligo, PAGE, Dry	50 nmol	A15605
DNA Oligo, PAGE, Liquid	50 nmol	A15604
Nucleic acid separation and analysis		
SYBR Safe DNA Gel Stain	400 µL	S33102
SYBR Gold Nucleic Acid Gel Stain	500 µL	S11494
UltraPure DNase/RNase-Free Distilled Water	500 mL	10977015
UltraPure Agarose	100 g	16500100
TrackIt 100 bp Plus DNA Ladder	100 applications	10488058
UltraPure TAE Buffer, 10X	4 L	15558026
E-Gel Agarose Gels with SYBR Safe stain	10 gels	A42135
E-Gel Double Comb Agarose Gels with SYBR Safe stain	10 gels	A42348
E-Gel EX Double Comb Agarose	10 gels	A42346
E-Gel CloneWell II Agarose Gels with SYBR Safe DNA Gel Stain, 0.8%	18 gels	G661818
E-Gel Agarose Gels with SYBR Safe DNA Gel Stain, 2%	18 gels	G521802
E-Gel EX Agarose Gels, 2% with SYBR Gold DNA stain	20 gels	G402002
E-Gel 1 Kb Plus DNA Ladder	100 applications	10488090
E-Gel Sample Loading Buffer, 1X	4 x 1.25 mL	10482055
E-Gel Power Snap Plus Electrophoresis System		G9301*
		G9311**



	Quantity	Cat. No.
Nucleic acid separation and analysis (cont.)		
E-Gel Power Snap Plus Electrophoresis System Starter kit, 48-well, 1%	1 kit	G9341* G9331**
E-Gel Power Snap Plus Electrophoresis System Starter kit, 48-well, 2%	1 kit	G9342* G9332**
E-Gel Power Snap Plus Electrophoresis System Starter kit, 96-well, 1%	1 kit	G9391* G9381**
E-Gel Power Snap Plus Electrophoresis System Starter kit, 96-well, 2%	1 kit	G9392* G9382**
E-Gel Power Snap Electrophoresis System Starter Kit, EX 2%	1 kit	G8342ST
E-Gel 48 Agarose Gels with SYBR Safe DNA Gel Stain, 2%	8 gels	G820802
E-Gel 96 Agarose Gels with SYBR Safe DNA Gel Stain, 2%	8 gels	G720802
E-Gel 48 Agarose Gels with SYBR Safe DNA Gel Stain, 4%	8 gels	G820804
Cloning and gene synthesis		
FastDigest BamHI	800 reactions	FD0054
	2,500 reactions	FD0055
FastDigest BclI	20 reactions	FD1253
	50 reactions	FD1254
FastDigest BshTI	20 reactions	FD1464
FastDigest DpnI	50 reactions	FD1703
	100 reactions	FD1704
FastDigest EcoRI	800 reactions	FD0274
	2,500 reactions	FD0275
FastDigest KpnI	300 reactions	FD0524
	20 reactions	FD0593
FastDigest NotI	50 reactions	FD0594
	150 reactions	FD0595
	250 reactions	FD0596
FastDigest Sall	200 reactions	FD0644
FastDigest XbaI	300 reactions	FD0684
	750 reactions	FD0685
FastDigest XhoI	400 reactions	FD0694
	1,200 reactions	FD0695
FastDigest Esp3I (BsmBI) (IIS class)	20 reactions	FD0454
FastDigest BpiI (BbsI) (IIS class)	20 reactions	FD1014
FastDigest Eco31I (BsaI) (IIS class)	50 reactions	FD0293
	100 reactions	FD0294

* North America, Europe, Middle East, and Africa.

** Asia Pacific, Japan, Latin America, and greater China.

	Quantity	Cat. No.
Cloning and gene synthesis (cont.)		
TOPO TA Cloning Kit for Subcloning, without competent cells	25 reactions	450641
Zero Blunt TOPO PCR Cloning Kit, without competent cells	25 reactions	450245
pENTR/D-TOPO Cloning Kit, with One Shot TOP10 Chemically Competent <i>E. coli</i>	20 reactions	K240020
One Shot TOP10 Chemically Competent <i>E. coli</i>	20 reactions	C404003
One Shot Stbl3 Chemically Competent <i>E. coli</i>	20 x 50 µL	C737303
MAX Efficiency DH5α Competent Cells	200 µL	18258012
ElectroMAX DH10B Cells	100 µL	18290015
MAX Efficiency Stbl2 Competent Cells	5 x 200 µL	10268019
MultiShot TOP10 Chemically Competent <i>E. coli</i>	5 plates	C40005
MultiShot StripWell TOP10 Chemically Competent <i>E. coli</i>	1 rack	C409601
MultiShot StripWell BL21 Star (DE3) Chemically Competent <i>E. coli</i>	1 rack	C609601
MultiShot FlexPlate TOP10 Chemically Competent <i>E. coli</i>	1 plate	C4081201
MultiShot FlexPlate DH5α T1 ^R Chemically Competent <i>E. coli</i>	1 plate	C4481201
MultiShot FlexPlate Stbl3 Chemically Competent <i>E. coli</i>	1 plate	C7381201
GeneArt Gibson Assembly HiFi Master Mix	50 reactions	A46628
GeneArt Gibson Assembly EX Master Mix	50 reactions	A46636
GeneArt Seamless Cloning and Assembly Enzyme Mix	20 reactions	A14606
GeneArt Type IIS Assembly Kit, AarI	10 reactions	A15916
GeneArt Type IIS Assembly Kit, BsaI	10 reactions	A15917
GeneArt Type IIS Assembly Kit, BbsI	10 reactions	A15918
GeneArt High-Order Genetic Assembly System	10 reactions	A13285
Gateway BP Clonase II Enzyme Mix	20 reactions	11789020
Gateway LR Clonase II Enzyme Mix	20 reactions	11791020
MultiSite Gateway Pro Plus	20 reactions	12537100
LR Clonase II Plus Enzyme	20 reactions	12538120
Gateway Vector Conversion System with One Shot <i>ccdB</i> Survival Cells	1 kit	11828029
PCR Cloning System with Gateway Technology with pDONR 221 and OmniMAX 2 Competent Cells	20 reactions	12535029
PCR Cloning System with Gateway Technology with pDONR/Zeo and OmniMAX 2 Competent Cells	20 reactions	12535037
Gateway pDONR 221 Vector	6 µg	12536017
pENTR/D-TOPO Cloning Kit, with One Shot TOP10 Chemically Competent <i>E. coli</i>	20 reactions	K240020
pCR 8/GW/TOPO TA Cloning Kit with One Shot TOP10 <i>E. coli</i>	20 reactions	K250020

GeneArt Gene Synthesis thermofisher.com/genesynthesis

GeneArt Strings DNA Fragments thermofisher.com/strings




School of Molecular Biology

Whether you want to review the basics, gain more in-depth knowledge, or discover the latest research tools, tap into our free education resources that are designed for new and experienced molecular biologists alike.

Learn more at thermofisher.com/molbioschool.

School of Molecular Biology




Speaking of Mol Bio podcast series

In this podcast, we have bold conversations that help inspire confidence in the future of molecular biology to tackle complex real-world problems.

Learn more at thermofisher.com/speakingofmolbio.

Speaking of Mol Bio podcast



Find out more at thermofisher.com/pcrandcloning