



Thermo Scientific Qtegra
Intelligent Scientific Data
Solution Software

Delivering quality

Driving productivity

Thermo
SCIENTIFIC


**There can be economy only
where there is efficiency**





The Thermo Scientific™ Qtegra™ Intelligent Scientific Data Solution™ (ISDS) Software is a software platform built on four pillars: Workflow, Data Management, Scalability and Compliance. Each pillar is carefully optimized to deliver simplicity, productivity, efficiency and quality in your laboratory.

- **Workflow**
- **Data Management**
- **Scalability**
- **Compliance**



Minimize steps needed to perform any task

A simple workflow and uncluttered user interface places the relevant controls at your disposal when you need them.

Less training and expertise are required to deliver your results.

Make all steps easy to understand and use

One click instrument "Get Ready" reduces operator time, freeing this valuable resource for other tasks.

Minimize time taken to perform any task

Just five clicks to create an analytical LabBook. The LabBook starts your intelligent sequence, with full QA/QC protocol and processes and reports your results. Unchain your time at the instrument to maximize your output.


Qtegra ISDS Software drives your Efficiency, Productivity and Quality

How do I get from sample to results?


A daily challenge for the laboratory or analyst. This challenge has many obstacles, including a complex collection of protocols, documents, and references. Due to this complexity, preparing to run an analysis on conventional software can be slow, tedious, and error-prone.

Qtegra ISDS Software has an intelligent workflow that drives you from sample to result. This intelligent and dynamic approach minimizes operator tasks between sample receipt and generation of the final result, dramatically boosting laboratory productivity.





Qtegra Workflow to Results



"Get Ready" – Automated System Set-Up
(Part of the Qtegra ISDS Software Dashboard)

Dashboard

The Qtegra Dashboard displays your system status. At a glance the status of your system and peripherals are displayed in a simple uncluttered layout.

Dashboard includes an intelligent automated instrument initialization algorithm - "Get Ready".

The Get Ready function sets your system automatically, saving you time and ensuring consistent system performance. The Dashboard takes you from instrument overview to your LabBook.

Template

A template contains the complete predefined measurement protocol: analyte details (elements, isotopes, wavelengths), internal standards, interference corrections, calibration parameters, QA/QC tests and sequence rules.

It contains comprehensive quality control tests that can be fully automated into your workflow, each with pass/fail actions.

The intelligent template enables continuous data quality monitoring and rule based automatic interaction on failures, such as re-calibration, repeat analysis and automatic dilution of samples outside the calibration range.

The user experience in Qtegra ISDS Software is based around simplicity. Extensive automation of tasks reduces the requirements for manual intervention. The intelligent workflow saves you time at the instrument, saves you time training operators and saves you time working with your results.

Workflow

Data Management

Scalability

Compliance

LabBook

The LabBook is home to your data and is prepared in just 5 clicks.

The user only needs to know two pieces of information to generate a LabBook:

- Which template?
- How many samples?

The analysis begins as soon as you create and schedule a LabBook; data is processed on the fly so your results are ready immediately. Your LabBook contains sample lists and fully interactive display of results, allowing you to dive into relevant areas of your data.


Is there a direct route to relevant data?

Once an analysis has been made, the Qtegra ISDS Software LabBook provides a layered approach to results through a series of views, where data can be displayed, queried and instantly recalculated.

This layered approach quickly directs the user to relevant areas of data, maximizing your manual data review, to maximize your productivity.



- 1 The initial display focuses on the most important values – the final concentrations.
- 2 More detailed information is available by expanding desired analyses.
- 3 Detailed calibration information.
- 4 Sensitivity, regression information as well as detection limits.
- 5 Commenting on individual analyses.



Data, data and more data. It floods a productive laboratory from a range of instruments and systems. Organizations that excel are those that master the disciplines of data handling, gathering and diving to accurate and relevant data from a range of sources.

- Workflow
- Data Management
- Scalability
- Compliance

Enhanced LabBook Query

Fast searches across multiple LabBooks can be performed directly within Qtegra LabBook Query.

Can't remember when you analyzed a sample? Then search for results based on any metadata or data contexts. The data can be queried over multiple LabBooks simultaneously.

Want to combine data from multiple analytical runs? Data can be combined from several LabBooks; for example, to assess long term reproducibility of quality control results.

WYSIWYG Data Export

Once reviewed, you have several options where to send your data. A "What you see, is what you get" (WYSIWYG) data export tool enables simple reporting of relevant results, or transfer of data to third-party software (including Iolite & Glitter). All information – whether it is graphical or data – can be instantly exported as an image or as text.

Full Integration with LIMS

Sample lists can be imported from LIMS into pre-defined protocol templates, incorporating the full range of quality control checks. Comprehensive, user-definable reports allow for flexible export to external LIMS software packages.

How to expand without adding complexity?

Whether you are challenged with globalization or consolidation, Qtegra ISDS Software allows organizations to deploy in a familiar environment. ISDS Software Qtegra software provides simple integration into your various workflows, maximizes capital investments and minimizes administrative costs.

Qtegra software scalability provides a balance between mainstream capabilities and specialized functions: focusing on the immediate task and creating awareness of greater possibilities.



We have a full range of instrumentation for the analysis of elements and isotopes. Qtegra ISDS Software is the unifying platform that links our innovative and market leading hardware to your application. Current, future and legacy instrumentation will in time be supported by this enabling software.

- Workflow
- Data Management
- **Scalability**
- Compliance

Standardization

Qtegra ISDS Software is standardized from instrument to instrument. This minimizes the requirement for specialist operators and dramatically reduces the requirement for extensive cross training.

Qtegra ISDS Software promotes collaboration and provides capabilities that allow a corporation to standardize:

- Harmonized workflow
- Streamlined method development and validation
- Multi-vendor accessory control for LC, GC, lasers and other sampling devices
- Reduced training efforts and cost
- Standardized data format
- Secure user access and tracking

Rapid Deployment

Qtegra ISDS Software is simple to install and configure. An intuitive drag-and-drop instrument configurator provides a bespoke system configuration for your particular application. Each Thermo Scientific system comes pre-deployed with Qtegra ISDS Software. Your installation engineer is able to deploy or restore factory configurations from a Qtegra file server to ensure rapid system set-up and security of system settings.

Localization

Qtegra ISDS Software is designed to be deployable across global laboratories. Qtegra ISDS Software is localized, already supporting English, Japanese, Chinese, Russian and German languages. Further language support is in progress.

Simple Configuration

Qtegra ISDS Software is modular, so it provides a flexible framework in which instrument and accessory plug-ins are combined into a single workflow. In addition to the instrument plug-in, Qtegra ISDS Software provides integrated plug-ins for autosamplers, autodilutors as well as major chromatography and laser ablation systems.

Plug-ins are available for all major Thermo Scientific liquid, ion and gas chromatography systems as well as support for third-party instruments. We also provide a software plug-in development kit to allow third-party manufacturers to interface their technology seamlessly to the powerful Qtegra ISDS Software.

Will new technology create a compliance challenge?

Companies in pharmaceutical, biotech or medical device industries continually push the boundaries of innovation to develop new products. These industries are constantly being challenged to meet rising standards of quality and to comply with rigorous regulatory requirements.

Thermo Scientific Qtegra ISDS Software recognizes these challenges and includes of the necessary tools for you to rapidly meet your compliance goals.



Thermo
SCIENTIFIC

Compliance violations can cost a company severely in fines and lost revenue. However, achieving full regulatory compliance can be time consuming and complex with legacy software or third-party packages. Qtegra ISDS Software simplifies your route to compliance.

- Workflow
- Data Management
- Scalability
- **Compliance**

21 CFR Part 11 Environment Support

Qtegra ISDS Software provides a full range of features to comply with your regulatory requirements, including:

- Audit trails
- Open / closed system control
- Electronic signatures
- Identification codes and passwords
- Tools for integrated data management
- Full record change traceability

Qtegra ISDS Software security supports the administrator by creating a controlled environment to ensure that the system is configured to comply with GxP and 21 CFR Part 11 regulations.

Data Security and Access Control

Qtegra ISDS Software ensures the utmost integrity and traceability for all of your raw data, metadata, and reported results.

Access control in Qtegra ISDS Software is provided by the Microsoft Windows® User Management system that is flexible, secure and customizable. Users are defined with privileges, roles, and access groups.

Qtegra ISDS Software enables you to establish audit trails, assign access privileges, set password expiry, and automated log-outs due to inactivity. Multi-level security is provided for every task.

Compliance Management

Qtegra ISDS Software also helps you manage your compliance on a day-to-day basis with such capabilities as secure time-stamped audit trails. Records can be easily reviewed so you can apply the appropriate electronic signatures. All data is relationally stored, so full documentation is automatic and complete, saving you time during compliance audits.

Qtegra driven – the journey has begun

Qtegra ISDS Software is developed using a platform-based approach; it is a core control software across different analytical devices. The introduction of Qtegra software for Thermo Scientific iCAP Q ICP-MS was the first step on the implementation roadmap for Qtegra ISDS Software control of a variety of technologies for the analysis of elements and isotopes.


Over the coming months the list of techniques and supported peripherals will rapidly increase to allow you to fully exploit the benefits of Qtegra ISDS Software. Cross-device training and adoption of new instrumentation is faster and easier, allowing for increased flexibility in modern multi-technique laboratories.

The range of instruments running Qtegra ISDS Software is rapidly expanding. This is supported on:

- Thermo Scientific™ iCAP™ 7000 ICP-OES
- Thermo Scientific™ iCAP™ 6000 ICP-OES
- Thermo Scientific™ iCAP™ Q ICP-MS
- Thermo Scientific™ Argus VI™ Noble Gas MS
- Thermo Scientific™ Helix SFT™ Noble Gas MS
- Thermo Scientific™ Helix MC™ Noble Gas MS
- Thermo Scientific™ Delta Ray™ Isotope Ratio Infrared Spectrometer (IRIS)
- Thermo Scientific™ 253 Ultra™ High Resolution Isotope Ratio Mass Spectrometer

Support coming soon:

- Thermo Scientific™ iCAP™ 7000 Plus ICP-OES
- Thermo Scientific™ ELEMENT 2 HR-ICP-MS
- Thermo Scientific™ ELEMENT XR HR-ICP-MS



thermofisher.com/Qtegra-ISDS

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
 SCIENTIFIC

A Thermo Fisher Scientific Brand