

**Thermo Scientific Nunc Dense Storage**

- 1.0mL Nunc™ CryoBank™ Tubes
- Freezer Boxes
- 13 x 13 Dividers

**maximize your space**

dense low temperature storage

**Thermo**  
SCIENTIFIC


# double your storage

within the s


## Think inside the box.

The Dense Low Temperature Storage format uses standard CryoBox sizes that work with the freezer racks you currently own. By using standard boxes with 13 x 13 tube separators, 169 1.0mL tubes fit in a standard size CryoBox – effectively doubling your capacity, and cutting storage costs by 50%.\*

For a free sample kit, go to:  
[www.thermoscientific.com/storemore](http://www.thermoscientific.com/storemore)



Cryogenic storage box systems contain the standard 1-2.0mL cryogenic tube in either a 9 x 9 (81 tubes) format, or a 10 x 10 (100 tubes) format



Since many samples are only 1.0mL or less, this storage system uses 1.0mL tubes in a 13 x 13 (169 tubes) array, to double existing freezer space capacity

\*Compared with a standard CryoBox that stores 81 samples

# age

space you already have.

## CryoBank Tubes (1.0mL)

- > The space-efficient shape of our CryoBank 1.0mL tubes enables dense storage to work. Although the most commonly used tube in sample storage is 2.0mL, very often only 1.0mL of sample is actually stored per tube.

## Plastic CryoBox with Divider

- > These white polycarbonate freezer boxes (132 x 132 x 52mm) come with a 13 x 13 divider and a clear lid.

## Cardboard CryoBox with Divider

- > These white cardboard storage boxes (132 x 132 x 75mm) have a 13 x 13 divider to hold 169 tubes.

## Plastic CryoBox without Divider


- > These white polycarbonate freezer boxes (132 x 132 x 52mm) have a clear lid, but are provided without a divider. Can be used to -150°C.

## Replacement 13 x 13 Dividers

- > Use the boxes you already own, but add this 13 x 13 divider to fit 169 tubes.

## Picker Pen (included with all CryoBank Tubes)

- > A simple, inexpensive pen with a specially shaped tip that fits like a lock-and-key into the top of any CryoBank tube for easy removal from storage boxes. Just click the top of the pen, and the tube is released.



Catalog #	Description	Units per Pack/Case
374115	1.0mL CryoBank Tubes, Blank, with Caps, Sterile, Loose packed	169/845
5050-0169	CryoBox with 13 x 13 divider, 133 x 133 x 52mm, white PC with clear lid	4/24
369640	CryoBox with 13 x 13 divider, 133 x 133 x 52mm, white cardboard with lid	4/24
5050-0001	CryoBox with no divider, 133 x 133 x 52mm, white PC with clear lid	4/24
374001	13 x 13 divider to fit 133 x 133 x 52mm CryoBox	4/24
369678	Picker Pen for CryoBank Tubes	1/24



# dense storage workflow

products to maximize capacity


Our low temperature storage portfolio includes all the tools you'll ever need for the most efficient sample storage.

**Double your storage** – Use CryoBank Tubes with CryoBoxes with 13 x 13 dividers to double your storage capacity.

**No more finger fumbling** – Just spend a minute or two with the Picker Pen, and you'll see how convenient it is to maneuver around your sample collection. Uncap or move any CryoBank tube from one rack to another – all with the help of the pen you write with. How cool is that?

**Expand your capacity even more** – Thermo Scientific Revco PLUS -86°C Freezers feature ultra space efficiency, fast temperature recovery after door openings, a 15% reduction in power, and more.

Visit [www.thermoscientific.com/coldstorage](http://www.thermoscientific.com/coldstorage)



© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of ThermoFisher Scientific Inc. and its subsidiaries.

**Asia:** China Toll-free: 800-810-5118 or 400-650-5118; India: +91 22 6716 2200,  
India Toll-free: 1 800 22 8374; Japan: +81-3-5826-1616;  
Other Asian countries: 65 68729717  
**Europe:** Austria: +43 1 801 40 0; Belgium: +32 2 482 30 30; Denmark: +45 4631 2000;  
France: +33 2 2803 2180; Germany: +49 6184 90 6000, Germany Toll-free: 0800 1 536 376;  
Italy: +39 02 95059554 or 434-254-375; Netherlands: +31 76 571 4440;  
Nordic/Baltic countries: +358 9 329 10200; Russia/CIS: +7 (812) 703 42 15;  
Spain/Portugal: +34 93 223 09 18; Switzerland: +41 44 454 12 22;  
UK/Ireland: +44 870 609 9203  
**North America:** USA/Canada +1 585 586 8800; USA Toll-free: 800 625 4327  
**South America:** USA sales support: +1 585 899 7198  
**Countries not listed:** +49 6184 90 6000

[www.thermoscientific.com](http://www.thermoscientific.com)

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific