



Clinical specimen collection

Thermo Scientific Samco Clicktainer Vials

Providing ultimate sample protection and user safety with a new, audible multi-click seal



Samco Clicktainer Vials

Multi-click system to help ensure secure and reliable sample protection

The scientists at Thermo Fisher Scientific understand the value of patient samples and the importance of safe transport. The Thermo Scientific™ Samco™ Clicktainer™ Vials offer improved sample protection by providing a leak tight seal, easy to grip surfaces, and straight sided design for easy handling.



Clicktainer cap

- Visual, audible, and tactile indication of full cap closure
- Safe tab design to help reduce risk of glove tears and accidental breakage

Leak proof*

- Product is leak tested in accordance with BS EN 14254

Transport

- Validated for use in hospital pneumatic transport systems

Sterile and non-sterile

- Available in both gamma sterilized and non-sterile configurations

Volume

- Choice of volume with wide mouth design
- 120mL or 90mL volume suitable for most applications

Base design

- Excellent clarity with easy-to-see graduations for visibility of samples within vial
- Knurls on cap and base make vial easy to grip with gloved hands

Stackable

- Optimize sample storage space



*The term "leakproof" applies to Thermo Scientific Samco Clicktainer products, 90ml/53mm and 120ml/53mm, that meet the following criteria: a) a vacuum test that draws 95 kPa of pressure within a sealed chamber with no leakage; b) leak testing in accordance with BS EN 14254, which involves filling the container with a dye, placing for 2 minutes on a roller mixer followed by 1 hour inverted in a water bath with no leakage. Note: these tests, using other liquids, may not yield the same results. To ensure safe usage, customers are advised to test Thermo Scientific Samco Clicktainer Vials under conditions of their planned applications. Thermo Scientific Samco Clicktainer Vials are leakproof at ambient temperature and pressure when used with their Clicktainer closures.

Performance, design, color, and format

Superior leak tight performance

- Proven leak tight in accordance with BS EN 14254
- Validated as 95kPa compliant
- Rigorously tested for suitability in pneumatic tubes

Clearly designed labels

- One label multiple languages – English, Spanish and French
- Fields to capture all relative clinical patient information
- Optional temperature strip
- Also available unlabeled

Wide range of cap colors and packaging configurations

- Cap colors to differentiate and manage samples
- Available pre-assembled, sterile and non-sterile, for convenient usage
- Separate closure and vial bulk packs for fixative pre-filling applications
- Individual peel pouches for operating rooms, cleanrooms, and other sterile environments
- Dispenser boxes ideal for neat storage in smaller labs and doctors' offices



Samco Clicktainer Vial label with fields to capture all FDA required information.

Samco Clicktainer Vials

Unique multi-click seal system to ensure sample protection and user safety

Samco Clicktainer Vials ordering information

Description	Closure color	Label type	Product configuration	Quantity (pk x pks per case)	Catalog number
120mL/53mm sterile					
120mL/53mm CLK SPC LBL OR CAP ST	Orange	Tabbed	Pre-assembled	75x4	120ORG53-1000
120mL/53mm CLK O/R N/LBL OR CAP ST	Orange	None	Peel pouch	1x100	120ORG53-4000
120mL/53mm CLK SPC LBL WT CAP ST	White	Tabbed	Pre-assembled	75x4	120WHT53-1000
120mL/53mm CLK MSU, NO FUNNEL, ST	White	Tabbed	Mid-stream kit	1x100	120WHT53-3000
120mL/53mm CLK O/R POUCH WT CAP ST	White	None	Peel pouch	1x100	120WHT53-4000
120mL/53mm non-sterile					
120mL/53mm CLK SPC LBL OR CAP NS	Orange	Tabbed	Pre-assembled	75x4	120ORG53-2000
120mL/53mm CLK SPC N/LBL WT CAP NS	White	None	Bagged separate	300x1	120WHT53-2000
90mL/53mm sterile					
90mL/53mm CLK SPC LBL OR CAP ST	Orange	Tabbed	Pre-assembled	100x4	90ORG53-1000
90mL/53mm CLK SPC LBL OR CAP ST DISP	Orange	Tabbed	Dispenser box	20x16	90ORG53-1001
90mL/53mm CLK SPC LBL WT CAP ST	White	Tabbed	Pre-assembled	100x3	90WHT53-1000
90mL/53mm CLK SPC TSLBL WT CAP ST	White	Label w/temp strip	Pre-assembled	100x3	90WHT53-1001
90mL/53mm CLK SPC TS WT CAP ST	White	Temp strip only	Pre-assembled	100x4	90WHT53-1002
90mL/53mm CLK SPC LBL GN CAP ST	Green	Tabbed	Pre-assembled	100x4	90GRN53-1000
90mL/53mm non-sterile					
90mL/53mm CLK SPC LBL WT CAP NS	White	Tabbed	Pre-assembled	75x4	90WHT53-2000
90mL/53mm CLK SPC TSLBL WT CAP NS	White	Label w/temp strip	Pre-assembled	300x1	90WHT53-2001
90mL/53mm CLK SPC N/LBL WT CAP NS	White	None	Bagged separate	300x1	90WHT53-2002
90mL/53mm CLK SPC N/LBL WT CAP NS	White	None	Pre-assembled	300x1	90WHT53-2003
90mL/53mm CLK SPC TS WT CAP NS	White	Temp strip only	Pre-assembled	100x4	90WHT53-2004
Non-sterile closures and containers					
53mm CLK CAP ORANGE NS	Orange	None	Bulk closures only	1500x1	CAPORG53-1000
53mm CAP WHITE NS	White	None	Bulk closures only	1500x1	CAPWHT53-1000
53mm CLK CAP GREEN NS	Green	None	Bulk closures only	1500x1	CAPGRN53-1000
120mL/53mm CLK SPC N/LBL NO CAP NS	120mL	None	Bulk container only	300x1	12000-5300
90mL/53mm CLK SPC N/LBL NO CAP NS	90mL	None	Bulk container only	300x1	9000-5300



Common applications

- Collection and transport of liquids and solids.
- Ideal for valuable and often hazardous patient samples, such as urine, stool, and sputum.
- Suitable for use with formalin or fixatives for histology applications.
- Sample management in research, food, and beverage.

Thermo Scientific™ Clinical Labware

The Thermo Scientific portfolio includes a full offering of consumables for your laboratory applications. Please contact us for more details.

Find out more at thermofisher.com/clicktainer

thermo scientific