

RapID™ ANA II QUALITY CONTROL REPORT FORM

1. *Record results of biochemical reactions under the expected reactions.
 Corrective Action: See Reverse Page

Date QC performed: _____
 Technologist: _____
 ANA II lot #: _____
 Expiration Date: _____
 Approved by: _____

Organism	URE	BLTS	aARA	ONPG	aGLU	BGLU	aGAL	aFUC	NAG	PO4	LGY	GLY	PRO	PAL	ARG	SER	PYR	IND
<i>C. sordelli</i> ^a	+	-	-	-	-	-	-	V	-	-	-	V	+	V	V	V	V	+
ATCC® 9714																		
	URE	BLTS	aARA	ONPG	aGLU	BGLU	aGAL	AFUC	NAG	PO4	LGY	GLY	PRO	PAL	ARG	SER	PYR	IND
<i>B. distasonis</i> ^a	-	+	V	+	+	+	+	-	+	+	+	+	-	+	+	+	+	-
ATCC® 8503																		
	URE	BLTS	aARA	ONPG	aGLU	BGLU	aGAL	AFUC	NAG	PO4	LGY	GLY	PRO	PAL	ARG	SER	PYR	IND
<i>B. uniformis</i>	-	+	+	+	+	+	+	+	+	+	+	-	-	-	-	-	-	+
ATCC® 8492																		

+ = Positive - = Negative V = Variable

*Visual Interpretation Only

^aKey indicator strain

CORRECTIVE ACTION

1. List organism with biochemical reaction that is outside the acceptable criteria:

2. Date corrective action started: _____

3. Corrective action procedure performed: _____

4. Results after corrective action:

Acceptable Reaction

Outside Acceptable Reaction

5. Further corrective action results, if the biochemical reaction is still outside the acceptable criteria:

6. Date corrective results accepted: _____