

Manual microbial identification

A simple, smart solution for a time-efficient, productive lab

Reliable reporting. Large database. Rapid results.



Manual microbial identification is now quicker than ever. With Thermo Scientific™ RapID™ Systems and ERIC™ Software. Experience fast identification, reduced turnaround time compared to traditional biochemical identification panels, visible color reactions and precise reporting for more organisms with a range of nine panels.

Designed exclusively for use with RapID Systems, the latest version of ERIC Software features an updated taxonomy for excellent coverage and confidence. All labs can take advantage of the easy-to-use software with updated nomenclature - available free online or in a CD format.

A comprehensive compendium

ERIC databases provide excellent reporting coverage for thorough analysis. Databases are updated with current organism nomenclature to keep your laboratory informed and productive.

Minimal hands-on time with one step inoculation and no oil overlays with RapID and ERIC Software. To learn more about our complete range of premium products, [visit our website](#).



Manual microbial identification made easy

No oil. No pipetting.
No 24-hr incubation.

One-step inoculation

Fast preparation time improves lab productivity

Four-hour incubation*

Quick same day time-to-result for a faster response

Common procedure

Reduced material usage saves time and money

Visible color reactions

Easy to read decreasing repeat testing

*Two-hour incubation for RapID SS/u.

RapidID systems and product guide

Color guides are included in every kit

Product	Format	Product Code
RapID ANA II	20 panels/kit	R8311002
RapID Inoculation Fluid 1mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 3	Each	R20413
RapID Spot Indole 15mL [†]	Each	R8309002
RapID ANA II QC Set/3	5 loops/pack	R4653050
ERIC RapID Compendium CD	Each	R8323600
RapID CB Plus	20 panels/kit	R8311008
RapID Inoculation Fluid 2mL [†]	20 tubes/pack	R8325106
McFarland Standard No. 4	Each	R20414
RapID Nitrate A 15mL [†]	Each	R8309003
RapID Nitrate B 15mL [†]	Each	R8309004
RapID CB Plus QC Set/3	5 loops/pack	R4653048
ERIC RapID Compendium CD	Each	R8323600
RapID NH Plus	20 panels/kit	R8311001
RapID Inoculation Fluid 1mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 3	Each	R20413
RapID Spot Indole 15mL [†]	Each	R8309002
RapID Nitrate A 15mL [†]	Each	R8309003
RapID Nitrate B 15mL [†]	Each	R8309004
BactiDrop™ Oxidase 0.75mL	50 ampules/pack	R21540
RapID NH Plus QC Set/4	5 loops/pack	R4653051
ERIC RapID Compendium CD	Each	R8323600
RapID NF Plus	20 panels/kit	R8311005
RapID Inoculation Fluid 1mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 1	Each	R20411
McFarland Standard No. 3	Each	R20413
RapID Spot Indole 15mL [†]	Each	R8309002
RapID Nitrate A 15mL [†]	Each	R8309003
BactiDrop Oxidase 0.75mL	50 ampules/pack	R21540
RapID NF Plus QC Set/4	5 loops/pack	R4653054
ERIC RapID Compendium CD	Each	R8323600

[†] Accessory product required for use with RapID Systems.

Product	Format	Product Code
RapID ONE	20 panels/kit	R8311006
RapID Inoculation Fluid 2mL [†]	20 tubes/pack	R8325106
McFarland Standard No. 2	Each	R20412
RapID Spot Indole 15mL [†]	Each	R8309002
BactiDrop Oxidase 0.75mL	50 ampules/pack	R21540
RapID ONE QC Set/4	5 loops/pack	R4653056
ERIC RapID Compendium CD	Each	R8323600
RapID Staph Plus	20 panels/kit	R8311009
RapID Inoculation Fluid 2mL [†]	20 tubes/pack	R8325106
McFarland Standard No. 3	Each	R20413
RapID Nitrate A 15mL [†]	Each	R8309003
RapID Nitrate B 15mL [†]	Each	R8309004
RapID Staph Plus QC Set/4	5 loops/pack	R4653047
ERIC RapID Compendium CD	Each	R8323600
RapID SS/u	20 panels/kit	R8311004
RapID Inoculation Fluid 1mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 1	Each	R20411
RapID Spot Indole 15mL [†]	Each	R8309002
RapID SS/u QC Set/5	5 loops/pack	R4653053
ERIC RapID Compendium CD	Each	R8323600
RapID STR	20 panels/kit	R8311003
RapID Inoculation Fluid 1mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 1	Each	R20411
RapID STR QC Set/5	5 loops/pack	R4653052
ERIC RapID Compendium CD	Each	R8323600
RapID Yeast Plus	20 panels/kit	R8311007
RapID Inoculation Fluid 2mL [†]	20 tubes/pack	R8325106
RapID Yeast Plus QC Set/5	5 loops/pack	R4653060
ERIC RapID Compendium CD	Each	R8323600

To learn more about our complete range of RapID products and the **free** online ERIC software visit: thermofisher.com/RapID

thermo scientific