



Related products:

Electronic Multichannel Pipette, 25-1250 μ l
Universal Charger
Pipette Tips; 960 tips per box
Pipette Troughs

Part #:

V4009
V4009PS
E4109
E1032

Fast, accurate inoculation

Thermo Scientific Sensititre Electronic Multichannel Pipette:
a quick, accurate and versatile manual plate inoculation system

Enhanced ergonomics

- Keeps your hand in the neutral position, increasing comfort and accuracy
- Upright design prevents tips from ever touching the bench and risking contamination

Versatile dispensing options

- Delivers 25 to 1250 μ l in 5 μ l increment settings, using one to eight tips
- Multiple programmable settings and easy single delivery set up with touch pad ease
- Dispensing speed option allows dosing of a plate in less than 10 seconds

Simple operation

- Effortless tip ejection
- Only requires 90-minute charging every 600 cycles
- On-board calibration software eliminates any need for tools

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

USA and Canada +1 800 871 8909
All Other Inquiries +44 1256 841144

Technical Support
USA and Canada +1 800 642 7029
All Other Inquiries +44 1342 318777

Thermo
S C I E N T I F I C

Part of Thermo Fisher Scientific

LT2033A/CS/02/12