

Cryo-TEM Supply Checklist

Equip your lab with the items needed for your cryo-TEM workflow

- | | |
|--|---|
| <input type="checkbox"/> Cryo Gloves | <input type="checkbox"/> 4L LN ₂ Dewar (2) |
| <input type="checkbox"/> 8" or 10" tweezers | <input type="checkbox"/> 1L LN ₂ Dewar |
| <input type="checkbox"/> Sharp tweezers | <input type="checkbox"/> C-clips |
| <input type="checkbox"/> Carbon fiber tipped tweezers | <input type="checkbox"/> Clip Rings |
| <input type="checkbox"/> Heating block for 1.5 or 2 ml tubes | <input type="checkbox"/> Hair dryer |
| <input type="checkbox"/> 50 ml conical tubes | <input type="checkbox"/> Drill |
| <input type="checkbox"/> Colored lab tape | <input type="checkbox"/> 0.8 mm nylon cord/thread |
| <input type="checkbox"/> LN ₂ Dewar with large neck | <input type="checkbox"/> Small screwdriver set |
| <input type="checkbox"/> LED light with magnifying lens | <input type="checkbox"/> Cryo-grid boxes with lids (get some with diamond holes and some with square holes) |

Learn more at thermofisher.com/emservices

thermo scientific