



# Amira-Avizo 2D Software version 2023.2

Release notes

# Contents

Introduction	3
<hr/>	
Operating Systems	4
<hr/>	
Resolved Issues	4
<hr/>	

# Introduction

This document covers the most important new features, improvements, and changes in this version of Thermo Scientific™ Amira™ 2D Software and Avizo™ 2D Software. In addition, you will find a list of new Xtras, including video tutorials, recipes, and workflows, which have been published in the [Xtras Library](#) since the previous release.

Visit our [Customer Support and Service Center](#) to access services and resources that help you make the best use of your product.

We value your feedback. If you have any suggestions for improvement, do not hesitate to [contact us](#).

## Operating Systems

Version 2023.2 of the Amira-Avizo 2D Analyzer and Trainer Packages run on:

- Microsoft Windows 10 (64-bit).
- Microsoft Windows 11 (64-bit)

Version 2023.2 of the Amira-Avizo 2D Labeler Package runs on:

- Microsoft Windows 10 (64-bit).

Our products are tested on the following configurations:

Platform	GPU	Driver number
Windows 10	NVIDIA Tesla T4	537.13
Windows 11	NVIDIA RTX A4500	537.34
Ubuntu-20.04	NVIDIA T1000	535.113.01

## Resolved issues

Name	ID	Description
Image Recipe Designer	AA-29155	Some color issues with the labeling module have been fixed.
FFT Filter	AA-28891	It is now possible to use the FFT Filter on images with 8-bit signed data type.
Licensing	AA-28495	Some minor licensing issues are now managed as warnings and no longer prevent the application from functioning.

Learn more at [thermofisher.com/amira-avizo2d](https://thermofisher.com/amira-avizo2d)

**thermo** scientific